

INTERPRETATION TABLE

DR low resolution SSP typing

Amplification patterns of the DRB1*01:01 to DRB1*10:04 alleles

		Well																						
		65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87
Length of spec.		255	200	215	210	220	210	210	175	235	250	180		170	105	215	225		155	170	225	240	215	175
PCR product(s)		255	200	215	210	220	210	210	175	235	250	180		170	105	215	225		155	170	225	240	215	175
Length of int. pos. control ¹		515	430	430	430	430	430	430	430	430	515	430	430	430	430	430	430	430	430	430	430	430	430	430
5'-primer(s) ²		5'-gAA ³ 14(129)	5'-gAA ³ 14(129)	5'-AAg ³ 13(126)	5'-AAg ³ 13(126)	5'-gTT ³ 47(227)	5'-gTC ³ 13(125)	5'-gTC ³ 13(125)	5'-gTC ³ 13(125)	5'-gTC ³ 13(125)	5'-ATA ³ 14(127)	5'-gTT ³ 16(133)	5'-gAg ³ 58(261)	5'-gCg ³ 31(178)	5'-gTC ³ 16(133)	5'-TAT ³ 12(122)	5'-TAT ³ 12(122)	5'-gTC ³ 13(125)	5'-gTT ³ 37(197)	5'-CAG ³ 34(189)	5'-CAG ³ 34(189)	5'-gTC ³ 10(116)	5'-gAT ³ 28(170)	5'-gTA ³ 13(125)
3'-primer(s) ³		5'-gAT ³ 67(286)	5'-gAT ³ 67(286)	5'-CTg ³ 70(295)	5'-gAg ³ 67(286)	5'-ggC ³ 73(305)	5'-gCT ³ 71(299)	5'-CTg ³ 70(295)	5'-CgG ³ 58(260)	5'-ggC ³ 73(305)	5'-gCTC ³ 71(298)	5'-CCT ³ 58(260)	5'-gAg ³ 73(305)	5'-CgC ³ 74(307)	5'-CCT ³ 58(260)	5'-gCT ³ 71(299)	5'-gAA ³ 67(286)	5'-CAG ³ 58(260)	5'-TCA ³ 42(213)	5'-CTg ³ 70(295)	5'-gTg ³ 60(265)	5'-AAT ³ 77(317)	5'-CTT ³ 87(346)	5'-CCT ³ 58(260)
71(299)	5'-gAT ³ 67(286)	5'-gAg ³ 67(286)	5'-gAg ³ 67(286)	5'-CTg ³ 70(295)	5'-gAg ³ 67(286)	5'-ggC ³ 73(305)	5'-gCT ³ 71(299)	5'-CTg ³ 70(295)	5'-CgG ³ 58(260)	5'-ggC ³ 73(305)	5'-gCTC ³ 71(298)	5'-CCT ³ 58(260)	5'-gAg ³ 73(305)	5'-CgC ³ 74(307)	5'-CCT ³ 58(260)	5'-gCT ³ 71(299)	5'-gAA ³ 67(286)	5'-CAG ³ 58(260)	5'-TCA ³ 42(213)	5'-CTg ³ 70(295)	5'-gTg ³ 60(265)	5'-AAT ³ 77(317)	5'-CTT ³ 87(346)	5'-CCT ³ 58(260)
71(298)	5'-Tg ³ 70(295)	5'-gAg ³ 67(286)	5'-gAg ³ 67(286)	5'-CTg ³ 70(295)	5'-gAg ³ 67(286)	5'-ggC ³ 73(305)	5'-gCT ³ 71(299)	5'-CTg ³ 70(295)	5'-CgG ³ 58(260)	5'-ggC ³ 73(305)	5'-gCTC ³ 71(298)	5'-CCT ³ 58(260)	5'-gAg ³ 73(305)	5'-CgC ³ 74(307)	5'-CCT ³ 58(260)	5'-gCT ³ 71(299)	5'-gAA ³ 67(286)	5'-CAG ³ 58(260)	5'-TCA ³ 42(213)	5'-CTg ³ 70(295)	5'-gTg ³ 60(265)	5'-AAT ³ 77(317)	5'-CTT ³ 87(346)	5'-CCT ³ 58(260)
72(301)	5'-CTg ³ 70(297)	5'-gAg ³ 67(286)	5'-gAg ³ 67(286)	5'-CTg ³ 70(295)	5'-gAg ³ 67(286)	5'-ggC ³ 73(305)	5'-gCT ³ 71(299)	5'-CTg ³ 70(295)	5'-CgG ³ 58(260)	5'-ggC ³ 73(305)	5'-gCTC ³ 71(298)	5'-CCT ³ 58(260)	5'-gAg ³ 73(305)	5'-CgC ³ 74(307)	5'-CCT ³ 58(260)	5'-gCT ³ 71(299)	5'-gAA ³ 67(286)	5'-CAG ³ 58(260)	5'-TCA ³ 42(213)	5'-CTg ³ 70(295)	5'-gTg ³ 60(265)	5'-AAT ³ 77(317)	5'-CTT ³ 87(346)	5'-CCT ³ 58(260)
74(307)	5'-TCC ³ 70(296)	5'-gAg ³ 67(286)	5'-gAg ³ 67(286)	5'-CTg ³ 70(295)	5'-gAg ³ 67(286)	5'-ggC ³ 73(305)	5'-gCT ³ 71(299)	5'-CTg ³ 70(295)	5'-CgG ³ 58(260)	5'-ggC ³ 73(305)	5'-gCTC ³ 71(298)	5'-CCT ³ 58(260)	5'-gAg ³ 73(305)	5'-CgC ³ 74(307)	5'-CCT ³ 58(260)	5'-gCT ³ 71(299)	5'-gAA ³ 67(286)	5'-CAG ³ 58(260)	5'-TCA ³ 42(213)	5'-CTg ³ 70(295)	5'-gTg ³ 60(265)	5'-AAT ³ 77(317)	5'-CTT ³ 87(346)	5'-CCT ³ 58(260)

		5'-gCg ^{3'}		5'-gCg ^{3'}	5'-ggC ^{3'}					5'-CAC ^{3'}									5'-Cag ^{3'}					
		5'-CCA ^{3'} 86(344)		5'-gCT ^{3'} 71(299)																				
				5'-ggC ^{3'} 73(305)																				
Well No.	DR	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87
Well No.	DR	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87
DRB1 allele ⁴	ser ⁹																							
*01:01:01-01:02:08, 01:04-01:38, 01:40N-01:41, 01:43-01:45	DR1, Null, -	65																						
*01:03, 01:39N	DR103, Null		66																					
*01:42	-	65	66																					
*03:01:01-01:03:01:19, 03:04:01-03:06, 03:09, 03:11:01-03:16, 03:18-03:20, 03:22-03:23, 03:25-03:26, 03:28, 03:30-03:31, 03:33-03:34, 03:36-03:37, 03:43-03:45, 03:47-03:48, 03:50-03:52, 03:54-03:63, 03:66-03:68N, 03:70-03:73, 03:75, 03:77-03:80, 03:11:02	DR17(3), DR3, Null, -					69	70											81						
*03:02:01-03:03, 03:27, 03:29, 03:38, 03:53, 03:74	DR18(3), -					69		71										81						
*03:07, 03:17, 03:21, 03:24, 03:32, 03:35, 03:39-03:41, 03:49	DR3, -					69												81						
*03:08, 03:65	-					69	70				75		77											
*03:10	DR17(3)					69	70													83	84			
*03:42, 03:69	-					69																		
*03:46, 03:64	-					69	70																	
*03:76	-						70					76		79				81						
*04:01:01-04:61, 04:63-04:68, 04:70-04:72:02, 04:74-04:104, 04:106-04:107	DR4, Null, -								72															
*04:62, 04:69, 04:73, 04:105	-								72											82				
*07:01:01-01:07:01:04, 07:03-07:22	DR7, Null, -									73														
*08:01:01-08:02:04, 08:04:01-08:07, 08:11, 08:16-08:17, 08:22, 08:24, 08:26, 08:28, 08:39, 08:42-08:44	DR8, -										74						80							

Well No.	DR	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87
DRB1 allele ⁴	ser ⁹																							
<i>DRB3*01:01:02:01-01:15,</i> <i>DRB3*02:01-02:26,</i> <i>DRB3*02:28-02:29N,</i> <i>DRB3*03:01:01-03:03</i>	DR52, Null, -																						85	
<i>DRB3*02:27</i>	-															79						85		
<i>DRB4*01:01:01:03:01:01,</i> <i>DRB4*01:03:01:03-01:08</i>	DR53, -																						86	
<i>DRB4*01:03:01:02N</i>	Null																		82				86	
<i>DRB5*01:01:01-01:14,</i> <i>DRB5*02:02-02:05</i>	DR51, Null, -																							87
Well No.	ser ⁵	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87