

Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21
HLA-A allele ^{1,2,6}	ser ³																					
*01:01:01:01:01:01:22, 01:01:24-01:01:47, 01:01:49-01:01:55, 01:01:57-01:01:64, 01:01:67-01:01:78, 01:01:80-01:04N, 01:06, 01:09:01:01-01:11N, 01:15N-01:18N, 01:20, 01:22N-01:25, 01:27N, 01:29, 01:31N-01:33, 01:35-01:42, 01:44-01:50, 01:52:01N-01:78, 01:80-01:94, 01:96-01:98, 01:100-01:120, 01:122-01:125, 01:127-01:135, 01:137-01:144, 01:146, 01:148, 01:150-01:158, 01:160N-01:166, 01:168-01:175, 01:177, 01:179N-01:191, 01:193, 01:195-01:199, 01:201-01:204, 01:206-01:207, 01:209-01:213, 01:215-01:227, 01:230-01:235, 01:237-01:243, 36:04 ⁴	A1, A36, Null, -		1		4																	
*01:01:23, 01:01:48, 01:01:65-01:01:66, 01:01:79, 01:08, 01:14, 01:30, 01:79, 01:99, 01:145, 01:147Q, 01:149, 01:159, 01:167, 01:178N, 01:205, 01:208Q, 01:214, 01:228Q, 01:236, 36:01, 36:03, 36:05 ⁴	A1, A36, Null, -		1																			
*01:01:56	-		1		4						11											
*01:07, 01:95	A1, -		1		4												16					
*01:12, 01:19, 01:200, 01:244	A1, -		1		3																	
*01:13, 01:176	A1, -		1		4						10	12		14								
*01:21, 01:126	A1, -		1		3	4																
*01:26, 01:136, 01:192, 11:94, 11:112, 11:211, 11:226 ⁵	A1, -		1		4						10											
*01:28, 01:229	-		1		4						10			14								
*01:43	-		1		4					9												
*01:51	-		1		4			7		9												
*01:121	-		1		4															18		
*01:194	-		1		4								12	14								
*02:01:01:01:02:01:15, 02:01:17-02:01:19, 02:01:21-02:01:81, 02:01:83-02:16, 02:18-02:22:02, 02:24:01-02:33, 02:36-02:45, 02:47, 02:49-02:54, 02:57-02:61, 02:63-02:64:02, 02:66-02:69, 02:71-02:77, 02:79:01-02:80, 02:82N-02:97:02, 02:99, 02:101:01-02:102, 02:104-02:123, 02:125N-02:128, 02:130-02:134, 02:136-02:148, 02:150-02:168, 02:170-02:236, 02:238-02:242, 02:244-02:308, 02:310-02:337, 02:339-02:406, 02:408-02:448, 02:450-02:453, 02:455-02:526, 02:528:01-02:570, 02:572-02:579, 02:581, 02:583-02:643N, 02:645-02:650, 02:652-02:661, 02:663-02:670, 02:672Q-02:700-02:706	A2, Low A2, A203, A210, Null, -				2																	
*02:17:01-02:17:04	A2, -		2				w															
*02:34-02:35:01, 02:35:03, 02:56:01-02:56:02, 02:62, 02:103, 02:580	A2, -		2										12								19	
*02:35:02	A2												12								19	
*02:46, 02:70, 02:651	A2, -		2																		19	
*02:48, 02:129, 02:571	A2, -																				19	
*02:55	A2, A28		2					7														20 21
*02:65, 02:407, 02:449	A2, A19, -		2																		18	
*02:78	-		w		4																19	
*02:81, 02:124	A2		2						8													
*02:135	-		2										12									
*02:149	-		2																			20
*02:169	-		2		4																	
*02:237	-		2											13	15							21
*02:243:01-02:243:02	-		2																			20
*02:309	-		2									11		13								20
*02:338	-		2		3																	
*02:454	-		2									11		13								
*02:527, 02:582	-		2					7														
*02:644	-							7		9												20
*02:662	-		2								10											
*02:671	-		2																			21
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21

*68:84	-						7	9											19			
*68:130:01	-						7		11										19			
*69:01:01-01-69:01:02, 69:03	A69, -						7												20			
*69:02	-						7		10										20			
*74:01:01-74:01:02, 74:01:04-74:03, 74:05-74:06, 74:08-74:12N, 74:14N-74:18, 74:20, 74:22, 74:24-74:28	A74, Null, -																		18			
*74:01:03	-								11										18			
*74:04, 74:21	A74, -																	17	18			
*74:07	A74																16		18			
*74:13	-									12						15			18			
*74:19	-				4														18			
*74:23	-			3															18			
*80:01:01-01-80:01:01:02, 80:03	A80, -						5													21		
*80:02	-						5		10													
B*07:81, B*08:52, B*18:67, B*38:41, B*38:67, B*51:185, B*53:05:01-53:05:02, B*53:16, B*53:33, B*57:60	B53, -							8														
C*03:365, C*12:131	-				4																	
C*06:187	-																			21		
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21