

<b>INTERPRETATION</b>		
<b>HLA-B*27 SSP typ</b>		
<b>Amplification pattern of the B*27:0</b>		
	<b>Well</b>	
	<b>1</b>	<b>2</b>
<b>Length of spec.</b>	<b>145</b>	<b>95</b>
<b>PCR product</b>		
<b>Length of int.</b>	<b>430</b>	<b>515</b>
<b>pos. control<sup>2</sup></b>		
<b>5'-primer<sup>3</sup></b>	<b>167</b>	<b>363</b>
	5' -gCT 3'	5' -AAT 3'
<b>3'-primer<sup>4</sup></b>	<b>272</b>	<b>418</b>
	5' -TgC 3'	5' -gTC 3'
<b>Well No.</b>	<b>1</b>	<b>2</b>
<b>HLA-B allele<sup>5</sup></b>		
<b>*27:01-27:05:08, 27:05:10-27:05:15, 27:08, 27:10, 27:13, 27:15, 27:17, 27:25, 27:28, 27:36-27:40, 27:42, 27:44-27:45, 27:47-27:69, 27:71-27:73</b>	<b>1</b>	<b>2</b>
<b>*27:05:09, 27:12, 27:16, 27:18, 27:23, 27:26, 27:29, 27:31, 37:02, 47:04-47:05</b>		<b>2</b>
<b>*27:06-27:07, 27:09, 27:11, 27:14, 27:19-27:21, 27:24, 27:27, 27:30, 27:32-27:35, 27:41, 27:43, 27:46, 27:70, 44:97</b>	<b>1</b>	
<b>HLA-B allele<sup>5</sup></b>		
<b>Well No.</b>	<b>1</b>	<b>2</b>

TABLE
ping
1 to 27: 73 alleles <sup>1</sup>
Length of spec.
PCR product
Length of int.
pos. control
5'-primer <sup>3</sup>
3'-primer <sup>4</sup>
Well No.
HLA-B allele <sup>5</sup>
*27:01-27:05:08, 27:05:10-27:05:15, 27:08, 27:10, 27:13, 27:15, 27:17, 27:25, 27:28, 27:36-27:40, 27:42, 27:44- 27:45, 27:47-27:69, 27:71-27:73
*27:05:09, 27:12, 27:16, 27:18, 27:23, 27:26, 27:29, 27:31, 37:02, 47:04-47:05
*27:06-27:07, 27:09, 27:11, 27:14, 27:19-27:21, 27:24, 27:27, 27:30, 27:32- 27:35, 27:41, 27:43, 27:46, 27:70, 44:97
HLA-B allele <sup>5</sup>
Well No.