

<b>INTEPRETATION TABLE</b>								
<b>DPA1 SSP typing</b>								
<b>Amplification patterns of the DPA1 alleles</b>								
	<b>Well<sup>5</sup></b>							
	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>
<b>Length of spec.</b>	<b>85</b>	<b>250</b>	<b>200</b>	<b>115</b>	<b>105</b>	<b>155</b>	<b>100</b>	<b>100</b>
<b>PCR product(s)</b>						<b>195</b>		
<b>Length of int.</b>	<b>515</b>	<b>515</b>	<b>430</b>	<b>430</b>	<b>430</b>	<b>515</b>	<b>430</b>	<b>430</b>
<b>pos. control<sup>1</sup></b>								
<b>5'-primer(s)<sup>2</sup></b>	<b>15</b>	<b>11</b>	<b>28</b>	<b>4</b>	<b>84</b>	<b>31</b>	<b>11</b>	<b>11</b>
	5'-ACg <sup>3'</sup>	5'-C gC <sup>3'</sup>	5'-gAA <sup>3'</sup>	5'-Cg g <sup>3'</sup>	5'-AAT <sup>3'</sup>	5'-g CA <sup>3'</sup>	5'-C gC <sup>3'</sup>	5'-C AT <sup>3'</sup>
						<b>43</b>		
						5'-TgT <sup>3'</sup>		
<b>3'-primer(s)<sup>3</sup></b>	<b>31</b>	<b>83</b>	<b>83</b>	<b>28</b>	<b>intr<sup>6</sup></b>	<b>69</b>	<b>31</b>	<b>31</b>
	5'-CAT <sup>3'</sup>	5'-ggT <sup>3'</sup>	5'-ggT <sup>3'</sup>	5'-TC g <sup>3'</sup>	5'-ggC <sup>3'</sup>	5'-gTC <sup>3'</sup>	5'-CTg <sup>3'</sup>	5'-CTg <sup>3'</sup>
						<b>83</b>		
						5'-ggT <sup>3'</sup>		
<b>Well No.</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>
<b>DPA1 allele<sup>4</sup></b>								
<b>*010301-010302, 010304</b>	<b>1</b>	<b>2</b>	<b>3</b>					
<b>*010303</b>		<b>2</b>	<b>3</b>					
<b>*0104</b>	<b>1</b>	<b>2</b>		<b>4</b>				
<b>*0105</b>	<b>1</b>				<b>5</b>			
<b>*010601</b>		<b>2</b>	<b>3</b>			<b>6</b>	<b>7</b>	
<b>*0107</b>	<b>1</b>	<b>2</b>	<b>3</b>					
<b>*0108</b>	<b>1</b>	<b>2</b>		<b>4</b>				
<b>*0109</b>	<b>1</b>	<b>2</b>	<b>3</b>					
<b>*0110</b>	<b>1</b>	<b>2</b>	<b>3</b>			<b>6</b>		
<b>*020101-020106</b>					<b>5</b>		<b>7</b>	
<b>*020201-020203</b>					<b>5</b>			<b>8</b>
<b>*0203</b>					<b>5</b>			
<b>*0204</b>					<b>5</b>	<b>6</b>		<b>8</b>
<b>*0301</b>			<b>3</b>					
<b>*0302</b>			<b>3</b>					
<b>*0303</b>				<b>4</b>				
<b>*0401</b>	<b>1</b>				<b>5</b>			
<b>DPA1 allele<sup>4</sup></b>								
<b>Well No.</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>

<b>INTERPRETATION TABLE</b>								
<b>DPA1 SSP typing</b>								
<b>Amplification patterns of the DPA1 alleles</b>								
<b>Well<sup>5</sup></b>								
<b>9</b>	<b>10</b>	<b>11</b>	<b>12</b>	<b>13</b>	<b>14</b>	<b>15</b>	<b>16</b>	
<b>205</b>	<b>85</b>	<b>90</b>	<b>205</b>	<b>135</b>	<b>140</b>	<b>245</b>	<b>220</b>	<b>Length of spec.</b>
								<b>PCR product(s)</b>
<b>430</b>	<b>515</b>	<b>515</b>	<b>430</b>	<b>430</b>	<b>430</b>	<b>430</b>	<b>430</b>	<b>Length of int.</b>
								<b>pos. control<sup>1</sup></b>
<b>11</b>	<b>15</b>	<b>66</b>	<b>18</b>	<b>51</b>	<b>50</b>	<b>15</b>	<b>23</b>	<b>5'-primer(s)<sup>2</sup></b>
5'-C AT <sup>3'</sup>	5'-ACC <sup>3'</sup>	5'-A TC <sup>3'</sup>	5'-gA A <sup>3'</sup>	5'-CA T <sup>3'</sup>	5'-C Cg <sup>3'</sup>	5'-ACC <sup>3'</sup>	5'-T AC <sup>3'</sup>	
	<b>15</b>					<b>15</b>		
	5'-ACC <sup>3'</sup>					5'-ACC <sup>3'</sup>		
<b>66</b>	<b>31</b>	<b>83</b>	<b>73</b>	<b>83</b>	<b>83</b>	<b>83</b>	<b>83</b>	<b>3'-primer(s)<sup>3</sup></b>
5'-T CA <sup>3'</sup>	5'-CAT <sup>3'</sup>	5'-ggT <sup>3'</sup>	5'-AgC <sup>3'</sup>	5'-ggT <sup>3'</sup>	5'-ggT <sup>3'</sup>	5'-ggT <sup>3'</sup>	5'-ggT <sup>3'</sup>	
<b>9</b>	<b>10</b>	<b>11</b>	<b>12</b>	<b>13</b>	<b>14</b>	<b>15</b>	<b>16</b>	<b>Well No.</b>
								<b>DPA1 allele<sup>4</sup></b>
								<b>*010301-010302,</b>
								<b>010304</b>
	<b>10</b>					<b>15</b>		<b>*010303</b>
								<b>*0104</b>
								<b>*0105</b>
								<b>*010601</b>
				<b>13</b>				<b>*0107</b>
					<b>14</b>			<b>*0108</b>
							<b>16</b>	<b>*0109</b>
								<b>*0110</b>
								<b>*020101-020106</b>
<b>9</b>								<b>*020201-020203</b>
	<b>10</b>							<b>*0203</b>
<b>9</b>								<b>*0204</b>
	<b>10</b>	<b>11</b>				<b>15</b>		<b>*0301</b>
<b>9</b>	<b>10</b>					<b>15</b>		<b>*0302</b>
	<b>10</b>	<b>11</b>				<b>15</b>		<b>*0303</b>
			<b>12</b>					<b>*0401</b>
								<b>DPA1 allele<sup>4</sup></b>
<b>9</b>	<b>10</b>	<b>11</b>	<b>12</b>	<b>13</b>	<b>14</b>	<b>15</b>	<b>16</b>	<b>Well No.</b>