

Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15
DQB1 allele <sup>1,2</sup>	ser. <sup>3</sup>															
*02:01:01:01-02:01:09, 02:01:11-02:01:13, 02:01:15-02:02:19, 02:04-02:76, 02:78-02:179, 02:181-02:206N	DQ2, Null -			3												
*02:03:02	-		2	3												
*02:77, 02:180	-			3			6									
*03:01:01:01-03:01:01:12, 03:01:01:14-03:01:01:48, 03:01:03-03:01:06, 03:01:09-03:01:56, 03:04:01:01-03:04:04, 03:09, 03:10:03, 03:16, 03:19:01:01-03:19:06, 03:21-03:22:02, 03:24, 03:27-03:29, 03:35-03:36, 03:42, 03:44, 03:46-03:57, 03:59-03:60, 03:69, 03:71, 03:76-03:77, 03:80, 03:82-03:84N, 03:92-03:94, 03:101-03:103, 03:108, 03:114-03:116, 03:119-03:122, 03:127-03:131, 03:133-03:135, 03:138-03:140, 03:142-03:144, 03:148, 03:150, 03:152, 03:157-03:160, 03:163-03:167, 03:169-03:173, 03:180, 03:182, 03:186-03:188, 03:191-03:193, 03:196-03:198:02, 03:201-03:202, 03:206-03:208, 03:216, 03:218-03:219, 03:231, 03:235-03:236, 03:241-03:243, 03:246, 03:253-03:257, 03:260, 03:264, 03:266-03:268, 03:271, 03:275-03:276N, 03:281, 03:284-03:286, 03:288, 03:290, 03:292-03:294, 03:297, 03:302-03:303N, 03:305-03:307, 03:309:01-03:309:02, 03:311-03:312, 03:317:01-03:318, 03:326-03:331, 03:338N, 03:340N-03:342, 03:347, 03:350, 03:353, 03:358N, 03:360-03:361, 03:366, 03:370, 03:372-03:373, 03:376N-03:378, 03:380-03:381, 03:385N, 03:387, 03:389-03:390, 03:394, 03:396, 03:399N-03:401, 03:404, 03:407N-03:408, 03:417, 03:419-03:421, 03:423-03:428, 03:430-03:432, 03:434-03:436, 03:438-03:439, 03:448, 03:451, 03:454-03:455, 03:458, 03:460, 03:465, 03:467-03:470, 03:472-03:473N, 03:475-03:476, 03:480Q, 03:482-03:483, 03:485-03:486, 03:488N, 03:491-03:492, 03:496-03:497, 03:499N	DQ7, Null -				4			7				12				
*03:01:02, 03:01:08, 03:73, 03:75, 03:109, 03:118N, 03:154, 03:162, 03:252, 03:291, 03:335, 03:354N, 03:391, 03:443	DQ7, Null -				4			7								
*03:01:07, 03:13, 03:58, 03:194, 03:232	DQ7, -				4			w					12			
*03:02:01:01-03:02:02:02, 03:02:04-03:02:05, 03:02:07-03:02:12, 03:02:14-03:02:17, 03:02:19-03:02:29, 03:02:31-03:02:36, 03:05:04, 03:07, 03:11, 03:18, 03:32, 03:45:01-03:45:02, 03:63-03:64, 03:66N-03:68, 03:85, 03:106, 03:146, 03:153, 03:161, 03:174-03:175, 03:178, 03:184-03:185, 03:190, 03:199, 03:203, 03:205, 03:210-03:211, 03:213N-03:215, 03:220-03:221, 03:224, 03:229, 03:240, 03:245, 03:247, 03:251, 03:261, 03:263:01:01-03:263:01:02, 03:265, 03:269N, 03:273-03:274, 03:277-03:279, 03:287, 03:289, 03:295-03:296, 03:298-03:301, 03:310N, 03:315, 03:320-03:324, 03:333-03:334N, 03:339N, 03:343-03:345, 03:348-03:349, 03:362, 03:364, 03:367-03:369, 03:379, 03:383, 03:386, 03:388, 03:392, 03:403N, 03:409-03:410, 03:412, 03:415-03:416, 03:422N, 03:429, 03:433, 03:437, 03:440N-03:442, 03:444, 03:446-03:447, 03:450, 03:452, 03:456-03:457, 03:459, 03:462-03:464, 03:466, 03:471, 03:479, 03:481, 03:484, 03:490, 03:493, 03:495, 03:498, 03:500	DQ8, Null, -					5		7					12			
*03:02:03, 03:02:13, 03:05:02, 03:37, 03:189, 03:204, 03:223, 03:226, 03:233, 03:237N, 03:352, 03:371	DQ8, Null, -					5		7								
*03:02:06, 03:02:18, 03:02:30, 03:125, 03:228	-					5							12			





