

# INTERPRETATION TABLE

## DRB1\*04 SSP subtyping

Allele-specific amplification patterns of the DRB1\*04:01 to \*04:89 alleles

|                           | Well <sup>4</sup>               |                                 |                                 |                                 |                                 |                                 |                                 |                                 |                                 |                                 |                                 |                                 |                                 |                                 |                                 |                                 |                                 |                                 |                                 |                                 |                                 |                                 |                                 |                                 |
|---------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
|                           | 1                               | 2                               | 3                               | 4                               | 5                               | 6                               | 7                               | 8                               | 9                               | 10                              | 11                              | 12                              | 13                              | 14                              | 15                              | 16                              | 17                              | 18                              | 19                              | 20                              | 21                              | 22                              | 23                              | 24                              |
| Length of spec.           | 215                             | 215                             | 225                             | 225                             | 170                             | 115                             | 260                             | 260                             | 215                             | 170                             | 175                             | 180                             | 170                             | 200                             | 215                             | 150                             | 190                             | 145                             | 140                             | 100                             | 155                             | 220                             | 120                             | 85                              |
| PCR product               |                                 |                                 |                                 |                                 |                                 |                                 |                                 |                                 |                                 |                                 | 220                             | 265                             | 225                             | 220                             |                                 |                                 |                                 |                                 | 210                             | 175                             | 205                             | 235                             | 235                             |                                 |
| Length of int.            | 515                             | 430                             | 430                             | 515                             | 430                             | 430                             | 430                             | 430                             | 430                             | 515                             | 430                             | 515                             | 430                             | 430                             | 430                             | 515                             | 430                             | 430                             | 515                             | 430                             | 430                             | 430                             | 430                             | 430                             |
| pos. control <sup>1</sup> |                                 |                                 |                                 |                                 |                                 |                                 |                                 |                                 |                                 |                                 |                                 |                                 |                                 |                                 |                                 |                                 |                                 |                                 |                                 | 255                             |                                 |                                 |                                 |                                 |
| 5'-primer(s) <sup>2</sup> | 5'-ACA <sup>3'</sup><br>13(125) | 5'-ACA <sup>3'</sup><br>13(125) | 5'-ACA <sup>3'</sup><br>13(125) | 5'-ACA <sup>3'</sup><br>13(125) | 5'-ACA <sup>3'</sup><br>13(125) | 5'-ACA <sup>3'</sup><br>13(125) | 5'-ACA <sup>3'</sup><br>13(125) | 5'-ACA <sup>3'</sup><br>13(125) | 5'-ACA <sup>3'</sup><br>13(125) | 5'-ACA <sup>3'</sup><br>13(125) | 5'-ACA <sup>3'</sup><br>13(125) | 5'-ACA <sup>3'</sup><br>13(125) | 5'-ACA <sup>3'</sup><br>13(125) | 5'-ACA <sup>3'</sup><br>13(125) | 5'-ACC <sup>3'</sup><br>13(125) | 5'-gTA <sup>3'</sup><br>37(197) | 5'-ACA <sup>3'</sup><br>13(125) | 5'-ACA <sup>3'</sup><br>13(125) | 5'-ACA <sup>3'</sup><br>13(125) | 5'-ACA <sup>3'</sup><br>13(125) | 5'-ACA <sup>3'</sup><br>13(125) | 5'-ACA <sup>3'</sup><br>13(125) | 5'-ACA <sup>3'</sup><br>13(125) | 5'-gAA <sup>3'</sup><br>71(299) |
| 3'-primer(s) <sup>3</sup> | 5'-CTT <sup>3'</sup><br>71(298) | 5'-CTC <sup>3'</sup><br>71(298) | 5'-CCT <sup>3'</sup><br>74(308) | 5'-CCg <sup>3'</sup><br>74(308) | 5'-gCT <sup>3'</sup><br>57(256) | 5'-Cgg <sup>3'</sup><br>37(197) | 5'-CAC <sup>3'</sup><br>86(344) | 5'-CCA <sup>3'</sup><br>86(344) | 5'-TCT <sup>3'</sup><br>71(298) | 5'-ATC <sup>3'</sup><br>57(256) | 5'-CCT <sup>3'</sup><br>58(260) | 5'-CTg <sup>3'</sup><br>59(262) | 5'-AgC <sup>3'</sup><br>56(253) | 5'-gAT <sup>3'</sup><br>67(286) | 5'-CTT <sup>3'</sup><br>71(298) | 5'-CCT <sup>3'</sup><br>74(308) | 5'-gTT <sup>3'</sup><br>60(265) | 5'-ggA <sup>3'</sup><br>47(227) | 5'-TAG <sup>3'</sup><br>46(225) | 5'-gTT <sup>3'</sup><br>33(184) | 5'-CCA <sup>3'</sup><br>51(239) | 5'-CCT <sup>3'</sup><br>72(302) | 5'-AgT <sup>3'</sup><br>39(203) | 5'-CAC <sup>3'</sup><br>86(344) |
| Well No.                  | 1                               | 2                               | 3                               | 4                               | 5                               | 6                               | 7                               | 8                               | 9                               | 10                              | 11                              | 12                              | 13                              | 14                              | 15                              | 16                              | 17                              | 18                              | 19                              | 20                              | 21                              | 22                              | 23                              | 24                              |
| Length of spec.           | 215                             | 215                             | 225                             | 225                             | 170                             | 115                             | 260                             | 260                             | 215                             | 170                             | 175                             | 180                             | 170                             | 200                             | 215                             | 150                             | 190                             | 145                             | 140                             | 100                             | 155                             | 220                             | 120                             | 85                              |
| PCR product               |                                 |                                 |                                 |                                 |                                 |                                 |                                 |                                 |                                 |                                 | 220                             | 265                             | 225                             | 220                             |                                 |                                 |                                 |                                 | 210                             | 175                             | 205                             | 235                             | 235                             |                                 |
| Well No.                  | 1                               | 2                               | 3                               | 4                               | 5                               | 6                               | 7                               | 8                               | 9                               | 10                              | 11                              | 12                              | 13                              | 14                              | 15                              | 16                              | 17                              | 18                              | 19                              | 20                              | 21                              | 22                              | 23                              | 24                              |

|   |          |          |          |          |          |          |          |          |          |           |           |           |           |           |           |           |           |           |           |           |           |           |           |           |  |
|---|----------|----------|----------|----------|----------|----------|----------|----------|----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|--|
| *04:01:01-04:01:06                              | 1        |          | 4        |          | 7        |          | 10       |          |          |           |           |           |           |           |           |           |           |           |           |           |           |           | 24        |           |  |
| *04:02  |          | 2        | 4        |          |          |          | 8        |          | 10       |           |           |           | 14        |           |           |           |           |           |           |           |           |           |           |           |  |
| *04:03:01-04:03:05                              |          |          | 3        |          |          |          | 8        | 9        | 10       |           |           |           |           |           | 16        |           |           |           |           |           |           |           |           |           |  |
| *04:04:01-04:04:04                              |          |          | 4        |          |          |          | 8        | 9        | 10       |           |           |           |           |           |           |           |           |           |           |           |           |           |           |           |  |
| *04:04:05                                       |          |          | 4        |          |          |          | 8        | w        | 10       |           |           |           |           |           |           |           |           |           |           |           |           |           |           |           |  |
| *04:05:01, 04:05:03-04:05:07, 04:05:09-04:05:10 |          |          | 4        | 5        | 7        |          | 9        |          |          |           |           |           |           |           |           |           |           |           |           |           |           |           |           |           |  |
| *04:05:02                                       |          |          | 4        | 5        | 7        |          | 9        |          |          |           |           |           |           |           |           |           |           |           |           |           |           |           |           |           |  |
| *04:05:08                                       |          |          | 4        | 5        | 7        |          | w        |          |          |           |           |           |           |           |           |           |           |           |           |           |           |           |           |           |  |
| *04:06:01-04:06:03                              |          |          | 3        |          | 6        | 8        | 9        | 10       |          |           |           |           |           |           |           |           |           |           |           |           |           |           |           |           |  |
| *04:07:01-04:07:04                              |          |          | 3        |          | 7        |          | 9        | 10       |          |           |           |           |           |           | 16        |           |           |           |           |           |           |           |           |           |  |
| *04:08:01-04:08:02                              |          |          | 4        |          | 7        |          | 9        | 10       |          |           |           |           |           |           |           |           |           |           |           |           |           |           |           |           |  |
| *04:09  | 1        |          | 4        | 5        | 7        |          |          |          |          |           |           |           |           |           |           |           |           |           |           |           |           |           | 24        |           |  |
| *04:10  |          |          | 4        | 5        |          | 8        | 9        |          |          |           |           |           |           |           |           |           |           |           |           |           |           |           |           |           |  |
| *04:11  |          |          | 3        |          | 5        |          | 8        | 9        |          |           |           |           |           |           | 16        |           |           |           |           |           |           |           |           |           |  |
| *04:12  |          |          |          | 5        |          | 8        | w        |          |          |           |           | 13        | 14        |           |           |           |           |           |           |           |           |           |           |           |  |
| *04:13  | 1        |          | 4        |          |          | 8        |          | 10       |          |           |           |           |           |           |           |           |           |           |           |           |           |           |           |           |  |
| *04:14  |          | 2        | 4        |          | 7        |          |          | 10       |          |           |           |           | 14        |           |           |           |           |           |           |           |           |           |           |           |  |
| *04:15  |          |          | 4        |          |          | 8        | 9        | 11       |          |           |           |           |           |           |           |           | 17        |           |           |           |           |           |           |           |  |
| *04:16  | 1        |          | 4        |          | 7        |          |          | 10       | 12       |           |           |           |           |           |           |           |           |           |           |           |           |           | 24        |           |  |
| *04:17:01-04:17:02                              |          |          | 3        |          | 5        | 7        | 9        |          |          |           |           |           |           |           | 16        |           |           |           |           |           |           |           |           |           |  |
| *04:18  |          |          |          |          |          |          | 8        | w        | 10       |           |           | 13        | 14        |           |           |           |           |           |           |           |           |           |           |           |  |
| *04:19  |          |          | 4        |          | 6        | 7        |          | 9        | 10       |           |           |           |           |           |           |           |           |           |           |           |           |           |           |           |  |
| *04:20  |          |          | 3        |          | 6        | 7        |          | 9        | 10       |           |           |           |           |           |           |           |           |           |           |           |           |           |           |           |  |
| *04:21  | 1        |          | 4        |          | 6        | 7        |          | 10       |          |           |           |           |           |           |           |           |           |           |           |           |           |           | 24        |           |  |
| *04:22  |          |          |          |          |          |          | 8        |          | 10       | 11        |           |           |           |           |           |           |           |           |           |           |           |           |           |           |  |
| *04:23  |          |          | 4        |          |          |          | 8        | 9        | 10       | 12        |           |           |           |           |           |           |           |           |           |           |           |           |           |           |  |
| *04:24  |          |          | 4        | 5        | 7        |          | 9        |          |          |           |           |           |           |           |           |           | 19        |           |           |           |           |           |           |           |  |
| *04:25  |          |          |          |          |          |          | 8        | w        | 10       |           |           | 13        |           |           |           | 17        |           |           |           |           |           |           |           |           |  |
| *04:26  | 1        |          | 4        |          | 7        |          |          | w        |          |           |           |           |           |           |           |           |           |           |           |           |           |           | 24        |           |  |
| *04:27  |          |          | 3        |          |          |          | 8        | 9        | 10       |           |           |           |           |           | 16        |           |           |           | 20        |           |           |           |           |           |  |
| *04:28  |          |          | 4        | 5        | 7        |          | 9        |          |          |           |           |           |           |           |           | 18        |           |           |           |           |           |           |           |           |  |
| *04:29  |          |          | 4        | 5        | 7        |          | 9        |          |          |           |           |           |           |           |           |           |           |           | 21        |           |           |           |           |           |  |
| *04:30  |          |          | 4        | 5        | 7        |          | 9        |          |          |           |           |           |           |           |           |           |           |           |           |           |           |           |           |           |  |
| *04:31  |          |          |          |          | 7        |          | w        | 10       |          |           |           | 13        |           |           |           |           |           |           |           |           |           |           |           |           |  |
| *04:32  |          |          | 4        |          |          | 8        | 9        | 10       |          |           |           |           |           |           |           |           |           |           |           |           | 22        |           |           |           |  |
| *04:33  | 1        |          | 4        |          | 7        |          |          | 10       |          |           |           |           |           |           |           |           |           |           |           |           |           | 23        | 24        |           |  |
| *04:34  |          |          |          |          |          |          |          |          |          |           |           |           | 15        |           |           |           |           |           |           |           |           |           | 24        |           |  |
| *04:35  | 1        |          | 4        |          | 7        |          |          | 10       |          |           |           |           |           |           |           | 18        |           |           |           |           |           |           | 24        |           |  |
| *04:36  |          |          | 4        |          |          | 8        | 9        | 10       |          |           |           |           |           |           |           | 17        |           |           |           |           |           |           |           |           |  |
| <b>Well No.</b>                                 | <b>1</b> | <b>2</b> | <b>3</b> | <b>4</b> | <b>5</b> | <b>6</b> | <b>7</b> | <b>8</b> | <b>9</b> | <b>10</b> | <b>11</b> | <b>12</b> | <b>13</b> | <b>14</b> | <b>15</b> | <b>16</b> | <b>17</b> | <b>18</b> | <b>19</b> | <b>20</b> | <b>21</b> | <b>22</b> | <b>23</b> | <b>24</b> |  |
| <b>Length of spec.</b>                          | 215      | 215      | 225      | 225      | 170      | 115      | 260      | 260      | 215      | 170       | 175       | 180       | 170       | 200       | 215       | 150       | 190       | 145       | 140       | 100       | 155       | 220       | 120       | 85        |  |
| <b>PCR product</b>                              |          |          |          |          |          |          |          |          |          |           | 220       | 265       | 225       | 220       |           |           |           |           | 210       | 175       | 205       | 235       | 235       |           |  |
|   |          |          |          |          |          |          |          |          |          |           |           |           |           |           |           |           |           |           |           | 255       |           |           |           |           |  |
| <b>Well No.</b>                                 | <b>1</b> | <b>2</b> | <b>3</b> | <b>4</b> | <b>5</b> | <b>6</b> | <b>7</b> | <b>8</b> | <b>9</b> | <b>10</b> | <b>11</b> | <b>12</b> | <b>13</b> | <b>14</b> | <b>15</b> | <b>16</b> | <b>17</b> | <b>18</b> | <b>19</b> | <b>20</b> | <b>21</b> | <b>22</b> | <b>23</b> | <b>24</b> |  |
| *04:37  |          | 2        | 4        |          |          |          |          | 8        |          | 10        |           |           |           |           |           |           |           |           |           |           |           |           |           | 21        |  |
| *04:38  | 1        |          | 4        |          |          | 7        |          |          |          | 10        |           |           | 14        |           |           |           |           |           |           |           |           | 21        |           | w         |  |
| *04:39  |          |          | 3        |          |          |          | 8        | 9        | 10       |           |           |           |           |           | 16        |           |           |           | 19        |           |           |           |           |           |  |

















