

HLA-A low screening 2018 8G1 R01

Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
HLA-A allele ^{1,2}	ser ³																								
*01:01:01:01-01:01:22, 01:01:24-01:01:47, 01:01:49-01:01:55, 01:01:57-01:01:64, 01:01:67-01:01:78, 01:01:80-01:04N, 01:06, 01:09:01:01-01:11N, 01:15N-01:18N, 01:20, 01:22N-01:25, 01:27N, 01:29, 01:31N-01:33, 01:35-01:42, 01:44-01:50, 01:52:01N-01:78, 01:80-01:94, 01:96-01:98, 01:100-01:120, 01:122-01:125, 01:127-01:135, 01:137-01:144, 01:146, 01:148, 01:150-01:158, 01:160N-01:166, 01:168-01:175, 01:177, 01:179N-01:191, 01:193, 01:195-01:199, 01:201-01:204, 01:206-01:207, 01:209-01:213, 01:215-01:227, 01:230-01:235, 01:237-01:243	A1, Null, -	1			4																				
*01:01:23, 01:01:48, 01:01:65-01:01:66, 01:01:79, 01:08, 01:14, 01:30, 01:79, 01:99, 01:145, 01:147Q, 01:149, 01:159, 01:167, 01:178N, 01:205, 01:208Q, 01:214, 01:228Q, 01:236	A1, Null, -	1																							
*01:01:56	-	1			4								12												
*01:07, 01:95	A1, -	1			4													17							
*01:12, 01:19, 01:200, 01:244	A1, -	1		3																					
*01:13, 01:176	A1, -	1			4							11		13		15									
*01:21, 01:126	A1, -	1		3	4																				
*01:26, 01:136, 01:192, 11:94, 11:112, 11:211 ⁴	A1, -	1			4							11													
*01:28, 01:229	-	1			4							11				15									
*01:43	-	1			4							10													
*01:51	-	1			4				8			10													
*01:121	-	1			4															19					
*01:194	-	1			4									13		15									
*02:01:01:01-02:01:15, 02:01:17-02:01:19, 02:01:21-02:01:81, 02:01:83-02:16, 02:18-02:22:02, 02:24:01-02:33, 02:36-02:45, 02:47, 02:49-02:54, 02:57-02:61, 02:63-02:64:02, 02:66-02:69, 02:71-02:77, 02:79:01-02:80, 02:82N-02:97:02, 02:99, 02:101:01-02:102, 02:104-02:123, 02:125N-02:128, 02:130-02:134, 02:136-02:148, 02:150-02:168, 02:170-02:236, 02:238-02:242, 02:244-02:308, 02:310-02:337, 02:339-02:406, 02:408-02:448, 02:450-02:453, 02:455-02:526, 02:528:01-02:570, 02:572-02:575, 02:577-02:579, 02:581, 02:583-02:643N, 02:645-02:650, 02:652-02:661, 02:663-02:670, 02:672Q-02:681, 02:683-02:706	A2, Low A2, A203, A210, Null, -		2																						
*02:17:01-02:17:04	A2, -		2					w																	
*02:34-02:35:01, 02:35:03, 02:56:01-02:56:02, 02:62, 02:103, 02:580	A2, -		2											13							20		22		
*02:35:02	A2													13							20		22		
*02:46, 02:70	A2		2																		20		22		
*02:48, 02:129, 02:571	A2, -																				20		22		
*02:55	A2, A28		2						8													21		23	
*02:65, 02:407, 02:449	A2, A19, -		2																	19					
*02:78	-		w		4																20		22		
*02:81, 02:124	A2		2								9														
*02:135	-		2											13											
*02:149	-		2																			21			
*02:169	-		2		4																				
*02:237	-		2												14		16								23
*02:243:01-02:243:02	-		2																18			21			
*02:309	-		2										12		14							21			
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
*02:338	-		2	3																					
*02:454	-		2										12		14										
*02:527, 02:582	-		2						8																
*02:576, 02:682	-		2																				22		
*02:644	-								8		10												21		
*02:651	-		2																			20			

Negative Control

*02:662	-		2								11														
*02:671	-		2																			23			
*03:01:01:01-03:01:18, 03:01:20-03:01:22, 03:01:24-03:01:37, 03:01:39-03:07:02, 03:09- 03:11N, 03:13-03:14, 03:16-03:17:02, 03:20- 03:23:02, 03:25-03:29, 03:31, 03:33-03:35, 03:37- 03:40, 03:42, 03:44:01-03:49, 03:51-03:56, 03:58, 03:60-03:62, 03:64-03:71, 03:73-03:74, 03:76- 03:81, 03:83-03:87, 03:90-03:94, 03:96- 03:103:01, 03:104-03:106, 03:109-03:110, 03:112- 03:123:01, 03:124-03:134, 03:136-03:141, 03:143- 03:151, 03:153-03:170, 03:174-03:175, 03:179- 03:185, 03:188-03:193, 03:195-03:197N, 03:201- 03:202, 03:204, 03:206-03:210, 03:212-03:214, 03:216-03:218, 03:220-03:224, 03:226-03:230, 03:231-02-03:245, 03:247-03:248, 03:250-03:251, 03:253-03:259, 03:261-03:266N, 03:268-03:272, 03:274-03:281, 03:283N-03:289	A3, Null, -		3									13													
*03:01:19, 03:103:02	-		3								12	13													
*03:01:23, 03:08, 03:36N, 03:57, 03:59, 03:72, 03:89:01-03:89:02, 03:107-03:108, 03:111, 03:142, 03:172-03:173, 03:176, 03:178N, 03:198- 03:199, 03:203, 03:205, 03:211, 03:252, 03:267	A3, Null, -		3																						
*03:01:38, 03:123:02, 03:171, 03:231:01	A3, Null, -		3									13	14												
*03:12	A3		3	4																					
*03:15, 03:19, 03:30	A3		3	5								13													
*03:18, 03:135	-	1		4								13													
*03:24, 03:50	A3, -		3				8					13													
*03:32, 03:219	-		3												17										
*03:41	-		3																			23			
*03:43, 03:82, 03:186	-		3									13	15	17											
*03:63	-		3								11	13										23			
*03:75, 03:177	-		3									13										23			
*03:88	-		3	4							11											23			
*03:95	-											13	14	16								23			
*03:152, 03:273	-		3	5											17										
*03:187	-											13										22			
*03:215, 03:225, 03:249, 25:46 ⁴	-											13													
*03:246, 03:282	-		3									13			19										
*11:01:01:01-11:02:06, 11:04-11:09, 11:12-11:16, 11:18-11:19, 11:21N-11:24:02, 11:26-11:27, 11:29-11:39, 11:41-11:52Q, 11:54-11:59, 11:61- 11:93, 11:95-11:111, 11:114-11:115N, 11:117- 11:129, 11:131-11:138, 11:141-11:154, 11:156- 11:157, 11:159-11:161, 11:163-11:165, 11:167- 11:174, 11:176-11:182Q, 11:184-11:190, 11:192- 11:198, 11:200-11:208N, 11:210N, 11:212- 11:221, 11:224-11:225, 11:227-11:268	A11, Null, -			4								11													
*11:03, 11:20, 11:25:01-11:25:02, 11:60, 11:158, 11:175, 11:209	A11, -		3	4								11													
*11:10, 11:191, 11:269	A11, -			4			8					11													
*11:11	-			4								11	12												
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
*11:17, 11:40, 11:223	-			4								10	11												
*11:113, 11:162	-			4								11			15										
*11:116, 11:140, 11:199:01, 11:222, 66:23 ⁴	-			4								11		13											
*11:130	-			3	4							11		14										23	
*11:139	-			4	5	7						11													
*11:155	-			4								11												22	
*11:166	-			4	6							11													
*11:183	-			3	4				8			11													
*11:199:02	-			4								11		13											20
*11:226	-	1		4								11													22
*23:01:01:01-23:01:12, 23:01:14-23:13, 23:15- 23:56, 23:58-23:63, 23:65, 23:67-23:68, 23:70- 23:83	A23, Null, -				5	6																			
*23:01:13	-				5	6																		20	

Negative Control

*23:14:01-23:14:02, 24:24, 24:71, 24:315, 24:392 ⁴	A23, A24, A9, -				5	6	7																				
*23:57	-				5																						21
*23:64	-				5	6											17										
*23:66, 24:14:01:01-24:15, 24:51-24:53, 24:57, 24:64, 24:94, 24:114, 24:138, 24:188, 24:222N, 24:228, 24:291, 24:296, 24:304, 24:316, 24:324, C*04:01:03 ⁴	A24, Null, Cw4, -				5																						
*24:02:01:01-24:11N, 24:13:01-24:13:02, 24:17, 24:20:01:01-24:23, 24:25-24:43, 24:45N-24:50, 24:54-24:56, 24:58-24:63, 24:66-24:70, 24:72-24:81, 24:83N-24:91, 24:93, 24:95-24:103, 24:105-24:113, 24:115-24:137, 24:139-24:187, 24:189-24:203, 24:205-24:210, 24:212, 24:214-24:221, 24:223-24:227, 24:229-24:242, 24:244-24:290, 24:292-24:295, 24:297-24:303N, 24:305-24:314, 24:317-24:323N, 24:325-24:339, 24:341-24:391, 24:393-24:394	A24, Low A24, A2403, A9, Null, -				5		7																				
*24:18, 24:204, 24:213	A24, A3, -				5		7																				23
*24:19, 24:44	A9				5		7		11																		
*24:82	-				5		7												18								21
*24:92	-		3		5																						
*24:104, 24:243	-				5		7												17								
*24:340	-				5		7																				20
*25:01:01:01-25:02, 25:04, 25:07-25:16, 25:18-25:29, 25:31-25:33, 25:35, 25:37-25:42N, 25:44-25:45	A25, Null, -							8	9			13															
*25:03	A25							8	9			13					17										
*25:05	A25							8	9		12	13															20
*25:06	A25							8	9		12																
*25:17, 25:34	-								9			13															
*25:30	-							8	9			13					17										21
*25:36	-							8	9																		
*25:43	-			4				8	9																		
*26:01:01:01-26:01:44, 26:08, 26:10-26:15, 26:17-26:18, 26:20, 26:23-26:28, 26:31-26:43:02, 26:45-26:48, 26:50-26:53, 26:55-26:63, 26:65-26:71N, 26:73-26:74, 26:76-26:77, 26:79-26:88, 26:90, 26:93-26:110, 26:112-26:119, 26:121, 26:123-26:149	A26, A10, Null, -							8		10		13															
*26:02:01-26:02:02, 26:29, 26:49	A26, -							8		10		w															
*26:03:01, 26:06, 26:21, 26:111	A26, -							8		11		13															23
*26:04, 26:75	A26, -							8		10																	
*26:05, 26:30	A26							8				13															23
*26:07:01-26:07:02, 26:64, 26:89	A26, -									10		13															
*26:09, 26:91	A26, -							8		10	12																
*26:16	A26						7			10		13															
*26:19	-			4						10		13	14														
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24		
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24		
*26:22	A26								8	10			13	14													21
*26:54	-								8	10	12	13															20
*26:72	-				4					10		13															
*26:78	-								8		11																23
*26:92	-										11		13														
*26:120	-	1							8	10		13															
*26:122	-								8	10		13															20
*29:01:01:01-29:06, 29:08N-29:12, 29:15-29:27, 29:29-29:38, 29:40-29:47, 29:50-29:66, 29:68-29:79, 29:81-29:97, 29:99-29:104, 29:106-29:111	A29, Null, -													14													
*29:07, 29:49	A29, -					5	6						14														
*29:13, 29:39, 29:98	-												14			17											
*29:14	-												14		16												
*29:28	-							8					14														
*29:48	-												14				18										
*29:67	-				4								14														
*29:80	-												14														21
*29:105	-									10			14						18								

Negative Control

HLA-A low screening 2018 8G1 R01

*66:12, 66:15	-							8				13												Negative Control				
*66:19	-							8			11	12	13															
*68:01:01:01-68:04:02, 68:06-68:08:01, 68:09-68:10, 68:12, 68:13:02-68:14, 68:16-68:18N, 68:21:01-68:24, 68:26-68:28, 68:30-68:40, 68:42-68:44, 68:46-68:65, 68:67-68:70, 68:72-68:83, 68:85-68:102, 68:104:01-68:104:02, 68:106-68:108, 68:110-68:116, 68:118-68:129, 68:131-68:133, 68:135, 68:137-68:156, 68:158-68:174	A68, A28, Null, -							8												20								
*68:05, 68:15, 68:20, 68:109, 68:136	A68, -							8												20			23					
*68:08:02, 68:105, 68:157	A68, -							8												20		21						
*68:11N, 68:41	Null, -							8												20			22					
*68:13:01, 68:66:01-68:66:02, 68:134	A68, -				4			8												20								
*68:19, 68:130:02	A68, -							8					14							20								
*68:25	A68							8											19		20							
*68:29	A68							8										18		20	21							
*68:45, 68:117	-						7	8											20									
*68:71	-							8					13							20								
*68:84	-							8			10								20									
*68:103:01-68:103:02	-				3			8											20									
*68:130:01	-							8				12							20									
*69:01:01:01-69:01:02, 69:03	A69, -							8														21						
*69:02	-							8			11											21						
*74:01:01-74:01:02, 74:01:04-74:03, 74:05-74:06, 74:08-74:12N, 74:14N-74:18, 74:20, 74:22, 74:24-74:28	A74, Null, -																					19						
*74:01:03	-											12										19						
*74:04, 74:21	A74, -																				18	19						
*74:07	A74																				17	19						
*74:13	-												13		16							19						
*74:19	-				4																	19						
*74:23	-				3																	19						
*80:01:01:01-80:01:01:02, 80:03	A80, -							6															23					
*80:02	-							6				11																
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24			
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24			
B*07:81, B*08:52, B*18:67, B*38:41, B*38:67, B*51:185, B*53:05:01-53:05:02, B*53:16, B*53:33, B*57:60	B53, -									9																		Negative Control
B*18:27	-							6																				
B*40:359, B*57:65, C*04:31, C*06:137, C*07:569	-																								22			
C*03:365, C*12:131	-				4																							
C*06:187	-																									23		
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24			