

INTERPRETATION TABLE

HLA-A*11 SSP subtyping

Amplification patterns of the A*11:01 to A*11:139 alleles

	Well ¹⁵																																	
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32		
Length of spec.	195	270	190	235	210	180	195	215	100	85	120	215	235	195	100	215	100	95	205	235	175	175	100	150	80	295	220	305	210	240	80	235	230	205
PCR product(s)					210	550	215	100			280	215	235		215	205			265	235	175	175		230		470		370	260	280		270	285	295
Length of int. pos control ¹	800	1070	1070	1070	800	1070	1070	800	800	800	1070	800	1070	800	1070	800	1070	1070	1070	1070	800	1070	1070	1070	800	1070	1070	1070	800	1070	1070	1070	1070	1070
5'-primer(s) ²	5'-ggg ³ 127	5'-TCC ³ 28	5'-ATA ³ 363	5'-ATA ³ 363	5'-ATA ³ 363	5'-CTC ³ 292	5'-ATA ³ 363	5'-TgC ³ 553	5'-TgC ³ 527	5'-TCA ³ 97	5'-ATA ³ 363	5'-ATA ³ 363	5'-ACC ³ 102	5'-CAG ³ 282	5'-CAG ³ 282	5'-ACA ³ 391	5'-ATA ³ 363	5'-TCC ³ 97	5'-CCT ³ 103	5'-CCT ³ 103	5'-CAC ³ 282	5'-CCT ³ 103	5'-gCA ³ 489	5'-gga ³ 341	5'-TCC ³ 97	5'-gga ³ 341	5'-CCT ³ 103	5'-CCT ³ 103	5'-CCA ³ 103	5'-AAC ³ 261	5'-CCT ³ 103	5'-TCA ³ 97	5'-ATA ³ 363	
3'-primer(s) ³	5'-gAC ³ 282	5'-CTT ³ 127	5'-CTg ³ 502	5'-CgT ³ 559	5'-CTC ³ 502	5'-CCg ³ 559	5'-TCC ³ 539	5'-CgT ³ 616	5'-CCg ³ 570	5'-CgT ³ 341	5'-TCT ³ 539	5'-CCg ³ 559	5'-gCA ³ 257	5'-Cgt ³ 341	5'-Cgg ³ 341	5'-CCg ³ 559	5'-CCg ³ 526	5'-gAC ³ 282	5'-CCA ³ 236	5'-CCA ³ 236	5'-Cgt ³ 341	5'-CCA ³ 214	5'-CCT ³ 527	5'-TgA ³ 353	5'-gAC ³ 282	5'-TCA ³ 362	5'-CAT ³ 271	5'-Cgt ³ 341	5'-CCA ³ 299	5'-CCg ³ 299	5'-CgT ³ 341	5'-CCA ³ 527		
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32		
HLA-A allele ^{4,5}	*11:01:01-11:01:20, 11:01:22-11:01:43, 11:01:45, 11:82-11:85, 11:87-11:88, 11:91-11:93, 11:95, 11:97, 11:103-11:105, 11:107, 11:112, 11:114, 11:116-11:117, 11:119, 11:122-11:123, 11:125- 11:126 , 11:128, 11:131-11:133, 11:135-11:136, 11:138 ⁶																																	
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32		
Length of spec.	195	270	190	235	210	180	195	215	100	85	120	215	235	195	100	215	100	95	205	235	175	175	100	150	80	295	220	305	210	240	80	235	230	205
PCR product(s)					210	550	215	100			280	215	235		215	205			265	235	175	175		230		470		370	260	280		270	285	295

*02:56:01-02:56:02, 02:62, 02:78, 02:103, 03:01:30, 03:40, 03:75, 03:95, 03:98, 03:105, 68:10, 68:13:01-68:14, 68:66, 74:13	1				15																														
*02:65, 02:80, 02:114, 02:152, 02:246, 02:279, 02:298, 02:304, 29:48, 33:08-33:09, 74:04, B*57:04, B*57:32																																			
*02:99																				22 23															
*02:101:01					6																														
*02:117																				w															
*02:154					6															22															
*02:156, 02:338	4				15																			32											
*02:294, 32:54	5																																		
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32			
Length of spec.	195	270	190	235	210	180	195	215	100	85	120	215	235	195	100	100	95	205	235	175	175	100	150	80	295	220	305	260	210	280	240	80	235	230	205
PCR product(s)					210	550					280					215	205			265	235		230		470		370	260	280			270	285	295	
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32			
*03:01:01:01- 03:01:01:03, 03:01:03- 03:01:08, 03:01:10- 03:01:22, 03:01:24- 03:01:29, 03:01:31- 03:01:37, 03:03N- 03:05:01, 03:06, 03:09, 03:11N-03:17:02, 03:19- 03:20, 03:25-03:30, 03:33-03:35, 03:37- 03:39, 03:43-03:49, 03:51-03:56, 03:58- 03:62, 03:64, 03:67- 03:68N, 03:70-03:72, 03:74, 03:77, 03:79- 03:81, 03:83-03:86, 03:91N-03:94, 03:96, 03:99-03:104, 03:109- 03:110, 03:112, 03:114- 03:121, 03:124-03:132, 03:134, 03:136-03:141, 03:143-03:151, 03:155- 03:159, 03:162N	1				15																														
*03:01:02, 03:05:02	1	4																		15				25											
*03:01:09, 03:01:23, 03:08, 03:36N, 03:57, 03:107-03:108, 03:111, 03:142, 03:152, 03:161N, 32:04																				15															
*03:01:38	1	4																		15				23											

*30:01:01-30:01:06, 30:02:01-30:02:04, 30:02:06-30:03, 30:07- 30:08, 30:10-30:13, 30:15-30:16, 30:19- 30:20, 30:22, 30:24- 30:26, 30:30-30:34, 30:36-30:38, 30:40- 30:42, 30:44-30:45, 30:48-30:50, 30:53, 30:56-30:71			3	4											15	17																														
*30:01:07, 30:23, 30:27N			3	4											15																															
*30:02:05			3	4													17											25																		
*30:04:01-30:04:02, 30:06, 30:29, 30:46			3	4			7																																				32			
*30:09, 31:04			3	4											15																												32			
*30:14L, 30:43			3												15	17																														
*30:17			3	4			7								15																												32			
*30:18, 30:55	1		3	4											15	17																														
*30:35, 30:54				4											15	17																														
*30:39			3	4											15	17	18																													
*30:51			3	4											15	17																											32			
*31:03			3	4					10						15														21														31	32		
*31:65, B*44:36																																														
*33:13									10																				21														31			
*33:49			3	4											15																												29	32		
*34:01:01-34:01:02, 34:05	1								10					w																												29				
*34:02:01, 34:02:03, 34:04, 34:08, 34:10N	1		w	4					10					14	15																											29				
*34:02:02	1		w	4					10					14																												29				
*34:03	1		3	4					10					14	15																											29				
*34:06, 66:06-66:07, 66:10, 66:18	1								10					14																												29				
*34:07	1		w	4					10					15																												29				
*34:09	1		w	4										15																												29				
*36:01-36:02, 36:05				4			7																																					21		
*43:01													13																																20	21
*66:01, 66:04, 66:09, 66:13-66:14, 66:19	1								10					14																														29		
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32														
Length of spec.	195	270	190	235	210	180	195	215	100	85	120	215	235	195	100	100	95	205	235	175	175	100	150	80	295	220	305	210	240	80	270	235	230	295	205											
PCR product(s)																																														
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32														

													58					*11:115N															
														63				*11:120															
43																	*11:121																
													59					*11:124															
													60					*11:127N, 11:137N ¹⁴															
																	64																
																	*11:129																
																	*11:139																
																	*01:01:01:01-01:01:22, 01:01:24-01:01:27, 01:01:29-01:01:47, 01:01:49-01:02, 01:09, 01:11N, 01:16N, 01:18N, 01:22N-01:24, 01:27N- 01:33, 01:35-01:40, 01:42, 01:44, 01:46- 01:50, 01:52N-01:55, 01:57N-01:60, 01:62, 01:64, 01:67-01:70, 01:74-01:78, 01:80- 01:82, 01:84-01:88, 01:90-01:99, 01:103, 01:106-01:108, 01:110- 01:114, 01:116-01:124																
																	*01:01:23, 01:01:48, 01:08, 01:15N																
																	*01:01:28, 01:07, 01:71, 01:100, 01:105																
																	*01:03, 36:03, 74:10																
																	*01:04N																
																	*01:06																
																	*01:10																
41																	*01:12, 01:19																
																	*01:13																
																	*01:14																
																	*01:17																
																	*01:20, 01:63, 01:66, 01:73, 01:89, 02:346																
33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		
110	85	105	65	105	120	280	90	205	175	445	145	210	120	55	200	185	90	265	230	195	190	120	230	245	90	370	105	130	75	275	340		
280	180	290	380	290	225	280	185	560	215	210	210	220	310	265	545	230	195	265	235	230	245	90	370	265	130	75	275	340					
305	265																																
33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64		
41																	*01:21, 01:126																
																	*01:25																
																	*01:26																
																	*01:41																
																	*01:43																
													52					*01:45															
																	*01:51																
42																	*01:56N																
																	*01:61																

																																*66:02, 66:16, 68:01:01:01-68:01:17, 68:02:01:01-68:02:05, 68:06-68:09, 68:11N- 68:12, 68:16-68:19, 68:21:01-68:23, 68:25- 68:29, 68:32-68:35, 68:37-68:39, 68:41, 68:43-68:45, 68:47- 68:53, 68:55-68:56, 68:58-68:60, 68:62, 68:64-68:65, 68:67- 68:71, 68:73-68:74, 68:76-68:83, 68:86- 68:89, 68:91, 68:94N- 68:96, 68:98-68:99, 69:01						
																																62	*66:05, 66:15					
																																62	*66:08					
																																62	*66:11					
																																62	*66:12					
																																		52	62	*66:17		
																																	*68:01:18, 68:24, 68:40, 68:72, 68:85, 68:93, 68:97					
																																	*68:15					
																																		43				
																																	*68:30, 68:75					
																																	*68:42, 68:54					
																																	*68:57, 68:92					
																																	*68:61					
																																	*68:63					
																																	*80:01:01:01-80:02					
																																	<i>B*07:102, B*08:18, B*35:178, B*73:01- 73:02</i>					
																																	<i>B*07:117</i>					
																																	<i>B*39:75</i>					
																																	<i>C*06:72, C*06:90, C*07:277</i>					
																																	<i>C*07:221</i>					
33																																						
33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64							
																																	Well No.					
																																	HLA-A allele ^{4,5}					