

INTERPRETATION TABLE

HLA-A*26 SSP subtyping

Amplification patterns of the A*26:01 to 26:72 alleles

	Well ¹⁴																							
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
Length of spec.	180	80	140	260	180	80	150	135	145	85	140	135	110	145	100	430	175	125	190	210	115	100	130	305
PCR product(s)						160			190	260		240	255					205	245		205	220	165	
Length of int.	800	800	###	800	###	###	800	###	###	###	800	###	###	###	800	###	###	800	800	###	###	###	###	###
pos. control ¹																								
5'-primer(s) ²	418	261	418	74	423	261	423	423	423	78	243	98	228	453	282	28	423	257	355	392	395	160	652	341
	5'-Agg ^{3'}	5'-AAC ^{3'}	5'-AgA ^{3'}	5'-C ^{3'}	5'-gCT ^{3'}	5'-AAC ^{3'}	5'-gCT ^{3'}	5'-gCT ^{3'}	5'-gCT ^{3'}	5'-TCT ^{3'}	5'-CCT ^{3'}	5'-CTC ^{3'}	5'-ATg ^{3'}	5'-AAA ^{3'}	5'-CAG ^{3'}	5'-TCg ^{3'}	5'-gCT ^{3'}	5'-CCC ^{3'}	5'-CCg ^{3'}	5'-CgA ^{3'}	5'-gCA ^{3'}	5'-ACT ^{3'}	5'-CTg ^{3'}	5'-ggA ^{3'}
				78		423				257	416	423	346					423	410		652	280		
				5'-TCC ^{3'}		5'-gCT ^{3'}				5'-Cgg ^{3'}	5'-gCg ^{3'}	5'-gCT ^{3'}	5'-gTA ^{3'}					5'-gCT ^{3'}	5'-gTT ^{3'}		5'-CTg ^{3'}	5'-CCC ^{3'}		
3'-primer(s) ³	559	299	517	292	560	299	524	517	527	299	341	299	299	559	341	282	559	341	559	559	559	341	743	362
	5'-CCg ^{3'}	5'-TCg ^{3'}	5'-CgT ^{3'}	5'-gTg ^{3'}	5'-ACA ^{3'}	5'-TCT ^{3'}	5'-CAC ^{3'}	5'-CgC ^{3'}	5'-CCA ^{3'}	5'-TCg ^{3'}	5'-CgT ^{3'}	5'-TCg ^{3'}	5'-TCg ^{3'}	5'-CCg ^{3'}	5'-CgT ^{3'}	5'-gAC ^{3'}	5'-CgT ^{3'}	5'-CgT ^{3'}	5'-CCg ^{3'}	5'-CCg ^{3'}	5'-CCg ^{3'}	5'-CgT ^{3'}	5'-TCC ^{3'}	5'-TCA ^{3'}
	559					542	538	519	570		517	517	559	559		290		589			728		776	
	5'-CCT ^{3'}					5'-CTT ^{3'}	5'-CTg ^{3'}	5'-ggA ^{3'}	5'-CCg ^{3'}		5'-CgT ^{3'}	5'-CgC ^{3'}	5'-CCg ^{3'}	5'-CCg ^{3'}		5'-gAA ^{3'}		5'-CTT ^{3'}			5'-CCT ^{3'}		5'-CAA ^{3'}	
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
HLA-A allele ^{4,5}																								
*26:01:01-26:01:22, 26:27 ¹⁵	1	2																						
*26:02		2	3																					
*26:03:01-26:03:02	1			4																				
*26:04		2			5																			
*26:05	1					6																		
*26:06				4							11													
*26:07:01-26:07:02	1									10														
*26:08, 26:67	1	2					7																	
*26:09		2																17						
*26:10	1	2												14										
*26:11N	1	2						8																
*26:12	1	2							9															
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
Length of spec.	##	80	##	##	##	80	##	##	##	85	##	##	##	##	##	##	##	##	##	##	##	##	##	##
PCR product(s)						##			##	##		##	##					##	##		##	##	##	

Length of spec.	180	80	140	260	180	80	150	135	145	85	140	135	110	145	100	430	175	125	190	210	115	100	130	305
PCR product(s)						160			190	260		240	255					205	245		205	220	165	
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
*02:81, 02:87, 02:112, 02:124, 02:129, 02:136, 23:01:01-23:08N, 23:11N- 23:12, 23:14-23:50, 24:02:01:01-24:02:32, 24:02:34-24:03:02, 24:05- 24:06, 24:08-24:09N, 24:11N, 24:13:01-24:15, 24:17-24:18, 24:20-24:23, 24:25, 24:27, 24:29-24:43, 24:45N, 24:47-24:64, 24:66-24:88, 24:90N- 24:99, 24:101-24:108, 24:110-24:111, 24:113- 24:128, 24:130, 24:132N- 24:180, 24:182-24:183N, 24:185N-24:190																								
*02:135	1																							
*02:154														14										
*02:241, 33:31																					21			
*02:309									9								17							
*03:01:19								8				12					17							
*03:09, 03:108, 11:18, 29:19, 30:13, 30:16, 30:44, 30:46, 33:24, 68:05, 68:15, 68:20, 74:06				4																				
*11:06				4												15								
*11:26																15								24
*11:91																15								
*23:09, 24:181																								
*23:10, 24:10, 24:46														14										
*23:13, 24:07, 24:24, 24:112, 24:131, 68:36																	16							
*24:129												12												
*25:01:01-25:01:05, 25:03, 25:07, 25:10, 25:12N, 25:14-25:16	1																							
*25:02	1															15								
*25:04	1						7																	
*25:05	1																							
*25:06																		17						
*25:08	1																							
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
Length of spec.	180	80	140	260	180	80	150	135	145	85	140	135	110	145	100	430	175	125	190	210	115	100	130	305
PCR product(s)						160			190	260		240	255					205	245		205	220	165	
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24

INTERPRETATION TABLE

HLA-A*26 SSP subtyping

Amplification patterns of the A*26:01 to 26:72 alleles

Well¹⁴

		25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	
		360	150	275	90	330	125	210	90	105	135	110	125	245	190	220	90	140	460	135	105	365	155	425	Length of spec.
		###	###	###	###	###	###	800	###	###	###	###	###	###	###	###	###	###	###	###	###	###	###	###	PCR product(s)
		###	###	###	###	###	###	###	###	###	###	###	###	###	###	###	###	###	###	###	###	###	###	###	Length of int.
		###	###	###	###	###	###	###	###	###	###	###	###	###	###	###	###	###	###	###	###	###	###	###	pos. control ¹
		5'-ggc ^{3'}	5'-cta ^{3'}	5'-cca ^{3'}	5'-ctc ^{3'}	5'-tcg ^{3'}	5'-cta ^{3'}	5'-acc ^{3'}	5'-gct ^{3'}	5'-cca ^{3'}	5'-gct ^{3'}	5'-gcc ^{3'}	5'-gct ^{3'}	5'-tcg ^{3'}	5'-cga ^{3'}	5'-tcg ^{3'}	5'-aac ^{3'}	5'-cct ^{3'}	5'-tcg ^{3'}	5'-gct ^{3'}	5'-gct ^{3'}	5'-tcg ^{3'}	5'-gct ^{3'}	5'-tcg ^{3'}	5'-primer(s) ²
		423	98	106	292	28	98	391	423	401	423	136	423	28	191	28	261	2 nd	28	423	423	28	423	28	
		5'-gct ^{3'}	5'-cca ^{3'}	5'-cgt ^{3'}	5'-cgt ^{3'}	5'-tcc ^{3'}	5'-tcc ^{3'}	5'-ccg ^{3'}	5'-ggc ^{3'}	5'-ccg ^{3'}	5'-cca ^{3'}	5'-cgt ^{3'}	5'-tgt ^{3'}	5'-cca ^{3'}	5'-cgt ^{3'}	5'-ggc ^{3'}	5'-gga ^{3'}	5'-tcc ^{3'}	5'-gcg ^{3'}	5'-cgc ^{3'}	5'-gcg ^{3'}	5'-tct ^{3'}	5'-cag ^{3'}	5'-tgc ^{3'}	3'-primer(s) ³
		418	208	341	341	186	180	559	472	559	518	341	506	265	341	77	311	426	319	517	486	224	538	283	
		559	559	743	290	538	570	595	559	559	570	595	559	559	271	265	271	426	319	517	486	224	538	283	
		5'-cca ^{3'}	5'-cca ^{3'}	5'-cca ^{3'}	5'-gaa ^{3'}	5'-caa ^{3'}	5'-ccg ^{3'}	5'-cca ^{3'}	5'-cca ^{3'}	5'-cca ^{3'}	5'-ccg ^{3'}	5'-cca ^{3'}	5'-cca ^{3'}	5'-ccg ^{3'}	5'-cag ^{3'}	5'-cca ^{3'}	5'-cca ^{3'}	5'-tcc ^{3'}	5'-gcg ^{3'}	5'-cgc ^{3'}	5'-gcg ^{3'}	5'-tct ^{3'}	5'-cag ^{3'}	5'-tgc ^{3'}	
25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	Well No.		
		HLA-A allele ^{4,5}																							
		*26:01:01-26:01:22, 26:27 ¹⁵ *26:02																							
		*26:03:01-26:03:02 *26:04																							
		*26:05 *26:06																							
		*26:07:01-26:07:02 *26:08, 26:67																							
		*26:09																							
		*26:10																							
		*26:11N																							
		*26:12																							
25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	Well No.		
##	##	##	90	##	##	##	90	##	##	##	##	##	##	##	##	90	##	##	##	##	##	##	##	Length of spec.	
##	##	##	##	##	##	##	##	##	##	##	##	##	##	##	##	##	##	##	##	##	##	##	##	PCR product(s)	

360	150	275	90	330	125	210	90	105	135	110	125	245	215	125	410	245	190	220	90	140	460	135	105	365	155	425	Length of spec.										
	175				235		155	200	190	245	215	410	220					410										PCR product(s)									
25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	Well No.														
																	42											*02:81, 02:87, 02:112, 02:124, 02:129, 02:136, 23:01:01-23:08N, 23:11N-23:12, 23:14-23:50, 24:02:01:01-24:02:32, 24:02:34-24:03:02, 24:05-24:06, 24:08-24:09N, 24:11N, 24:13:01-24:15, 24:17-24:18, 24:20-24:23, 24:25, 24:27, 24:29-24:43, 24:45N, 24:47-24:64, 24:66-24:88, 24:90N-24:99, 24:101-24:108, 24:110-24:111, 24:113-24:128, 24:130, 24:132N-24:180, 24:182-24:183N, 24:185N-24:190									
25																												*02:135									
																												*02:154									
																												*02:241, 33:31									
25							32																					*02:309									
25							32											43										*03:01:19									
																												*03:09, 03:108, 11:18, 29:19, 30:13, 30:16, 30:44, 30:46, 33:24, 68:05, 68:15, 68:20, 74:06									
			28																									*11:06									
					30																							*11:26									
															39													*11:91									
																	42											*23:09, 24:181									
																	42											*23:10, 24:10, 24:46									
																	42											*23:13, 24:07, 24:24, 24:112, 24:131, 68:36									
															39													*24:129									
																	42											*25:01:01-25:01:05, 25:03, 25:07, 25:10, 25:12N, 25:14-25:16									
																	42											*25:02									
												36					42											*25:04									
																	42											*25:05									
																	42											*25:06									
										34							42											*25:08									
25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	Well No.														
360	150	275	90	330	125	210	90	105	135	110	125	245	215	125	410	245	190	220	90	140	460	135	105	365	155	425	Length of spec.										
	175				235		155	200	190	245	215	410	220					410										PCR product(s)									
25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	Well No.														

