

Length of spec.	235	225	205	215	185	275	225	240	230	65	225	190	145	105	225	275	215	135	230	60	130	250	95	470	155	190	95	495	100	130	200	150			
PCR product(s)										205				235				210		175	195	155	210	185			240	145			210	155	230	270	
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32			
HLA-A allele ⁴																																			
*03:01:01:01, 03:01:01:03, 03:01:03- 03:01:05, 03:01:07, 03:01:09-03:01:11, 03:01:13, 03:01:15- 03:01:22, 03:01:24- 03:01:29, 03:01:31- 03:01:38, 03:26, 03:77, 03:80, 03:83-03:85, 03:92, 03:94, 03:100-03:101, 03:110, 03:114-03:117, 03:119-03:121, 03:124- 03:126, 03:131, 03:136- 03:139, 03:141, 03:144- 03:148, 03:151, 03:154- 03:159 ¹⁹	1	2	3					8									17																29		
*03:01:01:02N	1	2	3					8					13				17																29		
*03:01:02, 03:01:06, 03:01:14, 03:140	1	2						8									17																29		
*03:01:08, 03:109	1	2	3														17																29		
*03:01:12, 03:01:23	1	w	3					8									17																29		
*03:01:30, 03:98, 03:105		2	3					8									17																29		
*03:02:01-03:02:03, 03:82, 03:90, 03:106, 03:160	1	2		4				8		10																								29	
*03:03N	1	2	3		5			8									17																29		
*03:04:01-03:04:03	1	2	3			6		8									17																29		
*03:05:01		2	3					8	9								17																29		
*03:05:02, 30:55 ⁵		2						8	9								17																29		
*03:06	1	2	3				7	8									17																29		
*03:07	1	2						8		10							17																29		
*03:08	1		3					8			11						17																		
*03:09	1	2	3					8				12					17																29		
*03:10		2		4				8	9	10																								29	
*03:11N	1	2	3					8						14			17																29		
*03:12	1		3												15		17																	29	
*03:13	1	2	3					8									17	18															29		
*03:14	1	2	3					8								16	17																	29	
*03:15	1	2	3					8									17		19														29		
*03:16, 03:61 ⁷	1	2	3					8									17			20													29		
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32			
Length of spec.	235	225	205	215	185	275	225	240	230	65	225	190	145	105	225	275	215	135	230	60	130	250	95	470	155	190	95	495	100	130	200	150			
PCR product(s)										205				235				210		175	195	155	210	185			240	145			210	155	230	270	

Length of spec.	235	225	205	215	185	275	225	240	230	65	225	190	145	105	225	275	215	135	230	60	130	130	250	95	470	155	190	95	495	100	130	200	150							
PCR product(s)										205				235				210		175	195	155			240	210	185		240	145		210	155	230	270	220	150			
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	32							
*02:01:01:01-02:01:08, 02:01:10-02:01:15, 02:01:17- 02:01:19, 02:01:21-02:01:48, 02:01:50-02:01:81, 02:01:83- 02:04, 02:07:01-02:07:02, 02:09, 02:11:01, 02:11:03- 02:13, 02:15N-02:18, 02:20:01-02:20:02, 02:22:01- 02:22:02, 02:24:01-02:27, 02:29-02:33, 02:38-02:40:02, 02:42-02:43N, 02:45-02:47, 02:49, 02:52-02:53N, 02:55, 02:58-02:60:02, 02:63-02:71, 02:73-02:76:01, 02:77, 02:80, 02:82N-02:83N, 02:85-02:86, 02:88N-02:89, 02:92- 02:97:02, 02:101:01-02:102, 02:104-02:105, 02:107, 02:109-02:111, 02:113N- 02:121, 02:123, 02:125N, 02:128, 02:130-02:135, 02:138-02:141, 02:145, 02:147-02:153, 02:155, 02:157-02:163, 02:165- 02:168, 02:171:01-02:171:02, 02:173-02:177, 02:181- 02:183, 02:185-02:228, 02:230-02:231, 02:233-02:236, 02:238- 02:243, 02:245-02:246, 02:249-02:254, 02:256- 02:258, 02:260-02:270, 02:272-02:277, 02:279- 02:285, 02:287-02:289, 02:291-02:294, 02:296- 02:299, 02:301N-02:309, 02:311-02:323, 02:325- 02:327, 02:329, 02:332, 02:334-02:336, 02:339- 02:343, 02:345-02:354, 02:356N-02:357, 02:360- 02:372, 02:374-02:375, 02:377-02:379, 02:381, 02:383-02:386, 02:388- 02:392, 31:01:02-31:01:05, 31:01:07-31:02, 31:06, 31:09, 31:11-31:18, 31:20, 31:22- 31:32, 31:34-31:68, 31:70, 33:01:01-33:01:07, 33:03:01- 33:03:05, 33:03:07-33:07, 33:09, 33:11-33:12, 33:14- 33:16, 33:18:01-33:18:02, 33:20-33:23, 33:25-33:37, 33:39-33:47, 33:50, 33:52- 33:64, 74:05								8																																
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	32							

Length of spec.	235	225	205	215	185	275	225	240	230	65	225	190	145	105	225	275	215	135	230	60	130	250	95	470	155	190	95	495	100	130	200	150				
PCR product(s)										205				235				210		175	155		240	210	185		240	145			210	155	230	220	150	
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32				
*02:01:09, 02:50, 02:76:02								8				12														25										
*02:05:05, 02:06:07, 02:122												12														25										
*02:11:02, 33:08								8																		25										
*02:19, 02:36-02:37, 02:255, 24:44, 31:05, 33:10, 33:19								8																												
*02:34-02:35:02, 02:56:01-02:56:02, 02:62, 02:103		w						8																									29			
*02:35:03		w						8																		25						29				
*02:48, 02:90, 74:01-74:03, 74:08-74:09, 74:11-74:12N, 74:14N-74:17								8			11																									
*02:54, 23:01:01-23:03:02, 23:05-23:06, 23:08N-23:12, 23:14:01-23:15, 23:17-23:30, 23:32-23:41, 23:44, 23:46-23:54, 24:02:01:01-24:02:04, 24:02:06-24:02:28, 24:02:30-24:02:63, 24:04-24:06, 24:09N, 24:13:01-24:14, 24:17, 24:20-24:21:02, 24:25-24:27, 24:29-24:32, 24:34-24:40N, 24:43, 24:45N-24:50, 24:52-24:64, 24:66-24:74, 24:76-24:88, 24:90N-24:91, 24:93, 24:95-24:103, 24:105-24:107, 24:109-24:111, 24:113-24:117, 24:119-24:124, 24:126-24:130, 24:133-24:137, 24:139-24:166, 24:168-24:203, 24:205-24:206, 24:208-24:209, 24:212, 24:214, 24:216-24:218, 24:220-24:222N																																				
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32				

*33:49									8	9	10								17																					
*34:02:01, 34:02:03- 34:04, 34:07, 34:10N			3							9						15			17																					
*34:02:02										9						15			17																					
*34:08		2	3						8	9									17																		26			
*34:09			3							9						15			17																				29	
*36:01, 36:05										9		11																												
*36:02	1											11							18								23									28		31		
*66:17																15																								
*68:26, 68:65																15																							29	
*68:71		2							8																												26		29	
*74:04									8			11																								25				
*80:01:01:01- 80:01:01:02												w								17																			22	
*80:02									8			11								17																			22	
<i>C*07:35, C*08:48</i>																																								
<i>C*07:86</i>																																								
HLA-A allele⁴																																								
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32								

INTERPRETATION TABLE														
HLA-A*03 SSP subtyping														
Amplification patterns of the A*03:01 to A*03:161N alleles														
Well¹⁸														
	5'-T 3'	343	5'-CCC 3'	265	5'-gCT 3'	239								
			5'-TCT 3'	539	5'-gTT 3'	181								
			5'-TgT 3'	330	5'-CAA 3'	290								
			5'-CCA 3'	265	5'-CgA 3'	175								
			5'-ACg 3'	281	5'-CCC 3'	175	5'-Cgg 3'	152						
			5'-gAA 3'	574	5'-CCA 3'	518	5'-Cgg 3'	487						
					5'-AgT 3'	845	5'-CTC 3'	259						
					5'-ACA 3'	899								
					5'-Tgg 3'	506	5'-CTA 3'	331						
					5'-CCA 3'	555								
					5'-ACA 3'	899	5'-gCT 3'	316						
					5'-Cag 3'	561	5'-CCA 3'	518	5'-gCg 3'	486				
					5'-CCT 3'	527	5'-TTT 3'	327						
							5'-CgT 3'	559						
							5'-TCg 3'	299						
							5'-CgT 3'	559						
							5'-gAC 3'	282						
							5'-CCA 3'	555						
							5'-gCA 3'	605	5'-CCA 3'	538				
									5'-CgT 3'	143				
							5'-gTg 3'	583	5'-AgA 3'	545				
									5'-CTA 3'	296				
									5'-gCA 3'	777				
									5'-ACC 3'	270				
									5'-ACA 3'	899				
									5'-AgA 3'	545	5'-CCg 3'	527		
									5'-Cag 3'	298				
									5'-ACA 3'	899				
									5'-AgT 3'	490				
									5'-ACA 3'	899				
									5'-CAC 3'	524	5'-Tgg 3'	497		
											5'-CCg 3'	570		
33														
34														
35														
36														
37														
38														
39														
40														
41														
42														
43														
44														
45														
46														
47														
48														
49														
50														
51														
52														
53														
54														
55														
56														
57														
58														
59														
60														
61														
62														
63														
64														

Length of spec.
PCR product(s)

Length of int. pos.
control¹

5'-primer(s)²

3'-primer(s)³

Well No.
HLA-A allele⁴

285	210	180	160	125	275	235	210	160	120	225	120	95	250	195	165	235	200	240	275	185	220	190	260	150	240	195	165	270	150	190	110	240	145	110	70	545	195	130	285	215	85	260	225	240	175	215	275	180	215	240	390	170	190	200	175	225	Length of spec. PCR product(s)	
33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	Well No.																										
																																	*02:01:01:01-02:01:08, 02:01:10-02:01:15, 02:01:17-02:01:19, 02:01:21-02:01:48, 02:01:50-02:01:81, 02:01:83-02:04, 02:07:01-02:07:02, 02:09, 02:11:01, 02:11:03-02:13, 02:15N-02:18, 02:20:01-02:20:02, 02:22:01-02:22:02, 02:24:01-02:27, 02:29-02:33, 02:38-02:40:02, 02:42-02:43N, 02:45-02:47, 02:49, 02:52-02:53N, 02:55, 02:58-02:60:02, 02:63-02:71, 02:73-02:76:01, 02:77, 02:80, 02:82N-02:83N, 02:85-02:86, 02:88N-02:89, 02:92-02:97:02, 02:101:01-02:102, 02:104-02:105, 02:107, 02:109-02:111, 02:113N-02:121, 02:123, 02:125N, 02:128, 02:130-02:135, 02:138-02:141, 02:145, 02:147-02:153, 02:155, 02:157-02:163, 02:165-02:168, 02:171:01-02:171:02, 02:173-02:177, 02:181-02:183, 02:185-02:228, 02:230-02:231, 02:233-02:236, 02:238-02:243, 02:245-02:246, 02:249-02:254, 02:256-02:258, 02:260-02:270, 02:272-02:277, 02:279-02:285, 02:287-02:289, 02:291-02:294, 02:296-02:299, 02:301N-02:309, 02:311-02:323, 02:325-02:327, 02:329, 02:332, 02:334-02:336, 02:339-02:343, 02:345-02:354, 02:356N-02:357, 02:360-02:372, 02:374-02:375, 02:377-02:379, 02:381, 02:383-02:386, 02:388-02:392, 31:01:02-31:01:05, 31:01:07-31:02, 31:06, 31:09, 31:11-31:18, 31:20, 31:22-31:32, 31:34-31:68, 31:70, 33:01:01-33:01:07, 33:03:01-33:03:05, 33:03:07-33:07, 33:09, 33:11-33:12, 33:14-33:16, 33:18:01-33:18:02, 33:20-33:23, 33:25-33:37, 33:39-33:47, 33:50, 33:52-33:64, 74:05																									
33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	Well No.																										

		285	210	180			Length of spec.																												
		160	125	275	235			PCR product(s)																											
		210	160	120																															
		225	120	95																															
		250	195	165																															
		235	200	240																															
		275	185	220	190																														
		260	150	240	195	165																													
		270	150	190	110	240																													
		145	110	70																															
		545	195	130	285	215	85																												
		260	225	240	175	215																													
		275	180	215	240	390	170	190	175	225																									
												Well No.																							
33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	*02:78, 66:02, 66:12, 66:16, 68:01:01:01-68:02:05, 68:06-68:14, 68:16-68:19, 68:21:01-68:25, 68:27-68:28, 68:30, 68:32-68:35, 68:37-68:39, 68:41-68:44, 68:47-68:56, 68:58-68:64, 68:66-68:70, 68:73-68:83, 68:86-68:89, 68:91-68:92, 68:94N-68:96, 69:01			
																																		*02:81, 02:124, 02:146, 26:20, 32:01:01-32:02, 32:06, 32:09-32:14, 32:16, 32:18-32:26, 32:28-32:36, 32:38-32:41, 32:43-32:51, 32:53-32:55	
													46																63					*02:156, 02:338 *02:164, 31:69	
												46		48																				*02:237 *02:380, 31:19	
																																		*11:01:01-11:02:01, 11:02:03-11:03, 11:05, 11:07, 11:09-11:11, 11:13-11:15:02, 11:20, 11:22, 11:29, 11:32-11:34, 11:36-11:37, 11:41-11:42, 11:46, 11:48-11:49, 11:51-11:52Q, 11:54-11:56, 11:58-11:59, 11:61-11:72, 11:74-11:77, 11:79-11:87, 11:89-11:91, 11:93, 11:95-11:97, 11:99N-11:100, 11:102-11:105, 11:107, 11:109N-11:110, 11:112-11:117, 11:120, 11:122-11:129, 11:131-11:132	
		285	210	180			Length of spec.																												
		160	125	275	235			PCR product(s)																											
		210	160	120																															
		225	120	95																															
		250	195	165																															
		235	200	240																															
		275	185	220	190																														
		260	150	240	195	165																													
		270	150	190	110	240																													
		145	110	70																															
		545	195	130	285	215	85																												
		260	225	240	175	215																													
		275	180	215	240	390	170	190	175	225																									
												Well No.																							

