







*01:01:01:01-01:01:51, 01:02-01:04N, 01:07- 01:11N, 01:14-01:18N, 01:21-01:24, 01:26-01:29, 01:31N-01:33, 01:35- 01:50, 01:52N, 01:54- 01:64, 01:67-01:69:01, 01:70-01:71, 01:74-01:78, 01:80-01:83:01, 01:84- 01:97, 01:99, 01:101, 01:103-01:122																		<b>18</b>			<b>w</b>			
*01:01:52-01:01:53, 01:06, 01:25, 01:53N, 01:65, 01:69:02, 01:83:02, 01:100, 01:123N-01:125, 23:53, 29:03, 80:01:01-01- 80:02																		<b>18</b>						
*01:12, 01:19, 11:38, 24:17, 24:41, 24:208, 29:33									<b>10</b>									<b>18</b>						
*01:13, 03:18	<b>1</b>																	<b>18</b>			<b>w</b>			
*01:20, 01:66, 01:73																		<b>18</b>	<b>20</b>		<b>w</b>			
*01:51	<b>2</b>																	<b>18</b>			<b>w</b>			
*01:72, 01:102, 36:04																					<b>w</b>			
*01:98																		<b>18</b>			<b>22</b>			
*01:126																		<b>18</b>						
<b>Well No.</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>	<b>9</b>	<b>10</b>	<b>11</b>	<b>12</b>	<b>13</b>	<b>14</b>	<b>15</b>	<b>16</b>	<b>17</b>	<b>18</b>	<b>19</b>	<b>20</b>	<b>21</b>	<b>22</b>	<b>23</b>	<b>24</b>
<b>Length of spec.</b>	<b>#</b>	<b>#</b>	<b>#</b>	<b>#</b>	<b>#</b>	<b>#</b>	<b>#</b>	<b>#</b>	<b>#</b>	<b>#</b>	<b>#</b>	<b>#</b>	<b>#</b>	<b>#</b>	<b>#</b>	<b>#</b>	<b>#</b>	<b>#</b>	<b>#</b>	<b>9</b>	<b>#</b>	<b>#</b>	<b>#</b>	<b>#</b>
<b>PCR product(s)</b>				<b>#</b>	<b>#</b>			<b>#</b>	<b>#</b>		<b>#</b>	<b>#</b>	<b>#</b>		<b>#</b>			<b>#</b>			<b>#</b>			
									<b>#</b>															
<b>Well No.</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>	<b>9</b>	<b>10</b>	<b>11</b>	<b>12</b>	<b>13</b>	<b>14</b>	<b>15</b>	<b>16</b>	<b>17</b>	<b>18</b>	<b>19</b>	<b>20</b>	<b>21</b>	<b>22</b>	<b>23</b>	<b>24</b>
*02:02, 02:05:01-02:05:05, 02:08, 02:47, 02:63, 02:102, 02:115, 02:154- 02:155, 02:172, 02:179, 02:186, 02:209, 02:232, 02:286, 02:320, 02:324, 02:337, 02:344, 02:359, 02:373N																		<b>w</b>						
*02:03:01-02:03:04, 02:148, 02:253, 02:258, 02:264, 02:281, 02:315, 02:323, 02:345, 02:355, 02:370, 02:402																		<b>17</b>						
*02:05:06, 02:22:01- 02:22:02, 02:50, 02:104, 02:136, 02:191, 02:230, 02:267, 02:376, 02:382, 23:02, 24:06, 24:138, 24:167, 31:30, 32:10, B*27:25, C*02:07																		<b>17</b>						



*03:01:01:01-03:01:08, 03:01:10-03:01:22, 03:01:24-03:01:29, 03:01:31-03:01:37, 03:02:01-03:03N, 03:04:02- 03:07, 03:09, 03:11N- 03:16, 03:17:02, 03:20- 03:23:02, 03:25-03:31, 03:33-03:35, 03:37-03:39, 03:42, 03:44-03:49, 03:51- 03:52, 03:54-03:56, 03:58- 03:65, 03:67-03:81, 03:83- 03:88, 03:90-03:94, 03:96, 03:99-03:102, 03:104, 03:106, 03:109-03:110, 03:112-03:121, 03:123- 03:134, 03:136-03:141, 03:143-03:151, 03:153- 03:160, 03:162N, 11:35, 11:130	1																							23	
<b>Well No.</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>	<b>9</b>	<b>10</b>	<b>11</b>	<b>12</b>	<b>13</b>	<b>14</b>	<b>15</b>	<b>16</b>	<b>17</b>	<b>18</b>	<b>19</b>	<b>20</b>	<b>21</b>	<b>22</b>	<b>23</b>	<b>24</b>	
<b>Length of spec.</b>	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	#	9	#	#	#	#
<b>PCR product(s)</b>				#	#			#	#		#	#	#		#			#			#				
									#																
<b>Well No.</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>	<b>9</b>	<b>10</b>	<b>11</b>	<b>12</b>	<b>13</b>	<b>14</b>	<b>15</b>	<b>16</b>	<b>17</b>	<b>18</b>	<b>19</b>	<b>20</b>	<b>21</b>	<b>22</b>	<b>23</b>	<b>24</b>	
*03:01:09, 03:01:23, 03:08, 03:32, 03:36N, 03:41, 03:57, 03:89, 03:107-03:108, 03:111, 03:142, 03:152, 03:161N, 29:01:01:01-29:02:01:02, 29:02:03-29:02:13, 29:04, 29:06-29:21, 29:23-29:27, 29:29-29:31, 29:34-29:50, 30:26, 31:03-31:04, 31:06, 32:04, 32:30, 32:32																								23	
*03:01:38	1																							23	
*03:04:01	1																	18						23	
*03:10, 11:04	1								10															23	
*03:17:01	1													15										23	
*03:24	1	2																						23	
*03:95	1																			20				23	
*03:97, 03:122, 03:135, 11:94	1																	18							





*34:02:03	1	2							9													23		
*34:03-34:04, 34:07, 34:09	1	2							?													23		
*34:05	1															17								
*34:06	1								?															
*34:08	1	2							?													23		
*34:10N	1																					23		
*36:01, 36:03, 36:05																					w	23		
*66:06	1								9															
*66:09	1															17						24		
*66:10	1															17								
*66:15					5											17								
*66:19	1															17				21				
*69:01	1																					24		
B*44:03:19																								
<b>Well No.</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>	<b>9</b>	<b>10</b>	<b>11</b>	<b>12</b>	<b>13</b>	<b>14</b>	<b>15</b>	<b>16</b>	<b>17</b>	<b>18</b>	<b>19</b>	<b>20</b>	<b>21</b>	<b>22</b>	<b>23</b>	<b>24</b>

## INTERPRETATION TABLE

### HLA-A\*68 SSP subtyping

Amplification patterns of the A\*68:01 to A\*68:99 and the A\*69:01 alleles

Well <sup>8</sup>																							
25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46		
240	100	185	200	80	145	220	175	150	165	140	195	285	145	180	195	200	130	365	185	150	110	Length of spec.	
	185	235	230	250	245																	PCR product(s)	
1070	1070	1070	1070	800	1070	1070	800	1070	1070	1070	1070	1070	1070	1070	1070	1070	1070	1070	1070	1070	1070	Length of int.	
																						pos. control <sup>1</sup>	
5'-TCg <sup>3'</sup>	5'-AAA <sup>3'</sup>	5'-ATA <sup>3'</sup>	5'-gCT <sup>3'</sup>	5'-ACT <sup>3'</sup>	5'-ACT <sup>3'</sup>	5'-CCT <sup>3'</sup>	5'-AAA <sup>3'</sup>	5'-CCg <sup>3'</sup>	5'-ACT <sup>3'</sup>	5'-CCg <sup>3'</sup>	5'-CCg <sup>3'</sup>	5'-CCT <sup>3'</sup>	5'-CCT <sup>3'</sup>	5'-gCg <sup>3'</sup>	5'-AAA <sup>3'</sup>	5'-CTC <sup>3'</sup>	5'-CgA <sup>3'</sup>	5'-TCg <sup>3'</sup>	5'-gAA <sup>3'</sup>	5'-CCg <sup>3'</sup>	5'-ggA <sup>3'</sup>	5'-primer(s) <sup>2</sup>	
28	453	87	121	102	102	2 <sup>nd</sup> I	453	413	160	413	106	2 <sup>nd</sup> I	2 <sup>nd</sup> I	397	453	98	194	28	666	413	383		
5'-gAA <sup>3'</sup>	5'-CCA <sup>3'</sup>	5'-CTg <sup>3'</sup>	5'-gCA <sup>3'</sup>	5'-CAG <sup>3'</sup>																			
376																							
5'-gCA <sup>3'</sup>																							
5'-ggT <sup>3'</sup>	5'-ACT <sup>3'</sup>	5'-gAC <sup>3'</sup>	5'-gAC <sup>3'</sup>	5'-TCT <sup>3'</sup>	5'-TCT <sup>3'</sup>	5'-TgC <sup>3'</sup>	5'-Tgg <sup>3'</sup>	5'-CAA <sup>3'</sup>	5'-gAC <sup>3'</sup>	5'-CCT <sup>3'</sup>	5'-gTT <sup>3'</sup>	5'-CCg <sup>3'</sup>	5'-gCg <sup>3'</sup>	5'-CCA <sup>3'</sup>	5'-gCC <sup>3'</sup>	5'-gTT <sup>3'</sup>	5'-gAC <sup>3'</sup>	5'-TCT <sup>3'</sup>	5'-CAT <sup>3'</sup>	5'-ggA <sup>3'</sup>	5'-TCT <sup>3'</sup>	3'-primer(s) <sup>3</sup>	
97	512	282	282	308	308	506	578	524	282	512	259	570	430	538	605	259	282	224	809	521	453		
	5'-AgA <sup>3'</sup>	5'-CAT <sup>3'</sup>	5'-CCA <sup>3'</sup>	5'-ggA <sup>3'</sup>	5'-CgT <sup>3'</sup>		5'-CCg <sup>3'</sup>																
				527																			
25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	Well No.	
240	100	185	200	80	145	220	175	150	165	140	195	285	145	180	195	200	130	365	185	150	110	Length of spec.	
	185	235	230	250	245																	PCR product(s)	
25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	Well No.	
																						HLA-A allele <sup>4</sup>	





																				*01:01:01:01-01:01:51, 01:02-01:04N, 01:07- 01:11N, 01:14-01:18N, 01:21-01:24, 01:26-01:29, 01:31N-01:33, 01:35- 01:50, 01:52N, 01:54- 01:64, 01:67-01:69:01, 01:70-01:71, 01:74-01:78, 01:80-01:83:01, 01:84- 01:97, 01:99, 01:101, 01:103-01:122			
																				*01:01:52-01:01:53, 01:06, 01:25, 01:53N, 01:65, 01:69:02, 01:83:02, 01:100, 01:123N-01:125, 23:53, 29:03, 80:01:01:01- 80:02			
																				*01:12, 01:19, 11:38, 24:17, 24:41, 24:208, 29:33 *01:13, 03:18 *01:20, 01:66, 01:73 *01:51 *01:72, 01:102, 36:04 *01:98 *01:126			
																			45				
<b>25</b>	<b>26</b>	<b>27</b>	<b>28</b>	<b>29</b>	<b>30</b>	<b>31</b>	<b>32</b>	<b>33</b>	<b>34</b>	<b>35</b>	<b>36</b>	<b>37</b>	<b>38</b>	<b>39</b>	<b>40</b>	<b>41</b>	<b>42</b>	<b>43</b>	<b>44</b>	<b>45</b>	<b>46</b>	<b>Well No.</b>	
##	##	##	##	80	##	##	##	##	##	##	##	##	##	##	##	##	##	##	##	##	##	##	<b>Length of spec.</b>
	##	##	##	##	##																		<b>PCR product(s)</b>
<b>25</b>	<b>26</b>	<b>27</b>	<b>28</b>	<b>29</b>	<b>30</b>	<b>31</b>	<b>32</b>	<b>33</b>	<b>34</b>	<b>35</b>	<b>36</b>	<b>37</b>	<b>38</b>	<b>39</b>	<b>40</b>	<b>41</b>	<b>42</b>	<b>43</b>	<b>44</b>	<b>45</b>	<b>46</b>	<b>Well No.</b>	
				29																			*02:02, 02:05:01-02:05:05, 02:08, 02:47, 02:63, 02:102, 02:115, 02:154- 02:155, 02:172, 02:179, 02:186, 02:209, 02:232, 02:286, 02:320, 02:324, 02:337, 02:344, 02:359, 02:373N *02:03:01-02:03:04, 02:148, 02:253, 02:258, 02:264, 02:281, 02:315, 02:323, 02:345, 02:355, 02:370, 02:402 *02:05:06, 02:22:01- 02:22:02, 02:50, 02:104, 02:136, 02:191, 02:230, 02:267, 02:376, 02:382, 23:02, 24:06, 24:138, 24:167, 31:30, 32:10, B*27:25, C*02:07









