

INTERPRETATION TABLE

HLA-B low resolution SSP

Amplification patterns of the B*07:02 to B*83:01 alleles¹⁵

		Well																								
		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	
Length of spec.		110	215	140	130	185	190	290	165	165	180	105	280	195	105	115	80	150	135	60	210	170	110	395	160	
PCR product(s)				235	265	235			330	220	390	190			435	395	150			400	245				425	
Length of int. pos. control ¹		800	1070	1070	800	800	800	1070	1070	1070	800	1070	1070	1070	1070	1070	1070	800	1070	1070	800	800	1070	1070	1070	
5'-primer(s) ²		5'-gAC ^{3'} 540	5'-TCg ^{3'} 97	5'-AgC ^{3'} 363	5'-CCg ^{3'} 103	5'-AAT ^{3'} 363	5'-CCg ^{3'} 103	5'-Agg ^{3'} 418	5'-Tgg ^{3'} 357	5'-AAA ^{3'} 435	5'-TTA ^{3'} 420	5'-AgA ^{3'} 206	5'-gTC ^{3'} 419	5'-TCT ^{3'} 142	5'-gAC ^{3'} 540	5'-Cgg ^{3'} 161	5'-gCT ^{3'} 167	5'-TCA ^{3'} 355	5'-gAC ^{3'} 206	5'-gTT ^{3'} 368	5'-TCT ^{3'} 142	5'-gTC ^{3'} 368	5'-TTA ^{3'} 420	5'-gAA ^{3'} 246	5'-Tgg ^{3'} 357	
				5'-AgC ^{3'} 363	5'-AgT ^{3'} 361					5'-ATA ^{3'} 412										5'-ggA ^{3'} 557						
3'-primer(s) ³		5'-gCT ^{3'} 605	5'-TgA ^{3'} 272	5'-gTg ^{3'} 309	5'-CgT ^{3'} 193	5'-TAT ^{3'} 246	5'-TAT ^{3'} 246	5'-gCC ^{3'} 572	5'-gTC ^{3'} 538	5'-CAG ^{3'} 559	5'-CTC ^{3'} 559	5'-Tgg ^{3'} 272	5'-TAT ^{3'} 3 ^d 1	5'-gTC ^{3'} 301	5'-ATC ^{3'} 309	5'-TgC ^{3'} 272	5'-TCT ^{3'} 234	5'-TCT ^{3'} 204	5'-gCT ^{3'} 463	5'-ggC ^{3'} 302	5'-gTT ^{3'} 259	5'-ggg ^{3'} 311	5'-ggT ^{3'} 544	5'-ggA ^{3'} 317	5'-TgC ^{3'} 272	5'-ggC ^{3'} 302
			5'-TgA ^{3'} 272	5'-CAG ^{3'} 559	5'-gTg ^{3'} 309	5'-CTC ^{3'} 559			5'-gTg ^{3'} 583	5'-gTC ^{3'} 570				5'-gCg ^{3'} 572	5'-ATC ^{3'} 309	5'-TgC ^{3'} 272	5'-TgA ^{3'} 272				5'-gTT ^{3'} 259	5'-ggg ^{3'} 311				5'-ggT ^{3'} 302
					5'-gTg ^{3'} 583	5'-CTC ^{3'} 559				5'-gTg ^{3'} 583				5'-gCg ^{3'} 572	5'-gTC ^{3'} 570	5'-gCT ^{3'} 605	5'-Tgg ^{3'} 272				5'-TCC ^{3'} 266	5'-gTC ^{3'} 538				5'-gCg ^{3'} 477
														5'-gCg ^{3'} 572							5'-TCC ^{3'} 266					
																					5'-Tgg ^{3'} 292	5'-gTA ^{3'} 292				
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	

Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
HLA-B allele ⁴	ser. ⁵																								
*07:02:01-07:02:29, 07:05:01-07:07, 07:10, 07:12-07:15, 07:18:01- 07:18:02, 07:21-07:23, 07:26, 07:29-07:31, 07:34- 07:35, 07:39, 07:41-07:46, 07:49N, 07:52, 07:55- 07:59, 07:61-07:64, 07:66- 07:68:03, 07:70-07:72, 07:74-07:77, 07:79-07:80, 07:82-07:83, 07:86-07:99, 07:101-07:121, 07:124- 07:130, 07:132, 07:135N- 07:138, 07:140-07:143	B7, B42, Null, -	1																							24
*07:03, 07:08, 07:16, 07:32, 07:37	B703, B7, -	1																							24
*07:04, 07:25	B7	1													14										24
*07:09, 07:11, 07:17	B7	1										12													24
*07:19, 07:33, 07:40, 07:48, 07:51, 07:53, 07:122, 07:139	B7, 42, -																								24
*07:20, 07:24, 07:131	B7, -	1	3																						24
*07:27	-	1									11									20					
*07:28	B7	1				6																			24
*07:36, 07:81	B7, -	1																							
*07:38	-	1																							
*07:47	B7	1																							24
*07:50	-	1									11														24
*07:54, 07:123	B7, -	1																							24
*07:60	-		3																						24
*07:65	-	1																	18						
*07:69, 07:85	-	1																					23	24	
*07:73	-	1														16									24
*07:78	-	1									10														24
*07:84	-	1																							24
*07:100	-		3																						24
*07:133	-	1																							24
*07:134	-	1																	18						24
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24

Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
*15:43	B15							7					12							19					
*15:45	B62			3				7												19					
*15:46	-				4			7	8				12							19					
*15:48, 15:150, 15:188	B62, -							7												19					
*15:53, 15:212	-				4				8				12							19					
*15:54	B72								8			11	12							19					
*15:55	B15			3				7					12							19					
*15:57	B62							7					12							19					
*15:62	-								8			11	12					17		19					
*15:64	B71								8			11	12							19					
*15:66	B62							7					12		14					19					
*15:68	-			3								11	12												
*15:69	B72								8			11	12							19		21			
*15:73	B62							7												19					
*15:74	B72								8				12							19					
*15:76	-							7					12							19					
*15:77	B62						6	7										17		19					
*15:80	B70								8				12		14			17		19					
*15:83	B62							7	8											19					
*15:85, 15:154, 15:165	B62, -							7					12					17		19					
*15:86	-							7												19					
*15:89	B77			w				7					12					17		19					
*15:91, 15:131	-								8			11	12	13						19					
*15:101	-							7					12							19					
*15:106	-				4			7	8				12					17		19					
*15:116	B62							7					12							19					
*15:124	B71								8				12		14					19					
*15:137	-							7												19					
*15:138	-	1						7					12							19					
*15:142	-		2					7					12							19					
*15:143	-							7	8				12							19					
*15:144	-							7					12					17		19					
*15:161	-								8				12	13	14					19					
*15:180	-		2										12												
*15:186	-								8				12		14					19		21			
*15:189	-						6	7					12		14					19					
*15:202	-							7					12					18		19					
*15:222	-									9			12		14					19					
*15:224	-							7												19					
*15:230	-	1						7		9					14					19					
*15:233	-						6	7												19					
*15:236	-							7					12												
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24

Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	
*40:21, 40:137	B40, B15, -																								?	
*40:23	B40	1			4																					24
*40:25, 40:69, 40:106	B40, B60, -				4																					?
*40:26	B21									9																?
*40:28	-									9								17								?
*40:31, 40:45, 40:80	B60, -				4																					?
*40:38, 40:52	B60				4								12													?
*40:40	B61																									
*40:44, 40:75, 40:86, 40:103, 40:127, 40:131, 40:162, 40:165, 40:167	B61, -																									24
*40:46	B60											11					16									?
*40:47	B40				4																					
*40:48	B60				4							10														?
*40:53, 40:72:02	B61, B60				4																					24
*40:58	B40				4																					?
*40:59-40:60	B40, -				4								12													24
*40:68	-																									?
*40:71	B61			3																						?
*40:74, 40:149	B60, -																									?
*40:76	B60				4																					
*40:83, 40:177	-																									24
*40:93	-											11					16									24
*40:95	-																									?
*40:96	-																									24
*40:98	-	1																								?
*40:105	-												12													?
*40:109	-																									24
*40:110	-				4																					24
*40:117	-				4																					
*40:132	-				4																					?
*40:135	-				4																					
*40:136	-	1			4				8						14											?
*40:138	-																									24
*40:148, 40:161	-																									24
*40:157	-																									
*40:158	-	1			4								12													?
*40:159	-												12													24
*40:166	-				4																					?
*40:170	-																									?
*40:174	-			3																						?
*41:01, 41:06-41:07, 41:12, 41:14, 41:16-41:17	B41, -				4				8																	24
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	
Well No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	

25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	ser. ⁵	Well No. HLA-B allele ⁴
			28																				47	B7, B42, Null, -	*07:02:01-07:02:29, 07:05:01-07:07, 07:10, 07:12-07:15, 07:18:01-07:18:02, 07:21-07:23, 07:26, 07:29-07:31, 07:34-07:35, 07:39, 07:41-07:46, 07:49N, 07:52, 07:55-07:59, 07:61-07:64, 07:66-07:68:03, 07:70-07:72, 07:74-07:77, 07:79-07:80, 07:82-07:83, 07:86-07:99, 07:101-07:121, 07:124-07:130, 07:132, 07:135N-07:138, 07:140-07:143
																							47	B703, B7, -	*07:03, 07:08, 07:16, 07:32, 07:37
		27	28																				47	B7	*07:04, 07:25
			28																				47	B7	*07:09, 07:11, 07:17
																							47	B7, 42, -	*07:19, 07:33, 07:40, 07:48, 07:51, 07:53, 07:122, 07:139
			28																				47	B7, -	*07:20, 07:24, 07:131
																						46	-	*07:27	
			28																				47	B7	*07:28
			28																			46	B7, -	*07:36, 07:81	
																						46	-	*07:38	
			w																				47	B7	*07:47
																							47	-	*07:50
			28										38										47	B7, -	*07:54, 07:123
			28																				47	-	*07:60
			28				32								40				w				47	-	*07:65
																							47	-	*07:69, 07:85
			28																				47	-	*07:73
			28							35													47	-	*07:78
			28							35													47	-	*07:84
			28																		45	47	-	*07:100	
			28									37											47	-	*07:133
			28																w				47	-	*07:134
25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	Well No.	
25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	Well No.	

25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48		Well No.
													38							45	46			B15	*15:43
													38							45		47		B62	*15:45
												37	38							45		47		-	*15:46
													38										47	B62, -	*15:48, 15:150, 15:188
												37	38							45		47		-	*15:53, 15:212
													38							45		47		B72	*15:54
								33															47	B15	*15:55
					w								38							45		47		B62	*15:57
													38										47	-	*15:62
																							47	B71	*15:64
													38							45		47		B62	*15:66
													38							45		47		-	*15:68
													38										47	B72	*15:69
													38		41	43				45		47		B62	*15:73
													38										47	B72	*15:74
			28					w												45		47		-	*15:76
													38							45		47		B62	*15:77
																							47	B70	*15:80
					30								38										47	B62	*15:83
													38							45		47		B62, -	*15:85, 15:154, 15:165
													38		40								47	-	*15:86
								33													46			B77	*15:89
													38										47	-	*15:91, 15:131
			28										38							45		47		-	*15:101
												37	38										47	-	*15:106
25													38										47	B62	*15:116
25																							47	B71	*15:124
													38							45		47		-	*15:137
													38										47	-	*15:138
													38							45		47		-	*15:142
												37								45		47		-	*15:143
								33												45		47		-	*15:144
																							47	-	*15:161
																				45		47		-	*15:180
																							47	-	*15:186
								33												45		47		-	*15:189
													38							45		47		-	*15:202
																					46			-	*15:222
													38		40			43					47	-	*15:224
																					46			-	*15:230
													38							45		47		-	*15:233
													38							45		47		-	*15:236
25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48		Well No.
25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48		Well No.

25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48		Well No.
25												37	38									47		B40, B15, -	*40:21, 40:137
25												37	38									47		B40	*40:23
25	?											37										47		B40, B60, -	*40:25, 40:69, 40:106
	?											37	38									47		B21	*40:26
	?				31							37	38									47		-	*40:28
	?							33				37	38									47		B60, -	*40:31, 40:45, 40:80
25	?											37	38									47		B60	*40:38, 40:52
												37	38									47		B61	*40:40
	?											37	38									47		B61, -	*40:44, 40:75, 40:86, 40:103, 40:127, 40:131, 40:162, 40:165, 40:167
25												38										47		B60	*40:46
25	?											38									46			B40	*40:47
25	?											37	38				43					47		B60	*40:48
	?											37	38									47		B61, B60	*40:53, 40:72:02
	?											37	38		40							47		B40	*40:58
25												37	38									47		B40, -	*40:59-40:60
	?											37										47		-	*40:68
	?											37	38			41		43				47		B61	*40:71
25	?											37	38									47		B60, -	*40:74, 40:149
25	?											38												B60	*40:76
												37	38									47		-	*40:83, 40:177
												38										47		-	*40:93
	?									35		37	38									47		-	*40:95
	?											38										46		-	*40:96
	?											37	38									47		-	*40:98
	?											37	38									47		-	*40:105
	?								34		36	38										46		-	*40:109
	?											38										46		-	*40:110
25	?							34		36	38											46		-	*40:117
25	?			29								37	38									47		-	*40:132
25	?											37	38									47		-	*40:135
	?	27										37	38									47		-	*40:136
25												37	38									47		-	*40:138
	?									35		37	38									47		-	*40:148, 40:161
	?											38										46		-	*40:157
	?											37	38									47		-	*40:158
	?											37	38									47		-	*40:159
25	?		28									37										47		-	*40:166
25												37	38									47		-	*40:170
	?											37	38									47		-	*40:174
	27			30								37	38									47		B41, -	*41:01, 41:06-41:07, 41:12, 41:14, 41:16-41:17
25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48		Well No.
25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48		Well No.

Negative Control

				30						35	36							44	46			-	*51:68	
				30							36								44	46			-	*51:81
										35	36								44	46			-	*51:93
				30						35	36								44	46			-	*51:97
				30						35	36								44	46			-	*51:101
				30						35	36										46		Null, -	*51:104, 51:118N
				30						35	36				41				44	46			-	*51:106
				30				34	35	36											46		-	*51:112
										35	36		38									46	B52, -	*52:01:01-52:01:01:02, 52:01:03, 52:01:06- 52:01:08, 52:04-52:05, 52:07-52:08, 52:10:01- 52:13, 52:15, 52:17-52:18, 52:22-52:24, 52:27
				30						35	36		38									46	B52, -	*52:01:02, 52:01:04, 52:01:09, 52:02-52:03, 52:06:01-52:06:02, 52:09
										35	36		38									46	B52, -	*52:01:05, 52:14
										35	36		38									46	-	*52:01:10
										35	36		38									46	-	*52:16
				30						35	36		38									46	-	*52:19
										35			38									46	-	*52:20
				30						35	36		38								45	46	-	*52:21
										35	36		38		41	43						46	-	*52:25
										35	36		38									46	-	*52:26
											36								44	46			B53, -	*53:01:01-53:01:05, 53:05, 53:08:01-53:08:02, 53:10, 53:15-53:16, 53:18, 53:20-53:21, 53:23-53:27
				31						36									44	46			B53	*53:02
																			44	46			B53, B35	*53:03, 53:09, 53:11- 53:13
										36									44	46			B53, -	*53:04, 53:07, 53:19
				30	31					36									44	46			B53, B51	*53:06
										36									44	46			-	*53:14
										36		38										46	-	*53:17:01-53:17:02
				29						36									44	46			-	*53:22
			28			32							39	40								47	B54, Null, -	*54:01:01, 54:02, 54:05N, 54:07-54:08N, 54:10, 54:13, 54:16-54:19, 54:21- 54:24
25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	Well No.
25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	Well No.

