

INTERPRETATION TABLE

HLA-B*14 SSP subtyping

Amplification patterns of the B*14:01 to B*14:19 alleles

	Well ⁵							
	1	2	3	4	5	6	7	8
Length of spec.	265	185	210	85	135	95	260	250
PCR product(s)					220	175		
Length of int.	800	1070	800	1070	1070	1070	1070	1070
pos. control ¹								
5'-primer(s) ²	361	103	103	539	91	113	363	103
	5'-AgT 3'	5'-CCT 3'	5'-CCg 3'	5'-gCg 3'	5'-ggC 3'	5'-CCA 3'	5'-Agg 3'	5'-CCT 3'
		103			564	363		103
		5'-CCT 3'			5'-TgA 3'	5'-AgC 3'		5'-CCT 3'
3'-primer(s) ³	583	246	272	583	269	246	583	310
	5'-gTg 3'	5'-TAT 3'	5'-TgC 3'	5'-gTg 3'	5'-AgT 3'	5'-TAT 3'	5'-gTg 3'	5'-gTT 3'
		246			3 rd I	419		
		5'-TAT 3'			5'-TAT 3'	5'-CgA 3'		
Well No.	1	2	3	4	5	6	7	8
HLA-B allele								
*14:01:01-14:01:02	1	2						8
*14:02:01-14:02:06	1		3					
*14:03	1		3	4				
*14:04	1		w		5			
*14:05			3			6		
*14:06:01-14:06:02			3				7	
*14:07N	1	2			5			8
*14:08, 39:43 ⁴		2					7	8
*14:09	1		3					
*14:10		2						8
*14:11	1		3					
*14:12	1	2						8
*14:13			3					
*14:14	1	2						8
*14:15	1		3					
*14:16	1		3					
*14:17	1		3			6		
*14:18	1		3					
*14:19	1	2						8
Well No.	1	2	3	4	5	6	7	8
Length of spec.	265	185	210	85	135	95	260	250
PCR product(s)					220	175		

Well No.	1	2	3	4	5	6	7	8
*07:02:01-07:02:15, 07:02:17-07:05:04, 07:05:06-07:15, 07:17-07:18:02, 07:20- 07:26, 07:28, 07:30, 07:32-07:33, 07:35, 07:39-07:42, 07:44-07:49N, 07:51-07:64, 07:66-07:69, 07:73-07:74, 07:76-07:78, 07:80, 07:82-07:85, 07:87-07:90, 07:92- 07:110, 07:112-07:114, 07:116-07:117, 15:71, 15:175, 27:08, 27:12, 27:18, 27:26, 27:33, 27:40, 27:42, 27:44, 35:50, 35:84, 40:02:01-40:06:03, 40:08-40:09, 40:11:01- 40:11:02, 40:18, 40:20, 40:24, 40:26-40:27, 40:29, 40:35, 40:37, 40:40, 40:44, 40:50, 40:56, 40:64, 40:68, 40:70-40:71, 40:74- 40:75, 40:78, 40:82-40:83, 40:85-40:86, 40:89-40:91, 40:93-40:95, 40:97-40:99, 40:103-40:105, 40:107, 40:111, 40:115, 40:119-40:120, 40:122, 40:127, 40:131, 40:133Q, 40:142N-40:145, 40:148, 48:02:02, 73:01-73:02, 81:01-81:05								8
*07:16, 07:37, 07:50, 07:75, 07:91, 15:68, 18:10, 18:21, 18:35, 40:149, 48:01:01- 48:02:01, 48:03:01-48:06, 48:08-48:17, 48:19-48:21, 48:23-48:24								8
*07:19, 07:31, 07:34, 07:43, 39:01:09, 39:02:02, 39:10, 39:16-39:17, 39:19:01, 39:23, 39:45, 40:39, 42:01:01-42:02, 42:04- 42:05:02, 42:07-42:10, 42:12-42:14, 54:02, 55:10, 67:01:02-67:02								8
*07:27, 07:38, 15:03:01-15:03:03, 15:29, 15:47, 15:49, 15:54, 15:62, 15:64, 15:74, 15:91, 15:98, 15:103, 15:123, 15:127, 15:131- 15:132, 15:151, 15:156, 15:158, 15:173, 35:25, 35:124-35:125, 35:142, 37:01:05, 37:02, 37:07, 40:12, 48:18, 48:22, 58:12								
*07:115						6		8
Well No.	1	2	3	4	5	6	7	8

Length of spec.	265	185	210	85	135	95	260	250
PCR product(s)					220	175		
Well No.	1	2	3	4	5	6	7	8
*08:01:01-08:02, 08:04-08:05, 08:07-08:10, 08:12:01-08:20, 08:22, 08:24, 08:26-08:39, 08:41-08:53, 08:55-08:66, 13:06, 15:42, 15:50, 15:83, 15:86, 15:121, 15:188, 15:199, 27:15, 27:28, 27:62, 27:71, 35:35, 35:60, 35:87, 35:93, 37:22, 39:12, 41:01-41:16, 44:04, 44:56N, 45:08, 46:03, 46:18, 46:21:01-46:21:02, 51:36, 54:01, 54:04-54:05N, 54:07-54:23, 55:01:01-55:05, 55:07, 55:11-55:21, 55:23, 55:25-55:48, 56:10, 56:12, 56:18, 56:23, 57:26, 58:05, 58:20, 59:01-59:03, 59:05								
*08:03, 08:23, 08:54, 15:69, 37:01:01-37:01:04, 37:01:06-37:01:07, 37:03N, 37:06, 37:08-37:10, 37:12-37:21, 37:23-37:25, 38:02:03, 38:04, 38:06, 38:25								
*13:15, 51:01:01-51:01:18, 51:01:20, 51:03, 51:07:01-51:09:02, 51:11N-51:14, 51:16, 51:19-51:20, 51:24:01-51:24:04, 51:26-51:28, 51:30-51:33, 51:35, 51:38-51:39, 51:41N, 51:43-51:44N, 51:48-51:53, 51:55, 51:57-51:58, 51:60-51:61, 51:63-51:80, 51:83-51:84, 51:86-51:89, 51:91, 51:95-51:100, 51:102-51:106, 52:01:01-52:02, 52:04-52:10, 52:13-52:15, 52:17-52:22, 56:05:01-56:06, 56:21, 58:08, 78:01-78:02:02, 78:05-78:06								
*15:09-15:10:02, 15:18:01-15:18:04, 15:23, 15:51-15:52, 15:72, 15:80, 15:90, 15:99, 15:108, 15:114-15:115, 15:119, 15:124, 15:133-15:134, 15:153, 15:161, 15:176, 15:197-15:198, 15:200			3					
*15:21, 15:66, 35:82, 35:85, 35:135			3					
*15:37			3				7	
*15:38:01-15:38:02, 15:185, 51:45							7	
*15:44			3					
*15:93, 15:186, 38:05			3					
*15:189, 35:26			3					
Well No.	1	2	3	4	5	6	7	8
Length of spec.	265	185	210	85	135	95	260	250
PCR product(s)					220	175		

Well No.	1	2	3	4	5	6	7	8
*18:01:01, 18:01:03-18:01:09, 18:05-18:08, 18:12, 18:15, 18:17N-18:20, 18:23N-18:24, 18:28, 18:31-18:34, 18:37-18:38, 18:40-18:43, 18:45-18:49, 18:51-18:52							7	8
*18:01:02, 18:03, 18:29							7	8
*18:02, 18:11, 18:14, 18:27, 18:39, 37:05, 37:11, 39:13:01-39:13:02, 39:20, 39:49, 42:11								8
*18:04							7	
*18:09, 18:25							7	
*18:13, 18:30				4			7	8
*18:22								8
*18:26							7	8
*18:36, 18:50							7	8
*18:44	1							8
*35:21, 35:24:01-35:24:02, 35:81, 35:109, 53:02, 53:06, 57:14, 58:09, C*15:39								
*35:96								
*37:04:01-37:04:02								
*38:01:01-38:02:02, 38:07, 38:09-38:12, 38:15, 38:17-38:18, 38:20-38:21, 38:23-38:24		2						
*38:03, 38:13-38:14, 38:16, 38:26		2						
*38:08, 38:22		2						
*38:19		2				6		
*39:01:01:01-39:01:01:02L, 39:01:03-39:01:08, 39:01:10-39:02:01, 39:06:01-39:06:02, 39:09, 39:14-39:15, 39:18, 39:19:02, 39:22, 39:25N-39:31, 39:34-39:36, 39:38Q-39:39, 39:41-39:42, 39:44, 39:46, 39:51-39:54, 39:57-39:62, 67:01:01, 67:03		2						8
*39:03, 39:24		2				6		8
*39:04			3					
*39:05:01-39:05:02, 39:07-39:08, 39:11, 39:55-39:56		2						8
*39:32		2					7	8
*39:33, 39:47, 39:50		2						8
*39:37		2				6		8
Well No.	1	2	3	4	5	6	7	8
Length of spec.	265	185	210	85	135	95	260	250
PCR product(s)					220	175		
Well No.	1	2	3	4	5	6	7	8

*39:40N		2			5			8
*39:48		2					7	8
*40:28								8
*40:51							7	
*40:113, 51:06, 51:62							7	
*42:06						6		8
*48:07						6		8
*51:04, 51:46, 51:56								
*51:21, 51:101, 78:07								
*51:22, 78:03			3					
*51:29, 51:82				4				
*52:16								
*57:11, 58:02, 58:06-58:07, 58:25, C*01:35, C*03:17, C*03:71, C*14:20								
C*01:30	1			4				
C*07:102								w
C*08:08						6		
HLA-B allele								
Well No.	1	2	3	4	5	6	7	8

INTERPRETATION TABLE

HLA-B*14 SSP subtyping

Amplification patterns of the B*14:01 to B*14:19 alleles

Well ⁵								
9	10	11	12	13	14	15	16	
160	180	135	230	210	190	100	175	Length of spec.
	270							PCR product(s)
1070	1070	1070	1070	1070	1070	800	1070	Length of int.
								pos. control ¹
361	106	560	361	412	142	103	361	5'-primer(s) ²
5' -AgT 3'	5' -CCA 3'	5' -CAC 3'	5' -AgT 3'	5' -ATA 3'	5' -TCT 3'	5' -CCg 3'	5' -AgT 3'	
	355							
	5' -TCA 3'							
479	246	3 rd I	541	583	292	164	494	3'-primer(s) ³
5' -CCC 3'	5' -TAT 3'	5' -TAT 3'	5' -TCC 3'	5' -gTg 3'	5' -gTA 3'	5' -gCA 3'	5' -TCC 3'	
	583		559		292			
	5' -gTg 3'		5' -CAg 3'		5' -gTA 3'			
					292			
					5' -gTA 3'			
9	10	11	12	13	14	15	16	Well No.
								HLA-B allele
		11						*14:01:01-14:01:02
		11						*14:02:01-14:02:06
		11						*14:03
		11						*14:04
		11						*14:05
		11						*14:06:01-14:06:02
		11						*14:07N
		11						*14:08, 39:43 ⁴
9		11						*14:09
	10	11						*14:10
	10	11						*14:11
		11	12					*14:12
		11		13				*14:13
		11			14			*14:14
		11				15		*14:15
		11					16	*14:16
		11						*14:17
		11			14			*14:18
			12					*14:19
9	10	11	12	13	14	15	16	Well No.
160	180	135	230	210	190	100	175	Length of spec.
	270							PCR product(s)

9	10	11	12	13	14	15	16	Well No.
								*07:02:01-07:02:15, 07:02:17-07:05:04, 07:05:06-07:15, 07:17-07:18:02, 07:20-07:26, 07:28, 07:30, 07:32-07:33, 07:35, 07:39-07:42, 07:44-07:49N, 07:51-07:64, 07:66-07:69, 07:73-07:74, 07:76-07:78, 07:80, 07:82-07:85, 07:87-07:90, 07:92-07:110, 07:112-07:114, 07:116-07:117, 15:71, 15:175, 27:08, 27:12, 27:18, 27:26, 27:33, 27:40, 27:42, 27:44, 35:50, 35:84, 40:02:01-40:06:03, 40:08-40:09, 40:11:01-40:11:02, 40:18, 40:20, 40:24, 40:26-40:27, 40:29, 40:35, 40:37, 40:40, 40:44, 40:50, 40:56, 40:64, 40:68, 40:70-40:71, 40:74-40:75, 40:78, 40:82-40:83, 40:85-40:86, 40:89-40:91, 40:93-40:95, 40:97-40:99, 40:103-40:105, 40:107, 40:111, 40:115, 40:119-40:120, 40:122, 40:127, 40:131, 40:133Q, 40:142N-40:145, 40:148, 48:02:02, 73:01-73:02, 81:01-81:05
					14			*07:16, 07:37, 07:50, 07:75, 07:91, 15:68, 18:10, 18:21, 18:35, 40:149, 48:01:01-48:02:01, 48:03:01-48:06, 48:08-48:17, 48:19-48:21, 48:23-48:24
		11						*07:19, 07:31, 07:34, 07:43, 39:01:09, 39:02:02, 39:10, 39:16-39:17, 39:19:01, 39:23, 39:45, 40:39, 42:01:01-42:02, 42:04-42:05:02, 42:07-42:10, 42:12-42:14, 54:02, 55:10, 67:01:02-67:02
					14			*07:27, 07:38, 15:03:01-15:03:03, 15:29, 15:47, 15:49, 15:54, 15:62, 15:64, 15:74, 15:91, 15:98, 15:103, 15:123, 15:127, 15:131-15:132, 15:151, 15:156, 15:158, 15:173, 35:25, 35:124-35:125, 35:142, 37:01:05, 37:02, 37:07, 40:12, 48:18, 48:22, 58:12
								*07:115
9	10	11	12	13	14	15	16	Well No.

160	180	135	230	210	190	100	175	Length of spec.
	270							PCR product(s)
9	10	11	12	13	14	15	16	Well No.
		11						*08:01:01-08:02, 08:04-08:05, 08:07-08:10, 08:12:01-08:20, 08:22, 08:24, 08:26-08:39, 08:41-08:53, 08:55-08:66, 13:06, 15:42, 15:50, 15:83, 15:86, 15:121, 15:188, 15:199, 27:15, 27:28, 27:62, 27:71, 35:35, 35:60, 35:87, 35:93, 37:22, 39:12, 41:01-41:16, 44:04, 44:56N, 45:08, 46:03, 46:18, 46:21:01-46:21:02, 51:36, 54:01, 54:04-54:05N, 54:07-54:23, 55:01:01-55:05, 55:07, 55:11-55:21, 55:23, 55:25-55:48, 56:10, 56:12, 56:18, 56:23, 57:26, 58:05, 58:20, 59:01-59:03, 59:05
		11			14			*08:03, 08:23, 08:54, 15:69, 37:01:01-37:01:04, 37:01:06-37:01:07, 37:03N, 37:06, 37:08-37:10, 37:12-37:21, 37:23-37:25, 38:02:03, 38:04, 38:06, 38:25
				13				*13:15, 51:01:01-51:01:18, 51:01:20, 51:03, 51:07:01-51:09:02, 51:11N-51:14, 51:16, 51:19-51:20, 51:24:01-51:24:04, 51:26-51:28, 51:30-51:33, 51:35, 51:38-51:39, 51:41N, 51:43-51:44N, 51:48-51:53, 51:55, 51:57-51:58, 51:60-51:61, 51:63-51:80, 51:83-51:84, 51:86-51:89, 51:91, 51:95-51:100, 51:102-51:106, 52:01:01-52:02, 52:04-52:10, 52:13-52:15, 52:17-52:22, 56:05:01-56:06, 56:21, 58:08, 78:01-78:02:02, 78:05-78:06
					14			*15:09-15:10:02, 15:18:01-15:18:04, 15:23, 15:51-15:52, 15:72, 15:80, 15:90, 15:99, 15:108, 15:114-15:115, 15:119, 15:124, 15:133-15:134, 15:153, 15:161, 15:176, 15:197-15:198, 15:200
					14			*15:21, 15:66, 35:82, 35:85, 35:135
					14			*15:37
		11						*15:38:01-15:38:02, 15:185, 51:45
		11			14			*15:44
	10							*15:93, 15:186, 38:05
	10							*15:189, 35:26
9	10	11	12	13	14	15	16	Well No.
160	180	135	230	210	190	100	175	Length of spec.
	270							PCR product(s)

9	10	11	12	13	14	15	16	Well No.
		11			14			*18:01:01, 18:01:03-18:01:09, 18:05-18:08, 18:12, 18:15, 18:17N-18:20, 18:23N-18:24, 18:28, 18:31-18:34, 18:37-18:38, 18:40-18:43, 18:45-18:49, 18:51-18:52
		11						*18:01:02, 18:03, 18:29
		11			14			*18:02, 18:11, 18:14, 18:27, 18:39, 37:05, 37:11, 39:13:01-39:13:02, 39:20, 39:49, 42:11
		11						*18:04
		11			14			*18:09, 18:25
		11			14			*18:13, 18:30
	10	11			14			*18:22
		11		13	14			*18:26
					14			*18:36, 18:50
		11			14			*18:44
	10							*35:21, 35:24:01-35:24:02, 35:81, 35:109, 53:02, 53:06, 57:14, 58:09, C*15:39
	10	11						*35:96
	10	11			14			*37:04:01-37:04:02
		11			14			*38:01:01-38:02:02, 38:07, 38:09-38:12, 38:15, 38:17-38:18, 38:20-38:21, 38:23-38:24
		11						*38:03, 38:13-38:14, 38:16, 38:26
					14			*38:08, 38:22
		11			14			*38:19
		11						*39:01:01:01-39:01:01:02L, 39:01:03-39:01:08, 39:01:10-39:02:01, 39:06:01-39:06:02, 39:09, 39:14-39:15, 39:18, 39:19:02, 39:22, 39:25N-39:31, 39:34-39:36, 39:38Q-39:39, 39:41-39:42, 39:44, 39:46, 39:51-39:54, 39:57-39:62, 67:01:01, 67:03
		11						*39:03, 39:24
	10	11						*39:04
		11			14			*39:05:01-39:05:02, 39:07-39:08, 39:11, 39:55-39:56
				13				*39:32
								*39:33, 39:47, 39:50
		11			14			*39:37
9	10	11	12	13	14	15	16	Well No.
160	180	135	230	210	190	100	175	Length of spec.
	270							PCR product(s)
9	10	11	12	13	14	15	16	Well No.

		11							*39:40N
		11		13					*39:48
	10			13					*40:28
		11		13					*40:51
				13					*40:113, 51:06, 51:62
		11							*42:06
					14				*48:07
	10			13					*51:04, 51:46, 51:56
		11		13					*51:21, 51:101, 78:07
				13					*51:22, 78:03
				13					*51:29, 51:82
				13	14				*52:16
			12						*57:11, 58:02, 58:06-58:07, 58:25, C*01:35, C*03:17, C*03:71, C*14:20
									C*01:30
									C*07:102
									C*08:08
									HLA-B allele
9	10	11	12	13	14	15	16		Well No.