

## INTERPRETATION TABLE

### HLA-B\*51 SSP subtyping

Amplification patterns of the B\*51:01 to B\*51:95 alleles

	Well <sup>3</sup>																																
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	
Length of spec.	105	95	265	255	270	210	270	105	85	160	245	495	200	105	100	180	105	140	145	110	130	160	210	210	130	105	205	205	205	245	105	105	165
PCR product(s)									195	245					230			220			185				175	155	270	255			225		
Length of int. pos. control <sup>1</sup>	800	1070	1070	1070	1070	800	1070	800	1070	1070	1070	1070	1070	1070	800	1070	800	800	800	1070	1070	1070	1070	800	1070	800	1070	1070	1070	800	1070	1070	1070
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	
HLA-B allele <sup>2</sup>																																	
*51:01:01, 51:01:03-51:01:07, 51:01:10-51:01:14, 51:01:16-51:01:17	1	2																				22		24						29			
*51:01:02	1	2																				22		24						29			
*51:01:08, 51:01:15	1	2																				22		24						29			
*51:01:09	1	2								10																				29			
*51:02:01, 51:02:03-51:02:04	1		3																			22		24						29			
*51:02:02	1		3																			22		24						29			
*51:03	1	2		4																		22		24						29			
*51:04	1	2			5																	22		24			27						
*51:05	1		3			w																							29		31		
*51:06	1	2					7															22		24						30			
*51:07:01		2						8														22		24						29			
*51:07:02		2						8																						29			
*51:08	1								9								18					22		24						29			
*51:09:01	1					6																22		24						29			
*51:09:02	1					6																								29			
*51:10	1		3			6				10	11									19		22		24									
*51:11N	1	2										12										22		24						29			
*51:12	1	2											13									22		24						29			
*51:13:01-51:13:02	1	2												14								22		24		26				29		31	
*51:14	1	2													15							22		24						29			
*51:15	1		3			6										16														29			
*51:16	1	2									11											22		24								32	
*51:17	1	2							9													22		24						29			
*51:18	1	2														16						22		24						29			
*51:19	1					6									15							22		24						29			
*51:20	1								9								18					22		24						29			
*51:21	1	2								10												22		24									
*51:22		2															17					22		24						29			
*51:23	1		3																			22		24				28		29			
*51:24:01	1	2																	19			22		24						29			
*51:24:02-51:24:03	1	2								10									19			22		24						29			
*51:26	1	2																		20		22		24						29			
*51:27N	1	2																			21	22		24					29				
*51:28	1	2															18					22		24						29			
*51:29	1	2																				22		24						29		31	
*51:30	1	2																				22		24	25					29			
*51:31	1					6					11											22		24								32	
*51:32	1	2																				22		24		26				29			
*51:33	1	2																				22		24			27	28		29			
*51:34	1		3								11											22		24								32	
*51:35	1	2							9							16						22		24						29			
*51:36	1		3							10								18				22		24									
*51:37	1	2																				22		24		26				29		31	
*51:38	1	2																				21	22		24	25				29			
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	
*51:39	1	2																				22		24			27			29			
*51:40	1		3			6																22		24						29			

*51:41N	1	2																		22	24					29						
*51:42	1																			22	24											
*51:43	1	2																		22	w					29						
*51:44N	1							9								18				22	24					29						
*51:45	1	2					7													22	24											
*51:46	1	2			5															22	24			27								
*51:48	1	2																		22	24					29						
*51:49	1	2																		22	24					29						
*51:50	1	2																		22	24					29						
*51:51, 51:76 <sup>4</sup>	1	2																		22	24					29						
*51:52	1	2																		22	24					29						
*51:53	1																			22	24	25				29						
*51:54	1		3				w																			29		31				
*51:55	1	2																		22	24					29						
*51:56	1	2			5															22	24			27		30						
*51:57	1	2																		22	24					29						
*51:58	1	2																		22	24					29						
*51:59	1																			22	24					30						
*51:60	1	2																		22	24					29						
*51:61	1	2																		22	24					29		31				
*51:62	1	2					7								16					22	24											
*51:63	1	2																		22	24		26			29		31				
*51:64	1	2					7													22	24					30						
*51:65	1	2																		22	24					29						
*51:66-51:67 <sup>5</sup>	1	2																		22	24					29						
*51:68, 51:72 <sup>6</sup>	1	2																		22	24					29						
*51:69	1	2																		22	24					29						
*51:70-51:71 <sup>7</sup>	1	2																		22	24					29						
*51:73	1						6													22	24					29						
*51:74	1	2																		22	24					29						
*51:75, 51:77 <sup>8</sup>	1	2																		22	24					29						
*51:78	1	2																		22						29						
*51:79	1	2																		22	24					29						
*51:80	1	2																		22	24					29						
*51:81	1																			22	24											
*51:82	1	2							11											22	24						31	32				
*51:83, 51:88 <sup>9</sup>	1	2																		22	24					29						
*51:84, 51:94 <sup>10</sup>	1	2																		22	24					29						
*51:85	1	2																		22	24					29						
*51:86, 51:95 <sup>11</sup>	1	2																		22	24					29						
*51:87	1	2																		22	24					29						
*51:89	1	2																								29						
*51:90	1	2																		22	24					29						
*51:91	1	2																		22	24					29						
*51:92	1	2												14						22	24		26			29		31				
*51:93	1		3			6				11																		32				
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32









*57:14, 58:09					5																									27				
*57:24																																		
*59:02			3		6				10										19															
*73:01					6				10 11				14					19									26					31		
*78:01	1	2																					22 23									29		
*78:02:01	1	2																															29	
*78:02:02	1	2																					22 23									29		
*78:03		2																	17				22 23									29		
*78:04	1		3		6																												29	
*78:07	1	2							10														22 23											
C*08:24																																		
HLA-B allele <sup>4</sup>																																		
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32		















