

INTERPRETATION TABLE

HLA-C low resolution SSP typing

Amplification patterns of the C*01:02 to C*18:05 alleles

	Well ¹⁰															
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
Length of spec.	90	130	280	170	280	130	390	130	245	115	340	100	120	160	130	255
PCR product(s)	150	200		275		330	445	355	425	165		160	250	215	255	445
		270								265		220				
		300								390						
Length of int.	800	800	800	1070	800	800	1070	800	800	800	1070	800	800	800	800	1070
pos. control ¹																
5'-primer(s) ²	89	47	105	355	105	47	1 st	28	47	1 st	1 st	361	201	97	98	98
	5'-gAA ^{3'}	5'-Agg ^{3'}	5'-gCT ^{3'}	5'-TCA ^{3'}	5'-gCT ^{3'}	5'-Agg ^{3'}	5'-CgA ^{3'}	5'-TCA ^{3'}	5'-Agg ^{3'}	5'-CgA ^{3'}	5'-CgA ^{3'}	5'-AgT ^{3'}	5'-CCA ^{3'}	5'-TCg ^{3'}	5'-CTA ^{3'}	5'-CTC ^{3'}
	98	89		459		112	527	2 nd	648	176		419	218	361	368	527
	5'-CTT ^{3'}	5'-gAA ^{3'}		5'-gAT ^{3'}		5'-CCT ^{3'}	5'-TgA ^{3'}	5'-CCA ^{3'}	5'-CAC ^{3'}	5'-gCA ^{3'}		5'-gTC ^{3'}	5'-gga ^{3'}	5'-AgT ^{3'}	5'-gTg ^{3'}	5'-TgA ^{3'}
		2 nd				118				527			2 nd			
		5'-CCA ^{3'}				5'-CCA ^{3'}				5'-TAC ^{3'}			5'-CCA ^{3'}			
		703														
		5'-CTA ^{3'}														
3'-primer(s) ³	142	176	343	589	343	201	302	213	302	175	302	474	289	218	312	311
	5'-TgA ^{3'}	5'-ACT ^{3'}	5'-C ^{3'}	5'-CTT ^{3'}	5'-T ^{3'}	5'-CTT ^{3'}	5'-ggT ^{3'}	5'-Cgg ^{3'}	5'-ggC ^{3'}	5'-CCg ^{3'}	5'-ggT ^{3'}	5'-gCA ^{3'}	5'-AgC ^{3'}	5'-gCT ^{3'}	5'-AgT ^{3'}	5'-ggT ^{3'}
	201	559	343			218	3 rd	420	853	302	304	477	289	527	459	5 th
	5'-CTT ^{3'}	5'-CTC ^{3'}	5'-g ^{3'}			5'-gCT ^{3'}	5'-gCA ^{3'}	5'-gCT ^{3'}	5'-CAT ^{3'}	5'-ggC ^{3'}	5'-CAA ^{3'}	5'-gCA ^{3'}	5'-AgC ^{3'}	5'-CCg ^{3'}	5'-AgA ^{3'}	5'-gCA ^{3'}
	201	861								601		538	539	538		
	5'-CTC ^{3'}	5'-TCg ^{3'}								5'-CTT ^{3'}		5'-gCA ^{3'}	5'-TCT ^{3'}	5'-CCA ^{3'}		
														538		
														5'-gCA ^{3'}		

*03:130, 03:140			3	4								11							
*03:145			3	4															
*03:151					5														
*04:01:01:01-04:01:15, 04:01:17-04:01:22, 04:01:24-04:01:27, 04:01:29-04:01:41, 04:05, 04:07, 04:09N, 04:12, 04:15:02, 04:17-04:20, 04:24-04:28, 04:30-04:31, 04:35, 04:38-04:41, 04:43-04:53, 04:56-04:57, 04:59Q-04:67, 04:69-04:70, 04:73-04:76, 04:78-04:79, 04:81-04:102, 04:104-04:106, 04:109-04:110, 04:113, 04:115N, 04:117-04:121, 04:123N-04:127												6		11		13			
*04:01:16, 04:23, 04:54, 04:72												6		11		13			
*04:01:23, 04:33, 04:116												6		11		13			
*04:01:28														11		13			
*04:03, 04:80			3									6		11		13	15		
*04:04:01-04:04:02, 04:13, 04:111, 04:122												6		11		13			
*04:06			3									6		11		13	15		
*04:08												6		11		13			
*04:10												6		11					
*04:11												6					16		
*04:14, 04:68												6		11		13			
*04:15:01, 04:15:03												6		11		13			
*04:16			3									6		11		13	15		
*04:29												6				13	16		
*04:32, 04:77			2									6		11		13			
*04:34												6		11		13			
*04:36, 04:55												6					16		
*04:37												6		11		13	15		
*04:42												6		11		13	14		
*04:58												6		11		13			
*04:71	1											6		11		13			
*04:103												6		11		13	15		
*04:107			3									6		11		13	15		
*04:108												6		11		13			
*04:112												6		11		13			
*04:114												6		11		13	16		
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16			
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16			
HLA-C allele⁴																			
*05:01:01:01-05:01:21, 05:03, 05:05-05:10, 05:12-05:15, 05:18:01-05:19, 05:21, 05:24, 05:28-05:31, 05:33-05:35, 05:37-05:38, 05:40-05:41, 05:44:01-05:45, 05:47-05:50, 05:52-05:54, 05:56-05:57, 05:58:02-05:61, 05:63-05:67, 05:69, 05:71-05:77, 05:80-05:81																	7	11	15
*05:04, 05:22:01-05:22:02																	7	11	15
*05:11, 05:32, 05:51Q, 05:70																	7	11	15

*05:16							7				11				15	
*05:17, 05:68							7				11				15	
*05:20							7				11					
*05:23							7				11				15	
*05:25, 05:55							7				11				15	
*05:26						6	7				11				15	
*05:27							7				11				15	
*05:36							7				11				15	
*05:39							7				11				15	
*05:42							7				11		13		15	
*05:43							7				11					
*05:46							7				11		13		15	
*05:58:01			3				7				11				15	
*05:62							7				11				15	
*05:78							6	7			11					
*05:79, 08:10 ⁷							7			10	11				15	
*06:02:01:01-06:02:01:02, 06:02:03-06:02:11, 06:02:13-06:02:15, 06:02:17-06:02:21, 06:07, 06:10, 06:12-06:13, 06:15-06:16N, 06:18-06:30, 06:33-06:34, 06:36-06:39, 06:41-06:51, 06:53-06:66, 06:68-06:71, 06:73-06:75, 06:78, 06:80, 06:83-06:87								8			11			14		
*06:02:12								8						14		
*06:02:16, 06:81								8			11					
*06:03:01			3					8			11		14	15		
*06:03:02, 06:76:01								8			11		14	15		
*06:04, 06:40, 06:79N								8			11					
*06:05								8			11					
*06:06								8			11					
*06:08			2					8			11			14		
*06:09								8			11			14		
*06:11, 06:52								8						14		
*06:14, 06:35, 06:72								8			11					
*06:17											11			14		
*06:31								8			11			14		
*06:32											11			14		
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
HLA-C allele⁴																
*06:67								8			11			14		
*06:76:02											11	12	13	14	15	
*06:77								8			11			14		
*06:82								8			11			14		

*07:01:01-07:01:06, 07:01:08-07:02:07, 07:02:10-07:02:28, 07:02:30-07:04:02, 07:04:04-07:06, 07:08, 07:11-07:13, 07:15-07:19, 07:21-07:25, 07:27:01, 07:29-07:30, 07:32N, 07:35-07:38:02, 07:42, 07:44-07:48, 07:50-07:52, 07:54-07:63, 07:65-07:75, 07:77-07:91, 07:93-07:95, 07:97-07:100, 07:102-07:122, 07:124, 07:126-07:138, 07:140-07:147, 07:149-07:151, 07:153-07:155, 07:157-07:160, 07:162-07:172, 07:174-07:176, 07:178-07:183, 07:185-07:195, 07:197-07:199:01, 07:200-07:209, 07:211-07:215, 07:217-07:237, 07:239-07:241, 07:243-07:245, 07:248-07:262										9							
*07:01:07, 07:04:03, 07:14, 07:27:02, 07:53, 07:199:02, 07:216										9							
*07:02:08, 07:02:29, 07:26, 07:33N, 07:39-07:40, 07:92, 07:139, 07:152N, 07:156, 07:246										9							
*07:02:09					6					9				13			
*07:07, 07:09, 07:49, 07:76, 07:247										9		11					
*07:10, 07:28, 07:41, 07:43, 07:184, 07:196										9							
*07:20, 07:263										9							
*07:31										9					14		
*07:64										9							16
*07:96			3							9							
*07:101, 07:148, 07:161		2								9	10						
*07:123, 07:173										9							
*07:125						6				9				13	14		
*07:177										9					14		
*07:210, 07:238										9		11					
*07:242				4						9							
*08:01:01-08:01:08, 08:03:01-08:03:02, 08:06, 08:08, 08:14, 08:16:01-08:16:02, 08:20, 08:22, 08:24, 08:26N, 08:28, 08:36N, 08:38, 08:40, 08:42, 08:44, 08:46, 08:50, 08:56-08:61											10						
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	
HLA-C allele⁴																	
*08:02:01-08:02:07, 08:04:01-08:04:02, 08:13, 08:15:01-08:15:02, 08:17-08:19, 08:23, 08:27, 08:29-08:30, 08:32-08:35, 08:37, 08:39, 08:43, 08:45, 08:49, 08:51-08:55N, 08:62-08:64											10						
*08:05											10			13			
*08:07, 08:47											10						

*08:09, 08:11											10					
*08:12											10					
*08:21											10		13			
*08:25											10		13			
*08:31		2									10					
*08:41											10					
*08:48											10					
*12:02:01-12:02:02, 12:02:04-12:02:07, 12:02:09, 12:17-12:18:01, 12:22, 12:30, 12:36, 12:40, 12:56, 12:64, 12:67-12:69, 12:73-12:74, 12:80N												12	13			
*12:02:03, 12:02:08, 12:10:01-12:10:02, 12:20												12	13			
*12:03:01:01-12:03:01:02, 12:03:03, 12:03:05-12:03:07, 12:03:10-12:03:12, 12:03:14-12:03:15, 12:03:18-12:03:19, 12:03:21-12:03:22, 12:06-12:07, 12:11-12:13, 12:19, 12:23, 12:25, 12:29, 12:32, 12:34-12:35, 12:37-12:39N, 12:42Q-12:43, 12:45-12:48, 12:50-12:53, 12:57, 12:61-12:63, 12:65-12:66, 12:70-12:71, 12:75-12:79, 12:82												12	13	14		
*12:03:02, 12:03:08, 12:03:16, 12:03:20, 12:26, 12:59														13	14	
*12:03:04												12		14		
*12:03:09							8					12		14		
*12:03:13												12	w	14		
*12:03:17												12	13	14		
*12:04:01, 12:33											11			14	15	
*12:04:02, 12:60											11	12	13	14	15	
*12:05											11	12		14	15	
*12:08												12	13			
*12:09											11				15	
*12:14:01													13			
*12:14:02												12	13			
*12:15							8					12	13	14		
*12:16												12	13			
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
HLA-C allele⁴																
*12:18:02, 16:43 ⁸													13			
*12:21											11	12			15	
*12:24												12	13			
*12:27												12	13			
*12:28, 12:58												12	13	14		
*12:31												12	13	14		
*12:41											11	12	13		15	
*12:44												12	13			
*12:49												12				
*12:54											11	12	13	14		
*12:55												12	13	14		16

*12:72												12	13		15	
*12:81												12	13	14		
*14:02:01-14:02:02, 14:02:04-14:02:07, 14:02:09-14:02:11, 14:07N, 14:11, 14:13-14:14, 14:18-14:19, 14:21N, 14:23-14:24:02, 14:26-14:27, 14:29-14:34, 14:36-14:37, 14:39-14:40, 14:43																16
*14:02:03, 14:03, 14:22, 14:38																16
*14:02:08												12				16
*14:04											11					16
*14:05																16
*14:06, 14:09, 14:28, A*30:62 ⁹																16
*14:08																16
*14:10																16
*14:12											11					16
*14:15														14		16
*14:16														14		16
*14:17																16
*14:20																16
*14:25			2													16
*14:35N																16
*14:41																
*15:02:01-15:02:09, 15:13, 15:15, 15:18, 15:26, 15:28, 15:31-15:35, 15:38-15:39, 15:41, 15:44-15:45, 15:47-15:53, 15:56-15:57, 15:60			3					8				11				15
*15:02:10			3					8				11				15
*15:03			3					8				11		13		15
*15:04-15:06:03, 15:09, 15:19, 15:24, 15:29-15:30, 15:40, 15:46, 15:54, 15:59			3									11				15
*15:07			3					8								
*15:08			3					8				11				15
*15:10:01								8				11				15
*15:10:02, 15:17			3					8				11				15
*15:11			3			6		8				11				15
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
HLA-C allele⁴																
*15:12					5			8				11				15
*15:16			3					8				11		13		15
*15:21			3					8								15
*15:22, 15:55			3									11				15
*15:23												11				
*15:25			3											13		
*15:27			3									11				15
*15:36						6						11				
*15:37			3					8				11				15
*15:42		2	3					8				11				15
*15:43			3					8								
*15:58			3					8				11				15
*15:61												11				15

*16:01:01, 16:01:03-16:01:10, 16:06-16:08, 16:10-16:11, 16:13-16:15:01, 16:16Q-16:17, 16:21-16:24, 16:26-16:28, 16:30N-16:32, 16:36, 16:38-16:39, 16:41, 16:44, 16:49, 16:51-16:52													13			
*16:01:02													13			
*16:02:01-16:02:09, 16:09, 16:12, 16:19, 16:25, 16:46-16:47										11			13		15	
*16:04:01, 16:33, 16:42														14		
*16:15:02											12		13			
*16:18													13			
*16:20													13			
*16:29		2												14		
*16:34			3			6							13			
*16:35														14		
*16:37													13		15	
*16:40							7									16
*16:45													13			
*16:48											11			14	15	
*16:50		2											13			
*17:01:01:01-17:01:03, 17:01:05-17:04, 17:06, 17:08-17:13		2									11				15	
*17:01:04		2									11		13		15	
*17:05		2									11				15	
*17:07											11				15	
*18:01-18:02, 18:04-18:05											11					
*18:03		2									11					
A*24:96, A*24:106, A*24:174, B*46:25																
B*14:03																
B*35:08:02, B*35:08:05																
B*58:02																
B*67:02													13			
HLA-C allele ⁴																
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	Well No.
17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	Well No.
																HLA-C allele ⁴
			20						26							*01:02:01-01:02:16, 01:06-01:08, 01:11-01:13, 01:15-01:20, 01:23, 01:25-01:33, 01:37N-01:42, 01:44- 01:48, 01:51-01:54, 01:56N-01:58, 01:61-01:73
			20	21					26							*01:03, 01:24, 01:34
																*01:04
																*01:05, 01:22, 01:36, B*54:18 ⁵
																*01:09
			20						26							*01:10
			20						26					31		*01:14
										27			30			*01:21
												29				*01:35
			20						26							*01:43
				21												*01:49-01:50, 01:55
			20						26							*01:59
		19	20						26							*01:60
							24							31		*02:02:01
							24							31		*02:02:02-02:02:03, 02:02:06- 02:02:12, 02:02:15-02:02:24, 02:07- 02:09, 02:11, 02:15, 02:19-02:20, 02:23-02:25Q, 02:28-02:31, 02:34- 02:38N, 02:40, 02:42-02:44, 02:46, 02:48, 02:50, 02:52N-02:54, 02:56- 02:57, 02:59-02:61, 02:63
														31		*02:02:05, 02:02:14, 02:10, 02:32, 02:51
														31		*02:02:13
																*02:03, 02:16:01-02:16:02
							24							31		*02:04
														31		*02:05, 02:17
17														31		*02:06
										w			30	31		*02:12
	18						24							31		*02:13
														31		*02:14
	18															*02:18
							24							31		*02:21, 02:26:02-02:26:03, 02:39, 02:45
				21			24							31		*02:22, 02:62
							24							31		*02:26:01
							24			27						*02:27:01
										27						*02:27:02
	18													31		*02:33
							24							31		*02:47
	18						24						30	31		*02:49
							24						30	31		*02:55
												29		31		*02:58
17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	Well No.
17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	Well No.
																HLA-C allele ⁴

Negative Control

