

INTERPRETATION TABLE

HLA-C*02 SSP subtyping

Amplification patterns for the HLA-C*02:02 to 02:36 alleles

	Well ⁵											
	1	2	3	4	5	6	7	8	9	10	11	12
Length of spec.	250	95	95	150	145	160	130	70	200	125	85	150
PCR product			120		240			280			170	230
Length of int.	800	800	800	1070	1070	800	800	1070	1070	1070	800	1070
pos. control ¹												
5'-primer(s) ²	2 nd	486	486	92	361	419	2 nd	105	703	113	486	97
	5'-CCA ^{3'}	5'-ACA ^{3'}	5'-ACA ^{3'}	5'-gTg ^{3'}	5'-AgT ^{3'}	5'-gTA ^{3'}	5'-CCA ^{3'}	5'-gCT ^{3'}	5'-CTA ^{3'}	5'-CCA ^{3'}	5'-ACA ^{3'}	5'-TCg ^{3'}
					453	420				118		368
					5'-AAT ^{3'}	5'-TTA ^{3'}				5'-CCA ^{3'}		5'-gTT ^{3'}
												449
												5'-CCA ^{3'}
3'-primer(s) ³	538	538	538	201	559	538	418	134	861	201	527	201
	5'-CCA ^{3'}	5'-CCA ^{3'}	5'-CAG ^{3'}	5'-CTT ^{3'}	5'-CTC ^{3'}	5'-CCA ^{3'}	5'-gTC ^{3'}	5'-AgC ^{3'}	5'-TCg ^{3'}	5'-CTT ^{3'}	5'-CCg ^{3'}	5'-CTT ^{3'}
			555					343			538	559
			5'-CCg ^{3'}					5'-g ^{3'}			5'-CCg ^{3'}	5'-CTC ^{3'}
			578								613	
			5'-Tgt ^{3'}								5'-gCA ^{3'}	
Well No.	1	2	3	4	5	6	7	8	9	10	11	12
HLA-C allele ⁴												
*02:02:01-02:02:03, 02:02:06-02:02:11	1	2							9	10		
*02:02:05	1								9	10		
*02:03			3						9	10		
*02:04	1	2		4					9	10		
*02:05	1				5				9	10		
*02:06	1					6			9	10		
*02:07	1	2					7		9	10		
Well No.	1	2	3	4	5	6	7	8	9	10	11	12
Length of spec.	250	95	95	150	145	160	130	70	200	125	85	150
PCR product			120		240			280			170	230
Well No.	1	2	3	4	5	6	7	8	9	10	11	12
*02:08	1	2						8	9	10		
*02:09	1	2							9	10	11	

*02:10	1								9	10		
*02:11	1	2							9	10		12
*02:12	1								9	10		
*02:13	1	2							9	10		
*02:14	1									10		12
*02:15, 02:21 ⁶	1	2							9	10		
*02:16:01									9	10		
*02:16:02			3						9	10		
*02:17	1				5				9	10		12
*02:18			3						9	10	11	
*02:19	1	2			6				9	10		
*02:20	1	2	3						9	10		
*02:22	1	2			5				9	10		
*02:23	1	2							9	10		
*02:24	1	2							9	10		
*02:25Q, 02:30 ⁷	1	2							9	10		
*02:26:01-02:26:02	1	2							9			
*02:27	1	2							9	10		
*02:28	1	2							9	10		
*02:29	1	2								10		
*02:31	1	2	3						9	10		
*02:32									9	10	11	
*02:33	1							8	9	10		
*02:34	1	2							9	10		
*02:35	1	2							9	10		
*02:36	1	2							9	10		
*01:04, 01:21, 06:11, 06:18, 06:33, 08:27, 08:29, 12:02:01- 12:03:11, 12:06-12:08, 12:10:01- 12:13, 12:17, 12:22, 12:23, 12:25-12:32, 12:34, 12:35, 14:16, 16:04:01	1											
*01:06, 03:36, 03:77, 07:95												
*01:09	1				6							
*01:10, B*07:02:07					5							
Well No.	1	2	3	4	5	6	7	8	9	10	11	12

*15:02:01-15:03, 15:05:01-15:06:02, 15:09-15:10:01, 15:12, 15:13, 15:15-15:18, 15:20-15:22, 15:24, 15:26-15:28												
*15:07												
*15:08												
*15:10:02								8				
*15:11									10			
*15:19												
*15:23												12
*16:09												
*16:10							7					
B*27:34, B*40:02:07, 40:06:02			3									
B*57:03:02			w									
B*73:01												
Well No.	1	2	3	4	5	6	7	8	9	10	11	12

			16		18						24	*15:02:01-15:03, 15:05:01-15:06:02, 15:09-15:10:01, 15:12, 15:13, 15:15-15:18, 15:20-15:22, 15:24, 15:26-15:28
			16								24	*15:07
					18						24	*15:08
			16		18						24	*15:10:02
			16								24	*15:11
			16		18			21			24	*15:19
			16		18						24	*15:23
					18				22			*16:09
												*16:10
												<i>B*27:34, B*40:02:07, 40:06:02</i>
												<i>B*57:03:02</i>
						w						<i>B*73:01</i>
13	14	15	16	17	18	19	20	21	22	23	24	Well No.