


INTERPRETATION TABLE


HLA-C*02 SSP subtyping


Amplification patterns for the HLA-C*02:02 to 02:66 alleles


	Well ⁹															
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
Length of spec.	250	95	95	150	145	160	130	70	200	125	85	150	225	80	130	250
PCR product(s)			120		240	215		280			170	230		115	190	
Length of int. pos. control ¹	800	800	800	1070	1070	800	800	1070	1070	1070	800	1070	1070	800	1070	1070
5'-primer(s) ²	5'-CCA ^{3'} 2 nd I	5'-ACA ^{3'} 486	5'-ACA ^{3'} 486	5'-gTg ^{3'} 92	5'-AgT ^{3'} 361	5'-ggT ^{3'} 364	5'-CCA ^{3'} 2 nd I	5'-gCT ^{3'} 105	5'-CTA ^{3'} 703	5'-CCA ^{3'} 113	5'-ACA ^{3'} 486	5'-TCg ^{3'} 97	5'-CCA ^{3'} 118	5'-ACA ^{3'} 486	5'-CCA ^{3'} 113	5'-CCA ^{3'} 2 nd I
					5'-AAT ^{3'} 453	5'-TTA ^{3'} 420				5'-CCA ^{3'} 118		5'-gTT ^{3'} 368			5'-TAC ^{3'} 369	
												5'-CCA ^{3'} 449				
3'-primer(s) ³	5'-CCA ^{3'} 538	5'-CCA ^{3'} 538	5'-CAG ^{3'} 538	5'-CTT ^{3'} 201	5'-CTC ^{3'} 559	5'-CCA ^{3'} 538	5'-gTC ^{3'} 418	5'-AgC ^{3'} 134	5'-TCg ^{3'} 861	5'-CTT ^{3'} 201	5'-CCg ^{3'} 527	5'-CTT ^{3'} 201	5'-ggC ^{3'} 302	5'-CCg ^{3'} 527	5'-CTT ^{3'} 201	5'-CAG ^{3'} 538
			5'-CCg ^{3'} 555					5'-g ^{3'} 343			5'-CCg ^{3'} 538	5'-CTC ^{3'} 559	5'-ggC ^{3'} 302	5'-CgT ^{3'} 559	5'-CCA ^{3'} 518	
			5'-TgT ^{3'} 578								5'-gCA ^{3'} 613					
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
HLA-C allele ^{4,5}																
*02:02:01-02:02:03, 02:02:06-02:02:12, 02:02:15-02:02:25, 02:48, 02:50, 02:54, 02:56-02:57, 02:59-02:66	1	2							9	10						
*02:02:05, 02:02:14, 02:51	1								9	10						
*02:02:13	1								9							
*02:03			3						9	10						16
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
Length of spec.	250	95	95	150	145	160	130	70	200	125	85	150	225	80	130	250
PCR product(s)			120		240	215		280			170	230		115	190	
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
*02:04	1	2		4					9	10						
*02:05	1				5				9	10						
*02:06	1					6			9	10						

*02:07	1	2					7		9	10							
*02:08	1	2						8	9	10							
*02:09	1	2							9	10	11						
*02:10	1									10							
*02:11	1	2							9	10		12					
*02:12	1								9	10			13				
*02:13	1	2							9	10				14			
*02:14	1									10		12					
*02:15, 02:21 ⁶	1	2							9	10					15		
*02:16:01									9	10							16
*02:16:02			3						9	10							16
*02:17	1					5			9	10		12					
*02:18			3						9	10	11						16
*02:19	1	2							9	10							
*02:20	1	2	3						9	10							
*02:22	1	2				5			9	10							
*02:23	1	2							9	10							
*02:24	1	2							9	10							
*02:25Q, 02:30 ⁷	1	2							9	10							
*02:26:01-02:26:03	1	2							9								
*02:27:01	1	2							9	10			13				
*02:27:02	1								9	10			13				
*02:28	1	2							9	10							
*02:29	1	2								10							
*02:31	1	2	3						9	10							
*02:32									9	10	11						
*02:33	1							8	9	10							
*02:34	1	2							9	10							
*02:35	1	2							9	10							
*02:36	1	2							9	10							
*02:37, 02:52N ⁸	1	2							9	10							
*02:38N	1	2							9	10							
*02:39	1	2							9								
*02:40	1	2							9	10							
*02:42	1	2							9	10							
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	
Length of spec.	250	95	95	150	145	160	130	70	200	125	85	150	225	80	130	250	
PCR product(s)			120		240	215		280			170	230		115	190		
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	
*02:43	1	2							9	10				14			
*02:44	1	2							9	10							
*02:45	1	2							9								
*02:46	1	2				6			9	10							
*02:47	1	2				6			9	10							
*02:49, 02:55	1	2							9	10							
*02:53	1	2							9	10							
*02:58	1									10							












*02:36	1	2							9	10								
*02:37, 02:52N ⁸	1	2							9	10								
*02:38N	1	2							9	10								
*02:39	1	2							9									
*02:40	1	2							9	10								
*02:42	1	2							9	10								
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16		
Length of spec.	250	95	95	150	145	160	130	70	200	125	85	150	225	80	130	250		
PCR product(s)			120		240	215		280			170	230		115	190			
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16		
*02:43	1	2							9	10				14				
*02:44	1	2							9	10								
*02:45	1	2							9									
*02:46	1	2				6			9	10								
*02:47	1	2				6			9	10								
*02:49, 02:55	1	2							9	10								
*02:53	1	2							9	10								
*01:04, 04:94:01-04:94:02, 06:02:09, 06:11, 06:18, 06:33, 08:27, 08:29, 12:02:01-12:03:16, 12:03:18-12:03:20, 12:06-12:08, 12:10:01-12:13, 12:17, 12:22-12:23, 12:25-12:32, 12:34-12:40, 12:42Q-12:59, 12:61-12:71, 12:73-12:76, 14:16, 16:04, 16:33, 16:42	1																	
*01:09	1																	
*01:10, B*07:02:07, B*27:05:15					5													
*01:21	1																	
*01:22, 01:35, 04:04:01-04:04:02, 04:13, 04:34, 04:58, 08:01:01-08:01:05, 08:03:01-08:04:02, 08:06, 08:08-08:09, 08:11, 08:13, 08:16:01-08:16:02, 08:20-08:22, 08:24, 08:26N, 08:36N, 08:38-08:42, 08:44, 08:46, 08:50, 08:54, 08:56-08:60, 12:14:01-12:14:02, 12:18:01-12:18:02, 12:20, 14:06, 14:15, 16:35, 16:40. B*58:02																		16
*03:07, 03:10, 03:15, 03:29, 03:45, 04:16, 04:103, 05:01:01:01-05:01:18, 05:03-05:07N, 05:09:01-05:10, 05:12, 05:14-05:16, 05:18:01-05:19, 05:21-05:25, 05:28-05:29, 05:31-05:42, 05:44-05:51Q, 05:53-05:67, 05:69, 05:71-05:74, 06:06, 06:40, 07:07, 07:09, 07:76, 12:09, 16:02:01-16:02:07, 16:02:09, 16:12, 16:19, 16:25, 16:37, 18:01-18:02, 18:04-18:05																		
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16		


*16:10							7										
*16:34										10				13			
<i>B*07:113, B*08:64, B*40:192</i>																	
<i>B*27:34, B*40:02:07, 40:06:02</i>			3														
<i>B*35:01:30, B*57:03:02</i>			w														
<i>B*67:02</i>																	
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	







A blank sheet of graph paper with a grid pattern. The grid consists of 20 columns and 30 rows of small squares. The grid is centered on the page and is completely empty, with no text or markings inside the cells.



110	235	80	180	115	110	205	325	160	140	140	90	210	80	100	150	Length of spec.	
160				205		390					220					PCR product(s)	
17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	Well No.	
								25								*03:21, 03:80, 03:142, B*07:55, B*07:100, B*15:45, B*15:63	
																*03:64	
							24									*03:81, 04:108, 07:123, 07:173, 16:20	
												29				*04:03, 04:80, 04:107	
												29				*04:06	
	18															*04:42, 05:43	
																*05:08, 05:52, 06:02:01:01- 06:02:01:02, 06:02:03-06:02:08, 06:02:10-06:03:02, 06:07, 06:09- 06:10, 06:12-06:17, 06:19-06:32, 06:34-06:39, 06:41-06:76, 12:04:01- 12:05, 12:21, 12:33, 12:41, 12:60, 12:72, 18:03	
																*05:11, 05:17, 05:27, 05:68, 05:70, 06:04, 08:10, 15:04, 15:30, 15:40, 17:01:01:01-17:11	
																*05:26	
	18													31		*06:05	
	18															*06:08	
														31		*07:02:09	
														31		*08:14	
																*08:31, 16:29	
																*08:37, 16:21, B*27:84	
17																*12:03:17	
							24									*12:15	
														31		*12:16	
																*14:25	
																*15:02:01-15:03, 15:05:01-15:06:03, 15:09-15:10:01, 15:12-15:13, 15:15- 15:18, 15:20-15:22, 15:24, 15:26- 15:29, 15:31-15:35, 15:37-15:39, 15:41-15:42, 15:44-15:60	
							24									*15:07, 15:36, 15:43	
	18						24									*15:08	
	18						24									*15:10:02	
							24									*15:11	
17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	Well No.	
110	235	80	180	115	110	205	325	160	140	140	90	210	80	100	150	Length of spec.	
160				205		390					220					PCR product(s)	
17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	Well No.	
	18			21			24									*15:19	
	18						24									*15:23	
	18				22											*16:09	

