

INTERPRETATION TABLE

HLA-C*08 SSP subtyping

Amplification patterns of the HLA-C*08:01 to 08:60 alleles

	Well ⁷															
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
Length of spec.	250	115	115	110	155	270	100	225	385	110	280	165	105	170	265	100
PCR product(s)						310		290	505	140				280		545
Length of int.	800	1070	800	1070	1070	800	1070	1070	1070	800	800	800	800	1070	1070	1070
pos. control ¹																
5'-primer ²	5'-CCA ^{3'} 2 nd I	5'-TAC ^{3'} 527	5'-TgA ^{3'} 527	5'-TAC ^{3'} 527	5'-gCA ^{3'} 176	5'-CCA ^{3'} 2 nd I	5'-AAT ^{3'} 453	5'-CAA ^{3'} 351	5'-AAA ^{3'} 312	5'-gCA ^{3'} 176	5'-AgT ^{3'} 361	5'-gCA ^{3'} 176	5'-gCg ^{3'} 539	5'-CgC ^{3'} 173	5'-CgA ^{3'} 1 st I	5'-TCT ^{3'} 142
					485			419	736	527				363		972
								5'-gTC ^{3'}	5'-gCA ^{3'}	5'-TAC ^{3'}				5'-Agc ^{3'}		5'-CTA ^{3'}
3'-primer ³	5'-TCA ^{3'} 539	5'-CTT ^{3'} 601	5'-CTT ^{3'} 601	5'-CCT ^{3'} 595	5'-AgC ^{3'} 289	5'-CgC ^{3'} 559	5'-CCA ^{3'} 512	5'-CTT ^{3'} 601	5'-CgT ^{3'} 526	5'-gCA ^{3'} 277	5'-CTT ^{3'} 601	5'-ggC ^{3'} 302	5'-CTT ^{3'} 601	5'-ggC ^{3'} 302	5'-CCg ^{3'} 175	5'-CTT ^{3'} 201
					289	599			956	598				601		1034
					5'-AgC ^{3'}	5'-TCC ^{3'}			5'-CAG ^{3'}	5'-CTC ^{3'}				5'-CTT ^{3'}		5'-AgT ^{3'}
					601											
					5'-CTC ^{3'}											
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
HLA-C allele ⁴																
*08:01:01-08:01:05, 08:46, 08:50, 08:58, 08:60	1	2										12				
*08:02:01-08:02:05, 08:45, 08:48, 08:53			3									12	13			
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
Length of spec.	250	115	115	110	155	270	100	225	385	110	280	165	105	170	265	100
PCR product(s)						310		290	505	140				280		545
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
*08:02:06			3									12	13			
*08:02:07			3									12	13			
*08:03:01-08:03:02	1	2		4								12				
*08:04:01-08:04:02	1		3									12				
*08:05			3		5							12	13			
*08:06	1	2		4		6						12				
*08:07			3				7					12	13			

*05:18:01			3											13			
*05:18:02-05:18:03, 05:39			3											13			
*05:23			3				7							13			
*05:25			3											13			
*05:27	1		3														
*05:29			3											13		15	
*05:42			3		5									13			
*05:43			3											13			
*05:44			3											13			
*05:46			3		5									13			
*05:62							7										
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	
Length of spec.	250	115	115	110	155	270	100	225	385	110	280	165	105	170	265	100	
PCR product(s)						310		290	505	140				280		545	
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	
*05:64			3											13			
*05:70			3														
*06:05, 12:16																	16
*06:67					5												
*06:73																	
*07:02:09																	16
*07:04:01-07:04:07, 07:11-07:12, 07:63, 07:139, 07:142, 07:181, 07:199																	
*07:41			3											13			
*07:68								8									
*07:101					5							12					
*15:23	1																
*15:29	1																16
*16:33						6											
*17:05	1				5												
A*01:01:06																	
A*29:10																	
B*14:32			3								11						
B*44:148																	
HLA-C allele																	
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	

INTERPRETATION TABLE

HLA-C*08 SSP subtyping

Amplification patterns of the HLA-C*08:01 to 08:60 alleles

Well ⁷																																
17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32																	
																	Length of spec.															
430	200	250	110	195	80	95	145	115	125	100	185	205	405	500	125	PCR product(s)																
1070	1070	800	1070	1070	1070	1070	1070	1070	1070	1070	1070	1070	1070	1070	1070	Length of int.																
																pos. control ¹																
																5'-primer ²																
																3'-primer ³																
																Well No.																
																HLA-C allele ⁴																
				21												*08:01:01-08:01:05, 08:46, 08:50, 08:58, 08:60																
																*08:02:01-08:02:05, 08:45, 08:48, 08:53																
																Well No.																
																Length of spec.																
																PCR product(s)																
																Well No.																
																*08:02:06																
																*08:02:07																
				21												*08:03:01-08:03:02																
																*08:04:01-08:04:02																
																*08:05																
																*08:06																
																*08:07																
																Well No.																
																Length of spec.																
																PCR product(s)																
																Well No.																
																*08:02:06																
																*08:02:07																
																*08:03:01-08:03:02																
																*08:04:01-08:04:02																
																*08:05																
																*08:06																
																*08:07																

