

INTERPRETATION TABLE

HLA-C*12 SSP subtyping

Amplification patterns of the C*12:02 to 12:50 alleles

Well¹¹

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
Length of spec.	235	100	220	340	130	75	140	95	95	155	220	140	105	100	95	185
PCR product(s)						415 150		245 155					150 105	145 100	140 95	225 185
Length of int.	800	###	800	###	###	###	800	###	###	###	###	###	###	###	800	###
pos. control ¹																
5'-primer(s) ²	98 5'-CTA 3'	419 5'-gTC 3'	361 5'-AgT 3'	1 st I 5'-CgA 3'	201 5'-CCA 3'	28 5'-TCA 3'	201 5'-CCA 3'	98 5'-CTA 3'	289 5'-Agg 3'	361 5'-AgT 3'	361 5'-AgT 3'	201 5'-CCA 3'	368 5'-gTT 3'	142 5'-TCT 3'	201 5'-CCA 3'	2 nd I 5'-CCA 3'
						431 5'-CgT 3'		420 5'-TTA 3'					412 5'-ATA 3'	368 5'-gTC 3'	420 5'-TTA 3'	
						504 5'-CAT 3'										
3'-primer(s) ³	289 5'-AgC 3'	477 5'-gCA 3'	538 5'-CCA 3'	302 5'-ggT 3'	289 5'-AgT 3'	270 5'-TAG 3'	302 5'-ggT 3'	214 5'-CCA 3'	341 5'-Cgg 3'	474 5'-gCA 3'	538 5'-CCg 3'	302 5'-ggC 3'	474 5'-gCA 3'	201 5'-CTT 3'	299 5'-TCT 3'	473 5'-CAA 3'
	289 5'-AgC 3'				289 5'-AgT 3'	538 5'-CCA 3'		301 5'-gCC 3'					477 5'-gCA 3'	474 5'-gCA 3'	474 5'-gCA 3'	512 5'-CCA 3'
	295 5'-TCC 3'							474 5'-gCA 3'								
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
HLA-C allele ⁴																
*12:02:01-12:02:05	1	2							9			12				
*12:02:06	1								9			12				

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
Length of spec.	235	100	220	340	130	75	140	95	95	155	220	140	150	100	95	185
PCR product(s)						415	150						150	145	140	225
								245	155							
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
*12:03:01:01-12:03:01:02, 12:03:03, 12:03:05- 12:03:07, 12:03:10- 12:03:12, 12:03:14	1		3						9	10		12				
*12:03:02, 12:03:08	1		3						9			12				
*12:03:04, 12:03:09	1		3						9	10						
*12:03:13	w		3						9	10		12				
*12:04:01	1		3	4					9							
*12:04:02	1		3	4			7		9	10						
*12:05			3	4	5		7			10						
*12:06, 12:48 ⁶	1		3			6			9	10		12				
*12:07	1		3					8	9	10		12				
*12:08	1	2				6			9			12				
*12:09, 05:16 ⁵				4	5		7				11					
*12:10:01-12:10:02	1	2							9			12	13			
*12:11	1		3						9	10		12		14		
*12:12	1		3						9	10		12			15	
*12:13	1		3						9	10		12				16
*12:14:01	1								9			12				16
*12:14:02	1	2							9			12				16
*12:15	1		3					8	9	10		12			15	
*12:16, 01:21 ⁵		2							9			12		14		
*12:17	1	2							9			12				
*12:18	1	2							9			12				
*12:19	1								9			12				
*12:20	1								9	10		12				
*12:21				4	5		7									
*12:22	1	2							9			12				
*12:23	1		3						9	10		12				
*12:24	1								9	10	11	12				
*12:25	1		3						9	10		12				
*12:26	1		3						9			12				
*12:27	1	2							9							
*12:28	1		3							10		12				

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
*01:14, 06:06, 14:12				4							11					
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
Length of spec.	235	100	220	340	130	75	140	95	95	155	220	140	105	100	95	185
PCR product(s)						150		155					150	145	140	225
						415		245								
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
*01:17									9		11					
*01:22, 01:35, 14:06, 14:15, 15:43, B*58:02																
*01:32											11					
*02:02:01-02:02:03, 02:02:05-02:02:07, 02:02:09-02:02:11, 02:02:14-02:04, 02:06-02:11, 02:13-02:16:02, 02:18-02:26:02, 02:28-02:40, 02:43-02:46, 04:10, 05:01:01-05:01:11, 05:01:13-05:01:15, 05:03-05:10, 05:12-05:15, 05:18-05:22, 05:24-05:26, 05:28-05:35, 05:37-05:41, 05:43-05:45, 05:47-05:51Q, 15:08, 16:02:01-16:02:05, 16:09, 16:12, 16:19				4	5		7									
*02:02:08, 02:02:13, 03:07, 03:45, 04:01:10, 04:01:23, 05:01:12, 06:35, 06:40, 07:09, 07:49, 07:76, 18:01-18:04				4												
*02:02:12					5											
*02:12	1								9			w				
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
Length of spec.	235	100	220	340	130	75	140	95	95	155	220	140	105	100	95	185
PCR product(s)						150		155					150	145	140	225

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
						415		245								
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
*02:27:01-02:27:02, 04:11, 04:36, 04:55, 08:02:01-08:02:02, 08:02:04-08:02:05, 08:15, 08:17-08:19, 08:23, 08:27-08:35, 08:37, 08:43, 14:08, 16:01:01, 16:01:03-16:01:05, 16:06-16:08, 16:10-16:11, 16:16Q-16:18, 16:20-16:24, 16:26-16:28, 16:30N-16:32					5							12				
*02:42				4	5		7							14		
*03:08, 03:29, 03:31						6										
*03:15	1			4												
*03:27, B*07:13	1								9							
*03:38:01-03:38:02	1								9							
*03:53	1															
*03:69	1															
*03:102, 07:24, 07:81																
*04:01:01:01-04:01:04, 04:01:06-04:01:09, 04:01:11-04:01:22, 04:01:24-04:01:26, 04:04:01-04:05, 04:07-04:09N, 04:13-04:15:02, 04:17-04:20, 04:23-04:28, 04:30-04:32, 04:34-04:35, 04:38-04:54, 04:56-04:57, 04:59Q-04:79, 04:81					4		7									
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
Length of spec.	235	100	220	340	130	75	140	95	95	155	220	140	105	100	95	185
PCR product(s)						150		155					150	145	140	225
						415		245								
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
*06:44			3	4		6										
*07:02:09												12		14		
*07:07				4					9							
*07:10, 07:43																
*07:16, 07:51									9							
*07:26, 07:92, 07:96, B*07:15	1															
*08:01:01-08:01:03, 08:03:01-08:04, 08:06, 08:08-08:09, 08:11, 08:13, 08:16, 08:20, 08:22, 08:24, 08:26N, 08:36N, 08:38- 08:42, 08:44, 15:07					5							12				
*08:05, 08:25, 16:15:01	1								9			12				
*08:07					5							12				16
*08:12, 14:02:03, 14:03, 14:10, 14:22					5						11	12				
*08:14					5							12		14		
*08:21	1								9			12				
*14:04				4					9		11					16
*14:17											11					
*15:03, 15:16	1			4			7		9							
*15:21					5							w				
*15:25	1											12				
*15:45				4	w		7									
*16:04:01, 16:29			3		5							12				
*16:13					5							12				
*16:14					5			8				12				
*16:15:02	1								9	10		12				
*17:01:01:01-17:01:03, 17:01:05-17:08	1			4					9							16
*17:01:04	1			4			7		9							16
A*02:211																
B*35:08:02												12			15	

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
<i>B*40:29</i>																
<i>B*67:02</i>	1								9			12		14		
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	Well No.
												29				*12:29, 12:38 ⁷
										27						*12:30, 12:36 ⁸
17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	Well No.
565	140	100	105	250	110	150	185	80	355	100	275	125	90	120	115	Length of spec.
	165		170		580		425	145		205	350	185	230	185		PCR product(s)
	245		220									210				
	270															
17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	Well No.
																*12:31
													30			*12:32, 12:34 ⁹
																*12:33, 02:05, 02:17 ⁵
	18															*12:35
															32	*12:37
								25								*12:39N
																*12:40
																*12:41
												29		31		*12:42Q
							24									*12:43
									26							*12:44
											28					*12:45, 12:50 ¹⁰
			20													*12:46N
														31		*12:47
																*12:49
																*01:02:01-01:03, 01:06-01:08, 01:10-01:13, 01:15-01:16, 01:18-01:20, 01:23-01:31, 01:33-01:34, 01:37N-01:45, 03:58, 03:86, 03:94, 03:99, 14:02:01-14:02:02, 14:02:04-14:02:06, 14:05, 14:07N, 14:11, 14:13-14:14, 14:18-14:21N, 14:23-14:25, B*14:03
																*01:04
																*01:09, 06:11, 06:52, 14:16

17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	Well No.
																*01:14, 06:06, 14:12
17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	Well No.
565	140	100	105	250	110	150	185	80	355	100	275	125	90	120	115	Length of spec.
	165		170		580		425	145		205	350	185	230	185		PCR product(s)
	245		220									210				
	270															
17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	Well No.
																*01:17
				21												*01:22, 01:35, 14:06, 14:15, 15:43, B*58:02
										28						*01:32
																*02:02:01-02:02:03, 02:02:05-02:02:07, 02:02:09-02:02:11, 02:02:14-02:04, 02:06- 02:11, 02:13-02:16:02, 02:18-02:26:02, 02:28- 02:40, 02:43-02:46, 04:10, 05:01:01:01-05:01:11, 05:01:13-05:01:15, 05:03- 05:10, 05:12-05:15, 05:18- 05:22, 05:24-05:26, 05:28- 05:35, 05:37-05:41, 05:43- 05:45, 05:47-05:51Q, 15:08, 16:02:01-16:02:05, 16:09, 16:12, 16:19
																*02:02:08, 02:02:13, 03:07, 03:45, 04:01:10, 04:01:23, 05:01:12, 06:35, 06:40, 07:09, 07:49, 07:76, 18:01- 18:04
																*02:02:12
																*02:12
17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	Well No.
565	140	100	105	250	110	150	185	80	355	100	275	125	90	120	115	Length of spec.
	165		170		580		425	145		205	350	185	230	185		PCR product(s)

17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	Well No.
	245		220									210				
	270															
17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	Well No.
																*02:27:01-02:27:02, 04:11, 04:36, 04:55, 08:02:01- 08:02:02, 08:02:04- 08:02:05, 08:15, 08:17- 08:19, 08:23, 08:27-08:35, 08:37, 08:43, 14:08, 16:01:01, 16:01:03- 16:01:05, 16:06-16:08, 16:10-16:11, 16:16Q- 16:18, 16:20-16:24, 16:26- 16:28, 16:30N-16:32
																*02:42
																*03:08, 03:29, 03:31
																*03:15
																*03:27, B*07:13
									26							*03:38:01-03:38:02
	18															*03:53
									26							*03:69
											28					*03:102, 07:24, 07:81
																*04:01:01:01-04:01:04, 04:01:06-04:01:09, 04:01:11-04:01:22, 04:01:24-04:01:26, 04:04:01-04:05, 04:07- 04:09N, 04:13-04:15:02, 04:17-04:20, 04:23-04:28, 04:30-04:32, 04:34-04:35, 04:38-04:54, 04:56-04:57, 04:59Q-04:79, 04:81
17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	Well No.

17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	Well No.
565	140	100	105	250	110	150	185	80	355	100	275	125	90	120	115	Length of spec.
	165		170		580		425	145		205	350	185	230	185		PCR product(s)
	245		220									210				
	270															
17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	Well No.
							24									*04:01:05
																*04:03, 04:06, 04:16, 05:42
	18															*04:12
																*04:29, 08:02:03
																*04:33
																*04:37
																*04:58
								26								*04:80
				21												*05:11, 05:17, 05:27, 08:10, 15:02:01-15:02:07, 15:04-15:06:03, 15:09- 15:13, 15:15, 15:17-15:20, 15:22-15:24, 15:26-15:42, 15:44
																*05:23
																*05:36
																*05:46, 16:25
																*06:02:01:01-06:02:01:02, 06:02:03-06:02:11, 06:07- 06:10, 06:12-06:13, 06:15- 06:19, 06:21-06:29, 06:31- 06:34, 06:36-06:39, 06:42- 06:43, 06:45-06:51, 06:53- 06:55
																*06:03
				21												*06:04
																*06:05
								26								*06:14
										28						*06:20
																*06:30
													30			*06:41

17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	Well No.
17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	Well No.

17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	Well No.
		19														<i>B*40:29</i>
																<i>B*67:02</i>
17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	Well No.

