

INTERPRETATION TABLE								
HLA-Cw*01 SSP subtyping								
Amplification patterns of the Cw*0101 to Cw*0120 alleles								
	Well⁵							
	1	2	3	4	5	6	7	8
Length of spec.	90	90	150	210	150	195	150	210
PCR product(s)		270	265	240	200		195	
Length of int.	800	1070	800	1070	800	1070	1070	1070
pos. control¹								
5'-primer(s)²	368	368	89	89	368	406	368	368
	5'-gTg ^{3'}	5'-gTg ^{3'}	5'-gAA ^{3'}	5'-gAA ^{3'}	5'-gTg ^{3'}	5'-gCA ^{3'}	5'-gTg ^{3'}	5'-gTg ^{3'}
			368	368				
			5'-gTg ^{3'}	5'-gTg ^{3'}				
3'-primer(s)³	419	419	312	287	479	559	479	538
	5'-CgT ^{3'}	5'-CgA ^{3'}	5'-AgT ^{3'}	5'-TCg ^{3'}	5'-CCA ^{3'}	5'-CgT ^{3'}	5'-CCA ^{3'}	5'-CCA ^{3'}
		595	477	538	527		523	
		5'-Cct ^{3'}	5'-gCA ^{3'}	5'-CCg ^{3'}	5'-CCA ^{3'}		5'-ACA ^{3'}	
Well No.	1	2	3	4	5	6	7	8
HLA-Cw allele⁴								
*010201-010205	1							
*0103		2						
*0104			3					8
*0105				4				
*0106	1				5			
*0107	1					6		
*0108	1						7	
*0109	1							8
*0110	1							
*0111	1							
*0112								
*0113	1							
*0114	1		3					
*0115	1	2						
*0116	1			4				
*0117	1							
*0118	1							
*0119	1							
*0120	1				5		7	
*0509, 0517, 0815								
*0617, 120306, 1611								
*0737				4				
HLA-Cw allele⁴								
B*1533								
Well No.	1	2	3	4	5	6	7	8

INTERPRETATION TABLE								
HLA-Cw*01 SSP subtyping								
Amplification patterns of the Cw*0101 to Cw*0120 alleles								
Well ⁵								
9	10	11	12	13	14	15	16	
230	290	340	155	155	240	230	255	Length of spec. PCR product(s)
				370				
800	1070	800	1070	800	1070	800	1070	Length of int. pos. control ¹
368	368	368	485	89	89	89	89	5'-primer(s) ²
5'-gTg ^{3'}	5'-gTg ^{3'}	5'-gTA ^{3'}	5'-CAA ^{3'}	5'-gAA ^{3'}	5'-gAA ^{3'}	5'-gAA ^{3'}	5'-gAA ^{3'}	
				368				
				5'-gTg ^{3'}				
559	619	3 rd	601	201	289	278	301	3'-primer(s) ³
5'-CTC ^{3'}	5'-ACT ^{3'}	5'-ATg ^{3'}	5'-CTC ^{3'}	5'-CTC ^{3'}	5'-AgC ^{3'}	5'-ggT ^{3'}	5'-gCA ^{3'}	
				3 rd				
				5'-CTC ^{3'}				
9	10	11	12	13	14	15	16	Well No.
								HLA-Cw allele ⁴
				13				*010201-010205
				13				*0103
				13				*0104
				13				*0105
				13				*0106
				13				*0107
				13				*0108
				13				*0109
9				13				*0110
	10			13				*0111
		11		13				*0112
			12	13				*0113
				13				*0114
				13				*0115
				13				*0116
				13	14			*0117
				13		15		*0118
				13			16	*0119
				13				*0120
			12					*0509, 0517, 0815
		11						*0617, 120306, 1611
								*0737
								HLA-Cw allele ⁴
			12					B*1533
9	10	11	12	13	14	15	16	Well No.