

<b>INTERPRETATION TABLE</b>								
<b>HLA-Cw*06 SSP subtyping</b>								
<b>Amplification patterns of the Cw*0602 to Cw*0616N alleles</b>								
	<b>Tube<sup>5</sup></b>							
	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>
<b>Length of spec.</b>	<b>240</b>	<b>220</b>	<b>135</b>	<b>240</b>	<b>165</b>	<b>240</b>	<b>110</b>	<b>240</b>
<b>PCR product</b>								
<b>Length of int.</b>	<b>800</b>	<b>800</b>	<b>1070</b>	<b>1070</b>	<b>800</b>	<b>1070</b>	<b>1070</b>	<b>1070</b>
<b>pos. control<sup>1</sup></b>								
<b>5'-primer<sup>2</sup></b>	<b>28</b>	<b>361</b>	<b>105</b>	<b>2<sup>nd</sup> In</b>	<b>176</b>	<b>2<sup>nd</sup> In</b>	<b>232</b>	<b>361</b>
	5'-TCA <sup>3'</sup>	5'-AgT <sup>3'</sup>	5'-gCT <sup>3'</sup>	5'-CCA <sup>3'</sup>	5'-gCA <sup>3'</sup>	5'-CCA <sup>3'</sup>	5'-AgA <sup>3'</sup>	5'-AgT <sup>3'</sup>
<b>3'-primer<sup>3</sup></b>	<b>97</b>	<b>538</b>	<b>201</b>	<b>539</b>	<b>302</b>	<b>538</b>	<b>302</b>	<b>539</b>
	5'-gTC <sup>3'</sup>	5'-CCA <sup>3'</sup>	5'-CTC <sup>3'</sup>	5'-TCA <sup>3'</sup>	5'-ggT <sup>3'</sup>	5'-CCg <sup>3'</sup>	5'-ggT <sup>3'</sup>	5'-CTC <sup>3'</sup>
<b>Tube No.</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>
<b>HLA-Cw allele<sup>4</sup></b>								
<b>*06020101-060202</b>	<b>+</b>	<b>+</b>						
<b>*0603</b>		<b>+</b>	<b>+</b>					
<b>*0604</b>	<b>+</b>			<b>+</b>				
<b>*0605</b>	<b>+</b>				<b>+</b>	<b>+</b>		
<b>*0606</b>	<b>+</b>					<b>+</b>		
<b>*0607</b>	<b>+</b>	<b>+</b>					<b>+</b>	
<b>*0608</b>	<b>+</b>	<b>+</b>						<b>+</b>
<b>*0609</b>	<b>+</b>	<b>+</b>						
<b>*0610</b>	<b>+</b>	<b>+</b>						
<b>*0611</b>	<b>+</b>	<b>+</b>						
<b>*0612</b>	<b>+</b>	<b>+</b>						
<b>*0613</b>	<b>+</b>	<b>+</b>						
<b>*0614</b>	<b>+</b>							
<b>*0615</b>	<b>+</b>	<b>+</b>						
<b>*0616N</b>	<b>+</b>	<b>+</b>						
<b>*010201-0103, 0105-0108, 0111-0119, 04010101- 040104, 0403, 0405, 0407- 0412, 0415-0421, 0423-0427, 0802, 0805, 0807, 0812, 1209, 140201-1405, 1407N, 1408, 1508, 1801, 1802</b>						<b>+</b>		
<b>*0104, 0109, 12030101-1207, 1211-1213, 1215, 160401</b>		<b>+</b>						
<b>*0110</b>						<b>+</b>		<b>+</b>
<b>Tube No.</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>

Length of spec.	240	220	135	240	165	240	110	240
PCR product								
Tube No.	1	2	3	4	5	6	7	8
*0203, 0216, 0218, 040401-040402, 0406, 0413, 080101-080102, 0803, 0804, 0806, 0808, 0809, 0811, 0813, 0814, 1214, 1218, 1220, 1406, 150201-1507, 1509-1521, 1701-1704				+				
*0205, 0217		+						+
*030201-0317, 0319-0331, 0333-0338, 0340, 0341			+					
*0332			+					
*0339, 1216	+							
*050101-0507N, 0509, 0510, 0512-0516					+	+		
*0508					+			
*0511, 0517, 0810				+	+			
*070401-070402, 0711, 0712, 0745								
*0710, 0743								
*1221					+			
*1803								
HLA-Cw allele <sup>4</sup>								
Tube No.	1	2	3	4	5	6	7	8

# INTERPRETATION TABLE

## HLA-Cw\*06 SSP subtyping

Amplification patterns of the Cw\*0602 to 0616N alleles

Tube <sup>5</sup>								
9	10	11	12	13	14	15	16	
165	200	130	210	155	295	365	240	Length of spec. PCR product
1070	800	1070	1070	1070	1070	800	1070	Length of int. pos. control <sup>1</sup>
412	142	213	529	361	341	376	501	5'-primer <sup>2</sup>
5'-ATA <sup>3'</sup>	5'-TCC <sup>3'</sup>	5'-CCC <sup>3'</sup>	5'-AgA <sup>3'</sup>	5'-AgT <sup>3'</sup>	5'-ggA <sup>3'</sup>	5'-gCT <sup>3'</sup>	5'-ggC <sup>3'</sup>	
538	302	302	3 <sup>rd</sup> I	475	353	3 <sup>rd</sup> I	3 <sup>rd</sup> I	3'-primer <sup>3</sup>
5'-CCA <sup>3'</sup>	5'-ggT <sup>3'</sup>	5'-ggC <sup>3'</sup>	5'-CTC <sup>3'</sup>	5'-ggT <sup>3'</sup>	5'-TgA <sup>3'</sup>	5'-CTC <sup>3'</sup>	5'-CTC <sup>3'</sup>	
9	10	11	12	13	14	15	16	Tube No.
								HLA-Cw allele <sup>4</sup>
								*06020101-060202
								*0603
								*0604
								*0605
								*0606
								*0607
								*0608
+								*0609
	+							*0610
		+						*0611
			+					*0612
				+				*0613
					+			*0614
						+		*0615
							+	*0616N
								*010201-0103, 0105-0108, 0111-0119, 04010101- 040104, 0403, 0405, 0407- 0412, 0415-0421, 0423-0427, 0802, 0805, 0807, 0812, 1209, 140201-1405, 1407N, 1408, 1508, 1801, 1802
								*0104, 0109, 12030101-1207, 1211-1213, 1215, 160401
								*0110
9	10	11	12	13	14	15	16	Tube No.

165	200	130	210	155	295	365	240	Length of spec. PCR product
9	10	11	12	13	14	15	16	Tube No.
								*0203, 0216, 0218, 040401-040402, 0406, 0413, 080101-080102, 0803, 0804, 0806, 0808, 0809, 0811, 0813, 0814, 1214, 1218, 1220, 1406, 150201-1507, 1509-1521, 1701-1704
								*0205, 0217
								*030201-0317, 0319-0331, 0333-0338, 0340, 0341
					+			*0332
								*0339, 1216
								*050101-0507N, 0509, 0510, 0512-0516
+								*0508
								*0511, 0517, 0810
		+						*070401-070402, 0711, 0712, 0745
					+			*0710, 0743
								*1221
+								*1803
								HLA-Cw allele <sup>4</sup>
9	10	11	12	13	14	15	16	Tube No.