

INTERPRETATION TABLE					
KIR HLA Ligand					
	Well				
	1	2	3	4	5
Length of spec.	340	340	350	350	370
PCR product					
Length of int.	800	800	800	1070	1070
pos. control <sup>1</sup>					
5'-primer(s) <sup>2</sup>	1 <sup>st</sup> I	1 <sup>st</sup> I	1 <sup>st</sup> I	1 <sup>st</sup> I	1 <sup>st</sup> I
	5'-CgA <sup>3'</sup>	5'-CgA <sup>3'</sup>	5'-CAg <sup>3'</sup>	5'-CAg <sup>3'</sup>	5'-gCA <sup>3'</sup>
3'-primer <sup>3</sup>	302	302	309	309	317
	5'-ggC <sup>3'</sup>	5'-ggT <sup>3'</sup>	5'-gTg <sup>3'</sup>	5'-ATC <sup>3'</sup>	5'-ggA <sup>3'</sup>
Well No.	1	2	3	4	5
HLA allele					
C*01:02:01-01:13, 01:15-01:46, 01:48-01:57, 02:27:01-02:27:02, 03:02:01-03:03:14, 03:03:16-03:04:16, 03:04:18-03:06, 03:08-03:09, 03:11:01-03:14, 03:16-03:28, 03:30-03:44, 03:46-03:114, 03:116-03:128, 04:11, 04:29, 04:36, 04:55, 06:11, 07:01:01-07:06, 07:08, 07:10-07:33N, 07:35-07:48, 07:50-07:75, 07:77-07:114, 07:116-07:207, 08:01:01-08:09, 08:11-08:53, 12:02:01-12:03:17, 12:06-12:08, 12:10:01-12:20, 12:22-12:26, 12:28-12:32, 12:34-12:40, 12:42Q-12:53, 12:55-12:59, 12:61-12:64, 14:02:01-14:03, 14:05-14:11, 14:13-14:31, 15:07, 15:25, 15:43, 16:01:01-16:01:06, 16:04:01, 16:06-16:08, 16:10-16:11, 16:13-16:18, 16:20-16:24, 16:26-16:36, 16:38-16:40	1				
C*01:14, 02:02:01-02:02:03, 02:02:05-02:02:11, 02:02:13-02:11, 02:13-02:26:03, 02:28-02:40, 02:42-02:52N, 03:07, 03:15, 03:45, 04:01:01-04:01:28, 04:01:30-04:01:33, 04:03-04:10, 04:12-04:20, 04:23-04:28, 04:30-04:35, 04:37-04:54, 04:56-04:100, 05:01:01-05:01:17, 05:03-05:63, 06:02:01-06:02:01:02, 06:02:03-06:02:11, 06:02:13-06:10, 06:12-06:51, 06:53-06:68, 07:07, 07:09, 07:49, 07:76, 08:10, 12:04:01-12:05, 12:09, 12:21, 12:33, 12:41, 12:54, 14:04, 14:12, 15:02:01-15:06:03, 15:08-15:13, 15:15-15:20, 15:22-15:24, 15:26-15:42, 15:44-15:55, 16:02:01-16:02:08, 16:09, 16:12, 16:19, 16:25, 17:01:01-17:10, 18:01-18:05		2			
C*02:12, 03:03:15, 03:10, 03:29, 15:21, 16:37	w				
Well No.	1	2	3	4	5
Length of spec.	340	340	350	350	370
PCR product					

Well No.	1	2	3	4	5
B*07:36, 07:38, 07:81, 08:03, 08:52, 15:13:01-15:13:02, 15:16:01-15:17:02, 15:23-15:24, 15:67, 15:87, 15:95, 15:157, 15:162, 15:168, 15:177, 15:196, 15:208, 15:216, 15:222, 15:230, 27:02:01-27:02:02, 27:30, 27:53, 27:57, 27:62, 27:65N, 27:75, 27:77, 38:01:01-38:01:05, 38:05-38:07, 38:09-38:14, 38:16, 38:19-38:22, 38:24-38:28, 38:30-38:34N, 40:13, 40:19, 40:109, 40:117, 44:06, 44:18, 44:25, 44:50, 44:95, 48:18, 49:01:01-49:01:03, 49:03-49:20, 51:01:01-51:24:04, 51:26-51:46, 51:48-51:53, 51:55-51:77, 51:79-51:122, 52:01:01-52:19, 52:21-52:26, 53:01:01-53:02, 53:04-53:08:02, 53:10, 53:14-53:26, 54:12, 56:21, 57:01:01-57:11, 57:13-57:52, 58:01:01-58:02, 58:04-58:16, 58:18-58:34, 59:01:01-59:05				4	
B*08:02, 13:01:01-13:04, 13:06-13:08Q, 13:10-13:23, 13:25-13:38, 13:40-13:50, 15:36, 15:89, 15:115, 18:09, 27:01, 37:10, 38:02:01-38:04, 38:08, 38:15, 38:18, 38:23, 38:29, 38:35, 40:47, 40:96, 40:110, 40:157, 44:02:01:01-44:02:18, 44:02:20-44:05:02, 44:07-44:08, 44:10, 44:12-44:17, 44:19N-44:24, 44:26-44:45, 44:47-44:49, 44:51-44:74, 44:76-44:89, 44:91-44:94, 44:96-44:128, 44:130, 44:132-44:134, 47:04, 49:02, 51:54, 51:78, 52:20, 53:09, 53:11-53:13, 56:07	3				
A*01:95, 02:81, 02:87, 02:112, 02:124, 02:129, 02:136, 23:01:01-23:46, 24:02:01:01-24:03:02, 24:05-24:11N, 24:13:01-24:15, 24:17-24:18, 24:20-24:27, 24:29-24:43, 24:45N-24:64, 24:66-24:88, 24:90N-24:99, 24:101-24:108, 24:110-24:128, 24:130-24:182, 25:01:01-25:16, 29:13, 31:07-31:08, 31:10, 32:01:01-32:36, 68:36					5
HLA allele					
Well No.	1	2	3	4	5