

101.123-24/06 – including *Taq* polymerase, IFU-01  
101.123-24u/06u – without *Taq* polymerase, IFU-02

Visit [www.olerup.com](http://www.olerup.com) for  
“Instructions for Use” (IFU)

Lot No.: **5F2**

Lot-specific information  
**Olerup SSP<sup>®</sup> DRB5**

Product number:	101.123-24/06 – including <i>Taq</i> pol. 101.123-24u/06u – without <i>Taq</i> pol.
Lot number:	5F2
Expiry date:	2020-04-01
Number of tests:	24 tests – Product No. 101.123-24/24u 6 tests – Product No. 101.123-06/06u
Number of wells per test:	18+1
Storage - pre-aliquoted primers:	dark at -20°C
- PCR Master Mix:	-20°C
- Adhesive PCR seals	RT
- Product Insert	RT

**This Product Description is only valid for Lot No. 5F2.**

Complete product documentation consists of generic Instructions for Use (IFU), lot specific Product Insert, Worksheet and Certificate.

**CHANGES COMPARED TO THE PREVIOUS OLERUP SSP<sup>®</sup> DRB5 LOT (1E2)**

The DRB5 kit is updated to enable separation of:

- Confirmed DRB5 alleles as listed in the IMGT/HLA database<sup>1</sup>
- Polymorphisms in exons outside of the region encoding the peptide binding domain
- Null and Alternatively expressed alleles

Three wells have been added to DRB5, wells **17 to 19**.

The format of the Worksheet has been changed.

The DRB5 primer set, specificity and interpretation tables have been updated for the HLA-DRB alleles described since the previous *Olerup SSP<sup>®</sup> DRB5* lot was made (**Lot No. 1E2**). The kit design is based on IMGT/HLA database 3.28.0.

101.123-24/06 – including *Taq* polymerase, IFU-01

101.123-24u/06u – without *Taq* polymerase, IFU-02

Visit [www.olerup.com](http://www.olerup.com) for  
“Instructions for Use” (IFU)

Lot No.: **5F2**

Lot-specific information

The primers of the wells detailed below have been exchanged, added or modified compared to the previous lot.

Well	5'-primer	3'-primer	rationale
16	Added	Added	Negative control moved to well 19, primer pairs added for the DRB5*01:20, DRB5*01:21 and DRB5*02:10 alleles.
17	New	New	New primer pairs added for the DRB5*01:27N and DRB5*02:12 alleles.
18	New	New	New primer pair added for the DRB5*02:13 allele.
19	-	-	Negative control added from well 16.

101.123-24/06 – including *Taq* polymerase, IFU-01  
101.123-24u/06u – without *Taq* polymerase, IFU-02

Visit [www.olerup.com](http://www.olerup.com) for  
“Instructions for Use” (IFU)

Lot No.: **5F2**

Lot-specific information

Well **19** contains Negative Control primer pairs, that will amplify more than 95% of the Olerup SSP<sup>®</sup> HLA Class I, DRB, DQB1, DPB1 and DQA1 amplicons as well as all the amplicons generated by the control primer pairs matching the human growth hormone gene.

HLA-specific PCR product sizes range from 75 to 200 base pairs.

The PCR product generated by the positive control primer pair is 430 base pairs.

Length of PCR product	105	200	105	80	75	80	85
<b>5'-primer<sup>1</sup></b>	<b>164</b>	<b>340</b>	<b>440</b>	<b>45</b>	<b>45</b>	<b>43</b>	<b>36</b>
	5'-CAC <sup>3'</sup>	5'-Agg <sup>3'</sup>	5'-TTA <sup>3'</sup>	5'-Tgg <sup>3'</sup>	5'-Tgg <sup>3'</sup>	5'-Tgg <sup>3'</sup>	5'-TAC <sup>3'</sup>
							<b>36</b>
							5'-TAT <sup>3'</sup>
<b>3'-primer<sup>2</sup></b>	<b>231</b>	<b>2<sup>nd</sup> I</b>	<b>507</b>	<b>59</b>	<b>58</b>	<b>57</b>	<b>47</b>
	5'-TgC <sup>3'</sup>	5'-AAA <sup>3'</sup>	5'-TTg <sup>3'</sup>	5'-CTC <sup>3'</sup>	5'-ggC <sup>3'</sup>	5'-CTC <sup>3'</sup>	5'-ACA <sup>3'</sup>
							<b>48</b>
							5'-gCA <sup>3'</sup>
							<b>48</b>
							5'-gCC <sup>3'</sup>
							<b>52</b>
							5'-TgT <sup>3'</sup>
<b>A*</b>	+	+	+				
<b>B*</b>	+	+	+				
<b>C*</b>	+	+	+				
<b>DRB1</b>				+	+		
<b>DRB3</b>				+	+		
<b>DRB5</b>				+			
<b>DQB1</b>					+		
<b>DPB1</b>						+	
<b>DQA1</b>							+

<sup>1</sup>The nucleotide position for HLA class I genes and the codon for HLA class II genes, in the 2<sup>nd</sup> or 3<sup>rd</sup> exon, matching the specificity-determining 3'-end of the primer is given. Nucleotide and codon numbering as on the [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla) web site. The sequence of the 3 terminal nucleotides of the primer is given.

<sup>2</sup>The nucleotide position for HLA class I genes and the codon for HLA class II genes, in the 2<sup>nd</sup> or 3<sup>rd</sup> exon or the 2<sup>nd</sup> intron, matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide and codon numbering as on the [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla) web site. The sequence of the 3 terminal nucleotides of the primer is given.

101.123-24/06 – including *Taq* polymerase, IFU-01  
101.123-24u/06u – without *Taq* polymerase, IFU-02

Visit [www.olerup.com](http://www.olerup.com) for  
“Instructions for Use” (IFU)

Lot No.: **5F2**

Lot-specific information

## PRODUCT DESCRIPTION

### DRB5 SSP subtyping

#### CONTENT

The primer set contains 5'- and 3'-primers for identifying the DRB5\*01:01:01 to DRB5\*01:31 and the DRB5\*02:02 to DRB5\*02:16 alleles.

#### PLATE LAYOUT

Each test consists of 19 PCR reactions in a 24 well cut PCR plate. Wells 20 to 24 are empty.

1	2	3	4	5	6	7	8
9	10	11	12	13	14	15	16
17	18	NC	empty	empty	empty	empty	empty

The 24 well cut PCR plate is marked with 'DRB5' in silver/gray ink.

Well No. 1 is marked with the Lot No. '5F2'.

Wells 1 to 18 – DRB5 primers.

Well 19 – Negative Control (NC).

A faint row of numbers is seen between wells 1 and 2 or wells 7 and 8 of the PCR trays. These stem from the manufacture of the trays, and should be disregarded. The PCR plates are covered with a PCR-compatible foil.

**Please note:** When removing each 24 well PCR plate, make sure that the remaining plates stay covered. Use a scalpel or a similar instrument to carefully cut the foil between the plates.

#### INTERPRETATION

Due to the sharing of sequence motifs between DRB5 alleles, non-DRB5 alleles will be amplified by some primer mixes. For further details see Specificity Table.

#### UNIQUELY IDENTIFIED ALLELES

All the DRB5 alleles, i.e. **DRB5\*01:01:01 to DRB5\*01:31 and DRB5\*02:02 to DRB5\*02:16**, recognized by the HLA Nomenclature Committee in April 2016<sup>1,2</sup> will be amplified by the primers in the DRB5 subtyping kit.

The DRB5 kit enables separation of the confirmed DRB5 alleles as listed in the IMGT/HLA database. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources. Current allele confirmation status for DRB5 alleles is listed below.

The DRB5 kit also enables identification of polymorphisms in exons outside of the region encoding the peptide binding domain and of null and alternatively expressed alleles.

101.123-24/06 – including *Taq* polymerase, IFU-01  
101.123-24u/06u – without *Taq* polymerase, IFU-02

Visit [www.olerup.com](http://www.olerup.com) for  
“Instructions for Use” (IFU)

Lot No.: **5F2**

Lot-specific information

<sup>1</sup>DRB5 alleles listed on the IMGT/HLA web page 2017-April-13, release 3.28.0, [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla).

<sup>2</sup>Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

## ALLELE CONFIRMATION STATUS

Allele	Status <sup>1</sup>	Allele	Status <sup>1</sup>	Allele	Status <sup>1</sup>
DRB5*01:01:01	Confirmed	DRB5*01:19	Unconfirmed	DRB5*02:09	Unconfirmed
DRB5*01:01:02	Unconfirmed	DRB5*01:20	Confirmed	DRB5*02:10	Confirmed
DRB5*01:01:03	Unconfirmed	DRB5*01:21	Confirmed	DRB5*02:11	Unconfirmed
DRB5*01:02	Confirmed	DRB5*01:22	Unconfirmed	DRB5*02:12	Confirmed
DRB5*01:03	Confirmed	DRB5*01:23	Unconfirmed	DRB5*02:13	Confirmed
DRB5*01:04	Unconfirmed	DRB5*01:24	Unconfirmed	DRB5*02:14	Unconfirmed
DRB5*01:05	Unconfirmed	DRB5*01:25	Unconfirmed	DRB5*02:15	Unconfirmed
DRB5*01:06	Unconfirmed	DRB5*01:26	Unconfirmed	DRB5*02:16	Unconfirmed
DRB5*01:07	Unconfirmed	DRB5*01:27N	Unconfirmed		
DRB5*01:08N	Confirmed	DRB5*01:28	Unconfirmed		
DRB5*01:09	Unconfirmed	DRB5*01:29	Unconfirmed		
DRB5*01:10N	Confirmed	DRB5*01:30	Unconfirmed		
DRB5*01:11	Unconfirmed	DRB5*01:31	Unconfirmed		
DRB5*01:12	Unconfirmed	DRB5*02:02	Confirmed		
DRB5*01:13	Confirmed	DRB5*02:03	Confirmed		
DRB5*01:14	Unconfirmed	DRB5*02:04	Unconfirmed		
DRB5*01:15	Unconfirmed	DRB5*02:05	Unconfirmed		
DRB5*01:16	Unconfirmed	DRB5*02:06	Unconfirmed		
DRB5*01:17	Confirmed	DRB5*02:07	Confirmed		
DRB5*01:18	Confirmed	DRB5*02:08	Unconfirmed		

<sup>1</sup>Allele status “confirmed” or “unconfirmed” as listed on the IMGT/HLA web page 2017-April-13, release 3.28.0, [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla).

## RESOLUTION IN HOMO- AND HETEROZYGOTES

Results file with resolution in DRB5 homo- and heterozygotes is available upon request.

101.123-24/06 – including *Taq* polymerase, IFU-01  
101.123-24u/06u – without *Taq* polymerase, IFU-02

Visit [www.olerup.com](http://www.olerup.com) for  
“Instructions for Use” (IFU)

Lot No.: **5F2**

Lot-specific information  
**SPECIFICITY TABLE**

**DRB5 SSP subtyping**

**Specificities and sizes of the PCR products of the 18+1 primer mixes used for DRB5 SSP subtyping**

Primer Mix	Size of spec. PCR product <sup>1</sup>	Size of control band <sup>2</sup>	Amplified DRB5 alleles <sup>3</sup>	Other amplified DRB alleles
1	255 bp	515 bp	*01:01:01-01:05, 01:07-01:20, 01:22-01:30, 02:03	DRB1*09:07
2	210 bp	515 bp	*01:01:01-01:05, 01:07-01:10N, 01:12-01:20, 01:22-01:31, 02:04, 02:08, 02:12	DRB1*09:07
3	225 bp	430 bp	*01:01:01-01:02, 01:04-01:05, 01:07-01:10N, 01:12-01:20, 01:22-01:31, 02:05, 02:08, 02:12	DRB1*09:07
4 <sup>4</sup>	100 bp	515 bp	*01:01:01-01:01:03, 01:04, 01:06-01:07, 01:09, 01:11, 01:15-01:19, 01:21-01:24, 01:26, 01:29-01:31	
5	150 bp	515 bp	*02:06	
6	145 bp	430 bp	*01:01:01, 01:05, 01:07, 01:09, 01:13, 01:16-01:19, 01:22-01:24, 01:26, 01:29-01:31	
7	145 bp	430 bp	*01:02-01:03, 01:05, 01:08N, 01:10N, 01:20, 01:25, 01:27N-01:28, 02:08, 02:12	
8	215 bp	430 bp	*01:02-01:03, 01:08N, 01:10N, 01:17, 01:20, 01:25, 01:27N-01:28, 02:05, 02:08, 02:12	
9 <sup>4</sup>	85 bp	430 bp	*01:03, 01:06, 01:09, 01:11, 01:21, 02:02-02:04, 02:06-02:07, 02:09-02:11, 02:13-02:16	
10	175 bp	430 bp	*01:16	
11 <sup>4</sup>	225 bp	430 bp	*01:13	
12	130 bp	430 bp	*01:04	
13	160 bp	430 bp	*01:07	
14	110 bp	430 bp	*01:12, 01:15	
15	200 bp	430 bp	*01:14	
16		430 bp	*01:06, 01:11, 01:21, 02:02-02:03, 02:06-02:07, 02:09-02:11, 02:13-02:16	
17	185 bp	515 bp	*02:02, 02:04-02:12, 02:14-02:16	
18	150 bp	430 bp	*01:01:02 <sup>?</sup> , 01:03 <sup>?</sup> , 01:07 <sup>?</sup> , 01:09 <sup>?</sup> , 01:18, 02:04 <sup>?</sup>	DRB1*15:02:03 <sup>?</sup> , DRB1*15:86, DRB1*16:01:02 <sup>?</sup> , DRB1*16:02:02 <sup>?</sup> ,

101.123-24/06 – including *Taq* polymerase, IFU-01  
101.123-24u/06u – without *Taq* polymerase, IFU-02

Visit [www.olerup.com](http://www.olerup.com) for  
“Instructions for Use” (IFU)

Lot No.: **5F2**

Lot-specific information

				DRB1*16:05:01 <sup>?</sup> , DRB4*01:05 <sup>?</sup> , DRB1*01:07:01 <sup>?</sup>
	195 bp		*01:08N	
14	145 bp	430 bp	*02:07	
15	235 bp	430 bp	*01:10N, 01:12, 01:15	DRB1*09:07
16 <sup>4</sup>	125 bp	430 bp	*01:21, 02:10	
	225 bp		*01:20, 02:08	DRB1*08:59, DRB1*11:210
17	130 bp	430 bp	*02:12	
	180 bp		*01:27N	
18	185 bp	430 bp	*02:13	DRB1*03:112, DRB1*09:06
19 <sup>5</sup>	-	-	Negative Control	

<sup>1</sup>Alleles are assigned by the presence of specific PCR product(s). However, the sizes of the specific PCR products may be helpful in the interpretation of DRB5 SSP subtypings. When the primers in a primer mix can give rise to HLA-specific PCR products of more than one length this is indicated if the size difference is more than 20 base pairs. Size differences of 20 base pairs or less are not given. For high resolution SSP kits, the alleles listed are specified according to amplicon length.

Nonspecific amplifications, i.e. a ladder or a smear of bands, may sometimes be seen. GC-rich primers have a higher tendency of giving rise to nonspecific amplifications than other primers. PCR fragments longer than the control bands may sometimes be observed. Such bands should be disregarded and do not influence the interpretation of the SSP typings.

PCR fragments migrating faster than the control bands, but slower than a 400 bp fragment may be seen in some gel read-outs. Such bands can be disregarded and do not influence the interpretation of the SSP typings.

Some primers may give rise to primer oligomer artifacts. Sometimes this phenomenon is an inherent feature of the primer pair(s) of a primer mix. More often it is due to other factors such as too low amount of DNA in the PCR reactions, taking too long time in setting up the PCR reactions, working at elevated room temperature or using thermal cyclers that are not pre-heated.

<sup>2</sup>The internal positive control primer pairs amplify segments of the human growth hormone gene. The internal positive control bands are 430 or 515 base pairs respectively, well distribution as outlined in the table. Well number 1 contains the longer, 515 bp, internal positive control band. The well distribution of the internal controls can help in orientation of the kit on gel photo, as well as allow for kit identification. In the presence of a specific amplification the intensity of the control band often decreases. In the presence of a specific amplification the intensity of the control band often decreases.

<sup>3</sup>For several DRB alleles 1<sup>st</sup> and/or 3<sup>rd</sup> exon(s) and above, as well as intron nucleotide sequences, are not available. In these instances it is not known whether some of the primers of the SSP sets are completely matched with the target sequences or not. Assumption is made that unknown sequences in these regions are conserved within allelic groups and that unknown sequences of codons 87 to 92 are identical with the DRB1\*01:01 consensus sequence.

<sup>4</sup>HLA-specific PCR products shorter than 125 base pairs have a lower intensity and are less sharp than longer PCR products.

<sup>5</sup>Primer mix 19 contains a negative control, which will amplify more than 95% of HLA amplicons as well as the amplicons generated by the control primer pairs matching the human growth hormone gene. HLA-specific PCR product sizes range from 75 to 200 base pairs and the PCR product generated by the HGH positive control primer pair is 430 base pairs. ‘?’, nucleotide sequence information not available for the primer matching sequence.

101.123-24/06 – including *Taq* polymerase, IFU-01  
101.123-24u/06u – without *Taq* polymerase, IFU-02

Visit [www.olerup.com](http://www.olerup.com) for  
“Instructions for Use” (IFU)

Lot No.: **5F2**

Lot-specific information

## PRIMER SPECIFICATION

Well No.	1	2	3	4	5	6	7	8	9	10	11	12
Length of spec.	255	210	225	100	150	145	145	215	85	130	110	185
PCR product(s)				150					175	160	200	
									225			
Length of int.	515	515	430	515	515	430	430	430	430	430	430	515
pos. control <sup>1</sup>												
5'-primer(s) <sup>2</sup>	13(125) 5'-gTA 3'	13(125) 5'-gTA 3'	13(125) 5'-gTA 3'	38(199) 5'-ACT 3'	37(196) 5'-Agg 3'	38(199) 5'-ACg 3'	37(196) 5'-AgA 3'	13(125) 5'-gTA 3'	13(125) 5'-gTA 3'	38(199) 5'-ACT 3'	13(125) 5'-gTA 3'	37(196) 5'-AgA 3'
				98(379) 5'-CTg 3'		38(199) 5'-ACg 3'	41(209) 5'-Cgg 3'		120(446) 5'-gAC 3'			
3'-primer(s) <sup>3</sup>	85(341) 5'-CAA 3'	67(286) 5'-gAA 3'	71(299) 5'-gCC 3'	57(258) 5'-gCg 3'	72(303) 5'-gCg 3'	72(303) 5'-gCg 3'	70(295) 5'-CTg 3'	70(295) 5'-gTT 3'	58(260) 5'-CCT 3'	67(286) 5'-gAT 3'	37(196) 5'-gTA 3'	85(341) 5'-CAg 3'
		67(286) 5'-gAA 3'	74(307) 5'-CAg 3'	135(490) 5'-gCC 3'			72(303) 5'-gCg 3'	71(299) 5'-gCg 3'	74(307) 5'-CAg 3'	78(319) 5'-CAC 3'	67(286) 5'-gAT 3'	
		70(296) 5'-TCC 3'	78(319) 5'-CAC 3'					71(299) 5'-gCg 3'	135(490) 5'-gCT 3'			
		72(303) 5'-gCg 3'										
Well No.	1	2	3	4	5	6	7	8	9	10	11	12

Well No.	13	14	15	16	17	18
Length of spec.	150	145	235	125	130	185
PCR product(s)	195			225	180	
Length of int.	430	430	430	430	430	430
pos. control <sup>1</sup>						
5'-primer(s) <sup>2</sup>	57(258) 5'-gAC 3'	24(157) 5'-ggT 3'	13(125) 5'-gTA 3'	6(103) 5'-CAT 3'	37(196) 5'-AgA 3'	37(196) 5'-AgA 3'
3'-primer(s) <sup>3</sup>	108(409) 5'-AgA 3'					
	93(365) 5'-gCg 3'	58(261) 5'-TCA 3'	78(319) 5'-CAC 3'	30(176) 5'-TgT 3'	67(286) 5'-gAT 3'	85(341) 5'-CAA 3'
	160(565) 5'-CAT 3'		79(323) 5'-TgC 3'	38(199) 5'-CAA 3'	83(336) 5'-CCC 3'	
			67(286) 5'-gAA 3'			
Well No.	13	14	15	16	17	18

<sup>1</sup>The internal positive control primer pairs amplify segments of the human growth hormone gene. The internal positive control bands are 430 or 515 base pairs respectively, well distribution as outlined in the table. Well number 1 contains the longer, 515 bp, internal positive control band. The well distribution of the internal controls can help in orientation of the kit on gel photo, as well as allow for kit identification. In the presence of a specific amplification the intensity of the control band often decreases.

<sup>2</sup>The nucleotide position matching the specificity-determining 3'-end of the primer is given. Nucleotide numbering as on the [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla) web site. The sequence of the 3 terminal nucleotides of the primer is given.

<sup>3</sup>The nucleotide position matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide numbering as on the [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla) web site. The sequence of the 3 terminal nucleotides of the primer is given.



101.123-24/06 – including *Taq* polymerase, IFU-01  
101.123-24u/06u – without *Taq* polymerase, IFU-02

Visit [www.olerup.com](http://www.olerup.com) for  
“Instructions for Use” (IFU)

Lot No.: **5F2**

Lot-specific information

CELL LINE VALIDATION SHEET																				
DRB5 SSP subtyping kit																				
				Prod. No.:	Well <sup>2</sup>															
					1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
					201556401	201556402	201669803	201556404	201556405	201556406	201669807	201669808	201783309	201556410	201669811	201556412	201669813	201669814	201556415	201783316
	IHWC cell line <sup>1</sup>	DRB5																		
1	9001 SA				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
2	9280 LK707	*01:02			+	+	+	-	-	+	+	-	-	-	-	-	-	-	-	-
3	9011 E4181324	*01:02			+	+	+	-	-	+	+	-	-	-	-	-	-	-	-	-
4	9275 GU373				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
5	9009 KAS011	*02:02			-	-	-	-	-	-	-	+	-	-	+	+	-	-	-	-
6	9353 SM				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
7	9020 QBL				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
8	9025 DEU				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
9	9026 YAR				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
10	9107 LKT3				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
11	9051 PITOUT				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
12	9052 DBB				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
13	9004 JESTHOM				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
14	9071 OLGA				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
15	9075 DKB				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
16	9037 SWEIG007				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
17	9282 CTM3953540				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
18	9257 32367				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
19	9038 BM16				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
20	9059 SLE005				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
21	9064 AMALA				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
22	9056 KOSE				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
23	9124 IHL				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
24	9035 JBUSH				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
25	9049 IBW9				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
26	9285 WT49				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
27	9191 CH1007				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
28	9320 BEL5GB				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
29	9050 MOU				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
30	9021 RSH				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
31	9019 DUCAF				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
32	9297 HAG				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
33	9098 MT14B				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
34	9104 DHIF				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
35	9302 SSTO				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
36	9024 KT17				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
37	9065 HHKB				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
38	9099 LZL				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
39	9315 CML				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
40	9134 WHONP199				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
41	9055 H0301				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
42	9066 TAB089				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
43	9076 T7526				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
44	9057 TEM				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
45	9239 SHJO				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
46	9013 SCHU	*01:01			+	+	+	+	+	-	-	-	-	-	-	-	-	-	-	-
47	9045 TUBO				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
48	9303 TER-ND				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-

101.123-24/06 – including *Taq* polymerase, IFU-01  
101.123-24u/06u – without *Taq* polymerase, IFU-02

Visit [www.olerup.com](http://www.olerup.com) for  
“Instructions for Use” (IFU)

Lot No.: **5F2**

Lot-specific information

CELL LINE VALIDATION SHEET				
DRB5 SSP subtyping kit				
				Well <sup>2</sup>
				17 18
			Prod. No.:	201783317 201783318
	IHWC cell line <sup>1</sup>	DRB5		
1	9001 SA		-	-
2	9280 LK707	*01:02	-	-
3	9011 E4181324	*01:02	-	-
4	9275 GU373		-	-
5	9009 KAS011	*02:02	-	-
6	9353 SM		-	-
7	9020 QBL		-	-
8	9025 DEU		-	-
9	9026 YAR		-	-
10	9107 LKT3		-	-
11	9051 PITOUT		-	-
12	9052 DBB		-	-
13	9004 JESTHOM		-	-
14	9071 OLGA		-	-
15	9075 DKB		-	-
16	9037 SWEIG007		-	-
17	9282 CTM3953540		-	-
18	9257 32367		-	-
19	9038 BM16		-	-
20	9059 SLE005		-	-
21	9064 AMALA		-	-
22	9056 KOSE		-	-
23	9124 IHL		-	-
24	9035 JBUSH		-	-
25	9049 IBW9		-	-
26	9285 WT49		-	-
27	9191 CH1007		-	-
28	9320 BEL5GB		-	-
29	9050 MOU		-	-
30	9021 RSH		-	-
31	9019 DUCAF		-	-
32	9297 HAG		-	-
33	9098 MT14B		-	-
34	9104 DHIF		-	-
35	9302 SSTO		-	-
36	9024 KT17		-	-
37	9065 HHKB		-	-
38	9099 LZL		-	-
39	9315 CML		-	-
40	9134 WHONP199		-	-
41	9055 H0301		-	-
42	9066 TAB089		-	-
43	9076 T7526		-	-
44	9057 TEM		-	-
45	9239 SHJO		-	-
46	9013 SCHU	*01:01	-	-
47	9045 TUBO		-	-
48	9303 TER-ND		-	-

101.123-24/06 – including *Taq* polymerase, IFU-01  
101.123-24u/06u – without *Taq* polymerase, IFU-02

Visit [www.olerup.com](http://www.olerup.com) for  
“Instructions for Use” (IFU)

Lot No.: **5F2**

**Lot-specific information**

<sup>1</sup>The provided cell line HLA specificities are retrieved from the <http://www.ihwg.org/hla> web site. The specificity of an individual cell line may thus be subject to change.

<sup>2</sup>The specificity of each primer solution in the kit has been tested against 48 well characterized cell line DNAs and where applicable, additional cell line DNAs.

No DNAs carrying the alleles to be amplified by primer solutions 9, 10 and 14 to 18 were available. The specificities of the primers in primer solutions 9, 10, 15, 17 and 18 were tested by separately adding one, two or three additional 5'-primers, respectively one or two additional 3'-primers. In primer solutions 14 and 16 it was only possible to test the 3'-primers, the 5'-primer was not possible to test. In primer solutions 4, 7 and 9 one 5'-primer was not possible to test, and in primer solutions 2, 8, 13 and 15 to 17 one 3'-primer was not possible to test. Additional primers in primer solutions 1 to 4, 7, 8, 11 and 13 were tested by separately adding additional 5'-primers and/or 3'-primers.

101.123-24/06 – including *Taq* polymerase, IFU-01  
 101.123-24u/06u – without *Taq* polymerase, IFU-02

Visit [www.olerup.com](http://www.olerup.com) for  
 “Instructions for Use” (IFU)

Lot No.: **5F2**

Lot-specific information

**ADDRESSES:**

**Manufacturer:**

**Olerup SSP AB**, Franzengatan 5, SE-112 51 Stockholm, Sweden.

**Tel:** +46-8-717 88 27

**Fax:** +46-8-717 88 18

**E-mail:** [olerup-se@caredx.com](mailto:olerup-se@caredx.com)

**Web page:** <http://www.olerup.com>

**Distributed by:**

**Olerup GmbH**, Löwengasse 47 / 6, AT-1030 Vienna, Austria.

**Tel:** +43-1-710 15 00

**Fax:** +43-1-710 15 00 10

**E-mail:** [olerup-at@caredx.com](mailto:olerup-at@caredx.com)

**Web page:** <http://www.olerup.com>

**Olerup Inc.**, 901 S. Bolmar St., Suite R, West Chester, PA 19382

**Tel:** 1-877-OLERUP1

**Fax:** 610-344-7989

**E-mail:** [olerup-us@caredx.com](mailto:olerup-us@caredx.com)

**Web page:** <http://www.olerup.com>

For information on *Olerup* distributors worldwide, contact **Olerup GmbH**.