

## **Olerup SSP<sup>®</sup> Ethidium Bromide Dropper Bottle**

**Product number:** 103.301-10  
**Lot number:** 24G  
**Volume:** 10 ml  
**Concentration:** 625 µg/ml  
**Expiry date:** 2013-July-01  
**Storage:** Room temperature (25°C)

**PRODUCT USE** Adding ethidium has never been easier! For the recommended final concentration of 0.5µg/ml, simply add one drop from the bottle for every 50 ml of solution. For the most accurate and consistent results, hold the dropper bottle vertically and squeeze slowly until the forming drop falls from the tip.

$$\mathbf{0.040\ ml \times 625\ \mu g/ml \ / \ 50\ ml = 0.5\ \mu g/ml}$$

drop size	stock	final	final
	concentration	volume	concentration

Remember to allow agarose to cool to 55-60°C before addition of ethidium to prevent volatilization of the compound.

For post-run staining of gels, use a 1X buffer with 0.5 µg/ml ethidium bromide for 30 min at RT. This can be followed with an optional destaining soak in water or 1 mM MgSO<sub>4</sub> for 20 min at RT, to reduce background fluorescence.

**STORAGE** Store at room temperature (25°C) and shield from light. The storage vial supplied with your dropper bottle is specially designed to exclude damaging UV light. It also provides an extra level of protection from accidental spills.

**DISPOSAL** All ethidium containing solutions and gels must be disposed of properly. Follow your local recommendations.

**CAUTIONS** Ethidium bromide is a strong mutagen and a moderate acute toxin. Always use gloves and eye protection when handling ethidium containing solutions. For further information, consult the MSDS for this product on the reverse of this sheet. The UV light from sources used for visualization of ethidium stained DNA is damaging to your eyes. Always wear appropriate protective eyewear when using these instruments.

**Ethidium Bromide Dropper Bottle**

2

**103.301-10****Lot No.: 24G**

PRODUCT IDENTIFICATION	PRODUCT: Ethidium bromide solution, 0.625 mg/ml concentration. PRODUCT No.: 103.301-10 MANUFACTURE: Genesee Scientific Corporation, San Diego, CA 92126. CHEMTREC: 800-424-9300												
PRODUCT COMPOSITION	<table><thead><tr><th>CAS #</th><th>CHEMICAL NAME</th><th>%</th><th>EINECS #</th></tr></thead><tbody><tr><td>1239-45-8</td><td>ethidium bromide</td><td>1</td><td>214-984-6</td></tr><tr><td>7732-18-5</td><td>water</td><td>99</td><td>231-791-2</td></tr></tbody></table>	CAS #	CHEMICAL NAME	%	EINECS #	1239-45-8	ethidium bromide	1	214-984-6	7732-18-5	water	99	231-791-2
CAS #	CHEMICAL NAME	%	EINECS #										
1239-45-8	ethidium bromide	1	214-984-6										
7732-18-5	water	99	231-791-2										
PHYSICAL DATA	Aqueous solution, red-orange colour, no odour. Specific gravity: 1.0 pH @ 25°C: 6.9 +/- 1.0												
HAZARD DATA AND FIRST AID MEASURES	IRRITANT, SUSPECTED MUTAGEN EYE: irritant, in case of contact rinse thoroughly with water and consult physician. SKIN: irritant, may cause local dermatitis, in case of contact wash thoroughly with soap & water, do not handle material with bare hands. Avoid all contact with clothing, properly dispose of all contaminated clothing. INGESTION: toxic, irritant, causes nausea, vomiting. Seek medical aid immediately. INHALATION: irritant, avoid breathing fumes.												
FIRE HAZARD DATA	Non-flammable, wear MSHA/NIOSH approved self-contained breathing apparatus. Flashpoint: N/A NFPA RATING: Not published Hazardous decomposition products: nitrogen oxides, CO, CO2												
ACCIDENTAL RELEASE	Wear chemical safety goggles and gloves, collect material for disposal with absorbant material. May be disposed of by incineration in an approved device.												
EXPOSURE CONTROL	Chemical safety goggles as per 29 CFR 1910.133 Chemical resistant gloves Avoid contact with eyes, skin and clothing Wash thoroughly after exposure Use respiration protection devices that meet OSHA CFR 1910.134 and ANSI Z88.2 specifications												
REACTIVITY	Stable at room temperature under normal conditions. Avoid contact with strong oxidizing agents.												
TOXICOLOGY	RTECS: CAS# 1239-45-8: SF7950000 (ethidium bromide) LD50/LC50: UNLISTED CARCINOGENICITY: Not listed by ACGIH, IARC, NIOSH, NTP or OSHA. No other data available. MUTAGENICITY: structurally similar to compounds known to be potent mutagens.												
ENVIRONMENTAL CONCERNS	Dispose of in a manner consistent with federal, state and local regulations. RCRA: No information available. Not listed as a RCRA land disposed banned material.												
HANDLING & STORAGE	Store at room temperature in a UV light resistant container. Avoid eye and skin contact. Wear protective eye-wear, gloves and clothing. Follow directions for use.												
TRANSPORT INFORMATION	No specific USDOT, IMO, IATA, RID/ADR or TDG information available. Constituent materials not listed at the concentrations occurring in these products.												
REGULATIONS	Federal: TCSA: ingredients listed, no ingredients have a SNUR under TCSA. SARA: None of the ingredients have a RQ. None of the ingredients have a TPQ. None of the ingredients are listed under Section 313 CAA. This material does not contain any hazardous air pollutants or Class 1 or Class 2 ozone depleters OSHA. State RTK listings: none European: No hazards symbols. Risk not available. Safety S24/25.												
OTHER INFORMATION	The above information is believed correct as of the above date, but does not purport to be all inclusive and shall be used only as a guide. <i>Olerup</i> SSP AB shall not be held liable for any damage resulting from handling or from contact with the above product.												

July 2009  
Rev. No.: 00*Olerup* SSP AB, Hasselstigen 1,  
SE-133 33 Saltsjöbaden, Sweden.  
Tel: +46-8-717 88 27  
Fax: +46-8-717 88 18