

Lot No.: 17R

## **Olerup SSP® GelRed™ Dropper Bottle**

|                        |  |
|------------------------|--|
| <b>Product number:</b> | <b>103.302-05/25</b>   |
| <b>Lot number:</b>     | <b>17R</b>   |
| <b>Volume:</b>         | <b>Product No. 103.302-05 – 5 ml</b><br><b>Product No. 103.302-25 – 5x5 ml</b> |
| <b>Expiry date:</b>    | <b>2014-July-01</b>  |
| <b>Storage:</b>        | <b>Room temperature</b>  |

**PRODUCT USE** GelRed™ is a sensitive, stable and environmentally safe fluorescent nucleic acid dye designed to replace the toxic ethidium bromide for staining dsDNA, ssDNA or RNA in agarose gels or polyacrylamide gels.

For the recommended final concentration, simply add 4 drops à 0.030 ml from the bottle for every 100 to 120 ml of gel solution. For the most accurate and consistent results, hold the dropper bottle vertically and squeeze slowly until the forming drop falls from the tip.

**STORAGE** Store at room temperature.

The dropper bottle is specially designed to exclude damaging UV light. GelRed™ should be shielded from light and exposure to light should be avoided during long-term storage.

The Olerup SSP® GelRed™ dropper bottle should be stored at room temperature. At low temperatures, such as +4°C, dye precipitation may occur. Should this occur, heat the solution in a water bath at +45°C to +50°C for two minutes and/or vortex the solution.

**DISPOSAL** Not classified as hazardous waste, thus can be safely disposed of down the drain or as regular trash. Consult local, state or national regulations for proper disposal.

**CAUTIONS** The dye is noncytotoxic and nonmutagenic at concentrations well above the working concentrations used in gel staining. For further information, consult the MSDS for this product. GelRed™ is impermeable to latex gloves and cell membranes. A complete safety report on GelRed™ may be downloaded from [www.biotium.com](http://www.biotium.com).

The UV light from sources used for visualization of GelRed™ stained DNA is damaging to your eyes. Always wear appropriate protective eyewear when using these instruments.

Lot No.: 17R

## Material Safety Data Sheet Olerup SSP® GelRed™ Dropper Bottle

### Chemical Product and Company Identification

Date of Issue: 2010-06-30  
Product Name: Olerup SSP® GelRed™ Dropper Bottle  
Product number: 103.302-05  
Manufacturer: Olerup SSP AB, Franzengatan 5, SE-112 51 Stockholm, Sweden  
Tel: +46 8 717 88 27  
Fax: +46 8 717 88 18

### Composition and Information on Ingredients

| Component | Chemical Name    | Weight %    | CAS No.   | EINECS No. |
|-----------|------------------|-------------|-----------|------------|
| GelRed™   | Proprietary      | Proprietary | N.A.      | N.A.       |
| Water     | dihydrogen oxide | Proprietary | 7732-18-5 | 231-791-2  |

This product is not classified according to the EU regulations

### Hazards Identification

Emergency overview: Handle with care – substance has not been fully tested.

### First Aid Measures

Potentially harmful; avoid prolonged or repeated exposure. Wash hands thoroughly after handling. If eye or skin contact occurs, wash affected areas with water for 15 minutes and seek medical advice. If inhaled, move individual to fresh air and seek medical advice. If swallowed, seek medical advice.

### Fire Fighting Measures

Use carbon dioxide, dry chemical extinguishers, foam extinguishers or water. An approved self-contained breathing apparatus and protective clothing are recommended.

### Accidental Release Measures

Use appropriate protective equipment and methods to clean up spilled substance promptly. Absorb spill onto an appropriate material. Collect and dispose of all waste in accordance with applicable laws.

### Handling and Storage

Store at room temperature. The dropper bottle is specially designed to exclude damaging UV light. GelRed™ should be shielded from light and exposure to light should be avoided during long-term storage. The Olerup SSP® GelRed™ dropper bottle should be stored at room temperature. At low temperatures, such as +4°C, dye precipitation may occur. Should this occur, heat the solution in a water bath at +45°C to +50°C for two minutes and/or vortex the solution.

### Exposure Controls and Personal Protection

Gloves, protective clothing and eyewear should be worn and safe laboratory practices followed.

ACGIH/OSHA Permissible Exposure Limit Data: Not determined

### Physical and Chemical Properties

|        |                |
|--------|----------------|
| Form:  | Liquid         |
| Color: | Dark red       |
| Odour: | Not determined |

**Lot No.: 17R**

|  |                |
|--|----------------|
| Solubility in/ miscibility with water: | High           |
| pH in solution (20°C):                 | Not determined |
| Density (20°C):                        | Not determined |

**Stability and Reactivity**

**Thermal Decomposition:** No decomposition if used according to specifications

**Dangerous Reactions:** No dangerous reactions identified

**Dangerous Products of Decomposition:** No dangerous decomposition products identified

**Toxicological Information**

**RTECS Number:** None known

**Toxicity:** GelRed was subjected to a series of tests both by Biotium and by three independent testing services to assess the dye's safety for routine handling and disposal. These tests include: 1) glove penetration test; 2) cell membrane permeability and cytotoxicity test; 3) Ames test; and 4) environmental safety tests. Test results confirm that the dye is impenetrable to both latex gloves and cell membranes. The dye is noncytotoxic and nonmutagenic at concentrations well above the working concentrations used in gel staining. Furthermore, GelRed has successfully passed environmental safety tests in compliance with CCR Title 22 Hazardous Waste Characterization. As a result, GelRed is not classified as hazardous waste, thus can be safely disposed of down the drain or as regular trash, providing convenience and reducing cost in waste disposal. However, since these tests were not performed on human, we still advise that researchers exercise precautions when handling the dye or any other DNA-binding molecules by wearing protective gears. For detailed test results on GelRed, you may download a complete safety report at [www.biotium.com](http://www.biotium.com).

**Health Hazards:** See toxicity.

**Potential Hazards:** We are not aware of any potential hazards

**Carcinogenicity:** Not listed by NTP, IARC or OSHA

**Ecological and Disposal Information**

Not classified as hazardous waste, thus can be safely disposed of down the drain or as regular trash.

Consult local, state or national regulations for proper disposal.

**Transport Information**

**Hazards Class:** Not determined

**Identification Number:** Not determined

**Packing Group:** Not classified

**Proper Shipping Name (Technical Name):** Not determined

**Regulatory information**

This product is not classified according to the EU regulations.

**U.S Federal regulation TSCA:** Not listed.

**SARA:** Not regulated.

**Clean Water Act (CWA):** No products were found.

**Clean Air Act (CAA):** No products were found.

**WHMIS (Canada):** Not controlled under WHMIS (Canada), no products found.

**Other Information**

The information provided above is believed to be correct to our best knowledge, but does not purport to be all inclusive, and shall be used only as a guide. This material is sold for laboratory purposes only and is not required to appear on the TSCA inventory. It is not intended for food, drug, household, agricultural or cosmetic use. Its use must be supervised by a technically qualified individual experienced in handling potentially hazardous chemicals. Olerup SSP AB shall not be held liable for any damage resulting from handling or contact with the above product.

Lot No.: **17R**

**ADDRESSES:**

**Manufacturer:**

**Olerup SSP AB**, Franzengatan 5, SE-112 51 Stockholm, Sweden.

**Tel:** +46-8-717 88 27

**Fax:** +46-8-717 88 18

**E-mail:** [info-ssp@olerup.com](mailto:info-ssp@olerup.com)

**Web page:** <http://www.olerup-ssp.com>

**Distributed by:**

**Olerup GmbH**, Löwengasse 47 / 6, AT-1030 Vienna, Austria.

**Tel:** +43-1-710 15 00

**Fax:** +43-1-710 15 00 10

**E-mail:** [support-at@olerup.com](mailto:support-at@olerup.com)

**Web page:** <http://www.olerup.com>

**Olerup Inc.**, 901 S. Bolmar St., Suite R, West Chester, PA 19382

**Tel:** 1-877-OLERUP1

**Fax:** 610-344-7989

**E-mail:** [info.us@olerup.com](mailto:info.us@olerup.com)

**Web page:** <http://www.olerup.com>

For information on *Olerup* SSP distributors worldwide, contact **Olerup GmbH**.