

Lot No.: 2H7

Olerup SSP® GelRed™ Dropper Bottle

Product number:	103.302-05/25
Lot number:	2H7
Volume:	Product No. 103.302-05 – 5 ml Product No. 103.302-25 – 5x5 ml
Expiry date:	2020-11-01
Storage:	Room temperature

PRODUCT USE GelRed™ is a sensitive, stable and environmentally safe fluorescent nucleic acid dye designed to replace the toxic ethidium bromide for staining dsDNA, ssDNA or RNA in agarose gels or polyacrylamide gels.

For the recommended final concentration, simply add 4 drops à 0.030 ml from the bottle for every 100 to 120 ml of gel solution. For the most accurate and consistent results, hold the dropper bottle vertically and squeeze slowly until the forming drop falls from the tip.

STORAGE Store at room temperature.

The dropper bottle is specially designed to exclude damaging UV light. GelRed™ should be shielded from light and exposure to light should be avoided during long-term storage.

The Olerup SSP® GelRed™ dropper bottle should be stored at room temperature. At low temperatures, such as +4°C, dye precipitation may occur. Should this occur, heat the solution in a water bath at +45°C to +50°C for two minutes and/or vortex the solution.

DISPOSAL Not classified as hazardous waste, thus can be safely disposed of down the drain or as regular trash. Consult local, state or national regulations for proper disposal.

CAUTIONS The dye is noncytotoxic and nonmutagenic at concentrations well above the working concentrations used in gel staining. For further information, consult the MSDS for this product. GelRed™ is impermeable to latex gloves and cell membranes. A complete safety report on GelRed™ may be downloaded from www.biotium.com.

The UV light from sources used for visualization of GelRed™ stained DNA is damaging to your eyes. Always wear appropriate protective eyewear when using these instruments.

Lot No.: 2H7

Material Safety Data Sheet

Olerup SSP[®] GelRed™ Dropper Bottle

Chemical Product and Company Identification

Date of Issue: 2010-06-30
Product Name: Olerup SSP[®] GelRed™ Dropper Bottle
Product number: 103.302-05
Manufacturer: Olerup SSP AB, Franzengatan 5, SE-112 51 Stockholm, Sweden
Tel: +46 8 717 88 27
Fax: +46 8 717 88 18

Composition and Information on Ingredients

Component	Chemical Name	Weight %	CAS No.	EINECS No.
GelRed™	Proprietary	Proprietary	N.A.	N.A.
Water	dihydrogen oxide	Proprietary	7732-18-5	231-791-2

This product is not classified according to the EU regulations

Hazards Identification

Emergency overview: Handle with care – substance has not been fully tested.

First Aid Measures

Potentially harmful; avoid prolonged or repeated exposure. Wash hands thoroughly after handling. If eye or skin contact occurs, wash affected areas with water for 15 minutes and seek medical advice. If inhaled, move individual to fresh air and seek medical advice. If swallowed, seek medical advice.

Fire Fighting Measures

Use carbon dioxide, dry chemical extinguishers, foam extinguishers or water. An approved self-contained breathing apparatus and protective clothing are recommended.

Accidental Release Measures

Use appropriate protective equipment and methods to clean up spilled substance promptly. Absorb spill onto an appropriate material. Collect and dispose of all waste in accordance with applicable laws.

Handling and Storage

Store at room temperature. The dropper bottle is specially designed to exclude damaging UV light. GelRed™ should be shielded from light and exposure to light should be avoided during long-term storage. The Olerup SSP[®] GelRed™ dropper bottle should be stored at room temperature. At low temperatures, such as +4°C, dye precipitation may occur. Should this occur, heat the solution in a water bath at +45°C to +50°C for two minutes and/or vortex the solution.

Exposure Controls and Personal Protection

Gloves, protective clothing and eyewear should be worn and safe laboratory practices followed.

ACGIH/OSHA Permissible Exposure Limit Data: Not determined

Physical and Chemical Properties

Form:	Liquid
Color:	Dark red
Odour:	Not determined
Solubility in/ miscibility with water:	High
pH in solution (20°C):	Not determined

Lot No.: **2H7**

Density (20°C):

Not determined

Stability and Reactivity**Thermal Decomposition:** No decomposition if used according to specifications**Dangerous Reactions:** No dangerous reactions identified**Dangerous Products of Decomposition:** No dangerous decomposition products identified**Toxicological Information****RTECS Number:** None known

Toxicity: GelRed was subjected to a series of tests both by Biotium and by three independent testing services to assess the dye's safety for routine handling and disposal. These tests include: 1) glove penetration test; 2) cell membrane permeability and cytotoxicity test; 3) Ames test; and 4) environmental safety tests. Test results confirm that the dye is impenetrable to both latex gloves and cell membranes. The dye is noncytotoxic and nonmutagenic at concentrations well above the working concentrations used in gel staining. Furthermore, GelRed has successfully passed environmental safety tests in compliance with CCR Title 22 Hazardous Waste Characterization. As a result, GelRed is not classified as hazardous waste, thus can be safely disposed of down the drain or as regular trash, providing convenience and reducing cost in waste disposal. However, since these tests were not performed on human, we still advise that researchers exercise precautions when handling the dye or any other DNA-binding molecules by wearing protective gears. For detailed test results on GelRed, you may download a complete safety report at www.biotium.com.

Health Hazards: See toxicity.**Potential Hazards:** We are not aware of any potential hazards**Carcinogenicity:** Not listed by NTP, IARC or OSHA**Ecological and Disposal Information**

Not classified as hazardous waste, thus can be safely disposed of down the drain or as regular trash.

Consult local, state or national regulations for proper disposal.

Transport Information**Hazards Class:** Not determined**Identification Number:** Not determined**Packing Group:** Not classified**Proper Shipping Name (Technical Name):** Not determined**Regulatory information**

This product is not classified according to the EU regulations.

U.S Federal regulation TSCA:

Not listed.

SARA:

Not regulated.

Clean Water Act (CWA):

No products were found.

Clean Air Act (CAA):

No products were found.

WHMIS (Canada):

Not controlled under WHMIS (Canada), no products found.

Other Information

The information provided above is believed to be correct to our best knowledge, but does not purport to be all inclusive, and shall be used only as a guide. This material is sold for laboratory purposes only and is not required to appear on the TSCA inventory. It is not intended for food, drug, household, agricultural or cosmetic use. Its use must be supervised by a technically qualified individual experienced in handling potentially hazardous chemicals. Olerup SSP AB shall not be held liable for any damage resulting from handling or contact with the above product.

Lot No.: 2H7**ADDRESSES:****Manufacturer:****Olerup SSP AB**, Franzengatan 5, SE-112 51 Stockholm, Sweden.**Tel:** +46-8-717 88 27**Fax:** +46-8-717 88 18**E-mail:** olerup-se@caredx.com**Web page:** <http://www.olerup.com>**Distributed by:****Olerup GmbH**, Löwengasse 47 / 6, AT-1030 Vienna, Austria.**Tel:** +43-1-710 15 00**Fax:** +43-1-710 15 00 10**E-mail:** olerup-at@caredx.com**Web page:** <http://www.olerup.com>**Olerup Inc.**, 901 S. Bolmar St., Suite R, West Chester, PA 19382**Tel:** 1-877-OLERUP1**Fax:** 610-344-7989**E-mail:** olerup-us@caredx.com**Web page:** <http://www.olerup.com>For information on *Olerup* distributors worldwide, contact **Olerup GmbH**.