

## Olerup SSP® HLA-A\*29

Product number:	101.428-12u – without <i>Taq</i> polymerase
Lot number:	05M
Expiry date:	2013-November-01
Number of tests:	12
Number of wells per test:	23
Storage - pre-aliquoted primers:	dark at -20°C
- PCR Master Mix:	-20°C
- Adhesive PCR seals	RT
- Product Insert	RT

**This Product Description is only valid for Lot No. 05M.**

### CHANGES COMPARED TO THE PREVIOUS *OLERUP SSP*® HLA-A\*29 LOT

The HLA-A\*29 specificity and interpretation tables have been updated for the HLA-A alleles described since the previous *Olerup SSP*® HLA-A\*29 lot was made (Lot No. 19K).

Seven wells have been added to the HLA-A\*29 kit, well **17 to 23**.

The primers of the wells detailed below have been exchanged, added or modified compared to the previous lot.

Well	5'-primer	3'-primer	rationale
5	Moved	Moved	Primer pair moved to well 21, exchanged positive control primer pair.
8	Moved	Moved	Primer pair moved to well 22.
12	Moved	Moved	Primer pair moved to well 21.
17	New	New	New primer pairs for the A*29:24 and A*29:27 alleles.
18	New	New	New primer pair for the A*29:23 allele.
19	New	New	Primer added for the A*29:25 allele.
20	New	New	Primer added for the A*29:26 allele.
21	New	New	Primer pair from wells 5 and 12.
22	New	New	Primer pair from well 8.
23	New	New	Primer added for the A*29:29 allele.

Change in revision R01 compared to R00:

1. Primer mix 16 does not amplify the A\*29:21 allele. This has been corrected in the specificity and amplification tables. Thus, this lot of the HLA-A\*29 subtyping kit cannot distinguish the A\*29:02:01:01-29:02:03 and 29:02:05-29:02:08 and the A\*29:21 alleles.

## PRODUCT DESCRIPTION

### HLA-A\*29 SSP subtyping

#### CONTENT

The primer set contains 5'- and 3'-primers for identifying the A\*29:01 to A\*29:29 alleles.

#### PLATE LAYOUT

Each test consists of 23 PCR reactions in a 24 well cut PCR plate. Well 24 is empty

1	2	3	4	5	6	7	8
9	10	11	12	13	14	15	16
17	18	19	20	21	22	23	empty

The 24 well cut PCR plate is marked with 'HLA-A\*29' in silver/gray ink.

Well No. 1 is marked with the Lot No. '05M'.

A faint row of numbers is seen between wells 1 and 2 or wells 7 and 8 of the PCR trays. These stem from the manufacture of the trays, and should be disregarded.

The PCR plates are covered with a PCR-compatible foil.

**Please note:** When removing each 24 well PCR plate, make sure that the remaining plates stay covered. Use a scalpel or a similar instrument to carefully cut the foil between the plates.

#### INTERPRETATION

The interpretation of HLA-A\*29 SSP subtypings will be influenced by two A\*01, five A\*02, three A\*03, the A\*11:01:28, the A\*23:03:01, five A\*24, the A\*26:22, the A\*30:47, six A\*31, the A\*32, four A\*33, the A\*66\*09, the A\*68:58 and the A\*74 alleles when present on the other haplotype.

#### UNIQUELY IDENTIFIED ALLELES

All the HLA-A\*29 alleles, i.e. **A\*29:01 to A\*29:29 alleles**, recognized by the HLA Nomenclature Committee in April 2011<sup>1</sup> will give rise to unique amplification patterns by the primers in the HLA-A\*29 subtyping kit.

The HLA-A\*29 subtyping kit cannot distinguish the A\*29:01:01:01 and 29:01:02-29:01:03 alleles or the A\*29:02:01:01-29:02:03 and 29:02:05-29:02:08 alleles.

<sup>1</sup>HLA-A alleles listed on the IMGT/HLA web page 2011-April-08, release 3.4.0, [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla).

### RESOLUTION IN HOMO- AND HETEROZYGOTES

A total of 39 alleles generate 30 amplification patterns that can be combined in 465 homozygous and heterozygous combinations. 208 of these genotypes do not give rise to unique amplification patterns. The different lengths of the specific PCR products were not considered in these calculations.

+++-----	--+---+-	-----	*29:01:01:02N, *29:10 = *29:01:01:02N, *29:14
++-----	---+---+	-----	*29:01:01:02N, *29:02:04 = *29:01:01:02N, *29:18
++-----	-----++	-----	*29:01:01:01, *29:01:01:02N = *29:01:01:02N, *29:01:01:02N
+---+---+	-----+-	-----	*29:01:01:01, *29:19 = *29:06, *29:20 = *29:19, *29:20
+---+---+	+-----+	-----	*29:05, *29:16 = *29:08N, *29:20
+---+---+	-+-----	-----	*29:05, *29:17 = *29:09, *29:20
+---+---+	-----+-	-----	*29:10, *29:20 = *29:14, *29:20
+---+---+	-----++	-----	*29:05, *29:12 = *29:11, *29:20
+---+---+	-----+++	-----	*29:05, *29:15 = *29:20, *29:21
+---+---+	-----++	+-----	*29:05, *29:24 = *29:20, *29:27
+---+---+	-----+-	-----	*29:01:01:01, *29:05 = *29:02:01:01, *29:20 = *29:05, *29:20
+---+---+	+++-----	-----	*29:08N, *29:17 = *29:09, *29:16
+---+---+	+---+---+	-----	*29:10, *29:16 = *29:14, *29:16
+---+---+	+-----+	-----	*29:08N, *29:12 = *29:11, *29:16
+---+---+	-----+++	-----	*29:08N, *29:15 = *29:16, *29:21
+---+---+	+-----+	+-----	*29:08N, *29:24 = *29:16, *29:27
+---+---+	+-----+-	-----	*29:01:01:01, *29:08N = *29:02:01:01, *29:16 = *29:08N, *29:16
+---+---+	---+---+	-----	*29:10, *29:17 = *29:14, *29:17
+---+---+	---+---+	-----	*29:09, *29:12 = *29:11, *29:17
+---+---+	---+---+	-----	*29:09, *29:15 = *29:17, *29:21
+---+---+	---+---+	+-----	*29:09, *29:24 = *29:17, *29:27
+---+---+	---+---+	-----	*29:01:01:01, *29:09 = *29:02:01:01, *29:17 = *29:09, *29:17
+---+---+	---+---+	-----	*29:10, *29:12 = *29:12, *29:14
+---+---+	---+---+	-----	*29:10, *29:15 = *29:14, *29:15
+---+---+	---+---+	+-----	*29:10, *29:24 = *29:14, *29:24
+---+---+	---+---+	-----	*29:01:01:01, *29:10 = *29:01:01:01, *29:14
+---+---+	---+---+	-----	*29:11, *29:15 = *29:12, *29:21
+---+---+	---+---+	+-----	*29:11, *29:24 = *29:12, *29:27
+---+---+	---+---+	-----	*29:01:01:01, *29:11 = *29:02:01:01, *29:12 = *29:11, *29:12
+---+---+	---+---+	+-----	*29:15, *29:27 = *29:21, *29:24
+---+---+	---+---+	-----	*29:01:01:01, *29:21 = *29:02:01:01, *29:15 = *29:15, *29:21
+---+---+	---+---+	+-----	*29:01:01:01, *29:27 = *29:02:01:01, *29:24 = *29:24, *29:27
+-----+	---+---+	-----	*29:02:04, *29:20 = *29:18, *29:20
+-----+	---+---+	-----	*29:01:01:01, *29:20 = *29:20, *29:20
+-----+	+---+---+	-----	*29:02:04, *29:16 = *29:16, *29:18
+-----+	+-----+	-----	*29:01:01:01, *29:16 = *29:16, *29:16
+-----+	---+---+	-----	*29:02:04, *29:17 = *29:17, *29:18
+-----+	---+---+	-----	*29:01:01:01, *29:17 = *29:17, *29:17
+-----+	---+---+	-----	*29:02:04, *29:12 = *29:12, *29:18
+-----+	---+---+	-----	*29:02:04, *29:15 = *29:15, *29:18
+-----+	---+---+	+-----	*29:02:04, *29:24 = *29:18, *29:24
+-----+	---+---+	-----	*29:01:01:01, *29:02:04 = *29:01:01:01, *29:18
+-----+	---+---+	-----	*29:01:01:01, *29:12 = *29:12, *29:12
+-----+	---+---+	-----	*29:01:01:01, *29:15 = *29:15, *29:15
+-----+	---+---+	+-----	*29:01:01:01, *29:24 = *29:24, *29:24
--+-----	---+---+	-----	*29:03, *29:10 = *29:03, *29:14
--+-----	---+---+	-----	*29:02:04, *29:03 = *29:03, *29:18
--+-----	---+---+	-----	*29:02:01:01, *29:03 = *29:03, *29:03
--+-----	---+---+	-----	*29:04, *29:10 = *29:04, *29:14

Lot No.: **05M**

Lot-specific information

[www.olerup-ssp.com](http://www.olerup-ssp.com)

--+--+---	-----++-	-----	*29:02:04, *29:04 = *29:04, *29:18
--+--+---	-----++-	-----	*29:02:01:01, *29:04 = *29:04, *29:04
--+--+---	-----++-	-----	*29:10, *29:19 = *29:14, *29:19
--+--+---	-----++-	-----	*29:02:04, *29:19 = *29:18, *29:19
--+--+---	-----++-	-----	*29:02:01:01, *29:19 = *29:05, *29:06 = *29:05, *29:19 = *29:06, *29:19
--+--+---	-----++-	-----	*29:05, *29:10 = *29:05, *29:14
--+--+---	-----++-	-----	*29:02:04, *29:05 = *29:05, *29:18
--+--+---	-----++-	-----	*29:02:01:01, *29:05 = *29:05, *29:05
--+--+---	-----++-	-----	*29:06, *29:10 = *29:06, *29:14
--+--+---	-----++-	-----	*29:02:04, *29:06 = *29:06, *29:18
--+--+---	-----++-	-----	*29:02:01:01, *29:06 = *29:06, *29:06
--+--+---	-----++-	-----	*29:07, *29:10 = *29:07, *29:14
--+--+---	-----++-	-----	*29:02:04, *29:07 = *29:07, *29:18
--+--+---	-----++-	-----	*29:02:01:01, *29:07 = *29:07, *29:07
--+--+---	-----++-	-----	*29:08N, *29:10 = *29:08N, *29:14
--+--+---	-----++-	-----	*29:02:04, *29:08N = *29:08N, *29:18
--+--+---	-----++-	-----	*29:02:01:01, *29:08N = *29:08N, *29:08N
--+--+---	-----++-	-----	*29:09, *29:10 = *29:09, *29:14
--+--+---	-----++-	-----	*29:02:04, *29:09 = *29:09, *29:18
--+--+---	-----++-	-----	*29:02:01:01, *29:09 = *29:09, *29:09
--+--+---	-----++-	-----	*29:02:04, *29:10 = *29:02:04, *29:14 = *29:10, *29:18
--+--+---	-----++-	-----	*29:10, *29:11 = *29:11, *29:14
--+--+---	-----++-	-----	*29:10, *29:21 = *29:14, *29:21
--+--+---	-----++-	-----	*29:10, *29:27 = *29:14, *29:27
--+--+---	-----++-	-----	*29:10, *29:23 = *29:14, *29:23
--+--+---	-----++-	-----	*29:10, *29:25 = *29:14, *29:25
--+--+---	-----++-	-----	*29:10, *29:26 = *29:14, *29:26
--+--+---	-----++-	-----	*29:10, *29:22 = *29:14, *29:22
--+--+---	-----++-	-----	*29:10, *29:13 = *29:13, *29:14
--+--+---	-----++-	-----	*29:10, *29:29 = *29:14, *29:29
--+--+---	-----++-	-----	*29:02:01:01, *29:10 = *29:02:01:01, *29:14 = *29:10, *29:10 = *29:10, *29:14
--+--+---	-----++-	-----	*29:02:04, *29:11 = *29:11, *29:18
--+--+---	-----++-	-----	*29:02:04, *29:21 = *29:18, *29:21
--+--+---	-----++-	-----	*29:02:04, *29:27 = *29:18, *29:27
--+--+---	-----++-	-----	*29:02:04, *29:23 = *29:18, *29:23
--+--+---	-----++-	-----	*29:02:04, *29:26 = *29:18, *29:26
--+--+---	-----++-	-----	*29:02:04, *29:22 = *29:18, *29:22
--+--+---	-----++-	-----	*29:02:04, *29:13 = *29:13, *29:18
--+--+---	-----++-	-----	*29:02:04, *29:29 = *29:18, *29:29
--+--+---	-----++-	-----	*29:02:01:01, *29:02:04 = *29:02:01:01, *29:18
--+--+---	-----++-	-----	*29:02:01:01, *29:11 = *29:11, *29:11
--+--+---	-----++-	-----	*29:02:01:01, *29:21 = *29:21, *29:21
--+--+---	-----++-	-----	*29:02:01:01, *29:27 = *29:27, *29:27
--+--+---	-----++-	-----	*29:02:01:01, *29:23 = *29:23, *29:23
--+--+---	-----++-	-----	*29:02:01:01, *29:26 = *29:26, *29:26
--+--+---	-----++-	-----	*29:02:01:01, *29:22 = *29:22, *29:22
--+--+---	-----++-	-----	*29:02:01:01, *29:29 = *29:29, *29:29
-----	-----++-	-----	*29:02:04, *29:25 = *29:18, *29:25
-----	-----++-	-----	*29:02:04, *29:02:04 = *29:02:04, *29:18

\*29:01:01:01 = \*29:01:01:01 and 29:01:02-29:01:03  
\*29:02:01:01 = \*29:02:01:01-29:02:03 and 29:02:05-29:02:08

## SPECIFICITY TABLE

### HLA-A\*29 SSP subtyping

Specificities and sizes of the PCR products of the 23 primer mixes used for HLA-A\*29 SSP subtyping

Primer Mix	Size of spec. PCR product <sup>1</sup>	Size of control band <sup>2</sup>	Amplified HLA-A*29 alleles	Other amplified HLA-A alleles <sup>3</sup>
<b>1</b>	480 bp	<b>800 bp</b>	*29:01:01:01-29:01:03, 29:12, 29:15-29:17, 29:20, 29:24	
<b>2</b>	130 bp	1070 bp	*29:01:01:02N	
<b>3</b>	440 bp	<b>800 bp</b>	*29:02:01:01-29:02:03, 29:02:05-29:11, 29:13-29:14, 29:19, 29:21-29:23, 29:26-29:27, 29:29	
<b>4</b>	165 bp	<b>800 bp</b>	*29:03	*23:03:01, 31:05, 32:13, 33:10
<b>5</b>	130 bp	1070 bp	*29:04	
<b>6<sup>4,6</sup></b>	105 bp, 130 bp	<b>800 bp</b>	*29:05, 29:19-29:20	*11:01:28, 31:24, 32:02, 32:22
<b>7<sup>4,5,7</sup></b>	105 bp, 210 bp	1070 bp	*29:06, 29:19	*32:12
<b>8<sup>4</sup></b>	85 bp	<b>800 bp</b>	*29:07	*24:17, 24:41
<b>9<sup>4,5,8</sup></b>	80 bp, 170 bp	1070 bp	*29:08N, 29:16	*03:27
<b>10<sup>4,9</sup></b>	95 bp, 170 bp	<b>800 bp</b>	*29:09, 29:17	*03:01:18, 11:01:28, 31:24, 32:33, 33:34
<b>11<sup>4,10</sup></b>	100 bp, 195 bp	<b>800 bp</b>	*29:10, 29:14	
<b>12<sup>4,5</sup></b>	80 bp	<b>800 bp</b>	*29:02:04, 29:18	*32:01:01-32:01:07, 32:01:09-32:03, 32:05-32:33, 32:35, 74:01-74:13, 74:15
<b>13<sup>4,5,11</sup></b>	85 bp, 165 bp	1070 bp	*29:11-29:12	*31:16
<b>14</b>	200 bp	1070 bp	*29:01:01:01-29:13, 29:15-29:17, 29:19-29:27, 29:29	
<b>15</b>	240 bp	1070 bp	*29:01:01:01-29:12, 29:14-29:18, 29:20-29:27, 29:29	*33:13

Lot No.: **05M**

Lot-specific information

[www.olerup-ssp.com](http://www.olerup-ssp.com)

<b>16<sup>4</sup></b>	95 bp	1070 bp	*29:15	*02:221
<b>17<sup>12</sup></b>	140 bp, 190 bp	1070 bp	*29:24, 29:27	
<b>18</b>	185 bp	1070 bp	*29:23	
<b>19</b>	160 bp	1070 bp	*29:25	
<b>20</b>	505 bp	1070 bp	*29:26	
<b>21</b>	185 bp	1070 bp	*29:22	*01:20, 01:66, 02:19, 02:44, 02:149, 02:309, 03:95, 24:14, 24:93, 26:22, 30:47, 33:22, 66:09
<b>22</b>	260 bp	1070 bp	*29:13	*24:82, 31:07-31:08, 31:10
<b>23<sup>4</sup></b>	115 bp	<b>800 bp</b>	*29:29	*68:58

<sup>1</sup>Alleles are assigned by the presence of specific PCR product(s). However, the sizes of the specific PCR products may be helpful in the interpretation of HLA-A\*29 SSP typings.

When the primers in a primer mix can give rise to specific PCR products of more than one length this is indicated if the size difference is 20 base pairs or more. Size differences shorter than 20 base pairs are not given. For high resolution SSP kits the respective lengths of the specific PCR product(s) of the alleles amplified by these primer mixes are given.

Nonspecific amplifications, i.e. a ladder or a smear of bands, may sometimes be seen. GC-rich primers have a higher tendency of giving rise to nonspecific amplifications than other primers.

PCR fragments longer than the control bands may sometimes be observed. Such bands should be disregarded and do not influence the interpretation of the SSP typings.

PCR fragments migrating faster than the control bands, but slower than a 400 bp fragment may be seen in some gel read-outs. Such bands can be disregarded and do not influence the interpretation of the SSP typings.

Some primers may give rise to primer oligomer artifacts. Sometimes this phenomenon is an inherent feature of the primer pair(s) of a primer mix. More often it is due to other factors such as too low amount of DNA in the PCR reactions, taking too long time in setting up the PCR reactions, working at elevated room temperature or using thermal cyclers that are not pre-heated.

<sup>2</sup>The internal positive control primer pairs amplify segments of the human growth hormone gene. The two different control primer pairs give rise to either an internal positive control band of 1070 base pairs, for most wells, or a band of 800 base pairs, for some wells.

Well number 1 contains the primer pair giving rise to the shorter, 800 bp, internal positive control band in order to help in the correct orientation of the HLA-A\*29 subtyping.

In addition, wells number 3, 4, 6, 8, 10 to 12 and 23 contain the primer pair giving rise to the shorter, 800 bp, internal positive control band in order to allow kit identification.

In the presence of a specific amplification the intensity of the control band often decreases.

<sup>3</sup>Due to the sharing of sequence motifs between HLA-A alleles a few non-HLA-A\*29 alleles will be amplified by primer mixes 4, 6 to 10, 12, 13, 15, 16 and 21 to 23.

<sup>4</sup>Short specific PCR fragments are less intense and not as sharp as longer specific bands.

<sup>5</sup>Primer mixes 7, 9, 12 and 13 may give rise to nonspecific amplifications.

<sup>6</sup>Primer mix 6: Specific PCR fragment of 105 bp in the A\*29:19 and 29:20 alleles. Specific PCR fragment of 130 bp in the A\*29:05 and the A\*11:01:28, 31:24, 32:02 and 32:22 alleles.

<sup>7</sup>Primer mix 7: Specific PCR fragment of 105 bp in the A\*29:19 allele. Specific PCR fragment of 210 bp in the A\*29:06 and the A\*32:12 alleles.

<sup>8</sup>Primer mix 9: Specific PCR fragment of 80 bp in the A\*29:16 and in the A\*03:27 alleles. Specific PCR fragment of 170 bp in the A\*29:08N allele.

Lot No.: **05M**

Lot-specific information

[www.olerup-ssp.com](http://www.olerup-ssp.com)

<sup>9</sup>Primer mix 10: Specific PCR fragment of 95 bp in the A\*29:09 and the A\*03:01:18, 11:01:28, 31:24, 32:33 and 33:34 alleles. Specific PCR fragment of 170 bp in the A\*29:17 allele.

<sup>10</sup>Primer mix 11: Specific PCR fragment of 100 bp in the A\*29:14 allele. Specific PCR fragment of 195 bp in the A\*29:10 allele.

<sup>11</sup>Primer mix 13: Specific PCR fragment of 85 bp in the A\*29:11 allele. Specific PCR fragment of 165 bp in the A\*29:12 and in the A\*31:16 alleles.

<sup>12</sup>Primer mix 17: Specific PCR fragment of 140 bp in the A\*29:24 allele. Specific PCR fragment of 190 bp in the A\*29:27 allele.

<b>INTERPRETATION TABLE</b>												
<b>HLA-A*29 SSP subtyping</b>												
<b>Amplification patterns of the A*29:01 to 29:29 alleles</b>												
	<b>Well<sup>5</sup></b>											
	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>	<b>9</b>	<b>10</b>	<b>11</b>	<b>12</b>
<b>Length of spec. PCR product(s)</b>	480	130	440	165	130	105	105	85	80	95	100	80
<b>Length of int. pos. control<sup>1</sup></b>	800	1070	800	800	1070	800	1070	800	1070	800	800	800
<b>5'-primer(s)<sup>2</sup></b>	180	808	219	448	180	219	219	368	97	130	180	180
	5'-TTT <sup>3'</sup>	5'-CgT <sup>3'</sup>	5'-gCA <sup>3'</sup>	5'-CCT <sup>3'</sup>	5'-TTT <sup>3'</sup>	5'-gCA <sup>3'</sup>	5'-gCA <sup>3'</sup>	5'-gTT <sup>3'</sup>	5'-TCA <sup>3'</sup>	5'-AgA <sup>3'</sup>	5'-TTT <sup>3'</sup>	5'-TTT <sup>3'</sup>
						448	448		413	448	448	
						5'-CCT <sup>3'</sup>	5'-CCT <sup>3'</sup>		5'-CCg <sup>3'</sup>	5'-CCT <sup>3'</sup>	5'-CCT <sup>3'</sup>	
<b>3'-primer(s)<sup>3</sup></b>	376	895	376	570	268	282	282	413	224	257	238	218
	5'-gTg <sup>3'</sup>	5'-CTC <sup>3'</sup>	5'-gTC <sup>3'</sup>	5'-CCg <sup>3'</sup>	5'-ATg <sup>3'</sup>	5'-gAg <sup>3'</sup>	5'-gAg <sup>3'</sup>	5'-gCC <sup>3'</sup>	5'-TCT <sup>3'</sup>	5'-gCA <sup>3'</sup>	5'-CCT <sup>3'</sup>	5'-gCg <sup>3'</sup>
						282	616		454	502	601	
						5'-gAg <sup>3'</sup>	5'-CgC <sup>3'</sup>		5'-CTg <sup>3'</sup>	5'-CTT <sup>3'</sup>	5'-CTT <sup>3'</sup>	
						539						
						5'-TCT <sup>3'</sup>						
<b>Well No.</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>	<b>9</b>	<b>10</b>	<b>11</b>	<b>12</b>
<b>HLA-A allele</b>												
<b>*29:01:01:01, 29:01:02-29:01:03</b>	1											
<b>*29:01:01:02N</b>	1	2										
<b>*29:02:01:01-29:02:03, 29:02:05-29:02:08, *29:21<sup>4</sup></b>			3									
<b>*29:02:04</b>												12
<b>*29:03</b>			3	4								
<b>*29:04</b>			3		5							
<b>*29:05</b>			3			6						
<b>*29:06</b>			3				7					
<b>*29:07</b>			3					8				
<b>*29:08N</b>			3						9			
<b>*29:09</b>			3							10		
<b>Well No.</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>	<b>9</b>	<b>10</b>	<b>11</b>	<b>12</b>



INTERPRETATION TABLE											
HLA-A*29 SSP subtyping											
Amplification patterns of the A*29:01 to 29:29 alleles											
Well <sup>5</sup>											
13	14	15	16	17	18	19	20	21	22	23	
85	200	240	95	140	185	160	505	185	260	115	Length of spec. PCR product(s)
165				190							
1070	1070	1070	1070	1070	1070	1070	1070	1070	1070	800	Length of int. pos. control <sup>1</sup>
97	98	98	134	448	180	98	3 <sup>rd</sup> I	355	98	355	5'-primer(s) <sup>2</sup>
5'-TCA 3'	5'-CAC 3'	5'-CAC 3'	5'-CCT 3'	5'-CCT 3'	5'-TTT 3'	5'-CAC 3'	5'-ATA 3'	5'-CCg 3'	5'-CAC 3'	5'-CCA 3'	
448			484								
5'-CCT 3'			5'-ACg 3'								
221	257	299	257	545	326	217	667	497	317	430	3'-primer(s) <sup>3</sup>
5'-ACA 3'	5'-gCA 3'	5'-TCg 3'	5'-gCA 3'	5'-AgC 3'	5'-TgA 3'	5'-TgA 3'	5'-ggT 3'	5'-TgA 3'	5'-gga 3'	5'-gCT 3'	
494			538	595							
5'-TCg 3'			5'-CAA 3'	5'-CCA 3'							
13	14	15	16	17	18	19	20	21	22	23	Well No.
											HLA-A allele
	14	15									*29:01:01:01, 29:01:02- 29:01:03
	14	15									*29:01:01:02N
	14	15									*29:02:01:01-29:02:03, 29:02:05-29:02:08, *29:21 <sup>4</sup>
	14	15									*29:02:04
	14	15									*29:03
	14	15									*29:04
	14	15									*29:05
	14	15									*29:06
	14	15									*29:07
	14	15									*29:08N
	14	15									*29:09
13	14	15	16	17	18	19	20	21	22	23	Well No.

Lot No.: **05M**

Lot-specific information

[www.olerup-ssp.com](http://www.olerup-ssp.com)

Length of spec.	480	130	440	165	130	105	105	85	80	95	100	80
PCR product(s)						130	210		170	170	195	
Well No.	1	2	3	4	5	6	7	8	9	10	11	12
*29:10			3								11	
*29:11			3									
*29:12	1											
*29:13			3									
*29:14			3								11	
*29:15	1											
*29:16	1								9			
*29:17	1									10		
*29:18												12
*29:19			3			6	7					
*29:20	1					6						
*29:22			3									
*29:23			3									
*29:24	1											
*29:25												
*29:26			3									
*29:27			3									
*29:29			3									
*01:20, 01:66, 02:19, 02:44, 02:149, 02:309, 03:95, 24:14, 24:93, 26:22, 30:47, 33:22, 66:09												
*02:221												
*03:01:18, 33:34										10		
*03:27									9			
*11:01:28, 31:24						6				10		
*23:03:01, 31:05, 33:10				4								
*24:17, 24:41								8				
*24:82, 31:07-31:08, 31:10												
*31:16												
*32:01:01-32:01:07, 32:01:09-32:01:10, 32:03, 32:05-32:11Q, 32:14-32:21, 32:23-32:32, 32:35, 74:01- 74:13, 74:15												12
*32:02, 32:22						6						12
*32:12							7					12
Well No.	1	2	3	4	5	6	7	8	9	10	11	12

Lot No.: **05M**

Lot-specific information

www.olerup-ssp.com

85	200	240	95	140	185	160	505	185	260	115	Length of spec. PCR product(s)
165				190							Well No.
13	14	15	16	17	18	19	20	21	22	23	
	14	15									*29:10
13	14	15									*29:11
13	14	15									*29:12
	14								22		*29:13
		15									*29:14
	14	15	16								*29:15
	14	15									*29:16
	14	15									*29:17
		15									*29:18
	14										*29:19
	14	15									*29:20
	14	15						21			*29:22
	14	15			18						*29:23
	14	15		17							*29:24
	14	15				19					*29:25
	14	15					20				*29:26
	14	15		17							*29:27
	14	15								23	*29:29
								21			*01:20, 01:66, 02:19, 02:44, 02:149, 02:309, 03:95, 24:14, 24:93, 26:22, 30:47, 33:22, 66:09
			16								*02:221
											*03:01:18, 33:34
											*03:27
											*11:01:28, 31:24
											*23:03:01, 31:05, 33:10
											*24:17, 24:41
								22			*24:82, 31:07-31:08, 31:10
13											*31:16
											*32:01:01-32:01:07, 32:01:09-32:01:10, 32:03, 32:05-32:11Q, 32:14-32:21, 32:23-32:32, 32:35, 74:01- 74:13, 74:15
											*32:02, 32:22
											*32:12
13	14	15	16	17	18	19	20	21	22	23	Well No.

Lot No.: **05M**

Lot-specific information

[www.olerup-ssp.com](http://www.olerup-ssp.com)

Length of spec.	480	130	440	165	130	105	105	85	80	95	100	80
PCR product(s)						130	210		170	170	195	
Well No.	1	2	3	4	5	6	7	8	9	10	11	12
*32:13				4								12
*32:33										10		12
*33:13												
*68:58												
HLA-A allele												
Well No.	1	2	3	4	5	6	7	8	9	10	11	12

<sup>1</sup>The internal positive control primer pairs amplify segments of the human growth hormone gene. The two different control primer pairs give rise to either an internal positive control band of 1070 base pairs, for most wells, or a band of 800 base pairs, for some wells.

Well number 1 contains the primer pair giving rise to the shorter, 800 bp, internal positive control band in order to help in the correct orientation of the HLA-A\*29 subtyping. .

In addition, wells number 3, 4, 6, 8, 10 to 12 and 23 contain the primer pair giving rise to the shorter, 800 bp, internal positive control band in order to allow kit identification.

<sup>2</sup>The nucleotide position, in the 2<sup>nd</sup>, 3<sup>rd</sup> or 4<sup>th</sup> exons or in the 3<sup>rd</sup> intron, matching the specificity-determining 3'-end of the primer is given. Nucleotide numbering as on the [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla) web site. The sequence of the 3 terminal nucleotides of the primer is given.

<sup>3</sup>The nucleotide position, in the 2<sup>nd</sup>, 3<sup>rd</sup> or 4<sup>th</sup> exons, matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide numbering as on the [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla) web site. The sequence of the 3 terminal nucleotides of the primer is given.

Lot No.: **05M**

Lot-specific information

www.olerup-ssp.com

85	200	240	95	140	185	160	505	185	260	115	Length of spec. PCR product(s)
13	14	15	16	17	18	19	20	21	22	23	Well No.
											*32:13
											*32:33
		15									*33:13
										23	*68:58
											HLA-A allele
13	14	15	16	17	18	19	20	21	22	23	Well No.

<sup>4</sup>This lot of the HLA-A\*29 subtyping kit cannot distinguish the A\*29:02:01:01-29:02:03 and 29:02:05-29:02:08 and the A\*29:21 alleles.

<sup>5</sup>Primer mix 6: Specific PCR fragment of 105 bp in the A\*29:19 and 29:20 alleles. Specific PCR fragment of 130 bp in the A\*29:05 and the A\*11:01:28, 31:24, 32:02 and 32:22 alleles.

Primer mix 7: Specific PCR fragment of 105 bp in the A\*29:19 allele. Specific PCR fragment of 210 bp in the A\*29:06 and the A\*32:12 alleles.

Primer mix 9: Specific PCR fragment of 80 bp in the A\*29:16 and in the A\*03:27 alleles. Specific PCR fragment of 170 bp in the A\*29:08N allele.

Primer mix 10: Specific PCR fragment of 95 bp in the A\*29:09 and the A\*03:01:18, 11:01:28, 31:24, 32:33 and 33:34 alleles. Specific PCR fragment of 170 bp in the A\*29:17 allele.

Primer mix 11: Specific PCR fragment of 100 bp in the A\*29:14 allele. Specific PCR fragment of 195 bp in the A\*29:10 allele.

Primer mix 13: Specific PCR fragment of 85 bp in the A\*29:11 allele. Specific PCR fragment of 165 bp in the A\*29:12 and in the A\*31:16 alleles.

Primer mix 17: Specific PCR fragment of 140 bp in the A\*29:24 allele. Specific PCR fragment of 190 bp in the A\*29:27 allele.

CELL LINE VALIDATION SHEET																			
HLA-A*29 SSP subtyping kit																			
				Well															
				1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
				Lot No.:															
				200958401	200958402	200958403	200958404	201184705	201071406	200958407	201184708	200958409	200958410	201071411	201184712	200958413	201071414	201071415	201071416
	IHWC cell line	A*	A*																
1	9001 SA	*24:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
2	9280 LK707	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
3	9011 E4181324	*01:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
4	9275 GU373	*30:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
5	9009 KAS011	*01:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
6	9353 SM	*02:01	*26:03	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
7	9020 QBL	*26:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
8	9007 DEM	*31:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
9	9026 YAR	*26:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
10	9107 LKT3	*24:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
11	9051 PITOUT	*29:02		-	-	+	-	-	-	-	-	-	-	-	-	-	+	+	-
12	9052 DBB	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
13	9004 JESTHOM	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
14	9071 OLGA	*31:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
15	9075 DKB	*24:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
16	9037 SWEIG007	*29:02		-	-	+	-	-	-	-	-	-	-	-	-	-	+	+	-
17	9282 CTM3953540	*03:01	*80:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
18	9257 32367	*33:03	*74:01	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-
19	9038 BM16	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
20	9059 SLE005	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
21	9064 AMALA	*02:17		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
22	9056 KOSE	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
23	9124 IHL	*02:01	*34:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
24	9035 JBUSH	*32:01		-	-	-	-	-	-	-	-	-	-	+	-	-	-	-	-
25	9049 IBW9	*33:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
26	9285 WT49	*02:05		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
27	9191 CH1007	*24:10	*29:01	+	-	-	-	-	-	-	-	-	-	-	-	-	+	+	-
28	9320 BEL5GB	*02:01	*29:02	-	-	+	-	-	-	-	-	-	-	-	-	-	+	+	-
29	9050 MOU	*29:02		-	-	+	-	-	-	-	-	-	-	-	-	-	+	+	-
30	9021 RSH	*30:01	*68:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
31	9019 DUCAF	*30:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
32	9297 HAG	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
33	9098 MT14B	*31:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
34	9104 DHIF	*31:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
35	9302 SSTO	*32:01		-	-	-	-	-	-	-	-	-	-	-	+	-	-	-	-
36	9024 KT17	*02:06	*11:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
37	9065 HHKB	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
38	9099 LZL	*02:17		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
39	9315 CML	*01:01	*03:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
40	9134 WHONP199	*02:07	*30:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
41	9055 H0301	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
42	9066 TAB089	*02:07		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
43	9076 T7526	*02:06	*02:07	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
44	9057 TEM	*66:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
45	9239 SHJO	*23:01	*24:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
46	9013 SCHU	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
47	9045 TUBO	*02:16	*03:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
48	9303 TER-ND	*02:01	*11:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-

Lot No.: **05M**

Lot-specific information

www.olerup-ssp.com

CELL LINE VALIDATION SHEET										
HLA-A*29 SSP subtyping kit										
				Well						
				17	18	19	20	21	22	23
				201184717	201184718	201184719	201184720	201184721	201184722	201184723
	IHWC cell line	A*	A*	Lot No.:						
1	9001 SA	*24:02			-	-	-	-	-	-
2	9280 LK707	*02:01			-	-	-	-	-	-
3	9011 E4181324	*01:01			-	-	-	-	-	-
4	9275 GU373	*30:01			-	-	-	-	-	-
5	9009 KAS011	*01:01			-	-	-	-	-	-
6	9353 SM	*02:01	*26:03		-	-	-	-	-	-
7	9020 QBL	*26:01			-	-	-	-	-	-
8	9007 DEM	*31:01			-	-	-	-	-	-
9	9026 YAR	*26:01			-	-	-	-	-	-
10	9107 LKT3	*24:02			-	-	-	-	-	-
11	9051 PITOUT	*29:02			-	-	-	-	-	-
12	9052 DBB	*02:01			-	-	-	-	-	-
13	9004 JESTHOM	*02:01			-	-	-	-	-	-
14	9071 OLGA	*31:01			-	-	-	-	-	-
15	9075 DKB	*24:02			-	-	-	-	-	-
16	9037 SWEIG007	*29:02			-	-	-	-	-	-
17	9282 CTM3953540	*03:01	*80:01		-	-	-	-	-	-
18	9257 32367	*33:03	*74:01		-	-	-	-	-	-
19	9038 BM16	*02:01			-	-	-	-	-	-
20	9059 SLE005	*02:01			-	-	-	-	-	-
21	9064 AMALA	*02:17			-	-	-	-	-	-
22	9056 KOSE	*02:01			-	-	-	-	-	-
23	9124 IHL	*02:01	*34:01		-	-	-	-	-	-
24	9035 JBUSH	*32:01			-	-	-	-	-	-
25	9049 IBW9	*33:01			-	-	-	-	-	-
26	9285 WT49	*02:05			-	-	-	-	-	-
27	9191 CH1007	*24:10	*29:01		-	-	-	-	-	-
28	9320 BEL5GB	*02:01	*29:02		-	-	-	-	-	-
29	9050 MOU	*29:02			-	-	-	-	-	-
30	9021 RSH	*30:01	*68:02		-	-	-	-	-	-
31	9019 DUCAF	*30:02			-	-	-	-	-	-
32	9297 HAG	*02:01			-	-	-	-	-	-
33	9098 MT14B	*31:01			-	-	-	-	-	-
34	9104 DHIF	*31:01			-	-	-	-	-	-
35	9302 SSTO	*32:01			-	-	-	-	-	-
36	9024 KT17	*02:06	*11:01		-	-	-	-	-	-
37	9065 HHKB	*03:01			-	-	-	-	-	-
38	9099 LZL	*02:17			-	-	-	-	-	-
39	9315 CML	*01:01	*03:01		-	-	-	-	-	-
40	9134 WHONP199	*02:07	*30:01		-	-	-	-	-	-
41	9055 H0301	*03:01			-	-	-	-	-	-
42	9066 TAB089	*02:07			-	-	-	-	-	-
43	9076 T7526	*02:06	*02:07		-	-	-	-	-	-
44	9057 TEM	*66:01			-	-	-	-	-	-
45	9239 SHJO	*23:01	*24:02		-	-	-	-	-	-
46	9013 SCHU	*03:01			-	-	-	-	-	-
47	9045 TUBO	*02:16	*03:01		-	-	-	-	-	-
48	9303 TER-ND	*02:01	*11:01		-	-	-	-	-	-

## CERTIFICATE OF ANALYSIS

### Olerup SSP® HLA-A\*29 SSP

Product number: 101.428-12u – without *Taq* polymerase  
Lot number: 05M  
Expiry date: 2013-November-01  
Number of tests: 12  
Number of wells per test: 23

#### Well specifications:

Well No.	Production No.	Well No.	Production No.	Well No.	Production No.
1	2009-584-01	9	2009-584-09	17	2011-847-17
2	2009-584-02	10	2009-584-10	18	2011-847-18
3	2009-584-03	11	2010-714-11	19	2011-847-19
4	2009-584-04	12	2011-847-12	20	2011-847-20
5	2011-847-05	13	2009-584-13	21	2011-847-21
6	2010-714-06	14	2010-714-14	22	2011-847-22
7	2009-584-07	15	2010-714-15	23	2011-847-23
8	2011-847-08	16	2010-714-16		

The specificity of each primer solution of the kit has been tested against 48 well characterized IHWC cell line DNAs.

No DNAs carrying the alleles to be amplified by primer solutions 2, 4 to 7, 9 to 11, 13 and 16 to 23 were available. The specificities of the primers in primer solutions 4, 6, 7, 10, 11, 21 and 22 were tested by separately adding one additional 5'-primer, respectively one additional 3'-primer. In primer solutions 2, 5, 9, 13, 17 to 20 and 23 it was only possible to test the 5'-primer, the 3'-primer was not possible to test. In primer solution 16 it was only possible to the 3'-primer, the 5'-primer was not possible to test. In primer solution 6 and 7 one of the 3'-primers was not possible to test, and in primer solution 10 one of the 5'-primers was not possible to test.

**Results:** No false positive or false negative amplifications were obtained.

**Date of approval:** 2011-June-29

**Approved by:**

#### Production Quality Control



Lot No.: **05M**

Lot-specific information

[www.olerup-ssp.com](http://www.olerup-ssp.com)

## Declaration of Conformity

**Product name:** *Olerup* SSP® HLA-A\*29  
**Product number:** 101.428-12u  
**Lot number:** 05M

**Intended use:** HLA-A\*29 high resolution histocompatibility testing

**Manufacturer:** *Olerup* SSP AB  
Franzengatan 5  
SE-112 51 Stockholm, Sweden  
**Phone:** +46-8-717 88 27  
**Fax:** +46-8-717 88 18

We, *Olerup* SSP AB, hereby declare that this product, to which this Declaration of Conformity relates is in conformity with the following Standard(s) and other normative document(s) ISO 9001:2008 and ISO 13485:2003, following the provisions of the 98/79/EC Directive on *in vitro* diagnostic medical devices, Annex II List B, conformity assessed using Annex IV, as transposed into the national laws of the Member States of the European Union.

The Technical Documentation File is maintained at *Olerup* SSP AB, Franzengatan 5, SE-112 51 Stockholm, Sweden.

The Authorized Representative located within the Community is: *Olerup* SSP AB.

Notified Body: Lloyd’s Register Quality Assurance Limited, Hiramford, Middlemarch Office Village, Siskin Drive, Coventry CV3 4FJ, United Kingdom. (Notified Body number: 0088.)

Stockholm, Sweden  
2012-May-25

Ann-Cathrin Jareman  
Head of QA and Regulatory Affairs





Lot No.: **05M**

Lot-specific information

[www.olerup-ssp.com](http://www.olerup-ssp.com)

**ADDRESSES:**

**Manufacturer:**

**Olerup SSP AB**, Franzengatan 5, SE-112 51 Stockholm, Sweden.

**Tel:** +46-8-717 88 27

**Fax:** +46-8-717 88 18

**E-mail:** [info-ssp@olerup.com](mailto:info-ssp@olerup.com)

**Web page:** <http://www.olerup-ssp.com>

**Distributed by:**

**Olerup GmbH**, Löwengasse 47 / 6, AT-1030 Vienna, Austria.

**Tel:** +43-1-710 15 00

**Fax:** +43-1-710 15 00 10

**E-mail:** [support-at@olerup.com](mailto:support-at@olerup.com)

**Web page:** <http://www.olerup.com>

**Olerup Inc.**, 901 S. Bolmar St., Suite R, West Chester, PA 19382

**Tel:** 1-877-OLERUP1

**Fax:** 610-344-7989

**E-mail:** [info.us@olerup.com](mailto:info.us@olerup.com)

**Web page:** <http://www.olerup.com>

For information on *Olerup* SSP distributors worldwide, contact **Olerup GmbH**.