

101.429-12 – including *Taq* polymerase, IFU-01
101.429-12u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **38V**

Lot-specific Information
Olerup SSP® HLA-A*30

www.olerup-ssp.com

Product number:	101.429-12 – including <i>Taq</i> polymerase 101.429-12u – without <i>Taq</i> polymerase
Lot number:	38V
Expiry date:	2016-August-01
Number of tests:	12
Number of wells per test:	31+1
Storage - pre-aliquoted primers:	dark at -20°C
- PCR Master Mix:	-20°C
- Adhesive PCR seals	RT
- Product Insert	RT

This Product Description is only valid for Lot No. 38V.

**CHANGES COMPARED TO THE PREVIOUS OLERUP SSP®
HLA-A*30 LOT (21R)**

The HLA-A*30 kit is updated for new alleles to enable separation of:

- Confirmed¹ alleles as listed in the IMGT/HLA database
- Polymorphisms in exons outside of the region encoding the peptide binding domain
- Null and Alternatively expressed alleles

A well containing Negative Control primer pairs has been added.

The format of the Product Insert and Worksheet have been changed.

Two wells have been added to the HLA-A*30 kit, wells **31 to 32**.

¹As described in section Uniquely Identified Alleles.

The HLA-A*30 specificity and interpretation tables have been updated for the HLA-A alleles described since the previous *Olerup SSP®* HLA-A*30 lot was made (**Lot No. 21R**).

101.429-12 – including *Taq* polymerase, IFU-01
101.429-12u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **38V**

Lot-specific Information

www.olerup-ssp.com

The primers of the wells detailed below have been exchanged, added or modified compared to the previous lot.

Well	5'-primer	3'-primer	rationale
18	-	Added	3'-primer added for the A*30:71 allele.
19	Added	Added	Primer pair added for the A*30:77 allele, 3'-primer added for the A*30:70N allele.
22	-	Added	3'-primer added for the A*30:73N allele.
24	-	Added	3'-primer added for the A*30:78N allele.
25	Modified	-	Modified 5'-primer for improved HLA-specific amplification.
27	Added	Added	Primer pair added for the A*30:56 allele.
29	-	Added	3'-primer added for the A*30:73N allele.
31	New	New	New primer pair for the A*30:76N allele.
32			Negative Control.

Change in revision R01 compared to R00:

1. Primer mix 24 does not amplify the A*30:78N allele. This has been corrected in the Specificity and Interpretation tables. Thus, this lot of the HLA-A*30 subtyping kit cannot distinguish the A*30:01:01-30:01:07 alleles and the A*30:78N allele.

101.429-12 – including *Taq* polymerase, IFU-01
101.429-12u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **38V**

Lot-specific Information

www.olerup-ssp.com

Well **32** contains Negative Control primer pairs, that will amplify more than 95% of the *Olerup SSP®* HLA Class I, DRB, DQB1 and DPB1 amplicons as well as amplicons generated by a control primer pair.

PCR product sizes range from 75 to 430 base pairs.
The PCR product generated by the control primer pair is 430 base pairs.

Length of PCR product	105	200	105	80	75	80
5'-primer¹	164	340	440	45	45	43
	5'-CAC ^{3'}	5'-Agg ^{3'}	5'-TTA ^{3'}	5'-Tgg ^{3'}	5'-Tgg ^{3'}	5'-Tgg ^{3'}
3'-primer²	231	2nd I	507	59	58	57
	5'-TgC ^{3'}	5'-AAA ^{3'}	5'-TTg ^{3'}	5'-CTC ^{3'}	5'-ggC ^{3'}	5'-CTC ^{3'}
A*	+	+	+			
B*	+	+	+			
C*	+	+	+			
DRB1				+	+	
DRB3				+	+	
DRB5				+		
DQB1					+	
DPB1						+

¹The nucleotide position for HLA class I genes and the codon for HLA class II genes, in the 2nd or 3rd exon, matching the specificity-determining 3'-end of the primer is given. Nucleotide and codon numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

²The nucleotide position for HLA class I genes and the codon for HLA class II genes, in the 2nd or 3rd exon or the 2nd intron, matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide and codon numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

101.429-12 – including *Taq* polymerase, IFU-01
101.429-12u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **38V**

Lot-specific Information

www.olerup-ssp.com

PRODUCT DESCRIPTION

HLA-A*30 SSP subtyping

CONTENT

The primer set contains 5'- and 3'-primers for identifying the A*30:01 to A*30:78N alleles.

PLATE LAYOUT

Each test consists of 32 PCR reactions in a 32 well cut PCR plate

1	2	3	4	5	6	7	8
9	10	11	12	13	14	15	16
17	18	19	20	21	22	23	24
25	26	27	28	29	30	31	NC

The 32 well cut PCR plate is marked with 'HLA-A*30' in silver/gray ink.

Well No. 1 is marked with '38V'.

Wells 1 to 31 – HLA-A*30 high resolution primers.

Well 32 – Negative Control (NC).

A faint row of numbers is seen between wells 1 and 2 or wells 7 and 8 of the PCR trays. These stem from the manufacture of the trays, and should be disregarded.

The PCR plates are covered with a PCR-compatible foil.

Please note: When removing each 32 well PCR plate, make sure that the remaining plates stay covered. Use a scalpel or a similar instrument to carefully cut the foil between the plates.

INTERPRETATION

Due to the sharing of sequence motifs between HLA-A alleles non-HLA-A*30 alleles will be amplified by primer mixes 1, 3 to 5, 7, 9 to 12, 17, 19, 20, 28 and 29.

For further details see Specificity Table.

UNIQUELY IDENTIFIED ALLELES

All the HLA-A*30 alleles, i.e. **A*30:01 to A*30:78N alleles**, recognized by the HLA Nomenclature Committee in October 2013^{1,2} will be amplified by the primers in the HLA-A*30 subtyping kit.

The HLA-A*30 kit enables separation of the confirmed HLA-A*30 alleles as listed in the IMGT/HLA database. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources. Current allele confirmation status for HLA-A*30 alleles is listed below.

101.429-12 – including *Taq* polymerase, IFU-01
101.429-12u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **38V**

Lot-specific Information

www.olerup-ssp.com

The HLA-A*30 kit also enables identification of polymorphisms in exons outside of the region encoding the peptide binding domain and of null and alternatively expressed alleles

The HLA-A*30 kit cannot distinguish silent mutations in the A*30:01:01- A*30:01:07 alleles, the 30:02:01-30:02:04 and 30:02:06-30:02:12, the A*30:04:01-30:04:02 or the A*30:11:01-30:11:02 alleles.

¹HLA-A alleles listed on the IMGT/HLA web page 2013-October-11, release 3.14.0, www.ebi.ac.uk/imgt/hla.

²Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

ALLELE CONFIRMATION STATUS

Allele	Status ¹	Allele	Status ¹	Allele	Status ¹	Allele	Status ¹
A*30:01:01	Confirmed	A*30:13	Unconfirmed	A*30:44	Unconfirmed	A*30:74	Unconfirmed
A*30:01:02	Unconfirmed	A*30:14L	Unconfirmed	A*30:45	Confirmed	A*30:75	Unconfirmed
A*30:01:03	Unconfirmed	A*30:15	Unconfirmed	A*30:46	Unconfirmed	A*30:76N	Unconfirmed
A*30:01:04	Confirmed	A*30:16	Confirmed	A*30:47	Unconfirmed	A*30:77	Unconfirmed
A*30:01:05	Unconfirmed	A*30:17	Confirmed	A*30:48	Unconfirmed	A*30:78N	Unconfirmed
A*30:01:06	Unconfirmed	A*30:18	Unconfirmed	A*30:49	Unconfirmed		
A*30:01:07	Unconfirmed	A*30:19	Unconfirmed	A*30:50	Unconfirmed		
A*30:02:01	Confirmed	A*30:20	Confirmed	A*30:51	Confirmed		
A*30:02:02	Unconfirmed	A*30:22	Unconfirmed	A*30:52	Unconfirmed		
A*30:02:03	Confirmed	A*30:23	Unconfirmed	A*30:53	Unconfirmed		
A*30:02:04	Confirmed	A*30:24	Unconfirmed	A*30:54	Unconfirmed		
A*30:02:05	Confirmed	A*30:25	Unconfirmed	A*30:55	Unconfirmed		
A*30:02:06	Confirmed	A*30:26	Unconfirmed	A*30:56	Confirmed		
A*30:02:07	Unconfirmed	A*30:27N	Unconfirmed	A*30:57	Unconfirmed		
A*30:02:08	Unconfirmed	A*30:28	Unconfirmed	A*30:58	Unconfirmed		
A*30:02:09	Unconfirmed	A*30:29	Confirmed	A*30:59N	Confirmed		
A*30:02:10	Unconfirmed	A*30:30	Confirmed	A*30:60	Unconfirmed		
A*30:02:11	Unconfirmed	A*30:31	Confirmed	A*30:61	Confirmed		
A*30:02:12	Unconfirmed	A*30:32	Confirmed	A*30:62	Unconfirmed		
A*30:03	Confirmed	A*30:33	Unconfirmed	A*30:63	Unconfirmed		
A*30:04:01	Confirmed	A*30:34	Unconfirmed	A*30:64	Unconfirmed		
A*30:04:02	Confirmed	A*30:35	Unconfirmed	A*30:65	Unconfirmed		
A*30:06	Unconfirmed	A*30:36	Confirmed	A*30:66	Unconfirmed		
A*30:07	Confirmed	A*30:37	Unconfirmed	A*30:67	Unconfirmed		
A*30:08	Confirmed	A*30:38	Unconfirmed	A*30:68	Unconfirmed		
A*30:09	Confirmed	A*30:39	Unconfirmed	A*30:69	Unconfirmed		
A*30:10	Confirmed	A*30:40	Confirmed	A*30:70N	Unconfirmed		
A*30:11:01	Confirmed	A*30:41	Unconfirmed	A*30:71	Confirmed		
A*30:11:02	Unconfirmed	A*30:42	Confirmed	A*30:72	Unconfirmed		
A*30:12	Unconfirmed	A*30:43	Unconfirmed	A*30:73N	Unconfirmed		

¹Allele status “confirmed” or “unconfirmed” as listed on the IMGT/HLA web page 2013-October-11, release 3.14.0, www.ebi.ac.uk/imgt/hla.

RESOLUTION IN HOMO- AND HETEROZYGOTES

Results file with resolution in HLA-A*30 homo- and heterozygotes is available upon request.

101.429-12 – including *Taq* polymerase, IFU-01
101.429-12u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **38V**

Lot-specific Information

www.olerup-ssp.com

SPECIFICITY TABLE

HLA-A*30 SSP subtyping

Specificities and sizes of the PCR products of the 31+1 primer mixes used for HLA-A*30 SSP subtyping

Primer Mix	Size of spec. PCR product ¹	Size of control band ²	Amplified HLA-A*30 alleles ³	Other amplified HLA-A alleles ⁴
1	165 bp	800 bp	*30:01:01-30:02:04, 30:02:06-30:03, 30:07-30:16, 30:18-30:20, 30:22-30:25, 30:27N-30:28, 30:31-30:45, 30:47-30:76N, 30:78N	*02:52, 03:43, 23:46
2	205 bp	800 bp	*30:01:01-30:01:07, 30:08, 30:11:01-30:11:02, 30:14L-30:16, 30:18-30:20, 30:23-30:24, 30:26, 30:30-30:31, 30:35-30:44, 30:48-30:49, 30:53-30:56, 30:58-30:60, 30:62-30:63, 30:65, 30:71-30:75, 30:78N	
3 ⁸	210 bp	800 bp	*30:02:01-30:03, 30:07, 30:10, 30:12-30:13, 30:22, 30:25, 30:27N, 30:32-30:34, 30:45, 30:50-30:51, 30:57, 30:61, 30:64, 30:66-30:70N, 30:76N	*03:05:02, 11:24:01, 34:02:02, 80:01:01:01-80:03
4 ⁶	155 bp	1070 bp	*30:03, 30:11:01-30:11:02, 30:71	*01:02, 01:20
5	150 bp	1070 bp	*30:04:01-30:04:02, 30:06, 30:17, 30:29, 30:46, 30:77	*02:52, 03:82, 24:66, 68:06
6 ⁵	245 bp 80 bp 185 bp	1070 bp	*30:19 *30:06 *30:07	
7	200 bp	1070 bp	*30:08, 30:28	*02:185 ^W
8 ^{5,8}	85 bp	800 bp	*30:10	
9	180 bp	1070 bp	*30:01:01-30:02:12, 30:04:01-30:04:02, 30:06-30:07, 30:09-30:10, 30:12-30:20, 30:23-30:54, 30:56-30:70N, 30:72-30:78N	*24:124
10	150 bp	1070 bp	*30:12, 30:18, 30:55	*01:07, 02:185, 24:124, 26:19, 29:14, 31:01:02-31:01:06, 31:01:08-31:07, 31:09-31:43, 31:45-31:79
11	215 bp	800 bp	*30:01:01-30:04:02, 30:06, 30:09-30:20, 30:23-30:30, 30:32-30:54, 30:56-30:59N, 30:61-30:78N	*01:02 ^W , 01:20 ^W , 03:72, 11:88, 23:09 ^W , 23:51, 24:24, 24:67, 24:129 ^W , 24:145, 24:156, 24:191, 26:16, 29:37, 32:07, 68:45
12 ⁶	190 bp	1070 bp	*30:09, 30:35	*03:42, 03:131, 03:133, 31:03-31:04, 33:49

101.429-12 – including *Taq* polymerase, IFU-01
101.429-12u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: 38V		Lot-specific Information		www.olerup-ssp.com
13	155 bp	1070 bp	*30:01:01-30:02:12, 30:04:01-30:04:02, 30:06-30:10, 30:13-30:17, 30:19-30:20, 30:23-30:40, 30:42-30:54, 30:56-30:70N, 30:72-30:78N	
14 ⁶	210 bp	1070 bp	*30:13, 30:16, 30:44, 30:46	
15	245 bp	800 bp	*30:14L, 30:29	
16	265 bp	1070 bp	*30:15, 30:33	
17	225 bp	1070 bp	*30:01:01-30:01:07, 30:11:01-30:11:02, 30:14L-30:20, 30:23-30:24, 30:26, 30:30-30:31, 30:35-30:43, 30:48-30:49, 30:52-30:54, 30:56, 30:58-30:60, 30:62-30:63, 30:65, 30:72-30:75, 30:78N	*03:72, 11:88, 24:19, 29:37, 68:45
18	130 bp 210 bp	800 bp	*30:32 *30:20, 30:71	
19 ⁵	125 bp 160 bp 200 bp 235 bp	1070 bp	*30:70N *30:22 *30:31 *30:77	*03:45
20	210 bp 275 bp	800 bp	*30:23 *30:34	*03:04:02-03:04:03, 11:153:01
21	155 bp 185 bp	1070 bp	*30:30 *30:24	
22	180 bp 215 bp	800 bp	*30:36, 30:73N *30:27N	
23	155 bp	1070 bp	*30:26	
24 ⁷	200 bp	800 bp	*30:25	
25	170 bp	1070 bp	*30:40	
26 ⁵	80 bp	1070 bp	*30:42	
27	145 bp 195 bp	1070 bp	*30:56 *30:45	
28	270 bp	1070 bp	*30:51	*02:121, 02:425, 03:154, 23:47, 24:234, 29:06, 31:51, 68:14
29	130 bp	1070 bp	*30:59N, 30:61, 30:74	*01:57N, 02:156, 02:338, 03:17:01, 03:171, 11:119:01-11:119:02, 24:92, 68:103
30 ⁵	190 bp 80 bp	1070 bp	*30:73N *30:37	
31	150 bp	1070 bp	*30:76N	
32 ⁹	Negative Control			

¹Alleles are assigned by the presence of specific PCR product(s). However, the sizes of the specific PCR products may be helpful in the interpretation of HLA-A*30 SSP typings.

When the primers in a primer mix can give rise to HLA-specific PCR products of more than one length this is indicated if the size difference is more than 20 base pairs. Size differences of 20 base pairs or less are not given. For high resolution SSP kits, the alleles listed are specified according to amplicon length.

Nonspecific amplifications, i.e. a ladder or a smear of bands, may sometimes be seen. GC-rich primers have a higher tendency of giving rise to nonspecific amplifications than other primers.

101.429-12 – including *Taq* polymerase, IFU-01
101.429-12u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: 38V

Lot-specific Information

www.olerup-ssp.com

PCR fragments longer than the control bands may sometimes be observed. Such bands should be disregarded and do not influence the interpretation of the SSP typings.

PCR fragments migrating faster than the control bands, but slower than a 400 bp fragment may be seen in some gel read-outs. Such bands can be disregarded and do not influence the interpretation of the SSP typings.

Some primers may give rise to primer oligomer artifacts. Sometimes this phenomenon is an inherent feature of the primer pair(s) of a primer mix. More often it is due to other factors such as too low amount of DNA in the PCR reactions, taking too long time in setting up the PCR reactions, working at elevated room temperature or using thermal cyclers that are not pre-heated

²The internal positive control primer pairs amplify segments of the human growth hormone gene. The internal positive control bands are 1070 or 800 base pairs respectively, well distribution as outlined in the table. Well number 1 contains the shorter, 800 bp, internal positive control band. The well distribution of the internal controls can help in orientation of the kit on gel photo, as well as allow for kit identification. In the presence of a specific amplification the intensity of the control band often decreases.

³For several HLA Class I alleles 1st and/or 4th exon(s) and beyond, as well as intron nucleotide sequences, are not available. In these instances it is not known whether some of the primers of the SSP sets are completely matched with the target sequences or not. Assumption is made that unknown sequences in these regions are conserved within allelic groups.

⁴Due to the sharing of sequence motifs between HLA-A alleles non-HLA-A*30 alleles will be amplified by primer mixes 1, 3 to 5, 7, 9 to 12, 17, 19, 20, 28 and 29.

⁵HLA-specific PCR products shorter than 125 base pairs have a lower intensity and are less sharp than longer PCR products.

⁶Primer mixes 4, 12 and 14 may give rise to nonspecific amplifications, most pronounced in primer mixes 12 and 14.

⁷Primer mix 24 may have a tendency of giving rise to primer oligomer formation.

⁸Primer mixes 3 and 8 may give rise to a lower yield of HLA-specific PCR product than the other A*30 primer mixes.

⁹Primer mix 32 contains a negative control, which will amplify more than 95% of HLA amplicons as well as the amplicons generated by control primer pairs. PCR product sizes range from 75 to 200 base pairs. The PCR product generated by the control primer pair is 430 base pairs.

'w', might be weakly amplified.

101.429-12 – including *Taq* polymerase, IFU-01
101.429-12u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **38V**

Lot-specific Information

www.olerup-ssp.com

PRIMER SPECIFICATION

Well No.	1	2	3	4	5	6	7	8	9	10	11	12
Length of spec. PCR product	165	205	210	155	150	80	200	85	180	150	215	190
Length of int. pos. control ¹	800	800	800	1070	1070	1070	1070	800	1070	1070	800	1070
5'-primer(s) ²	414	363	363	123	123	123	78	367	98	127	98	363
	5'-gAA 3'	5'-ATA 3'	5'-ATA 3'	5'-AgT 3'	5'-AgT 3'	5'-AgT 3'	5'-TCT 3'	5'-TgC 3'	5'-CTC 3'	5'-ggg 3'	5'-CTC 3'	5'-ATA 3'
					414		362					
					5'-gAA 3'		5'-ggT 3'					
3'-primer(s) ³	538	526	526	238	325	163	238	411	238	238	270	505
	5'-CAA 3'	5'-CCA 3'	5'-CCg 3'	5'-CCC 3'	5'-gTg 3'	5'-CgC 3'	5'-CCT 3'	5'-TCA 3'	5'-CCT 3'	5'-CCT 3'	5'-ACA 3'	5'-gCT 3'
			535		524	265	526					524
			5'-CTA 3'		5'-CAT 3'	5'-CCC 3'	5'-CCg 3'					5'-CAC 3'
			538									
			5'-CAg 3'									
Well No.	1	2	3	4	5	6	7	8	9	10	11	12

Well No.	13	14	15	16	17	18	19	20	21	22	23	24
Length of spec. PCR product	155	210	245	265	225	130	125	210	155	180	155	105
						210	160	275	185	215		200
							200					
							235					
Length of int. pos. control ¹	1070	1070	800	1070	1070	800	1070	800	1070	800	1070	800
5'-primer(s) ²	123	123	363	363	98	123	123	363	410	363	413	123
	5'-AgT 3'	5'-AgT 3'	5'-ATA 3'	5'-ATA 3'	5'-CTC 3'	5'-AgT 3'	5'-AgT 3'	5'-ATA 3'	5'-gTg 3'	5'-ATA 3'	5'-CCg 3'	5'-AgT 3'
				784			650		634			
				5'-ggA 3'			5'-CCC 3'		5'-CAg 3'			
3'-primer(s) ³	238	292	563	586	282	212	206	530	526	494	526	187
	5'-CCT 3'	5'-gTg 3'	5'-Cgg 3'	5'-CAC 3'	5'-gAC 3'	5'-gCC 3'	5'-CA 3'	5'-CCT 3'	5'-CCA 3'	5'-TCC 3'	5'-CCA 3'	5'-g.C 3'
			572	899		289	240	595	777	513		282
			5'-gCg 3'	5'-ACA 3'		5'-AgC 3'	5'-ggA 3'	5'-CCT 3'	5'-gCA 3'	5'-TCC 3'		5'-gAC 3'
						294	281			535		
						5'-CgT 3'	5'-ACC 3'			5'-CTA 3'		
							845					
							5'-AgT 3'					
Well No.	13	14	15	16	17	18	19	20	21	22	23	24

101.429-12 – including *Taq* polymerase, IFU-01
101.429-12u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **38V**

Lot-specific Information

www.olerup-ssp.com

Well No.	25	26	27	28	29	30	31
Length of spec.	170	80	145	270	130	80	150
PCR product			195		190		
Length of int. pos. control ¹	1070	1070	1070	1070	1070	1070	1070
5'-primer(s) ²	397	363	176	385	363	485	127
	5' -gCC 3'	5' -ATA 3'	5' -gCA 3'	5' -ggT 3'	5' -ATA 3'	5' -CAC 3'	5' -gAT 3'
			414				
			5' -gAA 3'				
3'-primer(s) ³	526	403	282	616	453	526	238
	5' -CCA 3'	5' -gCT 3'	5' -gAC 3'	5' -CgC 3'	5' -TCT 3'	5' -CCA 3'	5' -CCT 3'
			570		454		
			5' -CAg 3'		5' -CTA 3'		
					513		
					5' -TCC 3'		
Well No.	25	26	27	28	29	30	31

¹The internal positive control primer pairs amplify segments of the human growth hormone gene. The internal positive control bands are 1070 or 800 base pairs respectively, well distribution as outlined in the table. Well number 1 contains the shorter, 800 bp, internal positive control band. The well distribution of the internal controls can help in orientation of the kit on gel photo, as well as allow for kit identification. In the presence of a specific amplification the intensity of the control band often decreases.

²The nucleotide position matching the specificity-determining 3'-end of the primer is given. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

³The nucleotide position matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

101.429-12 – including *Taq* polymerase, IFU-01
101.429-12u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **38V**

Lot-specific Information

www.olerup-ssp.com

CELL LINE VALIDATION SHEET																					
HLA-A*30 SSP subtyping kit ²																					
					Lot No.:	Well															
						1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
						201206201	201206202	201206203	201206204	201206205	201206206	201206207	201206208	201206209	201206210	201206211	201206212	201206213	201206214	201206215	201206216
	IHWC cell line ¹	A*	A*																		
1	9001 SA	*24:02				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
2	9280 LK707	*02:01				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
3	9011 E4181324	*01:01				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
4	9275 GU373	*30:01				+	+	-	-	-	-	-	-	+	-	+	-	+	-	-	-
5	9009 KAS011	*01:01				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
6	9353 SM	*02:01	*26:03			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
7	9020 QBL	*26:01				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
8	9025 DEU	*31:01				-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-
9	9026 YAR	*26:01				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
10	9107 LKT3	*24:02				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
11	9051 PITOUT	*29:02				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
12	9052 DBB	*02:01				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
13	9004 JESTHOM	*02:01				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
14	9071 OLGA	*31:01				-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-
15	9075 DKB	*24:02				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
16	9037 SWEIG007	*29:02				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
17	9282 CTM3953540	*03:01	*80:01			-	-	+	-	-	-	-	-	-	-	-	-	-	-	-	-
18	9257 32367	*33:03	*74:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
19	9038 BM16	*02:01				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
20	9059 SLE005	*02:01				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
21	9064 AMALA	*02:17				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
22	9056 KOSE	*02:01				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
23	9124 IHL	*02:01	*34:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
24	9035 JBUSH	*32:01				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
25	9049 IBW9	*33:01				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
26	9285 WT49	*02:05				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
27	9191 CH1007	*24:10	*29:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
28	9320 BEL5GB	*02:01	*29:02			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
29	9050 MOU	*29:02				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
30	9021 RSH	*30:01	*68:02			+	+	-	-	-	-	-	-	+	-	+	-	+	-	-	-
31	9019 DUCAF	*30:02				+	-	+	-	-	-	-	-	+	-	+	-	+	-	-	-
32	9297 HAG	*02:01				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
33	9098 MT14B	*31:01				-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-
34	9104 DHIF	*31:01				-	-	-	-	-	-	-	-	-	+	-	-	-	-	-	-
35	9302 SSTO	*32:01				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
36	9024 KT17	*02:06	*11:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
37	9065 HHKB	*03:01				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
38	9099 LZL	*02:17				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
39	9315 CML	*01:01	*03:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
40	9134 WHONP199	*02:07	*30:01			+	+	-	-	-	-	-	-	+	-	+	-	+	-	-	-
41	9055 H0301	*03:01				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
42	9066 TAB089	*02:07				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
43	9076 T7526	*02:06	*02:07			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
44	9057 TEM	*66:01				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
45	9239 SHJO	*23:01	*24:02			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
46	9013 SCHU	*03:01				-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
47	9045 TUBO	*02:16	*03:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
48	9303 TER-ND	*02:01	*11:01			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-

101.429-12 – including *Taq* polymerase, IFU-01
101.429-12u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **38V**

Lot-specific Information

www.olerup-ssp.com

CELL LINE VALIDATION SHEET				Well															
HLA-A*30 SSP subtyping kit ²				17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
				Lot No.:	201206217	201430918	201430919	201206220	201206221	201430922	201430923	201430924	201430925	201206226	201430927	201206228	201430929	201206230	201430931
	IHWC cell line ¹	A*	A*																
1	9001 SA	*24:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
2	9280 LK707	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
3	9011 E4181324	*01:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
4	9275 GU373	*30:01		+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
5	9009 KAS011	*01:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
6	9353 SM	*02:01	*26:03	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
7	9020 QBL	*26:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
8	9025 DEU	*31:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
9	9026 YAR	*26:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
10	9107 LKT3	*24:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
11	9051 PITOUT	*29:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
12	9052 DBB	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
13	9004 JESTHOM	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
14	9071 OLGA	*31:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
15	9075 DKB	*24:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
16	9037 SWEIG007	*29:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
17	9282 CTM3953540	*03:01	*80:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
18	9257 32367	*33:03	*74:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
19	9038 BM16	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
20	9059 SLE005	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
21	9064 AMALA	*02:17		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
22	9056 KOSE	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
23	9124 IHL	*02:01	*34:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
24	9035 JBUSH	*32:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
25	9049 IBW9	*33:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
26	9285 WT49	*02:05		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
27	9191 CH1007	*24:10	*29:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
28	9320 BEL5GB	*02:01	*29:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
29	9050 MOU	*29:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
30	9021 RSH	*30:01	*68:02	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
31	9019 DUCAF	*30:02		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
32	9297 HAG	*02:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
33	9098 MT14B	*31:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
34	9104 DHIF	*31:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
35	9302 SSTO	*32:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
36	9024 KT17	*02:06	*11:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
37	9065 HHKB	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
38	9099 LZL	*02:17		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
39	9315 CML	*01:01	*03:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
40	9134 WHONP199	*02:07	*30:01	+	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
41	9055 H0301	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
42	9066 TAB089	*02:07		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
43	9076 T7526	*02:06	*02:07	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
44	9057 TEM	*66:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
45	9239 SHJO	*23:01	*24:02	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
46	9013 SCHU	*03:01		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
47	9045 TUBO	*02:16	*03:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
48	9303 TER-ND	*02:01	*11:01	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	

¹The provided cell line HLA specificities are retrieved from the <http://www.ihwg.org/hla> web site. The specificity of an individual cell line may thus be subject to change.

101.429-12 – including *Taq* polymerase, IFU-01
101.429-12u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: 38V

Lot-specific Information

www.olerup-ssp.com

²The specificity of each primer solution in the kit has been tested against 48 well characterized cell line DNAs and where applicable, additional cell line DNAs.

No DNAs carrying the alleles to be amplified by primer solutions 6 to 7, 14 to 16 and 18 to 31 were available. The specificities of the primers in primer 6, 7, 14 to 16, 19, 23, 24, 27, 29 and 30 were tested by separately adding one additional 5'-primer, respectively one additional 3'-primer. In primer solutions 18, 20, 22, 26 and 28 it was only possible to test the 5'-primer, the 3'-primers were not possible to test. In primer solution 21, 25 and 31 it was only possible to test the 3'-primers, the 5'-primers were not possible to test.

In primer solutions 7 and 16 one 5'-primer was not possible to test, and in primer solutions 2, 3, 5, 6, 12, 15, 16, 19, 24, 27 and 29 one or two 3'-primers were not possible to test.

In addition, one 5'-primer in primer solution 5 was tested by adding an additional 3'-primer.

101.429-12 – including *Taq* polymerase, IFU-01
101.429-12u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **38V**

Lot-specific Information

www.olerup-ssp.com

101.429-12 – including *Taq* polymerase, IFU-01
101.429-12u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **38V**

Lot-specific Information

www.olerup-ssp.com

101.429-12 – including *Taq* polymerase, IFU-01
101.429-12u – without *Taq* polymerase, IFU-02

Visit www.olerup-ssp.com for
“Instructions for Use” (IFU)

Lot No.: **38V**

Lot-specific Information

www.olerup-ssp.com

ADDRESSES:

Manufacturer:

Olerup SSP AB, Franzengatan 5, SE-112 51 Stockholm, Sweden.

Tel: +46-8-717 88 27

Fax: +46-8-717 88 18

E-mail: info-ssp@olerup.com

Web page: <http://www.olerup-ssp.com>

Distributed by:

Olerup GmbH, Löwengasse 47 / 6, AT-1030 Vienna, Austria.

Tel: +43-1-710 15 00

Fax: +43-1-710 15 00 10

E-mail: support-at@olerup.com

Web page: <http://www.olerup.com>

Olerup Inc., 901 S. Bolmar St., Suite R, West Chester, PA 19382

Tel: 1-877-OLERUP1

Fax: 610-344-7989

E-mail: info.us@olerup.com

Web page: <http://www.olerup.com>

For information on *Olerup* SSP distributors worldwide, contact **Olerup GmbH**.