Visit www.olerup-ssp.com for

"Instructions for Use" (IFU)

101.852-12- including *Taq* polymerase, IFU-01 101.852-12u- without *Taq* polymerase, IFU-02

Lot No.: **40S**

Lot-specific information Olerup SSP[®] HLA-B*07 Add-on

Product number:	101.852-12 – including <i>Taq</i> polymerase 101.852-12u– without <i>Taq</i> polymerase
Lot number:	40S
Expiry date:	2015-October-01
Number of tests:	12
Number of wells per test:	6
Storage - pre-aliquoted primers:	dark at -20°C
- PCR Master Mix:	-20°C
 Adhesive PCR seals 	RT
- Product Insert	RT

This Product Description is only valid for Lot No. 40S

CHANGES COMPARED TO THE PREVIOUS OLERUP SSP[®] HLA-B*07 ADD-ON LOT (23N)

The Lot-specific information for HLA- B*07 Add-on including and without *Taq* polymerase is described in one common Product Insert.

One well has been added to the HLA-B*07 Add-on kit, well 6.

The HLA-B*07 Add-on specificity and interpretation tables have been updated for the HLA-B alleles described since the previous *Olerup* SSP[®] HLA-B*07 Add-on lot was made (Lot No. 23N).

The primers of the wells detailed below have been exchanged, added or modified compared to the previous lot.

Well	5'-primer	3'-primer	rationale
6	New	New	New primer pair for the B*07:161N allele.

101.852-12- including *Taq* polymerase, IFU-01 101.852-12u- without *Taq* polymerase, IFU-02 Visit <u>www.olerup-ssp.com</u> for "Instructions for Use" (IFU)

Lot No.: **40S**

Lot-specific information PRODUCT DESCRIPTION

HLA-B*07 SSP subtyping

CONTENT

The primer set contains 5'- and 3'-primers to separate the HLA-B*07:61 and B*07:161N from the B*07:02 alleles and the HLA-B*07:05 from the B*07:06 alleles.

PLATE LAYOUT

Each test consists of 6 PCR reactions in an 8 well cut PCR plate. Wells 7 to 8 are empty.

1 2	3	4	5	6	empty	empty
-----	---	---	---	---	-------	-------

The 8 well cut PCR plate is marked with 'HLA-B*07' in silver/gray ink.

Well No. 1 is marked with the Lot Number '40S'.

A faint row of numbers is seen between wells 1 and 2 or wells 7 and 8 of the PCR trays. These stem from the manufacture of the trays, and should be disregarded.

The PCR plates are covered with a PCR-compatible foil.

Please note: When removing each 8 well PCR plate, make sure that the remaining plates stay covered. Use a scalpel or a similar instrument to carefully cut the foil between the plates.

INTERPRETATION

The interpretation of HLA-B*07 Add-on SSP subtypings will be influenced by the B*07, the B*08:79, the B*27:36, two B*35, the B*37:07, most B*40, most B*42, two B*44, the B*46:06, the B*48, two B*54, seven B*55, six B*56 and the B*81 alleles. In addition, the C*04:77 allele will be amplified by primer mix 1.

UNIQUELY IDENTIFIED ALLELES

The HLA-B*07:02, 07:05, 07:06, 07:61 and 07:161N alleles give different patterns in the B*07 Add-on subtyping kit¹.

The HLA-B*07 Add-on kit cannot distinguish these silent mutations: the B*07:02:01-07:02:05, B*07:02:07-07:02:25 and 07:02:27-07:02:38, the B*07:05:01-07:05:06, the B*07:18:01-07:18:02, the B*07:22:01-07:22:02, the B*07:56:01-07:56:02, the B*07:68:01-07:68:03 and the B*07:96:01-07:96:02 alleles.

¹Based on HLA-B alleles listed on the IMGT/HLA web page 2013-January-11, release 3.11.0, <u>www.ebi.ac.uk/imgt/hla</u>.

101.852-12- including *Taq* polymerase, IFU-01 101.852-12u- without *Taq* polymerase, IFU-02

Lot No.: 40S

Lot-specific information

Visit <u>www.olerup-ssp.com</u> for "Instructions for Use" (IFU)

SPECIFICITY TABLE

HLA-B*07 Add-on SSP subtyping

Specificities and sizes of the PCR products of the 6 primer mixes used for HLA-B*07 Add-on SSP subtyping

Primer Mix	Size of spec. PCR product ¹	Size of control band ²	Amplified HLA-B*07 alleles ³	Other amplified HLA Class I alleles⁴
1 ⁵	105 bp	800 bp	*07:02:01-07:02:25, 07:02:27-07:03, 07:05:01- 07:18:02, 07:20-07:24, 07:27-07:32, 07:35-07:39, 07:41-07:42, 07:44-07:47, 07:49N-07:52, 07:54, 07:56:01-07:59, 07:61- 07:65, 07:67N-07:79, 07:81-07:99, 07:101- 07:121, 07:123-07:130, 07:132-07:138, 07:140- 07:145, 07:148-07:163, 07:166-07:181N	08:79, 35:66, 37:07, 40:15- 40:16, 40:23, 40:32, 40:98, 40:128, 40:198, 44:150, 48:05, 48:08, C*04:77
2 ⁶	355 bp	1070 bp	*07:05:01-07:05:06	35:46, 37:07
3	355 bp	1070 bp	*07:02:01-07:02:05, 07:02:07-07:04, 07:06- 07:99, 07:101-07:181N	08:79, 27:36, 40:01:01- 40:01:29, 40:07, 40:09 [?] , 40:10:01-40:10:02, 40:11:01 [?] -40:13 [?] , 40:14:01- 40:14:03, 40:15 [?] -40:16 [?] , 40:18 [?] , 40:20 [?] -40:21 [?] , 40:22N, 40:23 [?] -40:28 [?] , 40:30 [?] -40:38 [?] , 40:42 [?] -40:48 [?] , 40:49, 40:51-40:52, 40:54- 40:55, 40:59 [?] -40:60 [?] , 40:61, 40:62 [?] , 40:63, 40:64 [?] -40:65 [?] , 40:66-40:67, 40:69 [?] -40:71 [?] , 40:72:01-40:73, 40:74 [?] - 40:76 [?] , 40:77, 40:78 [?] , 40:79, 40:80 [?] , 40:81, 40:83 [?] -40:84 [?] , 40:86 [?] -40:87:01 [?] , 40:87:02- 40:88, 40:90 [?] , 40:92 [?] -40:93 [?] , 40:95 [?] -40:96 [?] , 40:98, 40:99 [?] - 40:121 [?] , 40:123 [?] , 40:124:01, 40:124:02 [?] -40:127 [?] , 40:128, 40:129 [?] -40:137 [?] , 40:138- 40:141, 40:145 [?] ,

101.852-12- including Taq polymerase, IFU-01	Visit <u>www.olerup-ssp.com</u> for
101.852-12u- without Taq polymerase, IFU-02	"Instructions for Use" (IFU)

Lot	: No.: 40S		Lot-specific information	
				40:146-40:147, 40:148 [?] - 40:149 [?] , 40:150-40:155N, 40:156 [?] -40:157 [?] , 40:158, 40:160:01 [?] -40:160:02 [?] , 40:162 [?] -40:163 [?] , 40:165 [?] - 40:166 [?] , 40:168 [?] -40:171 [?] , 40:172, 40:173 [?] , 40:175, 40:177 [?] -40:178 [?] , 40:179, 40:180 [?] -40:185 [?] , 40:186, 40:187 [?] -40:192 [?] , 40:193, 40:194 [?] -40:203 [?] , 40:204, 40:206 [?] , 40:207-40:208, 40:210, 40:212-40:213, 40:215-40:216N, 40:217 [?] - 40:219 [?] , 40:221-40:223, 40:224 [?] -40:228 [?] , 40:231 [?] - 40:233 [?] , 44:31, 46:06, 48:01:01-48:01:03, 48:03:01- 48:12, 48:14-48:28, 81:01- 81:02, 81:04N-81:06
4	425 bp	1070 bp	*07:05:01-07:06, 07:34, 07:40, 07:53, 07:80, 07:90, 07:97, 07:105, 07:112, 07:137-07:138, 07:140, 07:176	40:166, 42:01:01-42:02, 42:04-42:05:02, 42:07-42:08, 42:10, 42:12-42:15, 42:17- 42:18, 54:04, 54:11, 55:04, 55:08, 55:17, 55:20, 55:28, 55:49, 55:51, 56:05:01-56:06, 56:12, 56:15, 56:21, 81:01- 81:04N, 81:06
5	350 bp	1070 bp	*07:61	
6	405 bp	1070 bp	*07:161N	

¹Alleles are assigned by the presence of specific PCR product(s). However, the sizes of the specific PCR products may be helpful in the interpretation of HLA-B*07 Add-on SSP typings.

When the primers in a primer mix can give rise to specific PCR products of more than one length this is indicated if the size difference is 20 base pairs or more. Size differences shorter than 20 base pairs are not given. For high resolution SSP kits the respective lengths of the specific PCR product(s) of the alleles amplified by these primer mixes are given.

Nonspecific amplifications, i.e. a ladder or a smear of bands, may sometimes be seen. GC-rich primers have a higher tendency of giving rise to nonspecific amplifications than other primers.

PCR fragments longer than the control bands may sometimes be observed. Such bands should be disregarded and do not influence the interpretation of the SSP typings.

PCR fragments migrating faster than the control bands, but slower than a 400 bp fragment may be seen in some gel read-outs. Such bands can be disregarded and do not influence the interpretation of the SSP typings.

Some primers may give rise to primer oligomer artifacts. Sometimes this phenomenon is an inherit feature of the primer pair(s) of a primer mix. More often it is due to other factors such as too low amount of DNA in the PCR reactions, taking too long time in setting up the PCR reactions, working at elevated room temperature or using thermal cyclers that are not pre-heated.

²The internal positive control primer pairs amplify segments of the human growth hormone gene. The two different control primer pairs give rise to either an internal positive control band of 1070 base



Product Insert

101.852-12– including *Taq* polymerase, IFU-01 101.852-12u– without *Taq* polymerase, IFU-02

Visit <u>www.olerup-ssp.com</u> for "Instructions for Use" (IFU)

Lot No.: 40S

Lot-specific information

pairs, for most wells, or a band of 800 base pairs, for some wells.

Well number 1 contains the primer pair giving rise to the shorter, 800 bp, internal positive control band in order to help in the correct orientation of the HLA-B*07 Add-on subtyping.

In the presence of a specific amplification the intensity of the control band often decreases.

³For several HLA Class I alleles 1st and/or 4th exon(s) and beyond, as well as intron nucleotide sequences, are not available. In these instances it is not known whether some of the primers of the SSP sets are completely matched with the target sequences or not. We assume that unknown sequences in these regions are conserved within allelic groups.

⁴Due to the sharing of sequence motifs between HLA-B alleles non-HLA-B*07 alleles will be amplified by primer mixes 1 to 4. In addition, the C*04:77 allele will be amplified by primer mix 1.

⁵HLA-specific PCR products shorter than 125 base pairs have a lower intensity and are less sharp than longer PCR products.

⁶Primer solution 2 may give rise to a lower yield of HLA-specific PCR product than the other HLA-B*07 Add-on primer mixes.

'w', might be weakly amplified.

"?', nucleotide sequence information not available for the primer matching sequence.



Visit <u>www.olerup-ssp.com</u> for "Instructions for Use" (IFU)

101.852-12– including *Taq* polymerase, IFU-01 101.852-12u– without *Taq* polymerase, IFU-02

Lot No.: **40S**

Lot-specific information

INTERPRETATION TABLE						
HLA-B*07 Add-on SSP	typing	a				
	, , ,		W	ell		
	1	2	3	4	5	6
Length of spec.	105	355	355	425	350	405
PCR product						
Length of int.	800	1070	1070	1070	1070	1070
pos. control ¹						
5'-primer(s) ²	540	693	693	272	756	693
	^{5'} -Cgg ^{3'}	^{5′} -ggT ^{3′}	^{5'} -ggT ^{3'}	^{5'} -CTA ^{3'}	^{5′} -ACT ^{3′}	^{5′} -ggT ^{3′}
3'-nrimer(s) ³	605	916	916	412	970	962
	^{5'} -aCT ^{3'}	^{5'} -aAT ^{3'}	^{5'} -aAC ^{3'}	^{5'} -aTT ^{3'}	5' -TCT 3'	^{5'} -CAA ^{3'}
Well No.	1	2	3	4	5	6
HI A-B allele ⁴	-	_		-		
*07·02·01_07·02·05_07·02·07_07·02·25_07·02·27_07·03						
07.07.07.18.02 07.20.07.24 07.27.07.32 07.35.07.39						
07:41-07:42 07:44-07:47 07:49N-07:52 07:54 07:56:01-						
07:59 07:62-07:65 07:67N-07:79 07:81-07:89 07:91-						
07.96.02 07.98-07.99 07.101-07.104 07.106-07.111N	1		3			
07.113-07.121, 07.123-07.130, 07.132-07.136, 07.141-						
07:145.07:148-07:160.07:162-07:163.07:166-07:175.						
07:177-07:181N, 08:79, 40:98, 40:128, 48:05, 48:08						
*07:02:06. 35:66. 44:150. <i>C*04:77</i>	1					
*07:02:26, 07:04, 07:19, 07:25-07:26, 07:33, 07:43, 07:48,						
07:55, 07:60, 07:66, 07:122, 07:131, 07:139, 07:146-						
07:147, 07:164-07:165, 27:36, 40:01:01-40:01:29, 40:07,						
40:10:01-40:10:02, 40:14:01-40:14:03, 40:22N, 40:49,						
40:51-40:52, 40:54-40:55, 40:61, 40:63, 40:66-40:67,						
40:72:01-40:73, 40:77, 40:79, 40:81, 40:87:02-40:88,			•			
40:124:01, 40:138-40:141, 40:146-40:147, 40:150-			3			
40:155N, 40:158, 40:172, 40:175, 40:179, 40:186, 40:193,						
40:204, 40:207-40:208, 40:210, 40:212-40:213, 40:215-						
40:216N, 40:221-40:223, 44:31, 46:06, 48:01:01-48:01:03,						
48:03:01-48:04, 48:06-48:07, 48:09-48:12, 48:14-48:28,						
81:05						
*07:05:01-07:05:06	1	2		4		
*07:06, 07:90, 07:97, 07:105, 07:112, 07:137-07:138,	1		3	Λ		
07:140, 07:176	•		J	-		
*07:34, 07:40, 07:53, 07:80, 81:01-81:02, 81:04N, 81:06			3	4		
*07:61	1		3		5	
*07:161N	1		3			6
*35:46		2				
*37:07	1	2				
Well No.	1	2	3	4	5	6

Lot No · 40S

101.852-12- including Taq polymerase, IFU-01 101.852-12u- without Tag polymerase, IFU-02

Visit www.olerup-ssp.com for "Instructions for Use" (IFU)

Lot No.: 40S Lot-specific information	n					
Length of spec.	105	355	355	425	350	405
PCR product						
Well No.	1	2	3	4	5	6
*40:09, 40:11:01-40:13, 40:18, 40:20-40:21, 40:24-40:28,						
40:30-40:31, 40:33-40:38, 40:42-40:48, 40:59-40:60,						
40:62, 40:64-40:65, 40:69-40:71, 40:74-40:76, 40:78,						
40:80, 40:83-40:84, 40:86-40:87:01, 40:90, 40:92-40:93,						
40:95-40:96, 40:99-40:121, 40:123, 40:124:02-40:127,			2			
40:129-40:137, 40:145, 40:148-40:149, 40:156-40:157,			ſ			
40:160:01-40:160:02, 40:162-40:163, 40:165, 40:168-						
40:171, 40:173, 40:177-40:178, 40:180-40:185, 40:187-						
40:192, 40:194-40:197, 40:199-40:203, 40:206, 40:217-						
40:219, 40:224-40:228, 40:231-40:233						
*40:15-40:16, 40:23, 40:32, 40:198	1		?			
*40:166			?	4		
*42:01:01-42:02, 42:04-42:05:02, 42:07-42:08, 42:10,						
42:12-42:15, 42:17-42:18, 54:04, 54:11, 55:04, 55:08,				4		
55:17, 55:20, 55:28, 55:49, 55:51, 56:05:01-56:06, 56:12,				4		
56:15, 56:21, 81:03						
HLA-B allele						
Well No.	1	2	3	4	5	6

¹The internal positive control primer pairs amplify segments of the human growth hormone gene. The two different control primer pairs give rise to either an internal positive control band of 1070 base pairs, for most wells, or a band of 800 base pairs, for some wells.

Well number 1 contains the primer pair giving rise to the shorter, 800 bp, internal positive control band in order to help in the correct orientation of the HLA-B*07 Add-on subtyping.

In the presence of a specific amplification the intensity of the control band often decreases. ²The nucleotide position, in the 2nd, 3rd or 4th exon, matching the specificity-determining 3'-end of the primer is given. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

³The nucleotide position, in 3rd or 5th exon, matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

⁴The sequence of B*0701 has been shown to be in error.

'w', might be weakly amplified.

??, nucleotide sequence information not available for the primer matching sequence.



Lot-specific information

Visit <u>www.olerup-ssp.com</u> for "Instructions for Use" (IFU)

CELL LINE VALIDATION SHEET										
	HLA-	B*07 add-	on SS	6P sul	bty	piı	ng	kit		
						V	Vel	I		
					1	2	3	4	5	6
				Prod. No.:	201196101	201196102	201196103	201196104	201196105	201319106
	IHW	C cell line	E	3*						
1	9001	SA	*07:02		+	-	+	-	-	-
2	9280	LK707	*52:01	*73:01	-	-	-	-	-	-
3	9011	E4181324	*52:01	1	-	-	-	-	-	-
4	9275	GU373	*15:10	*53:01	-	-	-	-	-	-
5	9009	KAS011	*37:01	*= 4 04	-	-	-	-	-	-
6	9353	SM	*40:04	^51:01	-	-	-	-	-	-
/	9020		10:01		-	-	-	-	-	-
Ő	9025		35.01			-		-	-	-
9 10	9020 Q107	I KT3	*54.01		-	-	-	-	-	_
10	9107		*44.01		-	-	-	-	-	-
12	9052	DBB	*57.01		-	-	-	-	-	-
13	9025	JESTHOM	*27:05		-	-	-	-	-	-
14	9071	OLGA	*15:01	*15:20	-	-	-	-	-	-
15	9075	DKB	*40:01		-	-	+	-	-	-
16	9037	SWEIG007	*40:02		-	-	-	-	-	-
17	9282	CTM3953540	*08:01	*55:01	-	-	-	-	-	-
18	9257	32367	*14:01	*56:01	-	-	-	-	-	-
19	9038	BM16	*18:01		-	-	-	-	-	-
20	9059	SLE005	*40:01		-	-	+	-	-	-
21	9064	AMALA	*15:01		-	-	-	-	-	-
22	9056	KOSE	*35:03		-	-	-	-	-	-
23	9124	IHL	*40:02	*56:02	-	-	-	-	-	-
24	9035	JBUSH	*38:01		-	-	•	-	-	-
25	9049	IBW9	*14:02		-	-	-	-	-	-
26	9285	W149	*58:01	*= 4 0 4	-	-	-	-	-	-
27	9191	CH1007	^07:05	^51:01	+	+	-	+	-	-
28	9320	BELSGB	*44:02	44:03	-	-	-	-	-	-
29	9050		44:03			-	-	-	-	-
30	9021 0010		42.01 *18·01		-	-		+	-	_
32	9019	HAG	*41.02		_	-		-	-	-
33	9098	MT14B	*40.01		-	-	+	-	-	-
34	9104	DHIF	*38:01		-	-	-	-	-	-
35	9302	SSTO	*44:02		-	-	-	-	-	-
36	9024	KT17	*15:01	*35:01	-	-	-	-	-	-
37	9065	HHKB	*07:02		+	-	+	-	-	-
38	9099	LZL	*15:01		-	-	-	-	-	-
39	9315	CML	*08:01	*27:05	-	-	-	-	-	-
40	9134	WHONP199	*13:02	*46:01	-	-	-	-	-	-
41	9055	H0301	*14:02		-	-	-	-	-	-
42	9066	TAB089	*46:01		-	-	-	-	-	-
43	9076	T7526	*46:01		-	-	-	-	-	-
44	9057	TEM	*38:01		-	-	-	-	-	-
45	9239	SHJO	*42:01	*50:01	-	-	-	+	-	-
46	9013	SCHU	*07:02		+	-	+	-	-	-
47	9045	TUBO	*51:01		-	-	-	-	-	-
48	9303	TER-ND	*35:01	*44:03	-	-	-	-	-	-



Visit <u>www.olerup-ssp.com</u> for

"Instructions for Use" (IFU)

101.852-12- including *Taq* polymerase, IFU-01 101.852-12u- without *Taq* polymerase, IFU-02

Lot No.: **40S**

Lot-specific information CERTIFICATE OF ANALYSIS

Olerup SSP[®] HLA-B*07 Add-on SSP

101.852-12 – including <i>Taq</i> polymerase
101.852-12u– without <i>Taq</i> polymerase
40S
2015-October-01
12
6

Well specifications:

Well No.	Production No.
1	2011-961-01
2	2011-961-02
3	2011-961-03
4	2011-961-04
5	2011-961-05
6	2013-191-06

The specificity of each primer solution of the kit has been tested against 48 well characterized cell line DNAs.

No DNAs carrying the alleles to be amplified by primer solution 5 and 6 were available.

In primer solution 5 and 6 the 5'-primers were tested by adding one additional 3'primer, the 3'-primers were not possible to test.

Results: No false positive or false negative amplifications were obtained.

Date of approval: 2013-April-24

Approved by:

Production Quality Control



101.852-12- including *Taq* polymerase, IFU-01 101.852-12u- without *Taq* polymerase, IFU-02 Visit <u>www.olerup-ssp.com</u> for "Instructions for Use" (IFU)

Lot No.: 40S Lot-specific information Declaration of Conformity

Product name: Product number: Lot number:	<i>Olerup</i> SSP [®] HLA-B*07 add-on 101.852-12/12u 40S
Intended use:	HLA-B*07 high resolution histocompatibility testing
Manufacturer:	<i>Olerup</i> SSP AB Franzengatan 5 SE-112 51 Stockholm, Sweden <i>Phone:</i> +46-8-717 88 27 <i>Fax:</i> +46-8-717 88 18

We, *Olerup* SSP AB, hereby declare that this product, to which this Declaration of Conformity relates is in conformity with the following Standard(s) and other normative document(s) ISO 9001:2008 and ISO 13485:2012, following the provisions of the 98/79/EC Directive on *in vitro* diagnostic medical devices, Annex II List B, conformity assessed using Annex IV, as transposed into the national laws of the Member States of the European Union.

The Technical Documentation File is maintained at *Olerup* SSP AB, Franzengatan 5, SE-112 51 Stockholm, Sweden.

Notified Body: Lloyd's Register Quality Assurance Limited, Hiramford, Middlemarch Office Village, Siskin Drive, Coventry CV3 4FJ, United Kingdom. (Notified Body number: 0088.)

Stockholm, Sweden 2013-April-24

Ann-Cathrin Jareman Head of QA and Regulatory Affairs



101.852-12– including *Taq* polymerase, IFU-01 101.852-12u– without *Taq* polymerase, IFU-02

Lot No.: **40S**

Lot-specific information

Visit <u>www.olerup-ssp.com</u> for "Instructions for Use" (IFU)



101.852-12– including *Taq* polymerase, IFU-01 101.852-12u– without *Taq* polymerase, IFU-02 Visit <u>www.olerup-ssp.com</u> for "Instructions for Use" (IFU)

Lot No.: 40S ADDRESSES: Lot-specific information

Manufacturer: *Olerup* SSP AB, Franzengatan 5, SE-112 51 Stockholm, Sweden. *Tel:* +46-8-717 88 27 *Fax:* +46-8-717 88 18 *E-mail:* info-ssp@olerup.com *Web page:* http://www.olerup-ssp.com

Distributed by: Olerup GmbH, Löwengasse 47 / 6, AT-1030 Vienna, Austria. *Tel:* +43-1-710 15 00 *Fax:* +43-1-710 15 00 10 *E-mail:* support-at@olerup.com *Web page:* http://www.olerup.com

Olerup Inc., 901 S. Bolmar St., Suite R, West Chester, PA 19382 *Tel:* 1-877-OLERUP1 *Fax:* 610-344-7989 *E-mail:* info.us@olerup.com *Web page:* http://www.olerup.com

For information on Olerup SSP distributors worldwide, contact Olerup GmbH.

