

Olerup SSP® HLA-Cw*01

| | |
|----------------------------------|---|
| Product number: | 101.621-12 – including <i>Taq</i> polymerase 101.621-12u – without <i>Taq</i> polymerase |
| Lot number: | 44E |
| Expiry date: | 2010-April-01 |
| Number of tests: | 12 |
| Number of wells per test: | 16 |
| Storage - pre-aliquoted primers: | dark at -20°C |
| - PCR Master Mix: | -20°C |
| - Adhesive PCR seals | RT |
| - Product Insert | RT |

This Product Description is only valid for Lot No. 44E.

CHANGES COMPARED TO THE PREVIOUS *OLERUP SSP*® HLA-Cw*01 LOT

The HLA-Cw*01 specificity and interpretation tables have been updated for the HLA-Cw alleles described since the previous *Olerup SSP*® HLA-Cw*01 lot was made (Lot No. X34).

Three wells have been added to the HLA-Cw*01 kit,
wells **14 to 16**.

The primers of the wells detailed below have been exchanged, added or modified compared to the previous lot.

| Well | 5'-primer | 3'-primer | rationale |
|------|-----------|-----------|--|
| 2 | - | Added | Primer added for the Cw*0115 allele. |
| 3 | Added | Added | Primer pair added for the Cw*0114 allele. |
| 4 | Added | Added | Primer pair added for the Cw*0116 allele. |
| 5 | - | Added | Primer added for the Cw*0120 allele. |
| 7 | - | Added | Primer added for the Cw*0120 allele. |
| 8 | Modified | - | Increased specificity of specific primer pair. |
| 14 | New | New | New primer pair for the Cw*0117 allele. |
| 15 | New | New | New primer pair for the Cw*0118 allele. |
| 16 | New | New | New primer pair for the Cw*0119 allele. |

Changes in revision R01 compared to R00:

1. Fragment size corrected in footnotes to primer mix 13 in the Specificity and Interpretation Tables.

PRODUCT DESCRIPTION

HLA-Cw*01 SSP typing

CONTENT

The primer set contains 5'- and 3'-primers for identifying the Cw*0102 to Cw*0120 alleles.

PLATE LAYOUT

Each HLA-Cw*01 test consists of 16 PCR reactions in a 16 well cut PCR plate.

| | | | | | | | |
|---|----|----|----|----|----|----|----|
| 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 |
| 9 | 10 | 11 | 12 | 13 | 14 | 15 | 16 |

The 16 well PCR plate is marked with 'Cw*01'.

Well No. 1 is marked with the Lot No. '44E'.

The PCR plates are heat-sealed with a PCR-compatible foil.

Please note: When removing each 16 well PCR plate, make sure that the remaining plates stay sealed. Use a scalpel or a similar instrument to carefully cut the foil between the plates.

INTERPRETATION

Non HLA-Cw*01 alleles will be amplified by primer mixes 4, 11 and 12. In addition, the B*1533 allele will be amplified by primer mix 12.

UNIQUELY IDENTIFIED ALLELES

All the HLA-Cw*01 alleles, i.e. **Cw*0102 to Cw*0120**, recognized by the HLA Nomenclature Committee in April 2008¹ will be amplified by the primers in the HLA-Cw*01 SSP kit.

The HLA-Cw*01 subtyping kit cannot distinguish the Cw*010201 to Cw*010205 alleles.

¹HLA-Cw alleles listed on the IMGT/HLA web page 2008-April-08, release 2.21.0, www.ebi.ac.uk/imgt/hla.

RESOLUTION IN HOMO- AND HETEROZYGOTES

The 19 HLA-Cw*01 alleles can be combined in 190 homozygous and heterozygous combinations. Ninety of these genotypes do not give rise to unique amplification patterns. In these calculations the different sizes of the PCR products generated by primer mixes 2, 3, 4, 5, 7, and 13 have not been considered.

| | | |
|-----------|---------|-----------------------------------|
| +++----- | -----+ | 0103,0114 = 0114,0115 |
| ++-+----- | -----+ | 0103,0116 = 0105,0115 = 0115,0116 |
| ++-++-+- | -----+ | 0103,0120 = 0115,0120 |
| ++-+----- | -----+ | 0103,0106 = 0106,0115 |
| ++-+----- | -----+ | 0103,0107 = 0107,0115 |
| ++-+----- | -----+ | 0103,0108 = 0108,0115 |
| ++-+----- | -----+ | 0103,0109 = 0109,0115 |
| ++-+----- | +----- | 0103,0110 = 0110,0115 |
| ++-+----- | -+----- | 0103,0111 = 0111,0115 |

Lot No.: **44E**

Lot-specific information

www.olerup.com

| | | |
|-----------|-----------|--|
| ++----- | ----+---- | 0103,0113 = 0113,0115 |
| ++----- | ----+---- | 0103,0117 = 0115,0117 |
| ++----- | ----+---- | 0103,0118 = 0115,0118 |
| ++----- | ----+---- | 0103,0119 = 0115,0119 |
| ++----- | ----+---- | 0102,0103 = 0102,0115 = 0103,0115 = 0115,0115 |
| +++----- | ----+---- | 0105,0114 = 0114,0116 |
| +--+----- | ----+---- | 0102,0104 = 0104,0109 = 0104,0114 = 0109,0114 |
| +--+----- | ----+---- | 0102,0114 = 0114,0114 |
| +--+----- | ----+---- | 0105,0120 = 0116,0120 |
| +--+----- | ----+---- | 0105,0106 = 0106,0116 |
| +--+----- | ----+---- | 0105,0107 = 0107,0116 |
| +--+----- | ----+---- | 0105,0108 = 0108,0116 |
| +--+----- | ----+---- | 0105,0109 = 0109,0116 |
| +--+----- | +--+----- | 0105,0110 = 0110,0116 |
| +--+----- | +--+----- | 0105,0111 = 0111,0116 |
| +--+----- | ----+---- | 0105,0113 = 0113,0116 |
| +--+----- | ----+---- | 0105,0117 = 0116,0117 |
| +--+----- | ----+---- | 0105,0118 = 0116,0118 |
| +--+----- | ----+---- | 0105,0119 = 0116,0119 |
| +--+----- | ----+---- | 0102,0105 = 0102,0116 = 0105,0116 = 0116,0116 |
| +--+----- | ----+---- | 0102,0120 = 0106,0108 = 0106,0120 = 0108,0120 = 0120,0120 |
| +--+----- | ----+---- | 0102,0106 = 0106,0106 |
| +--+----- | ----+---- | 0102,0107 = 0107,0107 |
| +--+----- | ----+---- | 0102,0108 = 0108,0108 |
| +--+----- | ----+---- | 0102,0109 = 0109,0109 |
| +--+----- | +--+----- | 0102,0110 = 0110,0110 |
| +--+----- | +--+----- | 0102,0111 = 0111,0111 |
| +--+----- | ----+---- | 0102,0113 = 0113,0113 |
| +--+----- | ----+---- | 0102,0117 = 0117,0117 |
| +--+----- | ----+---- | 0102,0118 = 0118,0118 |
| +--+----- | ----+---- | 0102,0119 = 0119,0119 |

SPECIFICITY TABLE

HLA-Cw*01 SSP subtyping

Specificities and sizes of the PCR products of the 16 primer mixes used for HLA-Cw*01 SSP subtyping

| Primer Mix | Size of spec. PCR product ¹ | Size of control band ² | Amplified HLA-Cw*01 alleles | Other amplified HLA Class I alleles ³ |
|------------------------|--|-----------------------------------|-------------------------------------|--|
| 1⁴ | 90 bp | 800 bp | 010201-010205, 0106-0111, 0113-0120 | |
| 2^{4,5} | 90, 270 bp | 1070 bp | 0103, 0115 | |
| 3⁶ | 150, 265 bp | 800 bp | 0104, 0114 | |
| 4⁷ | 210, 240 bp | 1070 bp | 0105, 0116 | 0737 |
| 5⁸ | 150, 200 bp | 800 bp | 0106, 0120 | |
| 6 | 195 bp | 1070 bp | 0107 | |
| 7⁹ | 150, 195 bp | 1070 bp | 0108, 0120 | |
| 8 | 210 bp | 1070 bp | 0104, 0109 | |
| 9 | 230 bp | 800 bp | 0110 | |
| 10 | 290 bp | 1070 bp | 0111 | |
| 11 | 340 bp | 800 bp | 0112 | 0617, 120306, 1611 |
| 12 | 155 bp | 1070 bp | 0113 | 0509, 0517, 0815, B*1533 |
| 13¹⁰ | 155, 370 bp | 800 bp | 010201-0120 | |
| 14 | 240 bp | 1070 bp | 0117 | |
| 15 | 230 bp | 800 bp | 0118 | |
| 16 | 255 bp | 1070 bp | 0119 | |

¹Alleles are assigned by the presence of specific PCR product(s). However, the sizes of the specific PCR products may be helpful in the interpretation of HLA-Cw*01 SSP typings.

When the primers in a primer mix can give rise to specific PCR products of more than one length this is indicated if the size difference is 20 base pairs or more. Size differences shorter than 20 base pairs are not given. For high resolution SSP kits the respective length of the specific PCR product(s) of the alleles amplified by these primer mixes are given.

Nonspecific amplifications, i.e. a ladder or a smear of bands, may sometimes be seen. GC-rich primers have a higher tendency of giving rise to nonspecific amplifications than other primers.

PCR fragments longer than the control bands may sometimes be observed. Such bands should be disregarded and do not influence the interpretation of the SSP typings.

PCR fragments migrating faster than the control bands, but slower than a 400 bp fragment may be seen in some gel read-outs. Such bands can be disregarded and do not influence the interpretation of the SSP typings.

Some primers may give rise to primer oligomer artifacts. Sometimes this phenomenon is an inherent feature of the primer pair(s) of a primer mix. More often it is due to other factors such as too low amount of DNA in the PCR reactions, taking too long time in setting up the PCR reactions, working at elevated room temperature or using thermal cyclers that are not pre-heated.

²The internal positive control primer pairs amplify segments of the human growth hormone gene. The two different control primer pairs give rise to either an internal positive control band of 1070

Lot No.: **44E**

Lot-specific information

www.olerup.com

base pairs, for most wells, or a band of 800 base pairs, for some wells.

Well number 1 contains the primer pair giving rise to the shorter, 800 bp, internal positive control band in order to help in the correct orientation of the HLA-Cw*01 SSP subtyping.

In addition, wells number 3, 5, 9, 11, 13 and 15 contain the primer pair giving rise to the shorter, 800 bp, internal positive control band in order to allow kit identification.

In the presence of a specific amplification the intensity of the control band often decreases.

³Due to the sharing of sequence motifs between HLA Class I alleles some non-HLA-Cw*01 alleles will be amplified by primer mixes 4, 11 and 12. In addition, the B*1533 allele will be amplified by primer mix 12.

⁴Short specific PCR fragments are less intense and not as sharp as longer specific bands.

⁵Primer mix 2: Specific PCR product of 90 bp in the Cw*0103 allele. Specific PCR product of 270 bp in the Cw*0115 allele.

⁶Primer mix 3: Specific PCR product of 150 bp in the Cw*0104 allele. Specific PCR product of 265 bp in the Cw*0114 allele.

⁷Primer mix 4: Specific PCR product of 210 bp in the Cw*0105 and *0737 alleles. Specific PCR product of 240 bp in the Cw*0116 allele.

⁸Primer mix 5: Specific PCR product of 150 bp in the Cw*0120 allele. Specific PCR product of 200 bp in the Cw*0106 allele.

⁹Primer mix 7: Specific PCR product of 150 bp in the Cw*0120 allele. Specific PCR product of 195 bp in the Cw*0108 allele.

¹⁰Primer mix 13: Specific PCR fragment of 155 bp in the Cw*0105 and *0112 alleles. Specific PCR fragments of 155 and 370 bp in the Cw*010201 to *0104, *0106 to *0111 and *0113 to *0120 alleles (both specific PCR fragments may not always be obtained).

| INTERPRETATION TABLE | | | | | | | | |
|---|-------------------------|----------------------|----------------------|----------------------|----------------------|----------------------|----------------------|----------------------|
| HLA-Cw*01 SSP subtyping | | | | | | | | |
| Amplification patterns of the Cw*0101 to Cw*0120 alleles | | | | | | | | |
| | Well⁵ | | | | | | | |
| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 |
| Length of spec. | 90 | 90 | 150 | 210 | 150 | 195 | 150 | 210 |
| PCR product(s) | | 270 | 265 | 240 | 200 | | 195 | |
| Length of int. | 800 | 1070 | 800 | 1070 | 800 | 1070 | 1070 | 1070 |
| pos. control¹ | | | | | | | | |
| 5'-primer(s)² | 368 | 368 | 89 | 89 | 368 | 406 | 368 | 368 |
| | 5'-gTg ^{3'} | 5'-gTg ^{3'} | 5'-gAA ^{3'} | 5'-gAA ^{3'} | 5'-gTg ^{3'} | 5'-gCA ^{3'} | 5'-gTg ^{3'} | 5'-gTg ^{3'} |
| | | | 368 | 368 | | | | |
| | | | 5'-gTg ^{3'} | 5'-gTg ^{3'} | | | | |
| 3'-primer(s)³ | 419 | 419 | 312 | 287 | 479 | 559 | 479 | 538 |
| | 5'-CgT ^{3'} | 5'-CgA ^{3'} | 5'-AgT ^{3'} | 5'-TCg ^{3'} | 5'-CCA ^{3'} | 5'-CgT ^{3'} | 5'-CCA ^{3'} | 5'-CCA ^{3'} |
| | | 595 | 477 | 538 | 527 | | 523 | |
| | | 5'-Cct ^{3'} | 5'-gCA ^{3'} | 5'-CCg ^{3'} | 5'-CCA ^{3'} | | 5'-ACA ^{3'} | |
| Well No. | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 |
| HLA-Cw allele⁴ | | | | | | | | |
| *010201-010205 | 1 | | | | | | | |
| *0103 | | 2 | | | | | | |
| *0104 | | | 3 | | | | | 8 |
| *0105 | | | | 4 | | | | |
| *0106 | 1 | | | | 5 | | | |
| *0107 | 1 | | | | | 6 | | |
| *0108 | 1 | | | | | | 7 | |
| *0109 | 1 | | | | | | | 8 |
| *0110 | 1 | | | | | | | |
| *0111 | 1 | | | | | | | |
| *0112 | | | | | | | | |
| *0113 | 1 | | | | | | | |
| *0114 | 1 | | 3 | | | | | |
| *0115 | 1 | 2 | | | | | | |
| *0116 | 1 | | | 4 | | | | |
| *0117 | 1 | | | | | | | |
| *0118 | 1 | | | | | | | |
| *0119 | 1 | | | | | | | |
| *0120 | 1 | | | | 5 | | 7 | |
| *0509, 0517, 0815 | | | | | | | | |
| *0617, 120306, 1611 | | | | | | | | |
| *0737 | | | | 4 | | | | |
| HLA-Cw allele⁴ | | | | | | | | |
| B*1533 | | | | | | | | |
| Well No. | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 |

| INTERPRETATION TABLE | | | | | | | | |
|--|----------------------|----------------------|----------------------|----------------------|----------------------|----------------------|----------------------|--|
| HLA-Cw*01 SSP subtyping | | | | | | | | |
| Amplification patterns of the Cw*0101 to Cw*0120 alleles | | | | | | | | |
| Well ⁵ | | | | | | | | |
| 9 | 10 | 11 | 12 | 13 | 14 | 15 | 16 | |
| 230 | 290 | 340 | 155 | 155 | 240 | 230 | 255 | Length of spec. PCR product(s) |
| | | | | 370 | | | | |
| 800 | 1070 | 800 | 1070 | 800 | 1070 | 800 | 1070 | Length of int. pos. control ¹ |
| 368 | 368 | 368 | 485 | 89 | 89 | 89 | 89 | 5'-primer(s) ² |
| 5'-gTg ^{3'} | 5'-gTg ^{3'} | 5'-gTA ^{3'} | 5'-CAA ^{3'} | 5'-gAA ^{3'} | 5'-gAA ^{3'} | 5'-gAA ^{3'} | 5'-gAA ^{3'} | |
| | | | | 368 | | | | |
| | | | | 5'-gTg ^{3'} | | | | |
| 559 | 619 | 3 rd I | 601 | 201 | 289 | 278 | 301 | 3'-primer(s) ³ |
| 5'-CTC ^{3'} | 5'-ACT ^{3'} | 5'-ATg ^{3'} | 5'-CTC ^{3'} | 5'-CTC ^{3'} | 5'-AgC ^{3'} | 5'-ggT ^{3'} | 5'-gCA ^{3'} | |
| | | | | 3 rd I | | | | |
| | | | | 5'-CTC ^{3'} | | | | |
| 9 | 10 | 11 | 12 | 13 | 14 | 15 | 16 | Well No. |
| | | | | | | | | HLA-Cw allele ⁴ |
| | | | | 13 | | | | *010201-010205 |
| | | | | 13 | | | | *0103 |
| | | | | 13 | | | | *0104 |
| | | | | 13 | | | | *0105 |
| | | | | 13 | | | | *0106 |
| | | | | 13 | | | | *0107 |
| | | | | 13 | | | | *0108 |
| | | | | 13 | | | | *0109 |
| 9 | | | | 13 | | | | *0110 |
| | 10 | | | 13 | | | | *0111 |
| | | 11 | | 13 | | | | *0112 |
| | | | 12 | 13 | | | | *0113 |
| | | | | 13 | | | | *0114 |
| | | | | 13 | | | | *0115 |
| | | | | 13 | | | | *0116 |
| | | | | 13 | 14 | | | *0117 |
| | | | | 13 | | 15 | | *0118 |
| | | | | 13 | | | 16 | *0119 |
| | | | | 13 | | | | *0120 |
| | | | 12 | | | | | *0509, 0517, 0815 |
| | | 11 | | | | | | *0617, 120306, 1611 |
| | | | | | | | | *0737 |
| | | | | | | | | HLA-Cw allele ⁴ |
| | | | 12 | | | | | B*1533 |
| 9 | 10 | 11 | 12 | 13 | 14 | 15 | 16 | Well No. |

Lot No.: **44E**

Lot-specific information

www.olerup.com

¹The internal positive control primer pairs amplify segments of the human growth hormone gene. The two different control primer pairs give rise to either an internal positive control band of 1070 base pairs, for most wells, or a band of 800 base pairs, for some wells.

Well number 1 contains the primer pair giving rise to the shorter, 800 bp, internal positive control band in order to help in the correct orientation of the HLA-Cw*01 SSP subtyping.

In addition, wells number 3, 5, 9, 11, 13 and 15 contain the primer pair giving rise to the shorter, 800 bp, internal positive control band in order to allow kit identification.

²The nucleotide position, in the 2nd or 3rd exon, matching the specificity-determining 3'-end of the primer is given. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

³The nucleotide position, in the 2nd or 3rd exon or the 3rd intron, matching the specificity-determining 3'-end of the primer is given in the anti-sense direction. Nucleotide numbering as on the www.ebi.ac.uk/imgt/hla web site. The sequence of the 3 terminal nucleotides of the primer is given.

⁴The HLA-Cw*0101 nucleotide sequence has been shown to be identical to Cw*0102.

⁵Primer mix 2: Specific PCR product of 90 bp in the Cw*0103 allele. Specific PCR product of 270 bp in the Cw*0115 allele.

Primer mix 3: Specific PCR product of 150 bp in the Cw*0104 allele. Specific PCR product of 265 bp in the Cw*0114 allele.

Primer mix 4: Specific PCR product of 210 bp in the Cw*0105 and *0737 alleles. Specific PCR product of 240 bp in the Cw*0116 allele.

Primer mix 5: Specific PCR product of 150 bp in the Cw*0120 allele. Specific PCR product of 200 bp in the Cw*0106 allele.

Primer mix 7: Specific PCR product of 150 bp in the Cw*0120 allele. Specific PCR product of 195 bp in the Cw*0108 allele.

Primer mix 13: Specific PCR fragment of 155 bp in the Cw*0105 and *0112 alleles. Specific PCR fragments of 155 and 370 bp in the Cw*010201 to *0104, *0106 to *0111 and *0113 to *0120 alleles (both specific PCR fragments may not always be obtained).

| CELL LINE VALIDATION SHEET | | | | | | | | | | | | | | | | | | | | |
|-----------------------------------|------|----------|------------|-------------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
| HLA-Cw*01 SSP primer set | | | | | | | | | | | | | | | | | | | | |
| | | | | Well | | | | | | | | | | | | | | | | |
| | | | | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 | 13 | 14 | 15 | 16 | |
| | | | | Prod. No.: | 200617401 | 200844802 | 200844803 | 200844804 | 200844805 | 200617406 | 200844807 | 200844808 | 200844809 | 200617410 | 200844811 | 200617412 | 200729813 | 200844814 | 200844815 | 200844816 |
| IHC cell line | | | Cw* | | | | | | | | | | | | | | | | | |
| 1 | 9001 | SA | *0702 | | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 2 | 9280 | LK707 | *0701 | *1505 | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 3 | 9011 | E4181324 | *1202 | | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 4 | 9275 | GU373 | *0304 | *0401 | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 5 | 9009 | KAS011 | *0602 | | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 6 | 9353 | SM | *0304 | *0702 | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 7 | 9020 | QBL | *0501 | | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 8 | 9007 | DEM | *0602 | | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 9 | 9026 | YAR | *1203 | | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 10 | 9107 | LKT3 | *0102 | | + | - | - | - | - | - | - | - | - | - | - | - | + | - | - | - |
| 11 | 9051 | PITOUT | *1601 | | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 12 | 9052 | DBB | *0602 | | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 13 | 9067 | BTB | *0102 | | + | - | - | - | - | - | - | - | - | - | - | - | + | - | - | - |
| 14 | 9071 | OLGA | *0102 | *0304 | + | - | - | - | - | - | - | - | - | - | - | - | + | - | - | - |
| 15 | 9075 | DKB | *0304 | | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 16 | 9037 | SWEIG007 | *0202 | | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 17 | 9008 | WILJON | *1203 | | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 18 | 9257 | 32367 | *0102 | *0705 | + | - | - | - | - | - | - | - | - | - | - | - | + | - | - | - |
| 19 | 9038 | BM16 | *0701 | | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 20 | 9059 | SLE005 | *0304 | | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 21 | 9064 | AMALA | *0303 | | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 22 | 9056 | KOSE | *1203 | | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 23 | 9124 | IHL | *0102 | *1502 | + | - | - | - | - | - | - | - | - | - | - | - | + | - | - | - |
| 24 | 9035 | JBUSH | *1203 | | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 25 | 9049 | IBW9 | *0802 | | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 26 | 9285 | WT49 | *0701 | | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 27 | 9191 | CH1007 | *0704 | *1505 | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 28 | 9320 | BEL5GB | *0501 | *1601 | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 29 | 9050 | MOU | *1601 | | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 30 | 9021 | RSH | *1701 | | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 31 | 9019 | DUCAF | *0501 | | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 32 | 9297 | HAG | *1701 | *1703 | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 33 | 9098 | MT14B | *0304 | | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 34 | 9104 | DHIF | *1203 | | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 35 | 9302 | SSTO | *0501 | | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 36 | 9024 | KT17 | *0303 | *0401 | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 37 | 9065 | HHKB | *0702 | | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 38 | 9099 | LZL | *0303 | | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 39 | 9315 | CML | *0202 | *0701 | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 40 | 9134 | WHONP199 | *0102 | *0602 | + | - | - | - | - | - | - | - | - | - | - | - | + | - | - | - |
| 41 | 9055 | H0301 | *0802 | | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 42 | 9066 | TAB089 | *0102 | | + | - | - | - | - | - | - | - | - | - | - | - | + | - | - | - |
| 43 | 9076 | T7526 | *0102 | *0801 | + | - | - | - | - | - | - | - | - | - | - | - | + | - | - | - |
| 44 | 9057 | TEM | *1203 | | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 45 | 9239 | SHJO | *0602 | *1701 | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 46 | 9013 | SCHU | *0702 | | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 47 | 9045 | TUBO | *0704 | *1502 | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |
| 48 | 9303 | TER-ND | *0401 | *1601 | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - | - |



CERTIFICATE OF ANALYSIS

Olerup SSP[®] HLA-Cw*01 SSP

Product number: 101.621-12 – including *Taq* polymerase
101.621-12u – without *Taq* polymerase
Lot number: 44E
Expiry date: 2010-April-01
Number of tests: 12
Number of wells per test: 16

Well specifications:

| Well No. | Production No. | Well No. | Production No. |
|----------|----------------|----------|----------------|
| 1 | 2006-174-01 | 9 | 2008-448-09 |
| 2 | 2008-448-02 | 10 | 2006-174-10 |
| 3 | 2008-448-03 | 11 | 2008-448-11 |
| 4 | 2008-448-04 | 12 | 2006-174-12 |
| 5 | 2008-448-05 | 13 | 2007-298-13 |
| 6 | 2006-174-06 | 14 | 2008-448-14 |
| 7 | 2008-448-07 | 15 | 2008-448-15 |
| 8 | 2008-448-08 | 16 | 2008-448-16 |

The specificity of each primer solution of the HLA-Cw*01 primer set has been tested against 48 well characterized IHWC cell line DNAs.

No DNAs carrying the alleles to be amplified by primer solutions 2 to 12 and 14 to 16 were available. The specificity of the primers in primer solutions 2, 3, 4, 5, 8, 9, 11, 12 and 14 were tested by separately adding one additional 5'-primer, respectively one additional 3'-primer. In primer solution 6 it was only possible to test the 3'-primer, the 5'-primer was not possible to test. In primer solutions 7, 10, 15 and 16 it was only possible to test the 5'-primers, the 3'-primers were not possible to test.

Results: No false positive or false negative amplifications were obtained.

Date of approval: 2008-May-16

Approved by:

Quality Control, Supervisor

Lot No.: **44E**

Lot-specific information

www.olerup.com

Declaration of Conformity

Product name: *Olerup* SSP[®] HLA-Cw*01
Product number: 101.621-12, 101.621-12u
Lot number: 44E

Intended use: HLA-Cw*01 high resolution histocompatibility testing

Manufacturer: *Olerup* SSP AB
Hasselstigen 1
SE-133 33 Saltsjöbaden, Sweden
Phone: +46-8-717 88 27
Fax: +46-8-717 88 18

We, *Olerup* SSP AB, hereby declare that this product, to which this Declaration of Conformity relates is in conformity with the following Standard(s) and other normative document(s) ISO 9001:2000, ISO 17025:1999 and ISO 13485:2000, following the provisions of the 98/79/EC Directive on *in vitro* diagnostic medical devices, Annex III.

The Technical Construction File is maintained at *Olerup* SSP AB, Hasselstigen 1, SE-133 33 Saltsjöbaden, Sweden.

The Authorized Representative located within the Community is: *Olerup* SSP AB.

Saltsjöbaden, Sweden
2008-May-16

Olle Olerup
Managing Director

HLA-Cw*01
101.621-12 – including *Taq* polymerase
101.621-12u – without *Taq* polymerase

Product Insert

Page 12 of 12

Lot No.: **44E**

Lot-specific information

www.olerup.com

ADDRESSES:

Manufacturer:

Olerup SSP AB, Hasselstigen 1, SE-133 33 Saltsjöbaden, Sweden.

Tel: +46-8-717 88 27

Fax: +46-8-717 88 18

E-mail: info-ssp@olerup.com

Web page: <http://www.olerup.com>

Distributed by:

Olerup GmbH, Löwengasse 47 / 6, AT-1030 Vienna, Austria.

Tel: +43-1-710 15 00

Fax: +43-1-710 15 00 10

E-mail: support-at@olerup.com

Web page: <http://www.olerup.com>

Olerup Inc., 901 S. Bolmar St., Suite R, West Chester, PA 19382

Tel: 1-877-OLERUP1

Fax: 610-344-7989

E-mail: info.us@olerup.com

Web page: <http://www.olerup.com>

For information on *Olerup* SSP distributors worldwide, contact **Olerup GmbH**.