

101.331-24/06 – including *Taq* polymerase

101.331-24u/06u – without *Taq* polymerase

Lot No.: 7N5

Olerup SSP[®] DPA1

Product number: 101.331-24/06 – including *Taq* pol.
101.331-24u/06u – without *Taq* pol.

Lot number: 7N5

Expiry date: 2026-02-01

Number of tests: 24 tests – Product No. 101.331-24/24u
6 tests – Product No. 101.331-06/06u

Number of wells per test: 22+1

CHANGES COMPARED TO THE PREVIOUS DPA1 LOT (1L3):

Well	5'-primer	3'-primer	rationale
22	Added, moved	Added, moved	New primer pair added for the DPA1*02:12:01-02:12:02 alleles, negative control moved to well 23.
23	Added	Added	Negative control added from well 22.

THE NUMBER OF WELLS is increased from 22 to 23.

ALLELE COVERAGE:

DPA1*01:03 to DPA1*01:94, DPA1*02:01 to DPA1*02:64, DPA1*03:01 to DPA1*03:09 and DPA1*04:01 to DPA1*04:03, i.e. all the currently recognized DPA1 alleles, give rise to unique amplification patterns¹; www.ebi.ac.uk/imgt/hla, 2022-January-13, release 3.47.0.

¹Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

RESOLUTION IN DPA1 HOMO- AND HETEROZYGOTES:

Excellent.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

Well 22 has been added for improved resolution of the DPA1*02:12:01-02:12:02 alleles.

