

Olerup SSP[®] DR low resolution

Product number: 101.101-48/12 – including *Taq* pol.
101.101-48u/12u – without *Taq* pol.
Lot number: 93S
Expiry date: 2016-April-01
Number of tests: 48 tests – Product No. 101.101-48
12 tests – Product No. 101.101-12
Number of wells per test: 31 + 1

CHANGES COMPARED TO THE PREVIOUS DR LOW RESOLUTION LOT (13R):

| Well | 5'-primer | 3'-primer | rationale |
|------|-----------|-----------|---|
| 1 | - | Added | 3'-primer added for the DRB1*01:55 allele. |
| 9 | - | Added | 3'-primer added for the DRB1*07:25 allele. |
| 20 | Added | - | 5'-primer added for the DRB1*08:53 allele. |
| 24 | Added | Added | Primer pair added for resolution of the DRB5*01:08N allele. |
| 27 | Added | Added | Primer pair added for the DRB1*14:141 allele. |

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

All the HLA-DRB1, -DRB3, -DRB4¹ and -DRB5 alleles, i.e. **DRB1*01:01:01 to 10:06, DRB3*01:01:02:01 to DRB3*03:03, DRB4*01:01:01:01 to DRB4*01:08, DRB5*01:01:01 to DRB5*02:06**, recognized by the HLA Nomenclature Committee in July 2013² will be amplified by the primers in the DR low resolution SSP kit. The HLA-DRB alleles will be grouped into their corresponding serological specificities.

¹The DRB4*02:13R and DRB4*03:13R null alleles will not be amplified by the DR low resolution primer set.

²DRB alleles listed on the IMGT/HLA web page 2013-July-25, release 3.13.1, www.ebi.ac.uk/imgt/hla.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.