

**101.111-24/06 – including *Taq* polymerase**

**101.111-24u/06u – without *Taq* polymerase**

**Lot No.: 1H1**

## **Olerup SSP® DRB1\*01**

**Product number:** 101.111-24/06 – including *Taq* pol.  
101.111-24u/06u – without *Taq* pol.

**Lot number:** 1H1

**Expiry date:** 2022-11-01

**Number of tests:** 24 test – Product No. 101.111-24/24u  
6 tests – Product No. 101.111-06/06u

**Number of wells per test:** 23+1

The DRB1\*01 primer set is unchanged compared to the previous Olerup SSP® DRB1\*01 (Lot No. 9E8).

**THE NUMBER OF WELLS** is unchanged.

### **ALLELE COVERAGE:**

DRB1\*01:01 to DRB1\*01:89, i.e. all the currently recognized DRB1\*01 alleles, will be amplified by the primers in the DRB1\*01<sup>1</sup> subtyping kit; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2018-July-11, release 3.33.0.

The DRB1\*01 kit enables separation of the confirmed DRB1\*01 alleles as listed in the IMGT/HLA database 3.26.0. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources.

The DRB1\*01 kit also enables identification of null and alternatively expressed alleles.

The following DRB1\*01 alleles can be distinguished by the different sizes of the HLA-specific PCR product:

<b>Alleles</b>	<b>Primer mix</b>
DRB1*01:40N, 01:67	19

<sup>1</sup>Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

### **INFLUENCE ON THE INTERPRETATION OF DRB1\*01 SUBTYPINGS BY NON-DRB1\*01 ALLELES:**

None.

### **MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No modifications made.

Changes in revision R01 compared to R00:

1. The expiration date has been altered due to extension of shelf-life.