

**101.815-12 – including *Taq* polymerase**

**101.815-12u – without *Taq* polymerase**

**Lot No.: 9D0**

## **Olerup SSP® DRB1\*12 Add-on**

|                                  |   |
|----------------------------------|---|
| <b>Product number:</b>           | <b>101.815-12 – including <i>Taq</i> polymerase</b><br><b>101.815-12u – without <i>Taq</i> polymerase</b> |
| <b>Lot number:</b>               | <b>9D0</b>  |
| <b>Expiry date:</b>              | <b>2018-09-01</b>   |
| <b>Number of tests:</b>          | <b>12</b>   |
| <b>Number of wells per test:</b> | <b>3+1</b>  |

### **CHANGES COMPARED TO THE PREVIOUS DRB1\*12 ADD-ON LOT (8D2):**

The DRB1\*12 Add-on primer set is unchanged compared to the previous lot.

**THE NUMBER OF WELLS** is unchanged.

### **ALLELE COVERAGE:**

The DRB1\*12:01, DRB1\*12:10 and DRB1\*12:17 alleles give different patterns in the DRB1\*12 Add-on subtyping kit<sup>1</sup>; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2016-January-19, release 3.23.0.

The DRB1\*12 Add-on subtyping kit cannot distinguish the following silent mutations: the DRB1\*12:01:01-12:01:09 alleles.

<sup>1</sup>Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

### **RESOLUTION IN HOMO- AND HETEROZYGOTES:**

Excellent.

### **INFLUENCE ON THE INTERPRETATION OF DRB1\*12 ADD-ON SUBTYPINGS BY NON-DRB1\*12 ALLELES:**

None.