

## Olerup SSP<sup>®</sup> DRB1\*13

<b>Product number:</b>	<b>101.116-24/03 – including <i>Taq</i> pol. 101.116-24u/03u – without <i>Taq</i> pol.</b>
<b>Lot number:</b>	<b>0L5</b>
<b>Expiry date:</b>	<b>2024-04-01</b>
<b>Number of tests:</b>	<b>24 tests – Product No. 101.116-24/24u 3 tests – Product No. 101.116-03/03u</b>
<b>Number of wells per test:</b>	<b>57+1</b>

### CHANGES COMPARED TO THE PREVIOUS DRB1\*13 LOT (1K0):

The DRB1\*13 primer set is unchanged compared to the previous *Olerup SSP<sup>®</sup>* DRB1\*13 lot (**Lot No. 1K0**).

**THE NUMBER OF WELLS** is unchanged.

### ALLELE COVERAGE:

DRB1\*13:01 to DRB1\*13:295N, i.e., all the currently recognized DRB1\*13 alleles, will be amplified by the primers in the DRB1\*13 subtyping kit<sup>1,2</sup>; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2019-October-17, release 3.38.0.

The DRB1\*13 kit enables separation of the confirmed DRB1\*13 alleles as listed in the IMGT/HLA database 3.26.0. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources.

The DRB1\*13 kit also enables identification of null and alternatively expressed alleles.

The following DRB1\*13 alleles can be distinguished by the different sizes of the HLA-specific PCR product:

Alleles	Primer mix
DRB1*13:28, 13:221	23, 56
DRB1*13:80, 13:242:01	34
DRB1*13:92, 13:98	36

<sup>1</sup>Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

<sup>2</sup>The following alleles cannot be separated by the HLA-DRB1\*13 primer set. These alleles can be distinguished by the DR low resolution kit and/or the HLA- DRB1\*11 and DRB1\*14 high resolution kit.

### Alleles

DRB1*13:17, DRB1*11:229
DRB1*13:165, 13:171:01, DRB1*11:20
DRB1*13:193, DRB1*14:95

101.116-24/03 – including *Taq* polymerase  
101.116-24u/03u – without *Taq* polymerase

Lot No.: **OL5**

**RESOLUTION IN DRB1\*13 HOMOZYGOTES:**

Good.

**INFLUENCE ON THE INTERPRETATION OF DRB1\*13 SUBTYPINGS BY NON-DRB1\*13 ALLELES:**

None of importance.

**MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.