

101.117-12/03 – including *Taq* polymerase  
101.117-12u/03u – without *Taq* polymerase  
Lot No.: **5L2**

## Olerup SSP<sup>®</sup> DRB1\*14

**Product number:** 101.117-12/03 – including *Taq* pol.  
101.117-12u/03u – without *Taq* pol.  
**Lot number:** 5L2  
**Expiry date:** 2024-09-01  
**Number of tests:** 12 tests – Product No. 101.117-12/12u  
3 tests – Product No. 101.117-03/03u  
**Number of wells per test:** 47+1

### CHANGES COMPARED TO THE PREVIOUS DRB1\*14 LOT (5K3):

Well	5'-primer	3'-primer	rationale
21	-	Added	3'-primer added for the DRB1*14:114 allele.
43	Added	-	5'-primer added for the DRB1*14:111 allele.
47	-	Added	3'-primer added for the DRB1*14:114 allele.

**THE NUMBER OF WELLS** is unchanged.

### ALLELE COVERAGE:

DRB1\*14:01 to DRB1\*14:221, i.e. all the currently recognized DRB1\*14 alleles, will be amplified by the primers in the DRB1\*14 subtyping kit<sup>1,2</sup>; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2020-April-20, release 3.40.0.

The DRB1\*14 kit enables separation of the confirmed DRB1\*14 alleles as listed in the IMGT/HLA database 3.27.0. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources.

The DRB1\*14 kit also enables identification of many null and alternatively expressed alleles.

The following DRB1\*14 alleles can be distinguished by the different sizes of the HLA-specific PCR product:

Alleles	Primer mix
DRB1*14:01:02-14:01:04, 14:86	26
DRB1*14:26, 14:110	24

<sup>1</sup>Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

<sup>2</sup>The DRB1\*14 primer set cannot separate the following alleles. These alleles can be distinguished by the DR low resolution and/or DRB1\*03 and DRB1\*13 kits:

#### Alleles

DRB1\*14:21, DRB1\*03:76, DRB1\*13:178

### RESOLUTION IN DRB1\*14 HOMOZYGOTES:

Good.

101.117-12/03 – including *Taq* polymerase

101.117-12u/03u – without *Taq* polymerase

Lot No.: **5L2**

**INFLUENCE ON THE INTERPRETATION OF DRB1\*14 SUBTYPINGS BY NON-DRB1\*14 ALLELES:**

Most frequently encountered DRB1\*03,14; DRB1\*11,14 and DRB1\*13,14 genotypes give rise to unique amplification patterns.

**MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.