

101.117-12/03 – including *Taq* polymerase  
101.117-12u/03u – without *Taq* polymerase

Lot No.: **4R8**

## **Olerup SSP® DRB1\*14**

**Product number:** 101.117-12/03 – including *Taq* pol.  
101.117-12u/03u – without *Taq* pol.  
**Lot number:** 4R8  
**Expiry date:** 2026-12-01  
**Number of tests:** 12 tests – Product No. 101.117-12/12u  
3 tests – Product No. 101.117-03/03u  
**Number of wells per test:** 47+1

### **Changes compared to the previous DRB1\*14 Lot (6N0):**

The primers of the wells detailed below have been exchanged, added or modified compared to the previous lot (**Lot No. 6N0**).

Well	5'-primer	3'-primer	rationale
11	-	Exchanged	3'-primer added for improved allelic resolution.

**THE NUMBER OF WELLS** is unchanged.

### **ALLELE COVERAGE:**

DRB1\*14:01 to DRB1\*14:247, i.e. all the currently recognized DRB1\*14 alleles, will be amplified by the primers in the DRB1\*14 subtyping kit<sup>1,2</sup>; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2022-October-12, release 3.50.0.

The DRB1\*14 kit enables separation of the confirmed DRB1\*14 alleles as listed in the IMGT/HLA database 3.27.0. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources.

The DRB1\*14 kit also enables identification of many null and alternatively expressed alleles.

The following DRB1\*14 alleles can be distinguished by the different sizes of the HLA-specific PCR product:

Alleles	Primer mix
DRB1*14:01:02-14:01:04, 14:86	26
DRB1*14:26, 14:110	24

<sup>1</sup>Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

<sup>2</sup>The DRB1\*14 primer set cannot separate the following alleles. These alleles can be distinguished by the DR low resolution and/or DRB1\*03 and DRB1\*13 kits:

#### **Alleles**

DRB1\*14:21, DRB1\*03:76, DRB1\*03:174N, DRB1\*13:178



**0197**

101.117-12/03 – including *Taq* polymerase

101.117-12u/03u – without *Taq* polymerase

Lot No.: **4R8**

**RESOLUTION IN DRB1\*14 HOMOZYGOTES:**

Good.

**INFLUENCE ON THE INTERPRETATION OF DRB1\*14 SUBTYPINGS BY NON-DRB1\*14 ALLELES:**

Most frequently encountered DRB1\*03,14; DRB1\*11,14 and DRB1\*13,14 genotypes give rise to unique amplification patterns.

**MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.



**0197**

For *In Vitro* Diagnostic Use

MA125 v01 SSP Release Note Template

Date: January 2023, Rev. No: 00