

101.125-24/06 – including *Taq* polymerase
101.125-24u/06u – without *Taq* polymerase

Lot No.: **27M**

Olerup SSP® DRB1*15

Product number: 101.125-24/06 – including *Taq* pol.
101.125-24u/06u – without *Taq* pol.
Lot number: 27M
Expiry date: 2013-November-01
Number of tests: 24 test – Product No. 101.125-24
6 tests – Product No. 101.125-06
Number of wells per test: 32

CHANGES COMPARED TO THE PREVIOUS DRB1*15 LOT (09K):

Well	5'-primer	3'-primer	rationale
18	Exchanged	-	Exchanged 5'-primer and positive control primer pair, to decrease primer oligomer formation.
30	-	Added	Primer added for the DRB1*15:45 allele.
31	New	New	New primer pair for the DRB1*15:46 allele.
32	New	New	New primer pairs for the DRB1*15:47 and DRB1*15:49 alleles.

THE NUMBER OF WELLS is increased from 30 to 32.

ALLELE COVERAGE:

DRB1*15:01 to DRB1*15:49, i.e. all the currently recognized DRB1*15 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2010-July-16, release 3.1.0.

RESOLUTION IN DRB1*15 HOMOZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF DRB1*15 SUBTYPINGS BY NON-DRB1*15 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.