

## Olerup SSP® DRB1\*16

Product number: 101.126-12 – including *Taq* polymerase  
101.126-12u – without *Taq* polymerase  
Lot number: 65N  
Expiry date: 2014-October-01  
Number of tests: 12  
Number of wells per test: 16

### CHANGES COMPARED TO THE PREVIOUS DRB1\*16 LOT (02M):

Well	5'-primer	3'-primer	rationale
13	-	Added	Primer added for the DRB1*16:19 allele.

THE NUMBER OF WELLS is unchanged.

### ALLELE COVERAGE:

DRB1\*16:01 to DRB1\*16:19, i.e. all the currently recognized DRB1\*16 alleles, will be amplified by the primers in the DRB1\*16 subtyping kit; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2012-January-12, release 3.7.0.

The DRB1\*16 kit enables separation of the confirmed DRB1\*16 alleles as listed in the IMGT/HLA database. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources. Current allele confirmation status for DRB1\*16 alleles is listed below.

The DRB1\*16 kit also enables identification of polymorphisms in exons outside of the region encoding the peptide binding domain and of null and alternatively expressed alleles.

The DRB1\*16 subtyping kit cannot distinguish the following silent mutations: the DRB1\*16:01:01-16:01:02, DRB1\*16:02:01-16:02:02 or the DRB1\*16:05:01-16:05:02 alleles.

### RESOLUTION IN DRB1\*16 HOMOZYGOTES:

Excellent.

### INFLUENCE ON THE INTERPRETATION OF DRB1\*16 SUBTYPINGS BY NON-DRB1\*16 ALLELES:

None frequently occurring.

### MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.