

101.127-12/04 – including *Taq* polymerase  
101.127-12u/04u – without *Taq* polymerase

Lot No.: **13X**

## Olerup SSP® DRB1\*08

**Product number:** 101.127-12/04 – including *Taq* pol.  
101.127-12u/04u – without *Taq* pol.  
**Lot number:** 13X  
**Expiry date:** 2017-March-01  
**Number of tests:** 12 tests – Product No. 101.127-12/12u  
4 tests – Product No. 101.127-04/04u  
**Number of wells per test:** 23+1

### CHANGES COMPARED TO THE PREVIOUS DRB1\*08 LOT (79R):

Well	5'-primer	3'-primer	rationale
1	Added	-	5'-primer added for the DRB1*08:01:06 allele.
5	Added	-	5'-primer added for the DRB1*08:01:06 allele.
7	Added	-	5'-primer added for the DRB1*08:01:06 allele.
9	Moved, added	Added	5'-primer moved to well 18, primer pair added from well 24.
12	-	Added	3'-primer added for the DRB1*08:43 allele.
18	Added	-	5'-primer added from well 18.
23	-	Added	3'-primer added from well DRB1*08:43 allele.
24	Moved	Moved	Primer pair moved to well 9, negative Control.

**THE NUMBER OF WELLS** is unchanged.

### ALLELE COVERAGE:

DRB1\*08:01 to DRB1\*08:56, i.e. all the currently recognized DRB1\*08 alleles, will be amplified by the primers in the DRB1\*08 subtyping kit<sup>1,2</sup>; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2014-April-14, release 3.16.0.

The DRB1\*08 kit enables separation of the confirmed DRB1\*08 alleles as listed in the IMGT/HLA database. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources. The DRB1\*08 kit also enables identification of polymorphisms in exons outside of the region encoding the peptide binding domain and of null and alternatively expressed alleles.

The DRB1\*08 subtyping kit cannot distinguish the following silent mutations: DRB1\*08:01:01-08:01:06, the DRB1\*08:02:01-08:02:04, DRB1\*08:03:02-08:03:04, the DRB1\*08:04:01 and 08:04:04-08:04:07 alleles, the DRB1\*08:04:02-08:04:03 alleles, the DRB1\*08:30:01-08:30:03 or the DRB1\*08:36:01-08:36:02 alleles.

<sup>1</sup>Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

<sup>2</sup>The DRB1\*08:20 and the DRB1\*14:12:01-14:12:02 and 14:84 give rise to identical amplification patterns with the DRB1\*08 subtyping kit. These alleles can be distinguished by e.g. the DR low resolution kit and/or the DRB1\*14 subtyping kit.

101.127-12/04 – including *Taq* polymerase

101.127-12u/04u – without *Taq* polymerase

Lot No.: **13X**

**RESOLUTION IN DRB1\*08 HOMOZYGOTES:**

The DRB1\*08:01,08:01, DRB1\*08:02,08:02 and DRB1\*08:04,08:04 genotypes give rise to unique amplification patterns.

**INFLUENCE ON THE INTERPRETATION OF DRB1\*08 SUBTYPINGS BY NON-DRB1\*08 ALLELES:**

None frequently occurring.

**MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.