

101.122-24/06 – including *Taq* polymerase

101.122-24u/06u – without *Taq* polymerase

Lot No.: **1V5**

Olerup SSP® DRB4

Product number: 101.122-24/06 – including *Taq* pol.
101.122-24u/06u – without *Taq* pol.

Lot number: 1V5

Expiry date: 2029-10-01

Number of tests: 24 tests – Product No. 101.122-24/24u
6 tests – Product No. 101.122-06/06u

Number of wells per test: 31+1

The primers of the wells detailed below have been exchanged, added or modified compared to the previous lot (8R9).

Well	5'-primer	3'-primer	rationale
7	-	Modified	3'-primer modified due to unspecific amplification.
18	Modified	-	5'-primer modified due to unspecific amplification.
19	Modified	-	5'-primer modified due to unspecific amplification.
20	Modified	-	5'-primer modified due to unspecific amplification.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

DRB4*01:01:01:01 to DRB4*03:01N, i.e. all the currently recognized DRB4 alleles, will be amplified by the primers in the DRB4 subtyping kit¹; www.ebi.ac.uk/imgt/hla, 2025-January-15, release 3.59.0.

The DRB4 kit enables separation of the confirmed DRB4 alleles as listed in the IMGT/HLA database 3.29.0. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources.

The DRB4 kit also enables identification of many null and alternatively expressed alleles.

The following DRB4 alleles can be distinguished by the different sizes of the HLA-specific PCR product:

Alleles	Primer mix	Alleles	Primer mix
DRB4*01:08, 01:39	13	DRB4*01:32, 01:43, 01:54:01N-01:54:02N	24
DRB4*01:20, 01:36	21	DRB4*01:33, 01:38N, 01:40	22
DRB4*01:24, 01:37	23		

¹Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

RESOLUTION IN DRB4 HOMOZYGOTES:

Good.



0197

101.122-24/06 – including *Taq* polymerase

101.122-24u/06u – without *Taq* polymerase

Lot No.: 1V5

INFLUENCE ON THE INTERPRETATION OF DRB4 SUBTYPINGS BY NON-DRB4 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.



0197

For *In Vitro* Diagnostic Use

MA125 v01 SSP Release Note Template

Date: October 2025, Rev. No: 00