

101.412-24/04 – including *Taq* polymerase

101.412-24u/04u – without *Taq* polymerase

Lot No.: **15V**

## **Olerup SSP® HLA-A\*02**

**Product number:** 101.412-24/04 – including *Taq* pol.  
101.412-24u/04u – without *Taq* pol.

**Lot number:** 15V

**Expiry date:** 2016-July-01

**Number of tests:** 24 tests – Product No. 101.412-24/24u  
4 tests – Product No. 101.412-04/04u

**Number of wells per test:** 95+1

### **CHANGES COMPARED TO THE PREVIOUS HLA-A\*02 LOT (29S):**

| <b>Well</b> | <b>5'-primer</b> | <b>3'-primer</b> | <b>rationale</b>  |
|-------------|------------------|------------------|---|
| 5           | -                | Added            | 3'-primer added for the A*02:01:96 allele.  |
| 20          | -                | Added            | 3'-primer added for improved allelic resolution of the A*02:33 and 02:198 alleles.  |
| 24          | -                | Added            | 3'-primer added for the A*02:393 allele.  |
| 25          | -                | Added            | 3'-primer added for the A*02:384 allele.  |
| 27          | Added            | -                | 5'-primer added for the A*02:407 allele.  |
| 39          | -                | Added            | 3'-primer added from well 96, 3'-primer added for the A*02:439N allele.   |
| 40          | -                | Added            | 3'-primer added from well 96.   |
| 44          | Added            | -                | 5'-primer added for the A*02:342 allele.  |
| 65          | Added            | -                | 5'-primer added for the A*02:342 and A*02:440Q alleles.   |
| 66          | Added            | -                | 5'-primer added for improved allelic resolution of the A*02:30 and A*02:01:83 alleles.  |
| 67          | Added            | -                | 5'-primer added for improved allelic resolution of the A*02:07 and A*02:322 alleles, 5'-primer added for the A*02:407 allele. |
| 68          | -                | Added            | 3'-primer added for improved allelic resolution of the A*02:119 and A*02:158 alleles.   |
| 70          | -                | Added            | 3'-primer added for improved allelic resolution of the A*02:133 and A*02:250 alleles.   |
| 71          | -                | Added            | 3'-primer added from well 96.   |
| 75          | -                | Added            | 3'-primer added for the A*02:384 allele.  |
| 76          | Added            | -                | 5'-primer added for the A*02:342 allele.  |
| 79          | -                | Added            | 3'-primer added for the A*02:419 allele.  |
| 83          | -                | Added            | 3'-primer added for the A*02:440Q allele.   |
| 88          | -                | Added            | 3'-primer added from well 96.   |
| 89          | -                | Added            | 3'-primer added for the A*02:426 allele.  |
| 90          | Added            | -                | 5'-primer added for the A*02:377 allele.  |
| 92          | -                | Added            | 3'-primer added from well 96.   |
| 94          | Added            | -                | 5'-primer added for the A*02:377 allele.  |
| 96          | -                | Moved            | Negative Control, 3'-primers moved to wells 39, 40, 71, 88 and 92.  |

101.412-24/04 – including *Taq* polymerase  
101.412-24u/04u – without *Taq* polymerase

Lot No.: **15V**

**THE NUMBER OF WELLS** is unchanged.

#### ALLELE COVERAGE:

All HLA-A\*02 alleles recognized in July 2012, i.e. A\*02:01 to A\*02:456, will be amplified by the primers in the HLA-A\*02 subtyping kit; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2013-October-11, release 3.14.0.

The HLA-A\*02 kit enables separation of the confirmed HLA-A\*02 alleles as listed in the IMGT/HLA database. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources.

The HLA-A\*02 kit also enables identification of polymorphisms in exons outside of the region encoding the peptide binding domain and of null and alternatively expressed alleles

The following HLA-A\*02 alleles can be distinguished by the different sizes of the HLA-specific PCR product:

| Alleles                      | Primer mix | Alleles                  | Primer mix |
|------------------------------|------------|--------------------------|------------|
| A*02:28, 02:409              | 27         | A*02:134, 02:314N        | 85         |
| A*02:97:01-02:97:02, 02:305N | 59         | A*02:138, 02:284N        | 87         |
| A*02:107, 02:202, 02:251     | 65         | A*02:139, 02:235         | 88         |
| A*02:111, 02:350N            | 67         | A*02:141, 02:275         | 90         |
| A*02:116, 02:196             | 71         | A*02:150, 02:197, 02:325 | 94         |
| A*02:119, 02:263             | 73         | A*02:180, 02:358         | 29         |
| A*02:120, 02:223N            | 74         | A*02:193, 02:213         | 78         |
| A*02:132, 02:215, 02:237     | 83         | A*02:269, 02:367         | 70         |

The HLA-A\*02 subtyping kit cannot distinguish the silent mutations in the A\*02:01:01:01, 02:01:01:03-02:01:02, 02:01:04-02:01:15, 02:01:18-02:01:19, 02:01:21-02:01:51, 02:01:53-02:01:62, 02:01:64-02:01:73, 02:01:75-02:01:81, 02:01:84-02:01:86 and 02:01:89-02:01:90 and 02:01:92-02:01:99 alleles, the A\*02:243:01-02:243:02 alleles, the A\*02:01:17, 02:01:52, 02:01:63 and 02:01:74 alleles, the A\*02:01:87 and 02:01:91 alleles, the A\*02:03:01 and 02:03:03-02:03:06 alleles, the A\*02:05:01-02:05:06 alleles, the A\*02:06:01-02:06:11 and 02:06:13-02:06:15 alleles, the A\*02:07:01-02:07:03 alleles, the A\*02:11:01, 02:11:03-02:11:05 alleles, the A\*02:17:01-02:17:03 alleles, the 02:20:01-02:20:02 alleles, the A\*02:22:01-02:22:02 alleles, the A\*02:24:01-02:24:02 alleles, the A\*02:40:01-02:40:02 alleles, the A\*02:60:01-02:60:02 alleles, the A\*02:74:01-02:74:02 alleles, the A\*02:76:01-02:76:02 alleles, the A\*02:79:01-02:79:02 alleles, the A\*02:97:01-02:97:02 alleles, the A\*02:101:01-02:101:02 alleles and the A\*02:171:01-02:171:02 alleles.

<sup>1</sup>HLA-A alleles listed on the IMGT/HLA web page 2013-October-11, release 3.14.0, [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla).

<sup>2</sup>Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

**101.412-24/04 – including *Taq* polymerase**

**101.412-24u/04u – without *Taq* polymerase**

**Lot No.: 15V**

**RESOLUTION IN HLA-A\*02 HOMO- AND HETEROZYGOTES:**

Good.

**INFLUENCE ON THE INTERPRETATION OF HLA-A\*02 SUBTYPINGS BY NON-HLA-A\*02 ALLELES:**

A\*02:01,68:xx = A\*02:95,68:xx.

**MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.