

Olerup SSP[®] HLA-A*25

Product number: 101.423-06 – including *Taq* polymerase
 101.423-06u – without *Taq* polymerase

Lot number: 50X

Expiry date: 2017-July-01

Number of tests: 6

Number of wells per test: 15+1

CHANGES COMPARED TO THE PREVIOUS HLA-A*25 LOT (37S):

Well	5'-primer	3'-primer	rationale
4	Added	-	5'-primer added for the A*25:24 allele.
6	Added	Added	Primer pair added for the A*25:19 allele.
7	-	Added	3'-primer added from well 14, 3'-primer added for the A*25:19 allele.
8	-	Added	3'-primer added from well 14, 3'-primer added for the A*25:16 allele.
9	-	Added	3'-primers added for the A*25:16 and 25:19 alleles.
14	Added, moved	Added, moved	Primer pair moved to wells 7 and 8, primer pair added from well 16.
16	Moved	Moved	Primer pair moved to well 14, negative control.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

A*25:01 to A*25:29, i.e. all the currently recognized HLA-A*25 alleles, will be amplified by the primers in the HLA-A*25 subtyping kit^{1,2}; www.ebi.ac.uk/imgt/hla, 2014-October-10, release 3.18.0.

The HLA-A*25 kit enables separation of the confirmed HLA-A*25 alleles as listed in the IMGT/HLA database. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources.

The HLA-A*25 kit also enables identification of polymorphisms in exons outside of the region encoding the peptide binding domain and of null and alternatively expressed alleles.

The HLA-A*25 subtyping kit cannot distinguish the following silent mutations: the A*25:01:01 and 25:01:03-25:01:07 alleles.

¹Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

²The HLA-A*25 primer set cannot separate the A*25:18 and A*32:62 alleles. These alleles can be distinguished by the HLA-A low resolution kit and/or the HLA-A*32 high resolution kit.

101.423-06 – including *Taq* polymerase

101.423-06u – without *Taq* polymerase

Lot No.: 50X

RESOLUTION IN HLA-A*25 HOMO- AND HETEROZYGOTES:

Good.

INFLUENCE ON THE INTERPRETATION OF HLA-A*25 SUBTYPINGS BY NON-HLA-A*25 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No modifications made.