

Olerup SSP[®] HLA-A*25

Product number:	101.423-06 – including <i>Taq</i> polymerase 101.423-06u – without <i>Taq</i> polymerase
Lot number:	3N0
Expiry date:	2025-09-01
Number of tests:	6
Number of wells per test:	15+1

Changes compared to the previous HLA-A*25 Lot (6K4):

Well	5'-primer	3'-primer	rationale
1	Added	-	5'-primer added for the A*25:01:17 allele.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

A*25:01 to A*25:73, i.e. all the currently recognized HLA-A*25 alleles, will be amplified by the primers in the HLA-A*25 subtyping kit^{1,2}; www.ebi.ac.uk/imgt/hla, 2021-April-20, release 3.44.0.

The HLA-A*25 kit enables separation of the confirmed HLA-A*25 alleles as listed in the IMGT/HLA database 3.23.0. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources.

The HLA-A*25 kit also enables identification of many null and alternatively expressed alleles.

¹Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

² The HLA-A*25 primer set cannot separate the following alleles. These alleles can be distinguished by the HLA-A low resolution kit and/or the respective high-resolution kits:

Alleles

A*25:36, 25:43, 25:50, 25:52,
25:65, A*32:62
A*25:59, A*01:315
A*25:64, A*32:26:01-32:26:02

RESOLUTION IN HLA-A*25 HOMO- AND HETEROZYGOTES:

Good.

INFLUENCE ON THE INTERPRETATION OF HLA-A*25 SUBTYPINGS BY NON-HLA-A*25 ALLELES:

None of importance.



0197

HLA-A*25

Release Note

Page 2 of 2

101.423-06 – including *Taq* polymerase

101.423-06u – without *Taq* polymerase

Lot No.: **3N0**

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.



0197

For *In Vitro* Diagnostic Use

MA125 v01 SSP Release Note Template

Date: September 2021, Rev. No: 00