

Olerup SSP® HLA-A*26

Product number:	101.424-12 – including <i>Taq</i> polymerase 101.424-12u – without <i>Taq</i> polymerase
Lot number:	74Y
Expiry date:	2018-February-01
Number of tests:	12
Number of wells per test:	47+1

CHANGES COMPARED TO THE PREVIOUS HLA-A*26 LOT (01X):

Well	5'-primer	3'-primer	rationale
11	Added	-	5'-primer added for the allelic resolution of the A*26:106 allele.
22	Modified	-	5'-primer modified for increased yield and improved HLA-specific amplification.
33	Modified	-	5'-primer modified for increased yield and improved HLA-specific amplification.
43	Removed, added	Removed, added	Primer pair removed, primer pair added for the A*26:103 and A*26:107N alleles.
48	-	-	Updated negative control.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

A*26:01 to A*26:111, i.e. all the currently recognized HLA-A*26 alleles, will be amplified by the primers in the HLA-A*26 subtyping kit^{1,2}, www.ebi.ac.uk/imgt/hla, 2015-April-17, release 3.20.0.

The HLA-A*26 kit enables separation of the confirmed HLA-A*68 alleles as listed in the IMGT/HLA database. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources.

The HLA-A*26 kit also enables identification of polymorphisms in exons outside of the region encoding the peptide binding domain and of null and alternatively expressed alleles.

The HLA-A*26 subtyping kit cannot distinguish the following silent mutations: A*26:01:01-26:01:20 and 26:01:22-26:01:38 alleles, the A*26:02:01-26:02:02 alleles or the A*26:07:01-26:07:02 alleles.

101.424-12 – including *Taq* polymerase
 101.424-12u – without *Taq* polymerase

Lot No.: **74Y**

The following HLA-A*26 alleles can be distinguished by the different sizes of the HLA-specific PCR product:

Alleles	Primer mix	Alleles	Primer mix
A*26:17, 26:106	13	A*26:46, 26:53	33
A*26:23, 26:27	20	A*26:51, 26:75	38
A*26:24, 26:41	21	A*26:54, 26:55	36
A*26:25N, 26:38	22	A*26:59, 26:69	35
A*26:43:01, 26:61	32	A*26:62, 26:63	39

¹Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page

<http://hla.alleles.org/alleles/deleted.html>.

²The A*26:28 and A*26:52 alleles will give rise to identical amplification patterns with the HLA-A*26 subtyping kit.

RESOLUTION IN HLA-A*26 HOMO- AND HETEROZYGOTES:

The A*26:01,26:01, A*26:01,26:02 and A*26:02,26:02 genotypes give rise to unique amplification patterns with respect to confirmed alleles.

INFLUENCE ON THE INTERPRETATION OF HLA-A*26 SUBTYPINGS BY NON-HLA-A*26 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.