

## Olerup SSP<sup>®</sup> HLA-A\*29

**Product number:** 101.428-12 – including *Taq* polymerase  
 101.428-12u – without *Taq* polymerase  
**Lot number:** 5L0  
**Expiry date:** 2024-08-01  
**Number of tests:** 12  
**Number of wells per test:** 30+1

### CHANGES COMPARED TO THE PREVIOUS HLA-A\*29 LOT (3H6):

Well	5'-primer	3'-primer	rationale
30	Added, moved	Added, moved	Primer pair added for the A*29:88 allele. Negative control moved to primer mix 31.
31	Added	Added	Negative control added from primer mix 30.

**THE NUMBER OF WELLS** is increased from 30 to 31 wells.

### ALLELE COVERAGE:

A\*29:01 to A\*29:144, i.e. all the currently recognized HLA-A\*29 alleles, will be amplified by the primers in the HLA-A\*29 subtyping kit<sup>1</sup>; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2020-April-20, release 3.40.0.

The HLA-A\*29 kit enables separation of the confirmed HLA-A\*29 alleles as listed in the IMGT/HLA database 3.29.0. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources.

The HLA-A\*29 kit also enables identification of many null and alternatively expressed alleles.

The following HLA-A\*29 alleles can be distinguished by the different sizes of the HLA-specific PCR product:

Alleles	Primer mix
A*29:07, 29:46	8
A*29:08N, 29:139	9
A*29:11, 29:92	13
A*29:17, 29:144	10
A*29:20, 29:82	29

<sup>1</sup>Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

**RESOLUTION IN HLA-A\*29 HOMO- AND HETEROZYGOTES:**  
 Excellent.

**INFLUENCE ON THE INTERPRETATION OF HLA-A\*29 SUBTYPINGS BY NON-HLA-A\*29 ALLELES:**  
 None frequently occurring.

**HLA-A\*29**

**Release Note**

**Page 2 of 2**

**101.428-12 – including *Taq* polymerase**

**101.428-12u – without *Taq* polymerase**

**Lot No.: 5L0**

**MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.