

Olerup SSP[®] HLA-A*31

Product number: 101.430-12 – including *Taq* polymerase
 101.430-12u – without *Taq* polymerase
Lot number: 6G7
Expiry date: 2022-08-01
Number of tests: 12
Number of wells per test: 41+1

CHANGES COMPARED TO THE PREVIOUS HLA-A*31 LOT (4F2):

Well	5'-primer	3'-primer	rationale
40	Added	Added	Negative Control moved to well 42, new primer pair added for the A*31:126N allele.
41	New	New	New primer pair added for the A*31:131N allele.
42	-	-	Negative Control Added from well 40.

THE NUMBER OF WELLS is increased from 40 to 42 wells.

ALLELE COVERAGE:

A*31:01 to A*31:135, i.e. all the currently recognized HLA-A*31 alleles, will be amplified by the primers in the HLA-A*31 subtyping kit^{1,2}; www.ebi.ac.uk/imgt/hla, 2018-April-16, release 3.32.0.

The HLA-A*31 kit enables separation of the confirmed HLA- A*31 alleles as listed in the IMGT/HLA database 3.28.0. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources.

The HLA-A*31 kit also enables identification of polymorphisms in exons outside of the region encoding the peptide binding domain and of null and alternatively expressed alleles.

The following HLA-A*31 alleles can be distinguished by the different sizes of the HLA-specific PCR product:

Alleles	Primer mix	Alleles	Primer mix
A*31:01:02:03N, 31:35	21	A*31:29, 31:59	25
A*31:16, 31:46	16	A*31:36, 31:48	29

¹Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

²The A*31:44 and the A*33:15 alleles can be distinguished by the different sizes of the specific PCR products generated by primer mix 29.

The A*31 primer set cannot separate the A*31:89 from the A*33:125 and A*33:131 alleles. These alleles can be distinguished by the HLA-A low resolution and/or HLA-A*33 kits.

HLA-A*31

101.430-12 – including *Taq* polymerase

101.430-12u – without *Taq* polymerase

Lot No.: 6G7

RESOLUTION IN HLA-A*31 HOMO- AND HETEROZYGOTES:

Very good.

INFLUENCE ON THE INTERPRETATION OF HLA-A*31 SUBTYPINGS BY NON-HLA-A*31 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.

Changes in revision R01 compared to R00:

1. The expiration date has been altered due to extension of shelf-life.