

Olerup SSP[®] HLA-A*31

Product number: 101.430-12 – including *Taq* polymerase
 101.430-12u – without *Taq* polymerase
Lot number: 8K7
Expiry date: 2024-04-01
Number of tests: 12
Number of wells per test: 42+1

CHANGES COMPARED TO THE PREVIOUS HLA-A*31 LOT (6G7):

Well	5'-primer	3'-primer	rationale
1	Added	-	5'-primer added for the A*31:01:34 allele.
30	-	Modified	3'-primer modified in order to improve HLA-specific amplification.
33	Modified	-	5'-primer modified in order to improve HLA-specific amplification.
42	Added, moved	Added, moved	Primer pair added for the A*31:142 allele. Negative control moved to well 43.
43	Added	Added	Negative control moved from well 42.

THE NUMBER OF WELLS is increased from 42 to 43 wells.

ALLELE COVERAGE:

A*31:01 to A*31:167, i.e. all the currently recognized HLA-A*31 alleles, will be amplified by the primers in the HLA-A*31 subtyping kit^{1,2}; www.ebi.ac.uk/imgt/hla, 2019-October-17, release 3.38.0.

The HLA-A*31 kit enables separation of the confirmed HLA- A*31 alleles as listed in the IMGT/HLA database 3.28.0. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources.

The HLA-A*31 kit also enables identification of many null and alternatively expressed alleles.

The following HLA-A*31 alleles can be distinguished by the different sizes of the HLA-specific PCR product:

Alleles	Primer mix	Alleles	Primer mix
A*31:01:02:03N, 31:35	21	A*31:29, 31:59, 31:146	25
A*31:16, 31:46	16	A*31:36, 31:48	29

¹Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

HLA-A*31
101.430-12 – including *Taq* polymerase
101.430-12u – without *Taq* polymerase
Lot No.: 8K7

²The A*31 primer set cannot separate the following alleles. These alleles can be distinguished by the HLA-A low resolution and/or HLA-A*33 kits.

Alleles

A*31:89, A*33:125, 33:131

RESOLUTION IN HLA-A*31 HOMO- AND HETEROZYGOTES:

Very good.

INFLUENCE ON THE INTERPRETATION OF HLA-A*31 SUBTYPINGS BY NON-HLA-A*31 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.