

Olerup SSP[®] HLA-A*32

Product number: 101.431-12 – including *Taq* polymerase
101.431-12u – without *Taq* polymerase
Lot number: 99E
Expiry date: 2010-August-01
Number of tests: 12
Number of wells per test: 16

CHANGES COMPARED TO THE PREVIOUS HLA-A*32 LOT (Y13):

Well	5'-primer	3'-primer	rationale
3	-	Modified	Increased specificity of specific primer pair.
4	-	Modified	Increased specificity of specific primer pair.
12	-	New and modified	New primer for the A*3216 allele, and changed primer to increase specificity.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

A*3201 to A*3216, i.e. all the currently recognized HLA-A*32 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2008-July-11, release 2.22.0.

RESOLUTION IN HLA-A*32 HOMO- AND HETEROZYGOTES:

The A*3201,3201 genotype gives rise to a unique amplification pattern.

INFLUENCE ON THE INTERPRETATION OF HLA-A*32 SUBTYPINGS BY NON-HLA-A*32 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.