

Lot No.: **34G**

## **Olerup SSP<sup>®</sup> HLA-A\*66**

**Product number:** 101.427-06 – including *Taq* polymerase  
101.427-06u – without *Taq* polymerase  
**Lot number:** 34G  
**Expiry date:** 2011-October-01  
**Number of tests:** 6  
**Number of wells per test:** 12

### **CHANGES COMPARED TO THE PREVIOUS HLA-A\*66 LOT (Y03):**

<b>Well</b>	<b>5'-primer</b>	<b>3'-primer</b>	<b>rationale</b>
1	-	Added	Primer added for the A*6610 allele.
4	Moved	Moved	Primer pair moved to well 9.
9	Added	Added	Primer pair from well 4 for improved resolution of the A*6606 allele.
10	New	New	New primer pair for the A*6609 allele.
11	New	New	New primer pair for the A*6610 allele.
12	New	New	New primer pair for the A*6608 allele.

**THE NUMBER OF WELLS** has been increased from 8 to 12 wells.

#### **ALLELE COVERAGE:**

A\*6601 to A\*6610, i.e. all the currently recognized HLA-A\*66 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2009-July-17, release 2.26.0.

#### **RESOLUTION IN HLA-A\*66 HOMO- AND HETEROZYGOTES:**

Good.

#### **INFLUENCE ON THE INTERPRETATION OF HLA-A\*66 SUBTYPINGS BY NON-HLA-A\*66 ALLELES:**

None frequently occurring.

#### **MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No modifications made.